



Drug Analysis

Test No. 22-5001 Summary Report

Each sample pack consisted of two questioned samples containing powder substances, which participants were requested to analyze. Data were returned from 556 participants and are compiled into the following tables:

Report Contents:	<u>Page</u>
<u>Manufacturer's Information</u>	<u>2</u>
<u>Summary Comments</u>	<u>3</u>
<u>Table 1a: Drug Analysis (Item 1)</u>	<u>4</u>
<u>Table 1b: Examination Methods (Item 1)</u>	<u>14</u>
<u>Table 2a: Drug Analysis (Item 2)</u>	<u>29</u>
<u>Table 2b: Examination Methods (Item 2)</u>	<u>43</u>
<u>Table 3: Additional Comments</u>	<u>61</u>
<u>Appendix: Data Sheet</u>	

This report contains the data received from the participants in this test. Since these participants are located in many countries around the world, and it is their option how the samples are to be used (e.g., training exercise, known or blind proficiency testing, research and development of new techniques, etc.), the results compiled in the Summary Report are not intended to be an overview of the quality of work performed in the profession and cannot be interpreted as such. The Summary Comments are included for the benefit of participants to assist with maintaining or enhancing the quality of their results. These comments are not intended to reflect the general state of the art within the profession.

Participant results are reported using a randomly assigned "WebCode". This code maintains participant's anonymity, provides linking of the various report sections, and will change with every report.

Manufacturer's Information

Each sample pack consisted of two items. Item 1 contained 180 mg of powder consisting of 100% testosterone. Item 2 contained 300 mg of powder consisting of caffeine. Testosterone is a Schedule III controlled substance in the United States.

SAMPLE PREPARATION-

The testosterone in Item 1 was a pure powder. The caffeine in Item 2 was also a pure powder.

ITEM 1 (PREPARATION): Approximately 180 mg of the testosterone powder was weighed out and deposited into a glassine bag which was folded and secured with a label. The folded glassine bag was placed into a small zip top bag, which was heat sealed and then placed into a pre-labeled envelope.

ITEM 2 (PREPARATION): Approximately 300 mg of the caffeine powder was weighed out and deposited into a glassine bag which was folded and secured with a label. The folded glassine bag was placed into a small zip top bag, which was heat sealed and then placed into a pre-labeled envelope.

SAMPLE PACK ASSEMBLY: One of each of the Item 1 and Item 2 envelopes was placed into a larger pre-labeled sample pack envelope.

VERIFICATION: The above substances were identified by predistribution laboratories, who used the following combined list of techniques: color tests, FTIR, and GC/MS.

Summary Comments

This test was designed to allow participants to assess their proficiency in the qualitative analysis of questioned samples that may contain controlled substances. Item 1 consisted of approximately 180 mg of powder consisting of 100% testosterone. Item 2 consisted of approximately 300 mg of powder consisting of caffeine. In the United States, testosterone is a Schedule III controlled substance (See Manufacturer's Information for preparation details).

For Item 1, 548 of 556 participants (98.5%) reported testosterone. The remaining eight participants reported that no controlled substance was identified. For Item 2, 554 of 556 participants (99.6%) reported that no controlled substance was detected/identified. The two remaining participants reported testosterone.

The most common examination methods reported for both items include GC/MS, color tests, and FTIR.

Drug Analysis (Item 1)

TABLE 1a - Item 1 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
238KUR	testosterone	37PXCU	testosterone
23CH6B	Testosterone	37RMLX	Testosterone
24FTY6	testosterone	38JQ8D	Testosterone
24JECD	TESTOSTERONE	39THV8	testosterone
26NY8Z	testosterone (Schedule III)	3BT829	Testosterone
28HPDE	testosterone	3G66TU	Testosterone
28VEJC	Testosterone - Schedule III. Net weight of material - 0.18 (+/- 0.04) gram(s). (confidence level 99.7%)	3HYM3H	Testosterone
28XYBL	Testosterone	3KZPU6	Testosterone
29CDWA	Testosterone / Schedule III	3LEETB	testosterone
2BUPEC	Testosterone	3MCX3D	Testosterone.
2BXA7L	Testosterone	3MN9PB	Testosterone
2EGD74	Testosterone	3NNBP2	testosterone
2JCWHB	Testosterone	3Q7JLF	Testosterone
2KNPGJ	Testosterone	3Q8GWB	Testosterone
2NCGPN	Testosterone	3TRBUC	Testosterone
2NNXME	Testosterone	3TXDRW	No controlled substance detected/identified
2RA3GQ	Testosterone	3WWM3Y	Testosterone
2TG7ME	testosterone	3X8MRD	Testosterone
2TZRE6	Testosterone	3ZF27Y	Testosterone
2UJNB7	Testosterone	3ZGK8G	Testosterone
2UUM9T	Testosterone	43URZ9	TESTOSTERONE.
2WYCB3	Testosterone, a Schedule III controlled substance	462AY4	Testosterone
2Y2BB9	Testosterone / Schedule III	47XKYD	Testosterone
2YAH4D	Testosterone	48BQFV	Testosterone
2YPE76	Testosterone found CIII	4FALQT	testosterone
2Z3DQL	Testosterone	4GGBLP	testosterone
32RDBV	Testosterone	4GJZPL	Testosterone
36W96L	TESTOSTERONE	4M84WC	Testosterone
		4NGPBW	testosterone

TABLE 1a - Item 1 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
4PUHKL	Testosterone	7HBD8Y	Testosterone
4QQC7Q	testosterone	7HUFEZ	Testosterone
4RHQEZ	Testosterone	7KAQBN	Testosterone
4TBEXV	Testosterone	7MJRFV	Testosterone
4VKCKL	testosterone	7ND8KK	Testosterone
4XN6EJ	testosterone	7TCX87	Testosterone
4YWPBY	Testosterone	7W63LG	Testosterone
4YXCFR	Testosterone, Schedule III. Salt form not determined	7YDCT8	No controlled substance detected/identified
64W967	Testosterone	7Z4M47	Testosterone
67DVE6	testosterone	82VC88	Testosterone
67TLHA	Testosterone	83T3GM	Testosterone
6A2A3E	Testosterone	8443GN	Testosterone
6DV99W	Testosterone	84NF9L	Testosterone
6DVL2C	Testosterone	84ZJBC	Testosterone
6FHHA3	Testosterone	87WWDY	Testosterone
6MKAWE	Testosterone	882FBF	testosterone
6QG3L4	Testosterone	8D87F7	Testosterone
6QKH6V	Testosterone	8DMJ7N	Testosterone
6RF27W	testosterone	8L3VNV	Testosterone Indicated
6VJEP4	Testosterone	8MDWAU	Testosterone
6VJU62	Testosterone	8N8KRQ	Testosterone
72DKHN	Testosterone, Schedule III	8PZ7A3	Testosterone (a controlled substance)
72VJQF	Testosterone	8QDE4F	Testosterone
73MH28	Testosterone	8TGY3A	Testosterone
77TC8C	Testosterone	8UDEM L	testosterone
77TWYX	Testosterone	8UNJM3	Testosterone/Schedule III
78QZT2	Testosterone	8XQDPE	Testosterone
79T9A6	Testosterone	8YMTX6	Testosterone
7CQCP2	Testosterone	927Z2Y	testosterone
7CVW49	testosterone	94CME3	Testosterone

TABLE 1a - Item 1 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
967BWX	testosterone	A9RR9Y	Testosterone
97LDYZ	TESTOSTERONE	ADF97M	Testosterone (III)
98KPRM	containing testosterone	ADQJH4	No controlled substance detected/identified
9996WN	An Anabolic Steroid (Testosterone), a schedule III controlled non-narcotic substance.	AF3LWT	Testosterone
99TYGY	Testosterone	AFNNV2	testosterone
9EYF98	Testosterone	AFP33M	TESTOSTERONE
9EZARD	Testosterone	ANMNEG	testosterone
9GRUW6	Testosterone	APEK2P	Testosterone
9H39AX	Testosterone	AT82HN	Testosterone (C-III) Found
9HCA2Y	testosterone	AU8T2P	Item #1.1 contains Testosterone, an anabolic steroid, a schedule III controlled non-narcotic substance.
9HDAQQ	Testosterone	AX8UCK	Testosterone
9J9TTR	Testosterone	AZEC9V	Testosterone.
9JBH2V	testosterone	AZED83	Testosterone
9LGUT6	Testosterone	B3C9XN	testosterone
9LTRXA	Testosterone	B6KVZF	Testosterone (17beta-hydroxyandrost-4-en-3-one)
9MJE6E	TESTOSTERONE	B84J8L	testosterone
9MP33P	Testosterone	B9K2G8	Testosterone
9NKL4Q	Testosterone	B9KL88	Testosterone
9PJQKK	Testosterone	BEQ9RW	testosterone
9QRJAZ	testosterone	BGFGBX	Testosterone
9RRX6T	Contains Testosterone	BLPQBE	Testosterone
9YFYVB	Testosterone	BMGUXU	Testosterone
9ZCCH3	Testosterone	BPWUPQ	Contains Testosterone
A29ZQQ	Testosterone	BRXQLR	testosterone
A2QEAG	Testosterone	BYABMV	Testosterone
A7HCC6	Testosterone	BYU9H9	Testosterone
A8WQ8U	Testosterone	BYXFFV	testosterone
A93AQN	Testosterone was identified but was reported as an anabolic steroid, a schedule III controlled non-narcotic substance as per procedure.	BZ7NC9	Testosterone
		C2ZV4D	Testosterone

TABLE 1a - Item 1 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
C3CGUE	Testosterone	DPM2LY	Contains Testosterone
CAE8B7	Testosterone	DPZMXF	testosterone
CEKBWQ	Testosterone	DR6DRN	Testosterone
CF2EZJ	Testosterone	DTJAN2	Testosterone
CGTXWX	testosterone	DTZPV4	Testosterone
CHXEHM	Testosterone	DU8DQZ	Testosterone
CJX3CJ	Testosterone	DUA2AY	Testosterone
CQB9XP	testosterone	DZ277A	Testosterone
CQUZNZ	testosterone	DZ2YB9	Testosterone
CTFBHV	Testosterone	DZYLBU	Testosterone
CTZJMW	Testosterone	E2N7GP	Testosterone
CUXTV7	Testosterone	E4KUQU	Testosterone
CWL8KF	Testosterone	E4WEPA	Testosterone
CZ4CAY	Testosterone	E6NHCP	Testosterone
D2AG2R	testosterone	E6UCUC	testosterone
D386KV	Testosterone (a controlled substance)	E74F6F	Testosterone
D4K8GN	Testosterone	EC8EL7	TESTOSTERONE
D4MWKK	Testosterone	ECF9UM	Contains Testosterone
D73FD3	Testosterone	EDPEVN	Testosterone
DAMAA3	Testosterone (Schedule III)	EEJ3EJ	Testosterone III
DAVM6Y	Testosterone	EGL6V8	Testosterone
DCQA7T	Testosterone	EHJCJL	Testosterone
DCQZH3	Testosterone	EJRG3	Testosterone
DDQJNX	Testosterone	EK2L6C	Testosterone
DDQY2B	Testosterone	EL94Y2	Testosterone
DDY7RX	Testosterone	ELTPRR	Testosterone
DJTCEZ	Testosterone	ERYHLM	Testosterone
DK6R92	Testosterone	ETYEFQ	Testosterone
DKFBHT	Testosterone	EW22JG	Testosterone
DKNVG2	Testosterone	EWEQDT	Testosterone
DNNPQP	Testosterone	EYX4UA	Testosterone

TABLE 1a - Item 1 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
F4D8N7	Testosterone	GZCXXG	Testosterone
F6LYGZ	Testosterone	H2YWPK	testosterone
F99D7A	Testosterone	H4N8L6	Testosterone.
FA3KPU	Testosterone	H6BXGH	Testosterone
FBUHHD	testosterone	H9RAAR	Testosterone
FJ8R6U	Testosterone	HAAD3C	Testosterone
FJZFWG	Testosterone	HD3MQQ	testosterone
FLV7U4	Testosterone	HEX7TR	Testosterone.
FLZ2HN	Testosterone	HL9EV4	testosterone
FM6UV8	Testosterone	HNCERU	Testosterone
FPX2ER	Testosterone	HQGEWD	Testosterone
FREV9Y	testosterone	HQUZCG	Testosterone
FRKFNE	Testosterone	HVJNLT	Testosterone
FUP9XZ	Testosterone	HW8P7C	Testosterone
FV3636	Testosterone	HYLRBU	testosterone
FWDMLB	Testosterone	HYVLR4	Testosterone
FWRV6Q	Testosterone	J42AGH	Testosterone
FXFDZW	Testosterone, Boldenone	J4LYTC	Testosterone, Schedule III
FYM6PT	Testosterone	J6RK7Q	Testosterone
GCZ9YT	Testosterone in base form.	J8ZLGP	Testosterone
GEUTA8	Testosterone	J982NV	Testosterone
GGFR9M	Testosterone	JA8R6B	Testosterone
GML798	testosterone	JC3WX7	testosterone
GQG739	testosterone	JCBYUG	Testosterone was detected in Lab Item 1.
GRRTUV	Testosterone	JDMELW	Testosterone
GVTYBL	Testosterone	JE2C7T	Testosterone
GXDVRY	Testosterone	JE7PTP	Testosterone
GXHAFI	Testosterone	JEKPEJ	Testosterone
GYERDG	Testosterone (an anabolic steroid), a schedule III controlled non-narcotic substance.	JFEFHC	Testosterone
GYQLHV	Testosterone	JGCM8J	Testosterone
		JJ8PJX	testosterone

TABLE 1a - Item 1 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
JJZLFV	testosterone	L7DYTQ	Testosterone
JKNLCX	Testosterone	L7XJTW	Testosterone
JQ8TPW	Testosterone	L9488L	Testosterone
JQP67F	Testosterone	L98TKT	Testosterone
JTHVF2	Testosterone	LALGG	Testosterone.
JUCJUQ	Testosterone	LBD4Q6	Testosterone
JUP7C6	Testosterone	LBU9A4	Testosterone
JVWAGG	Testosterone	LJN42E	Testosterone
JXPKN	testosterone	LLTXXR	Testosterone
K3RPF	Testosterone	LM6V78	Testosterone
KC438N	Testosterone	M244JF	Testosterone
KCPJK7	Testosterone	M2UYKL	Testosterone
KCR9P3	Testosterone	M49J7L	testosterone
KDJBCH	Testosterone	M786MQ	Testosterone
KF9G6F	Testosterone	M8T8LY	Testosterone
KF9JQD	testosterone	MBAP4W	Testosterone
KFALCC	Testosterone, Schedule III	ME4WYQ	testosterone
KLF89L	Testosterone	MHCYCA	Testosterone
KLJZJV	Testosterone (17Beta-Hydroxyandrost-4-en-3-one)	MRCMXR	Testosterone
KMAPJV	Testosterone	MU88ZJ	Testosterone
KMLQXG	Testosterone	MULY73	testosterone
KNLMRK	testosterone	MVXBM	Testosterone
KVQAUM	Testosterone	MW6T8N	Testosterone
KXB4QM	Testosterone	MW9XYQ	testosterone
KYRGAK	Testosterone	MWRKVM	testosterone
L3HZKM	TESTOSTERONE	N7FX9R	Testosterone
L4W2R3	Testosterone	N7HZ2Q	Testosterone
L6ZFTU	Testosterone	NA6ZPP	Testosterone
L76LMB	testosterone	NB822W	Testosterone
L7BGKU	Testosterone	ND74JX	testosterone
		NFHR6H	Testosterone

TABLE 1a - Item 1 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
NFX3XV	Testosterone	PZDTGK	Testosterone
NGUCC8	Testosterone, Schedule III	PZY7J4	testosterone
NGXCCU	Testosterone	Q27ZLG	Testosterone
NL87AJ	No controlled substance detected/identified	Q2N2CQ	Testosterone
NQEUDT	No controlled substance detected/identified	Q4H4NT	Testosterone
NRNYVM	Testosterone	Q6NETN	Testosterone
NRQJMW	testosterone	Q9TTPU	Testosterone
NRYFJH	Testosterone	Q9WBW6	Testosterone
NT76WF	Testosterone	QAFLEH	Testosterone
NUHWU3	Testosterone	QAL6NM	Testosterone
NUTKMU	Testosterone	QFDWWW	testosterone
NYFFV2	Testosterone (CIII)	QG3FKF	Testosterone
NZ8AVQ	Testosterone	QGPL4F	testosterone
P2DQCE	Testosterone	QGZMRA	Testosterone
P46UXU	Testosterone	QHXYLG	Testosterone
P7UZWY	Testosterone	QMTKMR	Testosterone
P8KPXK	Testosterone	QRFEFJ	Testosterone
P8VECA	Testosterone	QUABWQ	Testosterone
PAR3FL	Testosterone	QURDNY	testosterone
PCVDG9	Testosterone	QVHJ24	Testosterone
PEYC6L	Testosterone	QWFRC9	Testosterone
PFR2HN	testosterone	QWWPLZ	Testosterone
PPW9XN	Testosterone	QXW3U	0.1779 ± 0.0018 grams of Testosterone (Schedule III)
PR66XJ	Yes, TESTOSTERONE which ia a controlled substance	QYHKX8	Testosterone
PT3NQD	Testosterone III	QYL6PG	Testosterone
PUCG2E	Testosterone	R2Q4NF	testosterone
PUTU4H	Testosterone	R6K3X6	Testosterone
PUURDT	Testosterone	R6PKW4	Testosterone
PVP9RB	Testosterone	R7KYR9	Testosterone
		R9QGQ6	Testosterone

TABLE 1a - Item 1 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
RB3Y7V	Testosterone	TM3DA7	Testosterone (C-III) found.
RBTHEF	Testosterone	TRF4PY	Testosterone
REY8BH	testosterone	TRW98W	Testosterone
RG9VPV	No controlled substance detected/identified	TTADVG	Testosterone
RJNJRD	Testosterone	TUKPG6	Testosterone
RKZ26Q	Testosterone	TYENUE	testosterone
RLFTBA	Testosterone	U7PGCN	Testosterone
RN3MUV	Testosterone	UAQD44	Testosterone
RP7NQF	Testosterone	UB636V	Testosterone
RPWBDR	Testosterone	UCYQMQ	Testosterone (salt form not determined)
RR3CEW	TESTOSTERONE	UD8MBG	Testosterone
RTBAMR	TESTOSTERONE	UFEC6C	Testosterone
RUJQCB	Testosterone	UFHBM9	testosterone
RWFCUP	Testosterone	UFV76M	Testosterone
RXM93R	Testosterone	UFX2JQ	Testosterone
T3QQLC	Testosterone	UJ973D	Testosterone
T8ACKY	testosterone	UMBWKM	testosterone
T8C2PV	Testosterone	UMJZBK	Testosterone
T9BH88	Testosterone	UQLTDW	Testosterone
T9CRHG	Testosterone	UQTF6G	Testosterone.
T9Q4XP	testosterone	UVC8UQ	Testosterone
TDRBN6	Testosterone	UX43YE	testosterone
TFRWAD	Testosterone	UYZGL8	Testosterone
TGV2YQ	TESTOSTERONE	UZ2CH8	Testosterone
TH42CK	Testosterone	UZ8KTF	testosterone
THAEP9	Testosterone	V2HJNP	Testosterone
THUJJN	Testosterone	V3HK2K	Testosterone
TJR87D	Testosterone	V3UWB7	Contains Testosterone
TLAX69	Testosterone (Anabolic Steroid)	V9CF6D	Testosterone
TLPVYD	Androstenedione*, Testosterone*, *Presumptively identified	VAJX2M	Testosterone
		VBN6UB	Testosterone

TABLE 1a - Item 1 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
VCJMVD	Testosterone	WWWWB9	Testosterone
VCQXF3	Testosterone	WZV889	Testosterone
VE3JGP	Testosterone	WZXQ6J	Testosterone
VE4YQN	Testosterone	X26JHH	Testosterone
VEFG4L	Testosterone	X4XVGB	Testosterone
VEJWKT	testosterone	X83WK3	Testosterone
VEKQJT	TESTOSTERONE	X9RUPF	Testosterone
VGLRD4	Testosterone	XBHVWN	Testosterone
VMFP2Z	testosterone	XGLGMJ	Testosterone
VU9UXG	testosterone	XLFWCD	Testosterone
VUCEQQ	Testosterone	XUMTVB	Testosterone
VURVKJ	Testosterone	Y3F9C3	Testosterone
VXV6YF	Testosterone	Y7CKP7	Testosterone found.
W6AR83	Testosterone, an anabolic steroid, a schedule III controlled non-narcotic substance.	Y7FJ82	Testosterone
W74W3H	Testosterone	Y7KBN7	Contains Testosterone
W7KUVM	Testosterone	Y7ZQGD	Testosterone
W8F9DK	Testosterone	Y8MKQ8	Testosterone
W8FAA6	Testosterone	YFEV3M	Testosterone
WBKAK	No controlled substance detected/identified	YGBBBD	Testosterone
WCBW8Y	testosterone	YK9P8B	No controlled substance detected/identified; Testosterone
WG2G33	Testosterone	YMKQM4	testosterone
WJ8VJF	Testosterone	YWQNDT	Testosterone
WLWVPB	Testosterone	YXGBP9	testosterone
WMAFWF	Testosterone	Z2CPNH	TESTOSTERONE
WMRNR3	Testosterone	Z2CTUE	Testosterone
WMUEZ7	Testosterone	Z2WCB8	Testosterone
WN87Q2	Testosterone (salt form not determined)	Z2WGMV	Testosterone
WQMPJJ	Testosterone	Z3BKGM	Testosterone
WUTX32	Testosterone	Z77LBN	Testosterone
		Z7BRZY	Testosterone

TABLE 1a - Item 1 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
Z7J8RC	Testosterone		
Z86Z7G	Testosterone		
ZACKPB	Testosterone		
ZB8J9V	No controlled substance detected/identified		
ZBDRE6	Testosterone		
ZEKFDX	testosterone		
ZEM7M2	Testosterone		
ZF3UY9	Testosterone		
ZKP9FF	Testosterone		
ZPHD8Y	Testosterone (preliminary ID)		
ZQT9VC	Testosterone		
ZUMNX2	Testosterone		
ZW2FE2	Testosterone		
ZWK27Y	testosterone		
ZXRQ39	Testosterone		
ZYGAKA	testosterone		
ZYNBY9	Testosterone		
ZYX2FY	Testosterone		
ZZJT2A	testosterone		

<i>Response Summary - Item 1</i>	
Participants: 556	
Testosterone:	548
No Controlled substance detected/identified:	8

Examination Methods (Item 1)

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
238KUR		✓							✓	
23CH6B									✓	TruNarc
24FTY6	✓								✓	
24JECD		✓							✓	
26NY8Z		✓					✓		✓	
28HPDE	✓	✓		✓					✓	
28VEJC	✓	✓							✓	
28XYBL		✓							✓	
29CDWA	✓	✓							✓	
2BUPEC			✓	✓				✓	✓	
2BXA7L		✓						✓	✓	
2EGD74	✓								✓	
2JCWHB	✓				✓				✓	
2KNPGJ	✓								✓	
2NCGPN	✓								✓	
2NNXME	✓				✓				✓	
2RA3GQ		✓					✓		✓	
2TG7ME	✓								✓	
2TZRE6	✓								✓	
2UJNB7	✓	✓							✓	
2UUM9T	✓							✓	✓	
2WYCB3	✓	✓							✓	
2Y2BB9	✓	✓							✓	
2YAH4D	✓	✓							✓	
2YPE76	✓	✓		✓					✓	
2Z3DQL	✓								✓	
32RDBV	✓								✓	
36W96L									✓	
37PXCW	✓								✓	
37RMLX	✓				✓				✓	
38JQ8D	✓								✓	GC/IR
39THV8		✓							✓	DART Agilent-TOF
3BT829	✓	✓						✓	✓	
3G66TU	✓				✓				✓	
3HYM3H	✓								✓	
3KZPU6	✓								✓	
3LEETB	✓								✓	FTIR with ATR accessory
3MCX3D	✓								✓	

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
3MN9PB	✓								✓	
3NNBP2	✓	✓							✓	
3Q7JLF	✓								✓	
3Q8GWB	✓				✓				✓	
3TRBUC		✓							✓	
3TXDRW	✓								✓	
3WWM3Y									✓	
3X8MRD	✓	✓		✓					✓	
3ZF27Y					✓				✓	
3ZGK8G	✓								✓	
43URZ9									✓	Two different GC/MS used
462AY4	✓								✓	
47XKYD		✓							✓	
48BQFV					✓				✓	
4FALQT					✓				✓	
4GGBLP		✓							✓	
4GJZPL	✓	✓			✓				✓	
4M84WC	✓								✓	GC/IR
4NGPBW									✓	
4PUHKL		✓							✓	XRD
4QQC7Q		✓							✓	
4RHQEZ	✓	✓							✓	
4TBEXV	✓				✓				✓	
4VKCKL	✓	✓							✓	
4XN6EJ	✓			✓					✓	
4YWBPY	✓	✓								
4YXCFR					✓				✓	
64W967					✓				✓	
67DVE6	✓								✓	
67TLHA	✓	✓			✓				✓	
6A2A3E		✓							✓	
6DV99W					✓				✓	
6DVL2C	✓	✓							✓	
6FHHA3	✓	✓							✓	
6MKAWE	✓								✓	
6QG3L4	✓	✓		✓		✓			✓	
6QKH6V		✓			✓				✓	
6RF27W	✓	✓							✓	
6VJEP4	✓								✓	

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
6VJU62	✓								✓	
72DKHN					✓				✓	
72VJQF		✓						✓	✓	
73MH28	✓								✓	
77TC8C				✓					✓	
77TWYX					✓				✓	
78QZT2	✓			✓					✓	
79T9A6									✓	
7CQCP2	✓			✓					✓	
7CVW49	✓	✓							✓	
7HBD8Y	✓	✓							✓	
7HUF EZ	✓	✓		✓					✓	
7KAQBN		✓							✓	
7MJRFV	✓	✓							✓	
7ND8KK		✓							✓	XRD
7TCX87	✓	✓		✓					✓	
7W63LG	✓	✓		✓					✓	
7YDCT8									✓	
7Z4M47	✓			✓					✓	
82VC88	✓								✓	
83T3GM	✓								✓	
8443GN	✓								✓	
84NF9L									✓	
84ZJBC	✓								✓	DART-TOF
87WWDY	✓			✓					✓	
882FBF		✓						✓	✓	
8D87F7					✓				✓	
8DMJ7N	✓								✓	
8L3VNV	✓							✓	✓	
8MDWAU	✓	✓							✓	
8N8KRQ	✓				✓				✓	
8PZ7A3		✓							✓	DART-Agilent TOF
8QDE4F					✓				✓	FID
8TGY3A		✓							✓	
8UDEML	✓	✓							✓	
8UNJM3	✓	✓							✓	
8XQDPE		✓						✓	✓	
8YMTX6	✓	✓							✓	
927Z2Y									✓	

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
94CME3	✓				✓				✓	
967BWX	✓								✓	
97LDYZ	✓								✓	
98KPRM									✓	
9996WN	✓								✓	
99TYGY	✓								✓	
9EYF98	✓			✓					✓	
9EZARD	✓								✓	
9GRUW6	✓								✓	
9H39AX	✓				✓				✓	
9HCA2Y	✓								✓	
9HDAQQ	✓	✓			✓				✓	
9J9TTR									✓	
9JBH2V	✓								✓	
9LGUT6	✓			✓					✓	
9LTRXA	✓								✓	
9MJE6E									✓	
9MP33P	✓	✓			✓				✓	
9NKL4Q	✓				✓				✓	
9PJQKK		✓							✓	
9QRJAZ									✓	DART-TOF
9RRX6T	✓	✓		✓					✓	
9YFYVB	✓	✓							✓	
9ZCCH3		✓							✓	
A29ZQQ	✓				✓				✓	
A2QEAG		✓							✓	
A7HCC6								✓	✓	
A8WQ8U					✓				✓	TruNarc
A93AQN	✓								✓	
A9RR9Y	✓	✓							✓	
ADF97M	✓				✓				✓	
ADQJH4									✓	
AF3LWT	✓								✓	
AFNNV2	✓								✓	
AFP33M					✓				✓	
ANMNEG					✓				✓	
APEK2P	✓				✓				✓	
AT82HN	✓	✓							✓	
AU8T2P	✓								✓	

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
AX8UCK	✓				✓				✓	
AZEC9V	✓								✓	
AZED83	✓			✓					✓	
B3C9XN	✓	✓						✓	✓	
B6KVZF		✓								X-Ray Diffractometry
B84J8L	✓							✓	✓	
B9K2G8									✓	
B9KL88	✓								✓	
BEQ9RW		✓							✓	
BGFGBX	✓								✓	
BLPQBE	✓								✓	
BMGUXU	✓	✓							✓	
BPWUPQ	✓			✓					✓	
BRXQLR					✓				✓	
BYABMV	✓	✓							✓	
BYU9H9	✓								✓	
BYXFFV	✓	✓							✓	
BZ7NC9	✓								✓	
C2ZV4D					✓				✓	
C3CGUE		✓			✓				✓	
CAE8B7					✓				✓	
CEKBWQ	✓	✓								
CF2EZJ	✓	✓							✓	
CGTXWX									✓	DART-Agilent ToF
CHXEHM		✓							✓	
CJX3CJ	✓				✓					GC-IRD
CQB9XP		✓							✓	weight determination
CQUZNZ	✓	✓						✓	✓	
CTFBHV	✓			✓					✓	
CTZJMW	✓								✓	
CUXTV7	✓								✓	
CWL8KF		✓							✓	
CZ4CAY	✓	✓							✓	
D2AG2R		✓							✓	
D386KV									✓	DART Agilent TOF
D4K8GN	✓								✓	
D4MWKK	✓	✓			✓				✓	
D73FD3	✓	✓							✓	
DAMAA3		✓						✓	✓	

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
DAVM6Y	✓								✓	TruNarc
DCQA7T	✓								✓	
DCQZH3	✓								✓	
DDQJNX	✓	✓					✓		✓	
DDQY2B	✓	✓		✓					✓	
DDY7RX	✓								✓	
DJTCEZ	✓								✓	FTIR with ATR accessory
DK6R92		✓							✓	
DKFBHT	✓	✓		✓					✓	
DKNVG2	✓				✓				✓	
DNNPQP								✓	✓	
DPM2LY	✓	✓		✓					✓	
DPZMXF	✓								✓	
DR6DRN									✓	Weight
DTJAN2	✓								✓	
DTZPV4	✓								✓	
DU8DQZ	✓								✓	
DUA2AY	✓								✓	net weight=0.176 (+/- 0.001) grams
DZ277A		✓							✓	
DZ2YB9	✓				✓				✓	FID, DART-MS
DZYLBU	✓								✓	
E2N7GP	✓	✓							✓	
E4KUQU	✓	✓							✓	
E4WEPA	✓								✓	
E6NHCP	✓	✓			✓				✓	
E6UCUC					✓				✓	
E74F6F	✓				✓				✓	
EC8EL7	✓	✓							✓	
ECF9UM	✓	✓		✓					✓	
EDPEVN	✓	✓							✓	
EEJ3EJ	✓	✓			✓				✓	
EGL6V8	✓								✓	
EHJCJL	✓	✓		✓					✓	
EJRG3		✓							✓	
EK2L6C	✓								✓	
EL94Y2		✓							✓	
ELTPRR	✓							✓	✓	
ERYHLM	✓	✓								
ETYEFQ	✓								✓	

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
EW22JG	✓	✓							✓	
EWEQDT	✓								✓	
EYX4UA	✓	✓							✓	
F4D8N7	✓	✓							✓	GC retention time comparison with standard
F6LYGZ		✓							✓	
F99D7A	✓	✓			✓				✓	
FA3KPU					✓					
FBUHHD	✓								✓	
FJ8R6U									✓	
FJZFWG							✓		✓	weight determination
FLV7U4		✓							✓	Microscope, XRF, XRD
FLZ2HN	✓			✓					✓	
FM6UV8	✓								✓	GC/IR
FPX2ER	✓	✓							✓	
FREV9Y	✓								✓	FTIR-ATR
FRKFNE	✓	✓							✓	
FUP9XZ	✓								✓	GC/IR
FV3636		✓			✓				✓	GC for retention time confirmation
FWDMLB		✓							✓	Electronic Balance
FWRV6Q	✓	✓							✓	
FXFDZW							✓			
FYM6PT	✓								✓	
GCZ9YT	✓								✓	
GEUTA8	✓								✓	
GGFR9M	✓	✓							✓	
GML798	✓	✓							✓	
GQG739		✓							✓	
GRRTUV		✓							✓	
GVTYBL	✓								✓	
GXDVRY			✓						✓	
GXHAFI	✓								✓	
GYERDG	✓								✓	
GYQLHV	✓							✓	✓	
GZCXXG	✓								✓	
H2YWPK		✓							✓	
H4N8L6	✓	✓		✓					✓	
H6BXGH		✓							✓	

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
H9RAAR	✓	✓							✓	
HAAD3C	✓	✓							✓	
HD3MQQ	✓							✓	✓	
HEX7TR		✓			✓				✓	
HL9EV4		✓							✓	weight, visual examination
HNCERU	✓							✓	✓	
HQGEWD	✓								✓	
HQUZCG		✓							✓	
HVJNLT	✓								✓	
HW8P7C	✓	✓							✓	
HYLRBU	✓								✓	
HYVLR4	✓	✓			✓				✓	Gas Chromatography retention time verification in comparison with a standard
J42AGH	✓	✓							✓	Atmospheric Sample Analysis Probe-MS
J4LYTC					✓				✓	
J6RK7Q	✓								✓	
J8ZLGP	✓								✓	
J982NV								✓	✓	
JA8R6B	✓	✓							✓	
JC3WX7	✓	✓							✓	
JCBYUG									✓	Weight
JDMELW	✓								✓	
JE2C7T	✓							✓	✓	
JE7PTP	✓								✓	
JEKPEJ		✓			✓				✓	
JFEFHC	✓				✓				✓	
JGCM8J	✓				✓				✓	
JJ8PJX	✓	✓						✓	✓	
JJZLFV	✓								✓	
JKNLCX	✓								✓	
JQ8TPW		✓							✓	
JQP67F		✓							✓	
JTHVF2	✓								✓	
JUCJUQ	✓	✓			✓				✓	
JUP7C6	✓								✓	
JVWAGG	✓								✓	GC/FID
JXPKVN	✓								✓	

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
K3RPF	✓	✓							✓	
KC438N	✓								✓	
KCPJK7		✓							✓	Polarized Light Microscopy
KCR9P3		✓							✓	
KDJBCH		✓							✓	
KF9G6F	✓				✓				✓	
KF9JQD	✓								✓	
KFALCC	✓				✓				✓	
KLF89L	✓				✓				✓	
KLJZJV		✓							✓	
KMAPJV					✓				✓	
KMLQXG	✓								✓	
KNLMRK	✓								✓	
KVQAUM	✓				✓				✓	
KXB4QM	✓								✓	
KYRGAK	✓	✓		✓					✓	
L3HZKM	✓								✓	
L4W2R3	✓								✓	
L6ZFTU	✓	✓							✓	
L76LMB					✓				✓	
L7BGKU	✓	✓							✓	
L7DYTQ	✓								✓	
L7XJTW	✓								✓	
L9488L	✓					✓			✓	
L98TKT	✓				✓				✓	
LAJLGG					✓	✓			✓	
LBD4Q6	✓								✓	
LBU9A4		✓							✓	
LJN42E	✓	✓				✓			✓	
LLTXR	✓	✓							✓	
LM6V78	✓								✓	
M244JF	✓				✓				✓	
M2UYKL	✓			✓					✓	
M49J7L	✓								✓	
M786MQ					✓				✓	
M8T8LY									✓	Gas Chromatography retention time verification in comparison with a standard

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
MBAP4W									✓	GC/IR
ME4WYQ	✓								✓	
MHCYCA				✓					✓	
MRCMXR		✓	✓	✓	✓		✓		✓	H-NMR
MU88ZJ	✓	✓			✓				✓	
MULY73	✓								✓	
MVXBMY	✓	✓							✓	
MW6T8N	✓								✓	
MW9XYQ	✓	✓							✓	
MWRKVM									✓	TruNarc
N7FX9R	✓	✓						✓	✓	
N7HZ2Q	✓								✓	
NA6ZPP								✓	✓	
NB822W		✓							✓	GC retention time
ND74JX		✓						✓	✓	
NFHR6H	✓								✓	
NFX3XV	✓	✓							✓	
NGUCC8					✓				✓	
NGXCCU				✓	✓				✓	
NL87AJ	✓			✓					✓	
NQEU DT									✓	
NRNYVM	✓								✓	
NRQJMW		✓							✓	
NRYFJH	✓	✓		✓					✓	
NT76WF					✓				✓	
NUHWU3		✓							✓	Balance, XRD
NUTKMU	✓								✓	
NYFFV2	✓								✓	
NZ8AVQ	✓	✓							✓	
P2DQCE	✓				✓				✓	
P46UXU	✓	✓							✓	
P7UZWY									✓	
P8KPXK									✓	
P8VECA		✓							✓	
PAR3FL	✓								✓	
PCVDG9					✓				✓	
PEYC6L	✓	✓							✓	
PFR2HN	✓								✓	
PPW9XN	✓								✓	

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
PR66XJ	✓								✓	
PT3NQD	✓			✓	✓				✓	
PUCG2E		✓			✓				✓	
PUTU4H	✓								✓	
PUURDT									✓	GCIR
PVP9RB									✓	
PZDTGK		✓							✓	Weight Determination, DART-TOFMS
PZY7J4		✓							✓	
Q27ZLG	✓								✓	
Q2N2CQ	✓								✓	
Q4H4NT	✓								✓	
Q6NETN								✓	✓	
Q9TTPU									✓	
Q9WBW6					✓				✓	
QAFLEH	✓								✓	
QAL6NM	✓								✓	
QFDWWW	✓	✓							✓	
QG3FKF	✓				✓				✓	
QGPL4F	✓	✓		✓					✓	
QGZMRA		✓							✓	
QHXYLG	✓								✓	
QMTKMR	✓								✓	
QRFEFJ	✓								✓	
QUABWQ									✓	GC/IR
QURDNY					✓				✓	
QVHJ24		✓			✓				✓	
QWFRC9	✓				✓				✓	
QWWPLZ	✓	✓							✓	
QXV3U	✓					✓				Weight Measurement, DART-TOF, GC-FID-MS
QYHKX8	✓								✓	Weight determination
QYL6PG	✓								✓	
R2Q4NF	✓								✓	
R6K3X6	✓	✓			✓				✓	
R6PKW4	✓				✓				✓	
R7KYR9	✓				✓				✓	
R9QGQ6					✓				✓	
RB3Y7V									✓	
RBTHEF	✓								✓	

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
REY8BH	✓				✓				✓	TruNarc
RG9VPV									✓	
RJNJRD	✓								✓	
RKZ26Q	✓								✓	
RLFTBA	✓	✓							✓	
RN3MUV	✓								✓	
RP7NQF	✓	✓							✓	
RPWBDR	✓								✓	
RR3CEW									✓	
RTBAMR									✓	
RUJQCB	✓								✓	
RWFCUP	✓								✓	
RXM93R	✓								✓	
T3QQLC	✓	✓		✓					✓	
T8ACKY		✓					✓		✓	
T8C2PV	✓								✓	
T9BH88	✓			✓					✓	
T9CRHG	✓	✓							✓	
T9Q4XP									✓	retention time
TDRBN6	✓								✓	
TFRWAD	✓	✓		✓					✓	
TGV2YQ	✓								✓	
TH42CK					✓				✓	
THAEP9		✓			✓				✓	
THUJJN	✓	✓							✓	
TJR87D	✓			✓					✓	
TLAX69	✓	✓							✓	
TLPVYD									✓	
TM3DA7	✓								✓	
TRF4PY	✓	✓			✓				✓	
TRW98W					✓				✓	
TTADVG									✓	DART Agilent-TOF
TUKPG6	✓								✓	
TYENUE	✓					✓		✓	✓	
U7PGCN	✓	✓							✓	
UAQD44	✓				✓				✓	
UB636V	✓				✓				✓	
UCYQMQ	✓	✓							✓	
UD8MBG	✓	✓							✓	

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
UFEC6C	✓								✓	
UFHBM9					✓				✓	
UFV76M	✓	✓						✓	✓	
UFX2JQ		✓							✓	
UJ973D	✓	✓							✓	
UMBWKM	✓								✓	
UMJZBK	✓								✓	GCIR
UQLTDW					✓				✓	
UQTF6G	✓							✓	✓	
UVC8UQ	✓								✓	
UX43YE	✓	✓		✓					✓	
UYZGL8	✓			✓					✓	
UZ2CH8					✓				✓	
UZ8KTF		✓							✓	DART-AccuTOF
V2HJNP	✓	✓							✓	GC retention time verification in comparison with a standard
V3HK2K	✓								✓	
V3UWB7	✓	✓		✓					✓	
V9CF6D		✓			✓				✓	
VAJ2M	✓								✓	
VBN6UB	✓	✓		✓					✓	
VCJMVD	✓	✓							✓	
VCQXF3									✓	weight determination
VE3JGP		✓							✓	
VE4YQN	✓								✓	
VEFG4L	✓	✓						✓	✓	
VEJWKT	✓	✓		✓					✓	
VEKQJT									✓	
VGLRD4	✓	✓							✓	
VMFP2Z	✓				✓				✓	
VU9UXG	✓								✓	
VUCEQQ		✓						✓	✓	
VURVKJ					✓				✓	
VXV6YF									✓	DART Agilent TOF
W6AR83	✓								✓	
W74W3H	✓								✓	
W7KUVM	✓								✓	
W8F9DK		✓							✓	

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
W8FAA6	✓	✓							✓	
WBEKAK									✓	
WCBW8Y	✓								✓	
WG2G33	✓				✓				✓	
WJ8VJF					✓				✓	
WLWVPB		✓							✓	
WMAWFW	✓								✓	
WMRNR3	✓								✓	
WMUEZ7	✓	✓							✓	
WN87Q2	✓				✓				✓	
WQMPJJ	✓								✓	
WUTX32	✓	✓			✓				✓	
WWWWB9	✓								✓	
WZV889	✓	✓							✓	
WZXQ6J	✓				✓				✓	
X26JHH	✓			✓					✓	
X4XVGB	✓								✓	
X83WK3	✓			✓					✓	
X9RUPF	✓								✓	
XBHVWN	✓	✓							✓	Gas Chromatography retention time verification in comparison with a standard
XGLGMJ	✓			✓					✓	
XLFWCD	✓								✓	
XUMTVB		✓							✓	
Y3F9C3	✓	✓				✓			✓	
Y7CKP7	✓	✓							✓	
Y7FJ82	✓	✓							✓	
Y7KBN7	✓	✓		✓					✓	
Y7ZQGD	✓	✓		✓					✓	
Y8MKQ8	✓								✓	
YFEV3M		✓							✓	
YGBBBD	✓								✓	
YK9P8B	✓	✓							✓	
YMKQM4					✓				✓	
YWQNDT		✓					✓		✓	
YXGBP9	✓								✓	
Z2CPNH	✓	✓							✓	
Z2CTUE	✓								✓	

TABLE 1b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
Z2WCB8	✓								✓	
Z2WGMV		✓							✓	
Z3BKGM	✓					✓				Weight Measurement/ DART-TOF MS/GC-FID-MS
Z77LBN	✓								✓	
Z7BRZY	✓								✓	
Z7J8RC	✓								✓	
Z86Z7G	✓	✓							✓	
ZACKPB	✓								✓	
ZB8J9V		✓							✓	
ZBDRE6	✓			✓					✓	
ZEKFDX	✓								✓	
ZEM7M2	✓								✓	
ZF3UY9	✓	✓							✓	
ZKP9FF								✓	✓	
ZPHD8Y									✓	
ZQT9VC	✓	✓		✓					✓	
ZUMNX2	✓			✓					✓	
ZW2FE2									✓	LC/UV/MS and Weight Determination
ZWK27Y					✓				✓	
ZXRQ39	✓								✓	
ZYGAKA	✓								✓	
ZYNBY9	✓								✓	
ZYX2FY	✓			✓					✓	
ZZJT2A		✓			✓				✓	

Response Summary										Total Participants: 556
	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	
Participants	381	214	3	54	103	8	8	31	547	
Percent	69%	38%	0.5%	10%	19%	1%	1%	6%	98%	

Drug Analysis (Item 2)

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
238KUR	No controlled substance detected/identified	2UUM9T	No controlled substance detected/identified
23CH6B	No controlled substance detected/identified	2WYCB3	No controlled substance detected/identified
24FTY6	No controlled substance detected/identified	2Y2BB9	No controlled substance detected/identified
24JECD	No controlled substance detected/identified	2YAH4D	No controlled substance detected/identified
26NY8Z	No controlled substance detected/identified	2YPE76	No controlled substance detected/identified
28HPDE	No controlled substance detected/identified	2Z3DQL	No controlled substance detected/identified
28VEJC	No controlled substance detected/identified	32RDBV	No controlled substance detected/identified
28XYBL	No controlled substance detected/identified	36W96L	No controlled substance detected/identified
29CDWA	No controlled substance detected/identified	37PXCW	No controlled substance detected/identified
2BUPEC	No controlled substance detected/identified	37RMLX	No controlled substance detected/identified
2BXA7L	No controlled substance detected/identified	38JQ8D	No controlled substance detected/identified
2EGD74	No controlled substance detected/identified	39THV8	No controlled substance detected/identified
2JCWHB	No controlled substance detected/identified	3BT829	No controlled substance detected/identified
2KNPGJ	No controlled substance detected/identified	3G66TU	No controlled substance detected/identified
2NCGPN	No controlled substance detected/identified	3HYM3H	No controlled substance detected/identified
2NNXME	No controlled substance detected/identified	3KZPU6	No controlled substance detected/identified
2RA3GQ	No controlled substance detected/identified	3LEETB	No controlled substance detected/identified
2TG7ME	No controlled substance detected/identified	3MCX3D	No controlled substance detected/identified
2TZRE6	No controlled substance detected/identified	3MN9PB	No controlled substance detected/identified
2UJNB7	No controlled substance detected/identified	3NNBP2	No controlled substance detected/identified

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
3Q7JLF	No controlled substance detected/identified	4VKCKL	No controlled substance detected/identified
3Q8GWB	No controlled substance detected/identified	4XN6EJ	No controlled substance detected/identified
3TRBUC	No controlled substance detected/identified	4YWBPY	No controlled substance detected/identified
3TXDRW	Testosterone	4YXCFR	No controlled substance detected/identified
3WWM3Y	No controlled substance detected/identified	64W967	No controlled substance detected/identified
3X8MRD	No controlled substance detected/identified	67DVE6	No controlled substance detected/identified
3ZF27Y	No controlled substance detected/identified	67TLHA	No controlled substance detected/identified
3ZGK8G	No controlled substance detected/identified	6A2A3E	No controlled substance detected/identified
43URZ9	No controlled substance detected/identified	6DV99W	No controlled substance detected/identified
462AY4	No controlled substance detected/identified	6DVL2C	No controlled substance detected/identified
47XKYD	No controlled substance detected/identified	6FHHA3	No controlled substance detected/identified
48BQFV	No controlled substance detected/identified	6MKAWA	No controlled substance detected/identified
4FALQT	No controlled substance detected/identified	6QG3L4	No controlled substance detected/identified; Caffeine
4GGBLP	No controlled substance detected/identified	6QKH6V	No controlled substance detected/identified
4GJZPL	No controlled substance detected/identified	6RF27W	No controlled substance detected/identified
4M84WC	No controlled substance detected/identified	6VJEP4	No controlled substance detected/identified
4NGPBW	No controlled substance detected/identified	6VJU62	No controlled substance detected/identified
4PUHKL	No controlled substance detected/identified	72DKHN	No controlled substance detected/identified
4QQC7Q	No controlled substance detected/identified	72VJQF	No controlled substance detected/identified
4RHQEZ	No controlled substance detected/identified	73MH28	No controlled substance detected/identified
4TBEXV	No controlled substance detected/identified	77TC8C	No controlled substance detected/identified

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
77TWYX	No controlled substance detected/identified	8D87F7	No controlled substance detected/identified
78QZT2	No controlled substance detected/identified	8DMJ7N	No controlled substance detected/identified
79T9A6	No controlled substance detected/identified	8L3VNV	No controlled substance detected/identified
7CQCP2	No controlled substance detected/identified	8MDWAU	No controlled substance detected/identified
7CVW49	No controlled substance detected/identified	8N8KRQ	No controlled substance detected/identified
7HBD8Y	No controlled substance detected/identified	8PZ7A3	No controlled substance detected/identified
7HUF EZ	No controlled substance detected/identified	8QDE4F	No controlled substance detected/identified
7KAQBN	No controlled substance detected/identified	8TGY3A	No controlled substance detected/identified
7MJRFV	No controlled substance detected/identified	8UDEML	No controlled substance detected/identified
7ND8KK	No controlled substance detected/identified	8UNJM3	No controlled substance detected/identified
7TCX87	No controlled substance detected/identified	8XQDPE	No controlled substance detected/identified
7W63LG	No controlled substance detected/identified	8YMTX6	No controlled substance detected/identified
7YDCT8	No controlled substance detected/identified	927Z2Y	No controlled substance detected/identified
7Z4M47	No controlled substance detected/identified	94CME3	No controlled substance detected/identified
82VC88	No controlled substance detected/identified	967BWX	No controlled substance detected/identified
83T3GM	No controlled substance detected/identified	97LDYZ	No controlled substance detected/identified
8443GN	No controlled substance detected/identified	98KPRM	No controlled substance detected/identified
84NF9L	No controlled substance detected/identified	9996WN	No controlled substance detected/identified
84ZJBC	No controlled substance detected/identified	99TYGY	No controlled substance detected/identified
87WWDY	No controlled substance detected/identified	9EYF98	No controlled substance detected/identified
882FBF	No controlled substance detected/identified	9EZARD	No controlled substance detected/identified

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
9GRUW6	No controlled substance detected/identified	A9RR9Y	No controlled substance detected/identified
9H39AX	No controlled substance detected/identified	ADF97M	No controlled substance detected/identified
9HCA2Y	No controlled substance detected/identified	ADQJH4	No controlled substance detected/identified
9HDAQQ	No controlled substance detected/identified	AF3LWT	No controlled substance detected/identified
9J9TTR	No controlled substance detected/identified	AFNNV2	No controlled substance detected/identified
9JBH2V	No controlled substance detected/identified	AFP33M	No controlled substance detected/identified
9LGUT6	No controlled substance detected/identified; Caffeine	ANMNEG	No controlled substance detected/identified
9LTRXA	No controlled substance detected/identified	APEK2P	No controlled substance detected/identified
9MJE6E	No controlled substance detected/identified	AT82HN	No controlled substance detected/identified
9MP33P	No controlled substance detected/identified	AU8T2P	No controlled substance detected/identified
9NKL4Q	No controlled substance detected/identified	AX8UCK	No controlled substance detected/identified
9PJQKK	No controlled substance detected/identified	AZEC9V	No controlled substance detected/identified
9QRJAZ	No controlled substance detected/identified	AZED83	No controlled substance detected/identified
9RRX6T	No controlled substance detected/identified	B3C9XN	No controlled substance detected/identified
9YFYVB	No controlled substance detected/identified	B6KVZF	No controlled substance detected/identified
9ZCCH3	No controlled substance detected/identified	B84J8L	No controlled substance detected/identified
A29ZQQ	No controlled substance detected/identified	B9K2G8	No controlled substance detected/identified
A2QEAG	No controlled substance detected/identified	B9KL88	No controlled substance detected/identified
A7HCC6	No controlled substance detected/identified	BEQ9RW	No controlled substance detected/identified
A8WQ8U	No controlled substance detected/identified	BGFGBX	No controlled substance detected/identified
A93AQN	No controlled substance detected/identified	BLPQBE	No controlled substance detected/identified

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
BMGUXU	No controlled substance detected/identified	CZ4CAY	No controlled substance detected/identified
BPWUPQ	No controlled substance detected/identified	D2AG2R	No controlled substance detected/identified
BRXQLR	No controlled substance detected/identified	D386KV	No controlled substance detected/identified
BYABMV	No controlled substance detected/identified	D4K8GN	No controlled substance detected/identified
BYU9H9	No controlled substance detected/identified	D4MWKK	No controlled substance detected/identified
BYXFFV	No controlled substance detected/identified	D73FD3	No controlled substance detected/identified
BZ7NC9	No controlled substance detected/identified	DAMAA3	No controlled substance detected/identified
C2ZV4D	No controlled substance detected/identified	DAVM6Y	No controlled substance detected/identified
C3CGUE	No controlled substance detected/identified	DCQA7T	No controlled substance detected/identified
CAE8B7	No controlled substance detected/identified	DCQZH3	No controlled substance detected/identified
CEKBWQ	No controlled substance detected/identified	DDQJNX	No controlled substance detected/identified
CF2EZJ	No controlled substance detected/identified	DDQY2B	No controlled substance detected/identified
CGTXWX	No controlled substance detected/identified	DDY7RX	No controlled substance detected/identified
CHXEHM	No controlled substance detected/identified	DJTCEZ	No controlled substance detected/identified
CJX3CJ	No controlled substance detected/identified	DK6R92	No controlled substance detected/identified
CQB9XP	No controlled substance detected/identified	DKFBHT	No controlled substance detected/identified
CQUZNZ	No controlled substance detected/identified	DKNVG2	No controlled substance detected/identified
CTFBHV	No controlled substance detected/identified	DNNPQP	No controlled substance detected/identified
CTZJMW	No controlled substance detected/identified	DPM2LY	No controlled substance detected/identified
CUXTV7	No controlled substance detected/identified	DPZMXF	No controlled substance detected/identified
CWL8KF	No controlled substance detected/identified	DR6DRN	No controlled substance detected/identified

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
DTJAN2	No controlled substance detected/identified	EL94Y2	No controlled substance detected/identified
DTZPV4	No controlled substance detected/identified	ELTPRR	No controlled substance detected/identified
DU8DQZ	No controlled substance detected/identified	ERYHLM	No controlled substance detected/identified
DUA2AY	No controlled substance detected/identified	ETYEFQ	No controlled substance detected/identified
DZ277A	No controlled substance detected/identified	EW22JG	No controlled substance detected/identified
DZ2YB9	No controlled substance detected/identified	EWEQDT	No controlled substance detected/identified
DZYLBU	No controlled substance detected/identified	EYX4UA	No controlled substance detected/identified
E2N7GP	No controlled substance detected/identified	F4D8N7	No controlled substance detected/identified
E4KUQU	No controlled substance detected/identified	F6LYGZ	No controlled substance detected/identified
E4WEPA	No controlled substance detected/identified	F99D7A	No controlled substance detected/identified
E6NHCP	No controlled substance detected/identified	FA3KPU	No controlled substance detected/identified
E6UCUC	No controlled substance detected/identified	FBUHHD	No controlled substance detected/identified
E74F6F	No controlled substance detected/identified	FJ8R6U	No controlled substance detected/identified
EC8EL7	No controlled substance detected/identified	FJZFWG	No controlled substance detected/identified
ECF9UM	No controlled substance detected/identified	FLV7U4	No controlled substance detected/identified
EDPEVN	No controlled substance detected/identified	FLZ2HN	No controlled substance detected/identified
EEJ3EJ	No controlled substance detected/identified	FM6UV8	No controlled substance detected/identified
EGL6V8	No controlled substance detected/identified	FPX2ER	No controlled substance detected/identified
EHJCJL	No controlled substance detected/identified	FREV9Y	No controlled substance detected/identified
EJRG3	No controlled substance detected/identified	FRKFNE	No controlled substance detected/identified
EK2L6C	No controlled substance detected/identified	FUP9XZ	No controlled substance detected/identified

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
FV3636	No controlled substance detected/identified	HAAD3C	No controlled substance detected/identified
FWDMLB	No controlled substance detected/identified	HD3MQQ	No controlled substance detected/identified
FWRV6Q	No controlled substance detected/identified	HEX7TR	No controlled substance detected/identified
FXFDZW	Caffeine	HL9EV4	No controlled substance detected/identified
FYM6PT	No controlled substance detected/identified	HNCERU	No controlled substance detected/identified
GCZ9YT	No controlled substance detected/identified	HQGEWD	No controlled substance detected/identified
GEUTA8	No controlled substance detected/identified	HQUZCG	No controlled substance detected/identified
GGFR9M	No controlled substance detected/identified	HVJNLT	No controlled substance detected/identified
GML798	No controlled substance detected/identified	HW8P7C	No controlled substance detected/identified
GQG739	No controlled substance detected/identified	HYLRBU	No controlled substance detected/identified
GRRTUV	Caffeine	HYVLR4	No controlled substance detected/identified
GVTYBL	No controlled substance detected/identified	J42AGH	No controlled substance detected/identified
GXDVRY	No controlled substance detected/identified	J4LYTC	No controlled substance detected/identified
GXHAFL	No controlled substance detected/identified	J6RK7Q	No controlled substance detected/identified
GYERDG	No controlled substance detected/identified	J8ZLGP	No controlled substance detected/identified
GYQLHV	No controlled substance detected/identified	J982NV	No controlled substance detected/identified
GZCXXG	No controlled substance detected/identified	JA8R6B	No controlled substance detected/identified
H2YWPK	No controlled substance detected/identified	JC3WX7	No controlled substance detected/identified
H4N8L6	No controlled substance detected/identified	JCBYUG	No controlled substance detected/identified
H6BXGH	No controlled substance detected/identified	JDMELW	No controlled substance detected/identified
H9RAAR	No controlled substance detected/identified	JE2C7T	No controlled substance detected/identified

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
JE7PTP	No controlled substance detected/identified	KFALCC	No controlled substance detected/identified
JEKPEJ	No controlled substance detected/identified	KLF89L	No controlled substance detected/identified
JFEFHC	No controlled substance detected/identified	KLJZJV	No controlled substance detected/identified
JGCM8J	No controlled substance detected/identified	KMAPJV	No controlled substance detected/identified
JJ8PJX	No controlled substance detected/identified	KMLQXG	No controlled substance detected/identified
JJZLFV	No controlled substance detected/identified	KNLMRK	No controlled substance detected/identified
JKNLCX	No controlled substance detected/identified	KVQAUM	No controlled substance detected/identified
JQ8TPW	No controlled substance detected/identified	KXB4QM	No controlled substance detected/identified
JQP67F	No controlled substance detected/identified	KYRGAK	No controlled substance detected/identified
JTHVF2	No controlled substance detected/identified	L3HZKM	No controlled substance detected/identified
JUCJUQ	No controlled substance detected/identified	L4W2R3	No controlled substance detected/identified
JUP7C6	No controlled substance detected/identified	L6ZFTU	No controlled substance detected/identified
JWAGG	No controlled substance detected/identified	L76LMB	No controlled substance detected/identified
JXPKVN	No controlled substance detected/identified	L7BGKU	No controlled substance detected/identified
K3RPF	No controlled substance detected/identified	L7DYTQ	No controlled substance detected/identified
KC438N	No controlled substance detected/identified	L7XJTW	No controlled substance detected/identified
KCPJK7	No controlled substance detected/identified	L9488L	No controlled substance detected/identified
KCR9P3	No controlled substance detected/identified	L98TKT	No controlled substance detected/identified
KDJBCH	No controlled substance detected/identified	LAJLGG	No controlled substance detected/identified
KF9G6F	No controlled substance detected/identified	LBD4Q6	No controlled substance detected/identified
KF9JQD	No controlled substance detected/identified	LBU9A4	No controlled substance detected/identified

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
LJN42E	No controlled substance detected/identified	NB822W	No controlled substance detected/identified
LLTXXR	No controlled substance detected/identified	ND74JX	No controlled substance detected/identified
LM6V78	No controlled substance detected/identified	NFHR6H	No controlled substance detected/identified
M244JF	No controlled substance detected/identified	NFX3XV	No controlled substance detected/identified
M2UYKL	No controlled substance detected/identified	NGUCC8	No controlled substance detected/identified
M49J7L	No controlled substance detected/identified	NGXCCU	No controlled substance detected/identified
M786MQ	No controlled substance detected/identified	NL87AJ	Testosterone
M8T8LY	No controlled substance detected/identified	NQEUDT	No controlled substance detected/identified
MBAP4W	No controlled substance detected/identified	NRNYVM	No controlled substance detected/identified
ME4WYQ	No controlled substance detected/identified	NRQJMW	No controlled substance detected/identified
MHCYCA	No controlled substance detected/identified	NRYFJH	No controlled substance detected/identified
MRCMXR	No controlled substance detected/identified	NT76WF	No controlled substance detected/identified
MU88ZJ	No controlled substance detected/identified	NUHWU3	No controlled substance detected/identified
MULY73	No controlled substance detected/identified	NUTKMU	No controlled substance detected/identified
MVXBMY	No controlled substance detected/identified	NYFFV2	No controlled substance detected/identified
MW6T8N	No controlled substance detected/identified	NZ8AVQ	No controlled substance detected/identified
MW9XYQ	No controlled substance detected/identified	P2DQCE	No controlled substance detected/identified
MWRKVM	No controlled substance detected/identified	P46UXU	No controlled substance detected/identified
N7FX9R	No controlled substance detected/identified	P7UZWY	No controlled substance detected/identified
N7HZ2Q	No controlled substance detected/identified	P8KPXK	No controlled substance detected/identified
NA6ZPP	No controlled substance detected/identified	P8VECA	No controlled substance detected/identified

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
PAR3FL	Caffeine	QFDWWW	No controlled substance detected/identified
PCVDG9	No controlled substance detected/identified	QG3FKF	No controlled substance detected/identified
PEYC6L	No controlled substance detected/identified	QGPL4F	No controlled substance detected/identified
PFR2HN	No controlled substance detected/identified	QGZMRA	No controlled substance detected/identified
PPW9XN	No controlled substance detected/identified	QHXYLG	No controlled substance detected/identified
PR66XJ	No controlled substance detected/identified	QMTKMR	No controlled substance detected/identified
PT3NQD	No controlled substance detected/identified	QRFEFJ	No controlled substance detected/identified
PUCG2E	No controlled substance detected/identified	QUABWQ	No controlled substance detected/identified
PUTU4H	No controlled substance detected/identified	QURDNY	No controlled substance detected/identified
PUURDT	No controlled substance detected/identified	QVHJ24	No controlled substance detected/identified
PVP9RB	No controlled substance detected/identified	QWFRC9	No controlled substance detected/identified
PZDTGK	No controlled substance detected/identified	QWWPLZ	No controlled substance detected/identified
PZY7J4	No controlled substance detected/identified	QXW3U	No controlled substance detected/identified
Q27ZLG	No controlled substance detected/identified	QYHKX8	No controlled substance detected/identified
Q2N2CQ	No controlled substance detected/identified	QYL6PG	No controlled substance detected/identified
Q4H4NT	No controlled substance detected/identified	R2Q4NF	No controlled substance detected/identified
Q6NETN	No controlled substance detected/identified	R6K3X6	No controlled substance detected/identified
Q9TTPU	No controlled substance detected/identified	R6PKW4	No controlled substance detected/identified
Q9WBW6	No controlled substance detected/identified	R7KYR9	No controlled substance detected/identified
QAFLEH	No controlled substance detected/identified	R9QGQ6	No controlled substance detected/identified
QAL6NM	No controlled substance detected/identified	RB3Y7V	No controlled substance detected/identified

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
RBTHEF	No controlled substance detected/identified	TFRWAD	No controlled substance detected/identified
REY8BH	No controlled substance detected/identified	TGV2YQ	No controlled substance detected/identified
RG9VPV	No controlled substance detected/identified	TH42CK	No controlled substance detected/identified
RJNJRD	No controlled substance detected/identified	THAEP9	No controlled substance detected/identified
RKZ26Q	No controlled substance detected/identified	THUJJN	No controlled substance detected/identified
RLFTBA	No controlled substance detected/identified	TJR87D	No controlled substance detected/identified
RN3MUV	No controlled substance detected/identified	TLAX69	No controlled substance detected/identified
RP7NQF	No controlled substance detected/identified	TLPVYD	No controlled substance detected/identified
RPWBDR	No controlled substance detected/identified	TM3DA7	No controlled substance detected/identified
RR3CEW	No controlled substance detected/identified	TRF4PY	No controlled substance detected/identified
RTBAMR	No controlled substance detected/identified	TRW98W	No controlled substance detected/identified
RUJQCB	No controlled substance detected/identified	TTADVG	No controlled substance detected/identified
RWFCUP	No controlled substance detected/identified	TUKPG6	No controlled substance detected/identified
RXM93R	No controlled substance detected/identified	TYENUE	No controlled substance detected/identified
T3QQLC	No controlled substance detected/identified	U7PGCN	No controlled substance detected/identified
T8ACKY	No controlled substance detected/identified	UAQD44	No controlled substance detected/identified
T8C2PV	No controlled substance detected/identified	UB636V	No controlled substance detected/identified
T9BH88	No controlled substance detected/identified	UCYQMQ	No controlled substance detected/identified
T9CRHG	No controlled substance detected/identified	UD8MBG	No controlled substance detected/identified
T9Q4XP	No controlled substance detected/identified	UFEC6C	No controlled substance detected/identified
TDRBN6	No controlled substance detected/identified	UFHBM9	No controlled substance detected/identified

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
UFV76M	No controlled substance detected/identified	VE4YQN	No controlled substance detected/identified
UFX2JQ	No controlled substance detected/identified	VEFG4L	No controlled substance detected/identified
UJ973D	No controlled substance detected/identified	VEJWKT	No controlled substance detected/identified
UMBWKM	No controlled substance detected/identified	VEKQJT	No controlled substance detected/identified
UMJZBK	No controlled substance detected/identified	VGLRD4	No controlled substance detected/identified
UQLTDW	No controlled substance detected/identified	VMFP2Z	No controlled substance detected/identified
UQTF6G	No controlled substance detected/identified	VU9UXG	caffeine
UVC8UQ	No controlled substance detected/identified	VUCEQQ	No controlled substance detected/identified
UX43YE	No controlled substance detected/identified	VURVKJ	No controlled substance detected/identified
UYZGL8	No controlled substance detected/identified	VXV6YF	No controlled substance detected/identified
UZ2CH8	No controlled substance detected/identified	W6AR83	No controlled substance detected/identified
UZ8KTF	No controlled substance detected/identified	W74W3H	No controlled substance detected/identified
V2HJNP	No controlled substance detected/identified	W7KUVM	No controlled substance detected/identified
V3HK2K	No controlled substance detected/identified	W8F9DK	No controlled substance detected/identified
V3UWB7	No controlled substance detected/identified	W8FAA6	No controlled substance detected/identified
V9CF6D	No controlled substance detected/identified	WBEKAK	No controlled substance detected/identified
VAJ2M	No controlled substance detected/identified	WCBW8Y	No controlled substance detected/identified; caffeine
VBN6UB	No controlled substance detected/identified	WG2G33	No controlled substance detected/identified
VCJMVD	No controlled substance detected/identified	WJ8VJF	No controlled substance detected/identified
VCQXF3	No controlled substance detected/identified	WLWVPB	No controlled substance detected/identified
VE3JGP	No controlled substance detected/identified	WMAFWF	No controlled substance detected/identified

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
WMRNR3	No controlled substance detected/identified	Y8MKQ8	No controlled substance detected/identified
WMUEZ7	No controlled substance detected/identified	YFEV3M	No controlled substance detected/identified
WN87Q2	No controlled substance detected/identified	YGBBBD	No controlled substance detected/identified
WQMPJJ	No controlled substance detected/identified	YK9P8B	No controlled substance detected/identified; Caffeine
WUTX32	No controlled substance detected/identified	YMKQM4	No controlled substance detected/identified
WWWWB9	No controlled substance detected/identified	YWQNDT	No controlled substance detected/identified
WZV889	No controlled substance detected/identified	YXGBP9	No controlled substance detected/identified
WZXQ6J	No controlled substance detected/identified	Z2CPNH	No controlled substance detected/identified
X26JHH	No controlled substance detected/identified	Z2CTUE	No controlled substance detected/identified
X4XVGB	No controlled substance detected/identified	Z2WCB8	No controlled substance detected/identified
X83WK3	No controlled substance detected/identified	Z2WGMV	No controlled substance detected/identified
X9RUPF	No controlled substance detected/identified	Z3BKGM	No controlled substance detected/identified
XBHVWN	No controlled substance detected/identified	Z77LBN	No controlled substance detected/identified
XGLGMJ	No controlled substance detected/identified	Z7BRZY	No controlled substance detected/identified
XLFWCD	No controlled substance detected/identified	Z7J8RC	No controlled substance detected/identified
XUMTVB	No controlled substance detected/identified	Z86Z7G	No controlled substance detected/identified
Y3F9C3	No controlled substance detected/identified	ZACKPB	No controlled substance detected/identified
Y7CKP7	No controlled substance detected/identified	ZB8J9V	No controlled substance detected/identified
Y7FJ82	No controlled substance detected/identified	ZBDRE6	No controlled substance detected/identified
Y7KBN7	No controlled substance detected/identified	ZEKFDX	No controlled substance detected/identified
Y7ZQGD	No controlled substance detected/identified	ZEM7M2	No controlled substance detected/identified

TABLE 2a - Item 2 Results

WebCode	Controlled Substance Detected	WebCode	Controlled Substance Detected
ZF3UY9	No controlled substance detected/identified		
ZKP9FF	No controlled substance detected/identified		
ZPHD8Y	No controlled substance detected/identified		
ZQT9VC	No controlled substance detected/identified		
ZUMNX2	No controlled substance detected/identified		
ZW2FE2	No controlled substance detected/identified		
ZWK27Y	No controlled substance detected/identified		
ZXRQ39	No controlled substance detected/identified		
ZYGAKA	No controlled substance detected/identified		
ZYNBY9	No controlled substance detected/identified		
ZYX2FY	No controlled substance detected/identified		
ZZJT2A	No controlled substance detected/identified		

<i>Response Summary - Item 2</i>	
Participants: 556	
No Controlled substance(s) detected/identified:	554
Other:	2

Examination Methods (Item 2)

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
238KUR		✓							✓	
23CH6B									✓	TruNarc
24FTY6	✓								✓	
24JECD		✓							✓	
26NY8Z		✓					✓		✓	
28HPDE	✓	✓		✓					✓	
28VEJC	✓	✓							✓	
28XYBL		✓							✓	
29CDWA	✓	✓							✓	
2BUPEC			✓	✓				✓	✓	
2BXA7L		✓						✓	✓	
2EGD74	✓								✓	
2JCWHB	✓								✓	
2KNPGJ	✓								✓	
2NCGPN	✓								✓	
2NNXME	✓								✓	
2RA3GQ		✓					✓		✓	
2TG7ME	✓								✓	
2TZRE6	✓								✓	
2UJNB7	✓	✓							✓	
2UUM9T	✓							✓	✓	
2WYCB3	✓	✓							✓	
2Y2BB9	✓	✓							✓	
2YAH4D	✓	✓							✓	
2YPE76	✓	✓		✓					✓	
2Z3DQL	✓								✓	
32RDBV	✓								✓	
36W96L									✓	
37PXCUI	✓								✓	
37RMLX	✓	✓			✓	✓			✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
38JQ8D	✓								✓	
39THV8		✓							✓	DART Agilent-TOF
3BT829	✓							✓	✓	
3G66TU	✓				✓				✓	
3HYM3H	✓	✓							✓	
3KZPU6	✓				✓				✓	
3LEETB	✓								✓	FTIR with ATR accessory
3MCX3D	✓								✓	
3MN9PB	✓								✓	
3NNBP2	✓	✓							✓	
3Q7JLF	✓								✓	
3Q8GWB	✓								✓	
3TRBUC		✓							✓	
3TXDRW	✓				✓				✓	
3WWM3Y									✓	
3X8MRD	✓	✓		✓					✓	
3ZF27Y		✓							✓	
3ZGK8G	✓								✓	
43URZ9									✓	Two different GC/MS used
462AY4	✓								✓	
47XKYD		✓							✓	
48BQFV					✓				✓	
4FALQT					✓				✓	
4GGBLP		✓							✓	
4GJZPL	✓	✓			✓				✓	
4M84WC	✓						✓		✓	
4NGPBW									✓	
4PUHKL		✓							✓	XRD
4QQC7Q		✓							✓	
4RHQEZ	✓	✓		✓					✓	
4TBEXV	✓				✓	✓			✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
4VKCKL	✓	✓		✓					✓	
4XN6EJ	✓	✓							✓	
4YWBPY	✓								✓	
4YXCFR					✓				✓	
64W967									✓	
67DVE6	✓								✓	
67TLHA	✓	✓							✓	
6A2A3E		✓							✓	XRD
6DV99W					✓				✓	
6DVL2C	✓	✓							✓	
6FHHA3	✓	✓							✓	
6MKawe	✓								✓	
6QG3L4	✓	✓		✓		✓			✓	
6QKH6V									✓	
6RF27W	✓	✓							✓	
6VJEP4	✓								✓	
6VJU62	✓								✓	
72DKHN					✓				✓	
72VJQF		✓						✓	✓	
73MH28	✓								✓	
77TC8C									✓	
77TWYX									✓	
78QZT2	✓			✓					✓	
79T9A6		✓							✓	
7CQCP2	✓			✓					✓	
7CVW49	✓	✓							✓	
7HBD8Y	✓	✓							✓	
7HUFEZ	✓	✓							✓	
7KAQBN		✓							✓	
7MJRFV	✓	✓							✓	
7ND8KK		✓							✓	XRD
7TCX87	✓	✓							✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
7W63LG	✓	✓		✓					✓	
7YDCT8									✓	
7Z4M47	✓			✓					✓	
82VC88	✓								✓	
83T3GM	✓								✓	
8443GN	✓								✓	
84NF9L		✓				✓			✓	
84ZJBC	✓									DART-TOF
87WWDY	✓	✓		✓					✓	
882FBF		✓						✓	✓	
8D87F7					✓				✓	
8DMJ7N	✓								✓	
8L3VNV								✓	✓	
8MDWAU	✓	✓							✓	
8N8KRQ	✓				✓				✓	
8PZ7A3		✓							✓	DART-Agilent TOF
8QDE4F					✓				✓	FID
8TGY3A		✓							✓	
8UDEML	✓								✓	
8UNJM3	✓	✓							✓	
8XQDPE		✓						✓	✓	
8YMTX6	✓	✓							✓	
927Z2Y									✓	
94CME3	✓								✓	
967BWX	✓				✓				✓	
97LDYZ	✓								✓	
98KPRM	✓								✓	
9996WN	✓								✓	
99TYGY	✓								✓	
9EYF98	✓								✓	
9EZARD	✓								✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
9GRUW6	✓								✓	
9H39AX	✓				✓	✓			✓	
9HCA2Y	✓								✓	
9HDAQQ	✓	✓							✓	
9J9TTR									✓	
9JBH2V	✓								✓	
9LGUT6	✓			✓					✓	
9LTRXA	✓								✓	
9MJE6E									✓	
9MP33P	✓	✓			✓				✓	
9NKL4Q	✓								✓	
9PJQKK		✓							✓	
9QRJAZ									✓	DART-TOF
9RRX6T	✓	✓		✓					✓	
9YFYVB	✓	✓							✓	
9ZCCH3		✓							✓	
A29ZQQ	✓				✓	✓			✓	
A2QEAG		✓							✓	
A7HCC6								✓	✓	
A8WQ8U									✓	TruNarc
A93AQN	✓								✓	
A9RR9Y	✓	✓							✓	
ADF97M	✓				✓				✓	
ADQJH4									✓	
AF3LWT	✓								✓	
AFNNV2	✓								✓	
AFP33M					✓				✓	
ANMNEG					✓				✓	
APEK2P	✓				✓	✓			✓	
AT82HN	✓	✓							✓	
AU8T2P	✓								✓	
AX8UCK	✓				✓				✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
AZEC9V	✓								✓	
AZED83	✓			✓					✓	
B3C9XN	✓							✓	✓	
B6KVZF		✓							✓	
B84J8L	✓	✓						✓	✓	
B9K2G8									✓	
B9KL88	✓								✓	
BEQ9RW									✓	
BGFBGX	✓								✓	
BLPQBE	✓								✓	
BMGUXU	✓	✓							✓	
BPWUPQ	✓			✓					✓	
BRXQLR					✓				✓	
BYABMV	✓	✓							✓	
BYU9H9	✓								✓	
BYXFFV	✓	✓							✓	
BZ7NC9	✓								✓	
C2ZV4D					✓				✓	
C3CGUE		✓			✓				✓	
CAE8B7									✓	
CEKBWQ	✓								✓	
CF2EZJ	✓								✓	
CGTXWX									✓	DART-Agilent ToF
CHXEHM		✓							✓	
CJX3CJ	✓				✓					GC-IRD
CQB9XP		✓							✓	weight determination
CQUZNZ	✓							✓	✓	
CTFBHV	✓			✓					✓	
CTZJMW	✓								✓	
CUXTV7	✓								✓	
CWL8KF		✓							✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
CZ4CAY	✓	✓							✓	
D2AG2R		✓							✓	
D386KV									✓	DART Agilent TOF
D4K8GN	✓								✓	
D4MWKK	✓				✓	✓			✓	
D73FD3	✓	✓							✓	
DAMAA3		✓						✓	✓	
DAVM6Y									✓	TruNarc
DCQA7T	✓								✓	
DCQZH3	✓								✓	
DDQJNX	✓	✓					✓		✓	
DDQY2B	✓	✓		✓					✓	
DDY7RX	✓								✓	
DJTCEZ	✓								✓	FTIR with ATR accessory
DK6R92		✓							✓	
DKFBHT	✓	✓		✓					✓	
DKNVG2	✓								✓	
DNNPQP								✓	✓	
DPM2LY	✓	✓		✓					✓	
DPZMXF	✓								✓	
DR6DRN									✓	Weight
DTJAN2	✓	✓							✓	
DTZPV4	✓								✓	
DU8DQZ	✓								✓	
DUA2AY	✓								✓	net weight=0.275 (+/- 0.001) grams
DZ277A		✓							✓	
DZ2YB9	✓				✓				✓	FID, DART-MS
DZYLBU	✓								✓	
E2N7GP	✓	✓							✓	
E4KUQU	✓	✓							✓	
E4WEPA	✓								✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
E6NHCP	✓	✓			✓				✓	
E6UCUC					✓				✓	
E74F6F	✓				✓				✓	
EC8EL7	✓	✓							✓	
ECF9UM	✓	✓		✓					✓	
EDPEVN	✓	✓							✓	
EEJ3EJ	✓	✓			✓				✓	
EGL6V8	✓								✓	
EHJCJL	✓	✓		✓					✓	
EJRGP3		✓							✓	
EK2L6C	✓								✓	
EL94Y2		✓							✓	
ELTPRR	✓							✓	✓	
ERYHLM	✓								✓	
ETYEFQ	✓				✓				✓	
EW22JG	✓	✓							✓	
EWEQDT	✓								✓	
EYX4UA	✓								✓	
F4D8N7	✓	✓							✓	
F6LYGZ		✓							✓	
F99D7A	✓	✓			✓				✓	
FA3KPU									✓	
FBUHHD	✓								✓	
FJ8R6U									✓	
FJZFWG							✓		✓	weight determination
FLV7U4		✓							✓	Microscope, XRF, XRD
FLZ2HN	✓								✓	
FM6UV8	✓								✓	
FPX2ER	✓	✓							✓	
FREV9Y	✓								✓	FTIR-ATR
FRKFNE	✓	✓							✓	
FUP9XZ	✓						✓		✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
FV3636		✓							✓	
FWDMLB		✓							✓	Electronic Balance
FWRV6Q	✓	✓							✓	
FXFDZW							✓			
FYM6PT	✓								✓	
GCZ9YT	✓								✓	
GEUTA8	✓								✓	
GGFR9M	✓	✓							✓	
GML798	✓								✓	
GQG739		✓							✓	
GRRTUV		✓							✓	
GVTYBL	✓								✓	
GXDVRY			✓						✓	
GXHAFI	✓								✓	
GYERDG	✓								✓	
GYQLHV	✓							✓	✓	
GZCXXG	✓				✓					
H2YWPK		✓							✓	
H4N8L6	✓	✓							✓	
H6BXGH		✓							✓	
H9RAAR	✓	✓							✓	
HAAD3C	✓	✓							✓	
HD3MQQ	✓							✓	✓	
HEX7TR		✓			✓				✓	
HL9EV4		✓							✓	weight, visual examination
HNCERU	✓							✓	✓	
HQGEWD	✓								✓	
HQUZCG		✓							✓	
HVJNLT	✓								✓	
HW8P7C	✓								✓	
HYLRBU	✓								✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
HYVLR4	✓	✓							✓	
J42AGH	✓	✓							✓	Atmospheric Sample Analysis Probe-MS
J4LYTC					✓				✓	
J6RK7Q	✓								✓	
J8ZLGP	✓								✓	
J982NV								✓	✓	
JA8R6B	✓	✓							✓	
JC3WX7	✓	✓							✓	
JCBYUG									✓	Weight
JDMELW	✓								✓	
JE2C7T	✓							✓	✓	
JE7PTP	✓	✓							✓	
JEKPEJ		✓							✓	
JFEFHC	✓				✓				✓	
JGCM8J	✓				✓	✓			✓	
JJ8PJX	✓	✓						✓	✓	
JJZLFV	✓								✓	
JKNLCX	✓								✓	
JQ8TPW		✓		✓						
JQP67F		✓							✓	
JTHVF2	✓								✓	
JUCJUQ	✓	✓							✓	
JUP7C6	✓								✓	
JWAGG	✓					✓			✓	GC/FID
JXPKVN	✓								✓	
K3RPFF	✓	✓							✓	
KC438N	✓								✓	
KCPJK7		✓							✓	Polarized Light Microscopy
KCR9P3		✓							✓	
KDJBCH		✓							✓	
KF9G6F	✓				✓	✓			✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
KF9JQD	✓								✓	
KFALCC	✓				✓				✓	
KLF89L	✓				✓	✓			✓	
KLJZJV		✓							✓	
KMAPJV		✓			✓				✓	
KMLQXG	✓								✓	
KNLMRK	✓				✓				✓	
KVQAUM	✓	✓			✓	✓			✓	
KXB4QM	✓								✓	
KYRGAK	✓	✓		✓					✓	
L3HZKM	✓								✓	
L4W2R3	✓								✓	
L6ZFTU	✓	✓							✓	
L76LMB					✓				✓	
L7BGKU	✓	✓							✓	
L7DYTQ	✓								✓	
L7XJTW	✓								✓	
L9488L	✓								✓	
L98TKT	✓								✓	
LAJLGG		✓							✓	
LBD4Q6	✓								✓	
LBU9A4		✓							✓	
LJN42E	✓	✓							✓	
LLTXR	✓	✓							✓	
LM6V78	✓								✓	
M244JF	✓								✓	
M2UYKL	✓			✓					✓	
M49J7L	✓				✓				✓	
M786MQ									✓	
M8T8LY									✓	Gas Chromatography retention time verification in comparison with a standard

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
MBAP4W									✓	
ME4WYQ	✓								✓	
MHCYCA				✓					✓	
MRCMXR		✓	✓	✓	✓		✓		✓	H-NMR
MU88ZJ	✓	✓							✓	
MULY73	✓								✓	
MVXBMV	✓	✓							✓	
MW6T8N	✓								✓	
MW9XYQ	✓	✓							✓	
MWRKVM									✓	TruNarc
N7FX9R	✓	✓						✓	✓	
N7HZ2Q	✓								✓	
NA6ZPP								✓	✓	
NB822W		✓							✓	
ND74JX		✓						✓	✓	
NFHR6H	✓								✓	
NFX3XV	✓	✓							✓	
NGUCC8					✓				✓	
NGXCCU					✓				✓	
NL87AJ	✓			✓					✓	
NQEU DT									✓	
NRNYVM	✓								✓	
NRQJMW		✓							✓	
NRYFJH	✓	✓		✓					✓	
NT76WF					✓				✓	
NUHWU3		✓							✓	Balance, XRD
NUTKMU	✓	✓							✓	
NYFFV2	✓								✓	
NZ8AVQ	✓	✓							✓	
P2DQCE	✓				✓				✓	
P46UXU	✓	✓							✓	
P7UZWY									✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
P8KPXK							✓		✓	
P8VECA		✓							✓	
PAR3FL	✓								✓	
PCVDG9					✓				✓	
PEYC6L	✓	✓							✓	
PFR2HN	✓							✓	✓	
PPW9XN	✓								✓	
PR66XJ	✓								✓	
PT3NQD	✓			✓	✓				✓	GC-IR
PUCG2E		✓							✓	
PUTU4H	✓								✓	
PUURDT									✓	
PVP9RB									✓	
PZDTGK		✓							✓	Weight Determination, DART-TOFMS
PZY7J4		✓							✓	
Q27ZLG	✓								✓	
Q2N2CQ	✓								✓	
Q4H4NT	✓								✓	
Q6NETN								✓	✓	
Q9TTPU									✓	
Q9WBW6					✓				✓	
QAFLEH	✓								✓	
QAL6NM	✓								✓	
QFDWWW	✓	✓							✓	
QG3FKF	✓				✓	✓			✓	
QGPL4F	✓	✓							✓	
QGZMRA		✓							✓	
QHXYLG	✓								✓	
QMTKMR	✓								✓	
QRFEFJ	✓								✓	
QUABWQ									✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
QURDNY					✓				✓	
QVHJ24		✓			✓				✓	
QWFRC9	✓				✓	✓			✓	
QWWPLZ	✓	✓							✓	
QXVV3U	✓					✓				Weight Measurement, DART-TOF
QYHKX8	✓								✓	Weight determination
QYL6PG	✓								✓	
R2Q4NF	✓								✓	
R6K3X6	✓				✓				✓	
R6PKW4	✓								✓	
R7KYR9	✓					✓			✓	
R9QGQ6					✓				✓	
RB3Y7V									✓	
RBTHEF	✓								✓	
REY8BH									✓	TruNarc
RG9VPV									✓	
RJNJRD	✓								✓	
RKZ26Q	✓								✓	
RLFTBA	✓	✓							✓	
RN3MUV	✓								✓	
RP7NQF	✓	✓							✓	
RPWBDR	✓								✓	
RR3CEW									✓	
RTBAMR	✓	✓							✓	ATR
RUJQCB	✓								✓	
RWFCUP	✓								✓	
RXM93R	✓								✓	
T3QQLC	✓	✓		✓					✓	
T8ACKY		✓					✓		✓	
T8C2PV	✓								✓	
T9BH88	✓			✓					✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
T9CRHG	✓	✓							✓	
T9Q4XP									✓	
TDRBN6	✓								✓	
TFRWAD	✓	✓		✓					✓	
TGV2YQ	✓								✓	
TH42CK					✓					
THAEP9		✓			✓				✓	
THUJJN	✓	✓							✓	
TJR87D	✓								✓	
TLAX69	✓								✓	
TLPVYD									✓	
TM3DA7	✓								✓	
TRF4PY	✓	✓							✓	
TRW98W					✓				✓	
TTADVG		✓							✓	DART Agilent-TOF
TUKPG6	✓								✓	
TYENUE	✓							✓	✓	
U7PGCN	✓	✓							✓	
UAQD44	✓								✓	
UB636V	✓								✓	
UCYQMQ	✓	✓							✓	
UD8MBG	✓	✓							✓	
UFEC6C	✓								✓	
UFHBM9					✓				✓	
UFV76M	✓	✓						✓	✓	
UFX2JQ		✓							✓	
UJ973D	✓	✓							✓	
UMBWKM	✓	✓							✓	
UMJZBK	✓								✓	
UQLTDW					✓				✓	
UQTF6G	✓							✓	✓	
UVC8UQ	✓								✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
UX43YE	✓	✓							✓	
UYZGL8	✓			✓					✓	
UZ2CH8					✓				✓	
UZ8KTF		✓							✓	DART-AccuTOF
V2HJNP	✓	✓							✓	
V3HK2K	✓								✓	
V3UWB7	✓	✓		✓					✓	
V9CF6D		✓							✓	
VAJX2M	✓	✓							✓	
VBN6UB	✓	✓		✓					✓	
VCJMVD	✓	✓							✓	
VCQXF3									✓	weight determination
VE3JGP		✓							✓	
VE4YQN	✓								✓	
VEFG4L	✓	✓						✓	✓	
VEJWKT	✓	✓		✓					✓	
VEKQJT									✓	
VGLRD4	✓	✓							✓	
VMFP2Z	✓								✓	
VU9UXG	✓								✓	
VUCEQQ		✓						✓	✓	
VURVKJ					✓				✓	
VXV6YF									✓	DART Agilent TOF
W6AR83	✓								✓	
W74W3H	✓								✓	
W7KUVM	✓								✓	
W8F9DK		✓							✓	
W8FAA6	✓	✓							✓	
WBEKAK									✓	
WCBW8Y	✓								✓	
WG2G33	✓								✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
WJ8VJF									✓	
WLWVPB		✓							✓	Agilent DART - TOF
WMAWFW	✓								✓	
WMRNR3	✓								✓	
WMUEZ7	✓	✓							✓	
WN87Q2	✓				✓	✓			✓	
WQMPJJ	✓								✓	
WUTX32	✓	✓			✓				✓	
WWWWB9	✓								✓	
WZV889	✓	✓							✓	
WZXQ6J	✓								✓	
X26JHH	✓								✓	
X4XVGB	✓								✓	
X83WK3	✓			✓					✓	
X9RUPF	✓								✓	
XBHVWN	✓	✓							✓	
XGLGMJ	✓								✓	
XLFWCD	✓								✓	
XUMTVB		✓							✓	
Y3F9C3	✓	✓				✓			✓	
Y7CKP7	✓	✓							✓	
Y7FJ82	✓	✓							✓	
Y7KBN7	✓	✓		✓					✓	
Y7ZQGD	✓	✓		✓					✓	
Y8MKQ8	✓								✓	
YFEV3M		✓							✓	
YGBBBD	✓								✓	
YK9P8B	✓	✓							✓	
YMKQM4					✓				✓	
YWQNDT		✓					✓		✓	
YXGBP9	✓								✓	
Z2CPNH	✓	✓							✓	

TABLE 2b

WebCode	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	Other
Z2CTUE	✓								✓	
Z2WCB8	✓								✓	solubility
Z2WGMV		✓							✓	
Z3BKGM	✓					✓				Weight Measurement/ DART-TOF MS/ GC-FID-MS
Z77LBN	✓								✓	
Z7BRZY	✓								✓	
Z7J8RC	✓								✓	
Z86Z7G	✓	✓							✓	
ZACKPB	✓								✓	
ZB8J9V		✓							✓	
ZBDRE6	✓			✓					✓	
ZEKFDX	✓								✓	
ZEM7M2	✓								✓	
ZF3UY9	✓	✓							✓	
ZKP9FF								✓	✓	
ZPHD8Y									✓	
ZQT9VC	✓	✓		✓					✓	
ZUMNX2	✓			✓					✓	
ZW2FE2									✓	LC/UV/MS and Weight Determination
ZWK27Y						✓			✓	
ZXRQ39	✓								✓	
ZYGAKA	✓								✓	
ZYNBY9	✓								✓	
ZYX2FY	✓			✓					✓	
ZZJT2A		✓							✓	

Response Summary										Total Participants: 556
	Color	FTIR	LC	UV	GC	TLC	LC/MS	Raman	GC/MS	
Participants	380	216	3	44	73	20	11	32	548	
Percent	68%	39%	0.5%	8%	13%	4%	2%	6%	99%	

Additional Comments

TABLE 3

WebCode	Additional Comments
238KUR	Item 2 was found to contain caffeine.
28XYBL	Caffeine detected in Item 2
2BUPEC	Item 2 was found to contain caffeine.
2RA3GQ	Exhibit 2: Caffeine identified
2UJNB7	Caffeine was identified in Item 2.
2YPE76	Caffeine was found in Item 2 but no controlled substances were detected.
2Z3DQL	Identified caffeine in item 2.
3BT829	For item 2 caffeine was indicated. The instrumental results were compared to reference spectra only.
3KZPU6	Item 2 contained caffeine.
3WWM3Y	Item #2) Caffeine is indicated.
3ZGK8G	Large caffeine peak in item #2.
43URZ9	ITEM 2 CONTAINED CAFFEINE (NOT A CONTROLLED SUBSTANCE), PRESENCE OF CAFFEINE WAS CONFIRMED.
462AY4	Caffeine was confirmed in Item 2
4GGBLP	Item 2 was found to contain caffeine.
4GJZPL	Caffeine was indicated in the powder from item 2.
4QQC7Q	Item 2 contains caffeine.
6A2A3E	Testosterone (C-III) was identified in Exhibit 1. Scheduling information reported is based on US Code of Federal Regulations Title 21 Section 1308. No controlled or related substances identified in Exhibit 2.
6RF27W	Item 2: Caffeine was identified.
73MH28	Item 2 contains caffeine.
78QZT2	Caffeine was detected in Item 2.
7CQCP2	Caffeine detected in Item 2
7CVW49	Caffeine was identified in Item 2.
7HBD8Y	Item 2: Caffeine was identified.
7KAQBN	Caffeine was identified in item 2.
7MJRFV	Item 1: Testosterone (Schedule III). Item 2: Caffeine (Noncontrolled).
7TCX87	Item 2 Caffeine
7W63LG	Item 2 caffeine present

TABLE 3

WebCode	Additional Comments
7YDCT8	Testosterone was presumptively identified in item #1. This compound is not included in the lab's confirmatory scope of testing. Testing for this compound would typically be sent out for completion by an outside reference laboratory.
7Z4M47	Caffeine was reported as the substance found in Item 2.
84NF9L	01: 176.3 MG +/- 1.4 MG OF FINE WHITE MATERIAL IN ONE FOLDED GLASSINE PACKET IN ONE ZIP-TOP BAG. THE FINE WHITE MATERIAL CONTAINS TESTOSTERONE. METHODS OF ANALYSIS: GC-MS X2. 02: 281.1 MG +/- 1.4 MG OF FINE WHITE MATERIAL IN ONE FOLDED GLASSINE PACKET IN ONE ZIP-TOP BAG. NO CONTROLLED SUBSTANCES DETECTED IN THE FINE WHITE MATERIAL. METHODS OF ANALYSIS: FT-IR, TLC, AND GC-MS X2.
87WWDY	Caffeine was detected in Item 2.
8L3VNV	Division does not have a reference material for this substance.
8MDWAU	Caffeine identified in item 2.
8N8KRQ	Item 2: Acetone wash was done on a second sampling of powder.
8TGY3A	Caffeine identified in item 2
8UDEML	item 2 contained caffeine
967BWX	Item 2 contained caffeine
9LGUT6	Though not controlled, caffeine was reported as the substance found in Item 2.
9MP33P	caffeine indicated in item 2
9YFYVB	Item 2 contained Caffeine
ADQJH4	Testosterone was presumptively identified in item #1. This compound is not included in the lab's confirmatory scope of testing. Testing for this compound would typically be sent out for completion by an outside reference laboratory.
AZED83	Item 2 contains caffeine.
B3C9XN	caffeine was indicated in item 2
B84J8L	Caffeine was identified in item 2
B9KL88	Caffeine was confirmed in Item 2
BGFGBX	Caffeine identified in Item 2.
BMGUXU	Caffeine indicated in item 2.
BYU9H9	Item 2 contains caffeine.
BZ7NC9	Caffeine was identified in item 2.
C2ZV4D	Item 2 was caffeine.
C3CGUE	Item 2 contained caffeine
CQUZNZ	I found Item 2 to contain caffeine. Comment re: packaging - since we need to repackage the items as we would casework, it would be great if the ziplock bags were heat sealed at the top of the bag and not right above the folded glassine/paper bag. Heat sealing the folded bag in a small compartment makes it more difficult to open without cutting the interior bag and also to repackage.

TABLE 3

WebCode	Additional Comments
CWL8KF	Item 2 contained caffeine.
DAMAA3	Item 2-Caffeine was identified.
DCQA7T	Caffeine was detected and reported in Item 2
DDQJNX	There were no controlled substances present in Item 2 but Caffeine was detected in Item 2.
DK6R92	Item 1 - Suspected androstenedione likely present as precursor/by-product of testosterone synthesis. Not confirmed due to low level and incomplete spectrum, and not further pursued. Peak just before 8 minutes shown with library match of testosterone acetate - there are significant spectral differences between the spectrum and the library match. No other good matches found, not further pursued. Item 2 - Caffeine identified. No controlled substances identified.
DKFBHT	item 2 result caffeine, non controlled substance
DPM2LY	Caffeine identified for Item 2.
DPZMXF	caffeine identified in Item 2.
DTZPV4	Analysis of Item 2 disclosed caffeine.
E2N7GP	Caffeine was indicated but not identified.
ECF9UM	Item 2 Contains Caffeine
EDPEVN	Caffeine (noncontrolled) was identified in item 2.
EGL6V8	Item 2 was found to contain caffeine.
EHJCJL	Item #2 contains Caffeine
EJRGP3	Item 2 contained caffeine
ETYEFQ	Item 2 contained caffeine
FJZFWG	Caffeine was detected in item 2.
FLZ2HN	Caffeine was matched via GCMS spectrum in powdered substance 2 despite containing no controlled substances.
GCZ9YT	Item 2 : caffein in base form (96.3 %).
GEUTA8	Possible trace of methamphetamine was detected in sample #2, but quantity was insufficient for identification.
GGFR9M	Item 2 indicated Caffeine, but was not confirmed.
GXDVRV	Item 1 and 2 for Test No. 22-5001 were distributed to Staff [Initials] of Toxicology Laboratory and were analysed by GCMS and HPLC.
GYQLHV	Caffeine was detected in Item 2.
H2YWPK	caffeine was identified in Item 2
H6BXGH	Item 2 is confirmed to contain caffeine.
HD3MQQ	Item 2 contains the non-controlled substance caffeine.
HL9EV4	Identification of caffeine was reported for item 2.

TABLE 3

WebCode	Additional Comments
HQUZCG	Caffeine was identified in item #2.
HVJNLT	Caffeine was identified in item 2.
J6RK7Q	Caffeine was confirmed in Item 2.
JA8R6B	Caffeine was identified in item 2.
JCBYUG	Caffeine was detected in Lab Item 2.
JQP67F	Item 2 contains caffeine.
JTHVF2	Item 2 contains caffeine
K3RPF	Confirmed Caffeine in Item 2.
KDJBCH	Item 2 contains caffeine.
KLJZJV	Item 2 was identified as caffeine.
KMAPJV	Item 2 was found to contain Caffeine.
KNLMRK	Item 2 contained caffeine.
KYRGAK	I found Caffeine to be in substance 2.
L7BGKU	Testosterone is not controlled substance in [Country]. Item 2 contained non-controlled substance - Caffeine.
LAJLGG	For Item #1 (2) solvent systems were used for TLC testing: MEOH: NH4OH (performed on 4/27/22) and MEOH (performed on 4/28/22). Both solvent systems displayed consistency between standard and sample however, both standard and sample traveled to the solvent line.
M2UYKL	Item 2 was identified as Containing Caffeine
M49J7L	Item 2 contained caffeine.
M8T8LY	The substance identified on item 2 was Caffeine.
ME4WYQ	Item 2 contains caffeine.
MRCMXR	Item 2 Caffeine
MVXBMY	Item 2: Caffeine
NL87AJ	Item 1 contains caffeine.
NQEUDT	Item 1 contains Testosterone, which is not a controlled drug according to the [Country] law of Narcotics.
NRYFJH	Item #2 contained Caffeine -not a controlled substance
NUHWU3	Caffeine identified in Item 2.
PFR2HN	Item 2 contains caffeine.
PT3NQD	Item 2 contains Caffeine
PZY7J4	Item 2 contained caffeine.

TABLE 3

WebCode	Additional Comments
Q2N2CQ	Work was performed at the Central Laboratory location in [Location].
Q4H4NT	Item 2 tentatively identified caffeine.
QHXYLG	Item 2 contained caffeine
QUABWQ	Testosterone Acetate indicated but not identified in GCMS
R6PKW4	Item 2 contained the non-controlled substance caffeine.
RBTHEF	Caffeine was confirmed in Item 2
RG9VPV	Testosterone was presumptively identified in item#1. This compound is not included in lab's confirmatory scope of testing. Testing for this compound would typically be sent out for completion by an outside reference laboratory.
RJNJRD	Caffeine was confirmed in Item 2
RKZ26Q	Caffeine was confirmed in Item 2
T9BH88	Item 2 was identified as Containing Caffeine.
TGV2YQ	Item 2 was presumptively identified to contain caffeine.
TLAX69	Unless stated otherwise, the isomer of the substance(s) reported was not determined.
TLPVYD	Caffeine (non-controlled substance) was presumptively identified in Item 2 as well.
UCYQMQ	Item 2 was found to contain caffeine.
UFV76M	Item 2 contains caffeine.
UFX2JQ	Item 2 contained caffeine.
UJ973D	Caffeine
UQLTDW	item 2 was found to contain caffeine.
UX43YE	Identified caffeine in Item 2.
UYZGL8	Caffeine was identified in item 2.
V3HK2K	Caffeine was confirmed in Item 2.
V3UWB7	Item 2 contains Caffeine.
V9CF6D	Caffeine was identified as a non-controlled substance for Item 2.
VBN6UB	Caffeine was detected by UV, FTIR/ATR, and GCMS in Item 2.
VCJMVD	Item 2 confirmed to contain caffeine.
VE4YQN	In Item 2, Caffeine was tentatively identified using the GC/MC method.
VURVKJ	Item 2 was found to contain Caffeine
W8F9DK	Item 2 contains Caffeine
W8FAA6	Caffeine was detected in sample 2.

TABLE 3

WebCode	Additional Comments
WG2G33	Test item 2 contained the non-controlled substance caffeine.
WMUEZ7	Caffeine (noncontrolled) was identified in Item 2.
WQMPJJ	In Item #2, caffeine was detected.
WWWWB9	Caffeine was confirmed in Item 2
X83WK3	Caffeine detected in item 2.
X9RUPF	Caffeine was detected in Item 2
XUMTVB	Item 2 contains caffeine.
Y3F9C3	Item 2 was found to contain Caffeine.
Y7FJ82	Item 2 was found to contain Caffeine.
YFEV3M	Item 1 contains testosterone. Item 2 contains caffeine.
YGBBBD	Caffeine was confirmed in Item 2.
YXGBP9	Item 2 contained caffeine.
Z2WGMV	Caffeine was identified in Item 2.
ZACKPB	Item 2 contained caffeine
ZB8J9V	In Item 1 Testosterone and Androst-5-en-17b-ol-3-one were identified (Testosterone was found only in the FTIR analysis). In Item 2 only Caffeine was identified in both FTIR and GC/MS analysis.
ZBDRE6	Item #2 contained Caffeine
ZF3UY9	Caffeine was detected in item 2
ZPHD8Y	In item 2 the non-controlled substance caffeine was preliminarily identified.
ZQT9VC	Item 2 contained Caffeine
ZYX2FY	Caffeine was detected in Item 2.
ZZJT2A	Caffeine was confirmed in Item 2.

**-End of Report-
(Appendix may follow)**

Test No. 22-5001: Drug Analysis

DATA MUST BE SUBMITTED BY **May 31, 2022, 11:59 p.m.** TO BE INCLUDED IN THE REPORT

Participant Code: U1234A

WebCode: Y3F7Q4

The Accreditation Release section can be accessed by using the "Continue to Final Submission" button above. This information can be entered at any time prior to submitting to CTS.

Scenario:

An investigator has submitted two questioned powdered substances for analysis. Please examine the samples and determine if any controlled substances are present.

Please note: The controlled status of substances is based on the list in section 1308 of Title 21 Code of Federal Regulations under the United States Controlled Substances Act.

Items Submitted (Sample Pack DA1):

Item 1: First powdered substance

Item 2: Second powdered substance

1a.) Were any controlled substances present in Item 1? If yes, list each substance.

- No controlled substance detected/identified Controlled Substance(s) detected/identified (please list below).
Indicate salt form if determined.

1b.) What methods were used to examine Item 1?

(If color and/or crystal tests are used, please simply check the appropriate method below. Specific color and crystal tests will not be listed in the summary report.)

Color Tests FTIR LC UV GC
Crystal Tests TLC LC/MS Raman GC/MS

Other (specify):

2a.) Were any controlled substances present in Item 2? If yes, list each substance.

- No controlled substance detected/identified Controlled Substance(s) detected/identified (please list below).
Indicate salt form if determined.

2b.) What methods were used to examine Item 2?

(If color and/or crystal tests are used, please simply check the appropriate method below. Specific color and crystal tests will not be listed in the summary report.)

- Color Tests FTIR LC UV GC
Crystal Tests TLC LC/MS Raman GC/MS

Other (specify):

3.) Additional Comments

Please note: Any additional formatting applied in the free form space below will not transfer to the Summary Report and may cause your information to be illegible. This includes additional spacing and returns that present your responses in lists and tabular formats.

RELEASE OF DATA TO ACCREDITATION BODIES

The Accreditation Release is accessed by pressing the "Continue to Final Submission" button online and can be completed at any time prior to submission to CTS.

CTS submits external proficiency test data directly to ASCLD/LAB, ANAB, and/or A2LA. Please select one of the following statements to ensure your data is handled appropriately.

- This participant's data is intended for submission to ASCLD/LAB, ANAB, and/or A2LA. (Accreditation Release section below must be completed.)
- This participant's data is **not** intended for submission to ASCLD/LAB, ANAB, and/or A2LA.

Have the laboratory's designated individual complete the following steps **only if your laboratory is accredited in this testing/calibration discipline** by one or more of the following Accreditation Bodies.

Step 1: Provide the applicable Accreditation Certificate Number(s) for your laboratory

ANAB Certificate No.
(Include ASCLD/LAB Certificate here)

A2LA Certificate No.

Step 2: Complete the Laboratory Identifying Information in its entirety.

Authorized Contact Person and Title

Laboratory Name

Location (City/State)