



## **Body Fluid Identification**

# **Test No. 23-5782 Summary Report**

---

Each participant received a sample set consisting of six items for examination. Participants were asked to examine each item and report the body fluid(s) identified. Data were returned from 178 participants and are compiled into the following tables:

<b>Report Contents:</b>	<b>Page:</b>
<a href="#"><u>Manufacturer's Information</u></a>	2
<a href="#"><u>Summary Comments</u></a>	3
<a href="#"><u>Key for Screening Tests Used</u></a>	4
<a href="#"><u>Table 1: Screening Conclusions &amp; Test(s) Performed</u></a>	5
<a href="#"><u>Table 2: Additional Comments</u></a>	114
<a href="#"><u>Appendix: Data Sheet</u></a>	

This report contains the data received from the participants in this test. Since these participants are located in many countries around the world, and it is their option how the samples are to be used (e.g., training exercise, known or blind proficiency testing, research and development of new techniques, etc.), the results compiled in the Summary Report are not intended to be an overview of the quality of work performed in the profession and cannot be interpreted as such. The Summary Comments are included for the benefit of participants to assist with maintaining or enhancing the quality of their results. These comments are not intended to reflect the general state of the art within the profession.

Participant results are reported using a randomly assigned "WebCode". This code maintains participant's anonymity, provides linking of the various report sections, and will change with every report.

# Manufacturer's Information

Each sample set consisted of six separate items for testing. Two items contained non-scenario related swabs (Items 1 and 2), two items contained swabs that related to the provided scenario (Items 3 and 4), and two items (Items 5 and 6) contained complete items of evidence that also related to the provided scenario. Participants were requested to examine Items 1 and 2 for body fluids that they typically encounter in routine casework examinations. For Items 3 through 6, the participants were requested to examine the items using the information provided by the scenario for guidance.

**SAMPLE PREPARATION:** Items were prepared as described in the table below. For the swab items where semen was used, the semen was mixed in a 1:1 ratio with TAE buffer. For the evidence items where semen was used, the semen was not mixed with TAE buffer. In all items where a mixture is present, the two body fluids were mixed in a 1:1 ratio. Stains from different sources were prepared at separate times and were packaged once thoroughly dried.


**SAMPLE SET ASSEMBLY:** For each sample set, all Items (1-6) were packaged into separate envelopes and then placed together in a pre-labeled sample set envelope and sealed. The sealed sample set envelopes were then packaged in pre-labeled heat seal envelopes and sealed. This process was repeated until all of the sample sets were prepared.

**VERIFICATION:** All predistribution laboratories confirmed the presence of the body fluid(s) that comprised each item.


<i>Sample Description</i>			
<u>Item</u>	<u>Substrate</u>	<u>Body Fluid</u>	<u>Volume</u>
1	4 Swabs: white plastic handle	Semen	100µL per swab
2	4 Swabs: pink plastic handle	Blood & Saliva Mixture	100µL per swab
3	4 Swabs: blue plastic handle	Semen	100µL per swab
4	4 Swabs: yellow plastic handle	Blood	100µL per swab
5	Lanyard	Blood & Saliva Mixture	200µL
6	Washcloth	Blood	200µL

The symbol "X" indicates the general location where the body fluids were deposited on Items 5 and 6.

**Item 5**



**Item 6**



Completed sample sets were stored at -20°C until shipment on September 18, 2023, following completion of the verification stage.

## **Summary Comments**

The Body Fluid Identification test was designed to allow participants to assess their proficiency in the identification of stains that may contain dried body fluids. Each participant received six items for testing; four items were swabs and two were whole evidence items. A scenario was not provided for Items 1 and 2; participants were requested to test for body fluids that they typically encounter in casework. A scenario was provided for Items 3 - 6, and participants were asked to treat these items as they would normal casework, testing for the fluids they deemed necessary after reviewing the scenario.

Item 1 consisted of swabs that were prepared using semen. The majority of responding participants reported the expected results: "Negative" for blood (98%), "Positive/Indicative" for semen (100%) and "Negative" for saliva (100%).

Item 2 consisted of swabs that were prepared using a mixture of blood and saliva. The majority of responding participants reported the expected results: "Positive/Indicative" for blood (100%), "Negative" for semen (99%), and "Positive/Indicative" for saliva (98%).

Item 3 consisted of swabs that were prepared using semen. The majority of responding participants reported the expected results: "Negative" for blood (100%), "Positive/Indicative" for semen (99%), and "Negative" for saliva (99%).

Item 4 consisted of swabs that were prepared using blood. The majority of responding participants reported the expected results: "Positive/Indicative" for blood (100%), "Negative" for semen (100%), and "Negative" for saliva (94%).

Item 5 consisted of a lanyard that was prepared using a mixture of blood and saliva. All responding participants reported the expected results: "Positive/Indicative" for blood (100%), "Negative" for semen (100%), and "Positive/Indicative" for saliva (99%).

Item 6 consisted of a washcloth that was prepared using blood. The majority of responding participants reported the expected results: "Positive/Indicative" for blood (100%), "Negative" for semen (99%), and "Negative" for saliva (94%).

## Key for Screening Tests Used

Participants were asked to use, where possible, the following chart of abbreviated screening test names. This was not an all inclusive list and was not designed to determine what tests should be performed. Participants were advised that tests not on this list may be used for screening.

<u>Test</u>	<u>Abbreviation</u>
Acid Phosphatase	AP
Alternate Light Source	ALS
Kastle Meyer	KM
Leucomalachite Green	LMG
Microscopic	Micro
Ortho-tolidine	O-tol
Phenolphthalein-Tetramethyl benzidine	PTMB
Prostate Specific Antigen	PSA
Rapid Stain Identification	RSID
Tetramethyl benzidine	TMB

# Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 1  
Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
2PK689	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
2PZNR7	Visual Examination <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID Semen <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID Saliva <i>Negative</i> <b>Conclusion:</b> Negative
2XF3XH	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
38Y4RX	LMG <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
3996DV	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	Phadebas <i>Negative</i> <b>Conclusion:</b> No Response
3H2GPC	PTMB <i>Negative</i> Bluestar <i>Negative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
3RNEGH	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Urine - Creatinine <i>Negative</i>		
3URBZK	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Spermatozoa Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Urine - Jaffe Test <i>Negative</i>		
3XDWT9	Visual examination did not detect red or brown stains therefore Blood testing was not performed <b>Conclusion:</b> Not Tested	Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	This laboratory does not have a current validated method for Saliva testing <b>Conclusion:</b> Not Tested

Table 1 - Item 1  
 Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
3YL8PC	PTMB	Negative	AP	Positive/Indicative	Phadebus Tablet	Negative
		<b>Conclusion:</b> Negative	PSA	Positive/Indicative		<b>Conclusion:</b> Negative
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
4DGTNC	TMB	Negative	AP	Positive/Indicative	SALigAE	Negative
	HemaTrace	Negative	Micro	Positive/Indicative		
		<b>Conclusion:</b> Negative	PSA	Positive/Indicative		<b>Conclusion:</b> Negative
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
4HDDVG	KM	Negative	AP	Positive/Indicative	Phadebas	Negative
		<b>Conclusion:</b> Negative	Micro	Positive/Indicative		<b>Conclusion:</b> Negative
	<u>Other Tests Performed</u>					
	Urine- Creatinine Negative					
4JLPBX	Phenolphthalein	Negative	Microscopic exam	Positive/Indicative	Phadebas	Negative
		<b>Conclusion:</b> Negative	p30	Positive/Indicative		<b>Conclusion:</b> Negative
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
4N3XXW	KM	Negative	Brentamine (AP)	Positive/Indicative	Human Amylase (alpha-amilosa)	Negative
		<b>Conclusion:</b> Negative	Human seminal protein (p30)	Positive/Indicative		<b>Conclusion:</b> Negative
		<b>Conclusion:</b> Negative	Microscopic	Positive/Indicative		<b>Conclusion:</b> Negative
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
4R39RC	TMB	Negative	AP	Positive/Indicative	SALigAE	Negative
	HemaTrace	Negative	p30	Positive/Indicative		
		<b>Conclusion:</b> Negative	Micro	Positive/Indicative		<b>Conclusion:</b> Negative
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
63V48H	KM	Negative	AP	Positive/Indicative	SERATEC AmylaseTest	Negative
	SERATEC HemDirect	Negative	SERATEC PSA	Positive/Indicative		<b>Conclusion:</b> No Response
		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response
69JRZ3	Visual	Negative	AP	Positive/Indicative	RSID	Negative
		<b>Conclusion:</b> Negative	RSID	Positive/Indicative		<b>Conclusion:</b> Negative
		<b>Conclusion:</b> Negative	Micro	Positive/Indicative		<b>Conclusion:</b> Negative
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
6BBMPW	benzidine	Negative	acid phosphotase solution	Positive/Indicative		<b>Conclusion:</b> Not Tested
		<b>Conclusion:</b> Negative	PSA kit	Positive/Indicative		<b>Conclusion:</b> Not Tested
		<b>Conclusion:</b> Negative	Microscope	Positive/Indicative		<b>Conclusion:</b> Not Tested
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested

Table 1 - Item 1  
 Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
6HUTL4	Visual Examination <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i>  <b>Conclusion:</b> Negative
6WTV3B	LMG <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Amylase Diffusion <i>Negative</i>  <b>Conclusion:</b> Negative
6YVZDB	Hemastix <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> ABACard p30 <i>Positive/Indicative</i> Sperm Search - Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	SALlgAE <i>Negative</i>  <b>Conclusion:</b> Negative
79FQD8	KM <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
7C4X2B	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
7FYPA8	KM <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i>  <b>Conclusion:</b> Negative
7G9YXY	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
7JZXCW	No brownish traces <b>Conclusion:</b> Not Tested	PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Amylase - Method CNPG3 <i>Negative</i> <b>Conclusion:</b> Negative
7LBQEC	Hemastix <i>Negative</i>  <b>Conclusion:</b> Negative	Microscopy <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID-Saliva <i>Negative</i>  <b>Conclusion:</b> Negative
7UGHWT	  <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RADIAL Diffusion <i>Negative</i>  <b>Conclusion:</b> Negative

Table 1 - Item 1  
 Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
83AAGT	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	   <b>Conclusion:</b> Not Tested
8AV3XY	KM <i>Negative</i>  <b>Conclusion:</b> No Response	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> No Response	Phadebas Paper <i>Negative</i> RSID <i>Negative</i>  <b>Conclusion:</b> No Response
8CE398	KM <i>Negative</i>  <b>Conclusion:</b> Negative	   <b>Conclusion:</b> Not Tested	   <b>Conclusion:</b> Not Tested
8GNHHW	TMB <i>Negative</i>  <b>Conclusion:</b> Negative	PSA <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Seratec Amylase <i>Negative</i>  <b>Conclusion:</b> Negative
8Z7GPE	No Blood-like stains observed  <b>Conclusion:</b> Not Tested	   <b>Conclusion:</b> Not Tested	   <b>Conclusion:</b> Not Tested
9BFE63	Visual <i>Negative</i> LMG <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID Semen <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	RSID Saliva <i>Negative</i>  <b>Conclusion:</b> Negative
9EJ86T	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Seratec <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i>  <b>Conclusion:</b> Negative
9K9PJ4	Micro <i>Negative</i> KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i>  <b>Conclusion:</b> Negative
9QA7Q4	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	RSID-Saliva <i>Negative</i>  <b>Conclusion:</b> Negative
9QVF4E	  <b>Conclusion:</b> Negative  <u>Other Tests Performed</u> ALS <i>Positive/Indicative</i>	   <b>Conclusion:</b> Not Tested	   <b>Conclusion:</b> Not Tested



Table 1 - Item 1  
 Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
9UVNAA	KM	Negative	AP	Positive/Indicative		
			PSA	Positive/Indicative		
			Micro	Positive/Indicative		
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
9VPA2F			RSID	Positive/Indicative	RSID	Negative
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
A7LW93	KM	Negative	AP	Positive/Indicative	RSID	Negative
			Micro	Positive/Indicative		
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
ACUZGP			PSA Semiquant	Positive/Indicative		
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
AEDFMB	LMG	Negative	AP	Positive/Indicative		
			Micro	Positive/Indicative		
			PSA	Positive/Indicative		
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
AF8ZDG	Visual appearance	Negative	ALS	Positive/Indicative		Our lab does not perform screening tests for Saliva
	Bluestar forensic tablets	Negative	Micro	Positive/Indicative		
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
ALCYQR	No Visual indications of Blood to test.		ALS, AP, Micro	Positive/Indicative	ALS, RSID Saliva	Negative
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
AQ64HB	No Blood-like discoloration noted.		AP	Positive/Indicative	Phadebas Amylase Test	Negative
			PSA	Positive/Indicative		
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
ARYUNP	KM	Negative	ALS	Positive/Indicative		
	HemaTrace	Negative	AP	Positive/Indicative		
			Micro	Positive/Indicative		
			PSA	Positive/Indicative		
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
B7XV2Q	No red brown staining		AP	Positive/Indicative	Radial Diffusion	Negative
			Micro	Positive/Indicative		
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
BLQQBW	Visual	Negative	AP	Positive/Indicative	RSID-Saliva	Negative
			RSID-Semen	Positive/Indicative		
			Micro	Positive/Indicative		
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative

Table 1 - Item 1  
Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
BNVB33	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i>  <b>Conclusion:</b> Negative
BQJU4Z	Visual <i>Negative</i> LMG <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i>  <b>Conclusion:</b> Negative
BTQ6Z4	Phenolphthalein <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
BU9Y2D	PEROXIDASES <i>Negative</i> HEXAGON OBTI <i>Negative</i> <b>Conclusion:</b> Negative	PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	TEST IS NOT PERFORMED  <b>Conclusion:</b> Not Tested
C84ZT6	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
CDC42R	RSID Blood <i>Negative</i> <b>Conclusion:</b> Negative	RSID Semen <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID Saliva <i>Negative</i> <b>Conclusion:</b> Negative
CGD64B	HEXAGON OBTI <i>Negative</i>  <b>Conclusion:</b> Negative	PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
CKDD9K	KM <i>Negative</i> <b>Conclusion:</b> No Response	AP <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	Phadebas <i>Negative</i> <b>Conclusion:</b> No Response
CNWLH7	  <b>Conclusion:</b> Not Tested	Acid phosphatase test <i>Positive/Indicative</i> Microscopy <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas paper press test <i>Negative</i> <b>Conclusion:</b> Negative
CQFMDR	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
CYTE6B	LMG <i>Negative</i> OBTI <i>Negative</i> <b>Conclusion:</b> No Response	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested

Table 1 - Item 1  
Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
DH7FGP	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
DKQAC4	FOB <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
DQFZL7	Visual appearance of Blood <i>Negative</i> KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro (Spermatozoa) <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas (press test) <i>Negative</i> <b>Conclusion:</b> Negative
DWKLQ3	KM <i>Negative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Sperm cells - Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Nucleated epithelial cells - Micro <i>Negative</i>		
DZQ388	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas paper press test <i>Negative</i> <b>Conclusion:</b> Negative
EBRRPL	LMG <i>Negative</i> RSID-Blood <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID-Saliva <i>Negative</i> <b>Conclusion:</b> Negative
ED299W	Visual <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID-Saliva <i>Negative</i> <b>Conclusion:</b> Negative
EGQD7A	HEXAGON OBTI IMMUNOCHROMATOGRAPHIC RAPID TEST <i>Negative</i> PEROXIDASE TEST <i>Negative</i> <b>Conclusion:</b> Negative	PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	TEST NOT PERFORMED <b>Conclusion:</b> Not Tested
EQF6UM	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
FBHAXM	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Urine - Creatinine <i>Negative</i>		

Table 1 - Item 1  
Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
FDLCGW	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID-Saliva <i>Negative</i>  <b>Conclusion:</b> Negative
FGHUXX	KM <i>Negative</i>  <b>Conclusion:</b> Negative	Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	[Laboratory] does not perform Saliva testing. <b>Conclusion:</b> Not Tested
FNJ9HY	PTMB <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas tablet <i>Negative</i> <b>Conclusion:</b> Negative
FV2LGX	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Sperm Hy-Liter PI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID-Saliva <i>Negative</i>  <b>Conclusion:</b> Negative
FWF8TK	kastle meyer KM <i>Negative</i>  <b>Conclusion:</b> No Response	acid phosphatase ap <i>Positive/Indicative</i> microscopy - christmas tree staining <i>Positive/Indicative</i> RSID Semen <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	RSID Saliva <i>Negative</i>  <b>Conclusion:</b> No Response
FXTH7R	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Urine - Jaffe Test <i>Negative</i>		
FZA2QL	O-tolidine <i>Negative</i> <b>Conclusion:</b> No Response	AP and p30 <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	Lab does not test for Saliva <b>Conclusion:</b> Not Tested
G6M6E3	LMG <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
GBR2ZN	No red/brown stain present  <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID Saliva <i>Negative</i> <b>Conclusion:</b> Negative
GCMLNM	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative

Table 1 - Item 1  
 Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
GFJDUY	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID-Saliva <i>Negative</i>  <b>Conclusion:</b> Negative
GTTNRM	Kastle Meyer <i>Negative</i>  <b>Conclusion:</b> Negative	Acid Phosphatase <i>Positive/Indicative</i> Microscopic <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas Amylase <i>Negative</i>  <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Urine - Creatinine <i>Negative</i>		
GVXEUV	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i>  <b>Conclusion:</b> Negative
H2GYEM	no reddish brown stains observed on swabs  <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID Saliva <i>Negative</i>  <b>Conclusion:</b> Negative
H3CJ3M	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
HFJHYZ	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
HQZC42	LMG <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	  <b>Conclusion:</b> Not Tested
HVCZ4X	Phenolphthalein <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	  <b>Conclusion:</b> Not Tested
J2ETAD	Visual exam <i>Negative</i> KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Sperm elution and Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas tube test <i>Negative</i>  <b>Conclusion:</b> Negative

Table 1 - Item 1  
Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
J6WQ4V	KM Negative Micro Negative	<b>Conclusion:</b> Negative	AP Positive/Indicative Micro Positive/Indicative	<b>Conclusion:</b> Positive/Indicative	RSID Negative	<b>Conclusion:</b> Negative
JQFVNY	O-tol Negative	<b>Conclusion:</b> Negative	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative	<b>Conclusion:</b> Positive/Indicative	RSID Negative	<b>Conclusion:</b> Negative
K4X8BR	KM Negative	<b>Conclusion:</b> Negative	AP Positive/Indicative RSID-Semen Positive/Indicative PSA Positive/Indicative	<b>Conclusion:</b> Positive/Indicative	RSID-Saliva Negative	<b>Conclusion:</b> Negative
KBH3F3	Negative	<b>Conclusion:</b> Negative	Positive/Indicative	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
KLZMZY	KM Negative	<b>Conclusion:</b> Negative	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
KMTD7D	LMG Negative	<b>Conclusion:</b> Negative	AP Positive/Indicative p30 Positive/Indicative	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
KQ99JP	KM Negative	<b>Conclusion:</b> Negative	Micro Positive/Indicative	<b>Conclusion:</b> Positive/Indicative	Our laboratory does not currently have a validated test for Saliva.	<b>Conclusion:</b> Not Tested
KTFL8J	KM Negative	<b>Conclusion:</b> Negative	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative	<b>Conclusion:</b> Positive/Indicative	RSID Negative	<b>Conclusion:</b> Negative
KXT4PZ		<b>Conclusion:</b> Not Tested	AP Positive/Indicative	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
	<u>Other Tests Performed</u> ALS Positive/Indicative					
KZXWGD	Molecular Serology Negative	<b>Conclusion:</b> Negative	Molecular Serology Positive/Indicative	<b>Conclusion:</b> Positive/Indicative	Molecular Serology Negative	<b>Conclusion:</b> Negative

Table 1 - Item 1  
Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
L2CKH6	ALS BLUE (420-470 nm) with yellow glasses <i>Negative</i> Hemastix (SIEMENS) <i>Negative</i>  <b>Conclusion:</b> Negative	ALS BLUE (420-470 nm) with yellow glasses <i>Positive/Indicative</i> ALS UV (350-380 nm) with UV glasses <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> AP <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	     <b>Conclusion:</b> Not Tested
L73VR3	PEROXIDASES <i>Negative</i> HEXAGON OBTI <i>Negative</i>  <b>Conclusion:</b> Negative	PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	TEST IS NOT PERFORMED.  <b>Conclusion:</b> Not Tested
LECJ72	  <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	  <b>Conclusion:</b> Not Tested
	<u>Other Tests Performed</u> Microscopic Slide <i>Positive/Indicative</i>		
LJAZHV	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	     <b>Conclusion:</b> Not Tested
M4VWUF	LMG (no Blood-like stains observed) <i>Negative</i>  <b>Conclusion:</b> Negative	     <b>Conclusion:</b> Not Tested	     <b>Conclusion:</b> Not Tested
MCLHXT	no red/brown staining  <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	SALigAE <i>Negative</i>     <b>Conclusion:</b> Negative
MFNJZB	No Visible Bloodlike staining  <b>Conclusion:</b> Not Tested	AP, Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
MUMKGH	reddish-brown staining not observed, so no KM test performed  <b>Conclusion:</b> Not Tested	PSA <i>Positive/Indicative</i> Microscopic <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i>     <b>Conclusion:</b> No Response
N44NZ2	Bluestar OBTI <i>Negative</i>  <b>Conclusion:</b> Negative	PSA <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	     <b>Conclusion:</b> Not Tested
NAT4TF	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>     <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Urine - Creatinine <i>Negative</i>		

Table 1 - Item 1  
Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
NBNQ6J	Visual <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i>  <b>Conclusion:</b> Negative
NDTENE	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Xmas/Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Urine-Creatinine <i>Negative</i>		
NDVZKV	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
NGE2EG	KM Test <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Spermatozoa Identification (Microscopic) <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
NNYWKR	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
NVFQJU	Visual examination <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas press test <i>Negative</i>  <b>Conclusion:</b> Negative
NXYREG	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
P7B7XX	KM <i>Negative</i> HEM DIRECT HEMOGLOBIN <i>Negative</i>  <b>Conclusion:</b> Negative	PSA <i>Positive/Indicative</i> Micro - CHRISTMAS TREE <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
P9FVYK	<b>Conclusion:</b> Not Tested	Micro <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
PFY2VR	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Urine - creatinine <i>Positive/Indicative</i>		



Table 1 - Item 1  
 Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
PGCTPU	LMG <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
PTJJQG	<b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Saliva testing is not performed at this laboratory.  <b>Conclusion:</b> Not Tested
Q9A4JA	No RBS present <b>Conclusion:</b> Not Tested	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
Q9UKZC	<b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
QCDPKT	No visible staining. Not tested per laboratory policy <b>Conclusion:</b> Not Tested	PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
QCVL2U	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro Exam <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
QD4XJW	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
QD6VPL	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> Sperm Hy-Liter PI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID-Saliva <i>Negative</i>  <b>Conclusion:</b> Negative
QDPC8A	No visible appearance of Blood or under ALS/IR <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> ALS <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	CSI does not test for Saliva  <b>Conclusion:</b> Not Tested
QEXQEF	No Visual characteristics of Blood  <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested

Table 1 - Item 1  
Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
QH2EXB	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Phadebas Amylase <i>Negative</i>  <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Urine - Creatinine <i>Negative</i>		
QMUJMA	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
QV428T	<b>Conclusion:</b> Not Tested	Microscopy <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
QZ2CEL	Visual, based on color of swabs testing was not warranted.  <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Saliva Tests are not performed at the [Laboratory]  <b>Conclusion:</b> Not Tested
T4XQ3M	<b>Conclusion:</b> Not Tested	Micro <i>Positive/Indicative</i> AP <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	RSID  <b>Conclusion:</b> Negative
TCMMGA	LMG <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Laboratory no longer performs this testing  <b>Conclusion:</b> Not Tested
THCZHx	Bluestar Obti <i>Negative</i>  <b>Conclusion:</b> Negative	Bluestar Identi-PSA <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> No Response
TLCCEW	Visual <i>Negative</i> fluorescence 445 nm to 510 nm <i>Negative</i> fluorescence 420 nm to 470 nm <i>Negative</i>  <b>Conclusion:</b> No Response	Visual <i>Negative</i> fluorescence 420 nm to 470 nm <i>Positive/Indicative</i> fluorescence 445 nm to 510 nm <i>Positive/Indicative</i> Bluestar PSA <i>Positive/Indicative</i>  <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested
TM7X8N	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
TMQJZE	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative

Table 1 - Item 1  
 Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
TQPWLH	Phenolphthalein HemaTrace	Negative Negative	Acid Phosphatase p30 Microscopic	Positive/Indicative Positive/Indicative Positive/Indicative	[Laboratory] does not conduct Saliva examinations	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
TRY22R	KM	Negative	AP PSA Microscopic	Positive/Indicative Positive/Indicative Positive/Indicative	Phadebas	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
TW9FKC	No Reddish/Brown Staining Observed		p30 Micro	Positive/Indicative Positive/Indicative	RSID	Negative
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
TXPCBQ	KM	Negative	AP PSA Micro	Positive/Indicative Positive/Indicative Positive/Indicative		
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
U7VQFP	KM	Negative	AP Micro	Positive/Indicative Positive/Indicative	Phadebas	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
U87FTA						
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
UJKFP4	Visual Kastle-Mayer test of supernatant	Negative Negative	Acid phosphatase test of supernatant Microscopy	Positive/Indicative Positive/Indicative	Phadebas tube test of supernatant (Amylase)	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
UMZ8EB	KM	Negative	AP Microscopic	Positive/Indicative Positive/Indicative	Phadebas	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
VGR42G	KM	Negative	AP RSID-Semen PSA	Positive/Indicative Positive/Indicative Positive/Indicative	Not competent / authorised to test for Saliva	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
VHJYQB	TMB	Negative	AP RSID	Positive/Indicative Positive/Indicative	RSID	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
	<u>Other Tests Performed</u>					
	Blood - White light Visual examination Negative					
	Semen - White light Visual examination Inconclusive					
	Saliva - White light Visual examination Inconclusive					

Table 1 - Item 1  
Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
VMDWF3	no red-brown staining observed <b>Conclusion:</b> Not Tested	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
VXWD89	<b>Conclusion:</b> No Response	brentamine spot test <i>Positive/Indicative</i> ABAcad p30 <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> No Response
W4ZMZM	O-tol <i>Negative</i> <b>Conclusion:</b> Negative	No current test in laboratory for Semen <b>Conclusion:</b> Not Tested	No current test in laboratory for Saliva <b>Conclusion:</b> Not Tested
WEGF6N	O-tol <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
WL47RJ	KM <i>Negative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Inconclusive	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Positive/Indicative</i> Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
WP2REF	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Sperm Hy-Liter PI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID-Saliva <i>Negative</i> <b>Conclusion:</b> Negative
WVQ6HP	LMG <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
WZYRRK	KM <i>Negative</i> HemaTrace <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
X32Q3D	Visual examination did not detect brown stains <b>Conclusion:</b> Not Tested	Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Currently, the lab does not have a validated method for testing for Saliva <b>Conclusion:</b> Not Tested
X7HDGH	KM <i>Negative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Inconclusive	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative

Table 1 - Item 1  
 Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
XBE2V4	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
XEYPW4	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
XPVZZ6	LMG <i>Negative</i> HemaTrace <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	No longer performed <b>Conclusion:</b> Not Tested
XTUBW6	LMG <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	No longer performed <b>Conclusion:</b> Not Tested
Y48XHB	No visible red/brown stains. <b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
Y7YU86	RSID Blood <i>Negative</i> <b>Conclusion:</b> Negative	RSID Semen <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID Saliva <i>Negative</i> <b>Conclusion:</b> Negative
Y94GNF	PTMB <i>Negative</i> Bluestar <i>Negative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
YDXBLG	KM <i>Negative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Prostate specific antigen <i>Positive/Indicative</i> ALS <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> ALS <i>Positive/Indicative</i> <b>Conclusion:</b> Negative
YEUW9F	TMB <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	SALigAE <i>Negative</i> <b>Conclusion:</b> Negative
YNE396	Bluestar Hexagon OBTI <i>Negative</i> <b>Conclusion:</b> Negative	Seratec PSA Semiquant <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Seratec Amylase Test <i>Negative</i> <b>Conclusion:</b> Negative
YNVC27	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative

Table 1 - Item 1  
 Evidence - Four white plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
YRECCF	O-tol <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i>  <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Urine - RSID <i>Negative</i>		
YTQ493	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Phadebas Amylase <i>Negative</i>  <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Urine - Creatinine <i>Negative</i>		
Z7EAZC	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
Z9YG2K	LMG <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	     <b>Conclusion:</b> Not Tested
ZADNKG	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> Sperm Hy-Liter PI <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	RSID-Saliva <i>Negative</i>     <b>Conclusion:</b> Negative
ZRC2V3	Combur test <i>Negative</i>  <b>Conclusion:</b> Negative	PSA-RSID <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Amylase <i>Negative</i>  <b>Conclusion:</b> Negative
ZYAT4W	Phenolphthalein <i>Negative</i>  <b>Conclusion:</b> Negative	     <b>Conclusion:</b> Not Tested	     <b>Conclusion:</b> Not Tested
ZZ3LYF	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Not competent / authorised to test for Saliva.     <b>Conclusion:</b> Not Tested

**Response Summary for Table 1 -  
Item 1: Four white plastic swabs containing Semen**

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	0	155	0
Negative	132	0	101
Inconclusive	2	0	0
Not Tested (NT)	34	15	69
No Response	10	8	8
<hr/>			
Total	178	178	178

# Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
2PK689	Img <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
2PZNR7	Visual Examination <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID Saliva <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
2XF3XH	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
38Y4RX	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
3996DV	KM <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	AP <i>Negative</i> <b>Conclusion:</b> No Response	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> No Response
3H2GPC	PTMB <i>Positive/Indicative</i> Bluestar <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
3RNEGH	KM <i>Positive/Indicative</i> Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
	<u>Other Tests Performed</u> Urine - Creatinine <i>Negative</i> Species origin - HemaTrace <i>Positive/Indicative</i>		
3URBZK	KM <i>Positive/Indicative</i> Seratec Hb <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
	<u>Other Tests Performed</u> Urine - Jaffe Test <i>Negative</i>		
3XDWT9	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	This laboratory does not have a current validated method for Saliva testing <b>Conclusion:</b> Not Tested



Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
3YL8PC	PTMB HemaTrace <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative	AP PSA <b>Conclusion:</b> Negative	Negative Negative	Phadebus Tablet <b>Conclusion:</b> Positive/Indicative	Positive/Indicative
4DGTNC	TMB HemaTrace <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative	AP Micro PSA <b>Conclusion:</b> Negative	Negative Negative Negative	SALigAE <b>Conclusion:</b> Positive/Indicative	Positive/Indicative
4HDDVG	KM Modified Takayama <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative	AP <b>Conclusion:</b> Negative	Negative	Phadebas <b>Conclusion:</b> Positive/Indicative	Positive/Indicative
	<u>Other Tests Performed</u> Species Origin- HemaTrace Positive/Indicative Urine- Creatinine Negative					
4JLPBX	Phenolphthalein <b>Conclusion:</b> Positive/Indicative	Positive/Indicative	Microscopic exam p30 <b>Conclusion:</b> Negative	Negative Negative	Phadebas <b>Conclusion:</b> Positive/Indicative	Positive/Indicative
4N3XXW	KM HemaTrace <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative	Bremtamine (AP) Human seminal protein (p30) Microscopic <b>Conclusion:</b> Negative	Negative Negative Negative	Human Amylase (alpha-amilosa) <b>Conclusion:</b> Positive/Indicative	Positive/Indicative
4R39RC	TMB HemaTrace <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative	AP p30 Micro <b>Conclusion:</b> Negative	Negative Negative Negative	SALigAE <b>Conclusion:</b> Positive/Indicative	Positive/Indicative
63V48H	KM SERATEC HemDirect <b>Conclusion:</b> No Response	Positive/Indicative Positive/Indicative	AP SERATEC PSA <b>Conclusion:</b> No Response	Negative Negative	SERATEC AmylaseTest <b>Conclusion:</b> Positive/Indicative	Positive/Indicative
69JRZ3	Visual LMG HemaTrace <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative Positive/Indicative	AP <b>Conclusion:</b> Negative	Negative	RSID <b>Conclusion:</b> Positive/Indicative	Positive/Indicative
6BBMPW	benzidine Hem Direct kit <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative	acid phosphotase solution <b>Conclusion:</b> Negative	Negative	<b>Conclusion:</b> Not Tested	

**Table 1 - Item 2**  
**Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture**

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
6HUTL4	Visual Examination LMG HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	AP	Negative	RSID	Positive/Indicative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Positive/Indicative
6WTV3B	LMG	Positive/Indicative	AP Micro PSA	Negative Negative Negative	Amylase Diffusion	Positive/Indicative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Positive/Indicative
6YVZDB	Hemastix ABACard HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative	SALigAE	Positive/Indicative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Positive/Indicative
79FQD8	KM HemaTrace	Positive/Indicative Positive/Indicative				
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Not Tested	<b>Conclusion:</b>	Not Tested
7C4X2B	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative		
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Not Tested
7FYPA8	KM Micro	Positive/Indicative Positive/Indicative	AP Micro	Negative Negative	RSID	Positive/Indicative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Positive/Indicative
7G9YXY	KM Seratec HB	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Positive/Indicative
7JZXCW	TMB Hexagon OBTI	Positive/Indicative Positive/Indicative	PSA	Negative	Amylase - Method CNPG3	Positive/Indicative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Positive/Indicative
7LBQEC	Hemastix HemaTrace	Positive/Indicative Positive/Indicative	Microscopy AP	Negative Negative	RSID-Saliva	Positive/Indicative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Positive/Indicative
7UGHWT	KM	Positive/Indicative	AP Micro	Negative Negative	RADIAL Diffusion	Positive/Indicative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Positive/Indicative

Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
83AAGT	KM HemaTrace	Positive/Indicative Positive/Indicative	AP Micro PSA	Negative Negative Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested	
8AV3XY	KM HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative	RSID	Positive/Indicative
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response	
8CE398	KM HemaTrace	Positive/Indicative Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
8GNHHW	OBTI	Positive/Indicative	PSA	Negative	Seratec Amylase	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
8Z7GPE	LMG	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
9BFE63	Visual LMG HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	AP	Negative	RSID Saliva	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
9EJ86T	KM HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative	RSID	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
9K9PJ4	Micro KM	Positive/Indicative Positive/Indicative	AP Micro	Negative Negative	RSID	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
9QA7Q4	KM Bluestar OBTI	Positive/Indicative Positive/Indicative	AP RSID-Semen PSA	Negative Negative Negative	RSID-Saliva	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
9QVF4E	LMG	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
9UVNAA	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	

Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
9VPA2F	<b>Conclusion:</b> Not Tested	RSID <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
A7LW93	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
ACUZGP	Phenolphthalein <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
	<u>Other Tests Performed</u> Hemedirect <i>Positive/Indicative</i>		
AEDFMB	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
AF8ZDG	Visual appearance <i>Positive/Indicative</i> Bluestar forensic tablets <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	Our lab does not perform screening tests for Saliva <b>Conclusion:</b> Not Tested
ALCYQR	Visual, KM, HemDirect <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS, AP <i>Negative</i> <b>Conclusion:</b> Negative	ALS, RSID Saliva <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
AQ64HB	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas Amylase Test <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
ARYUNP	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Positive/Indicative</i> Micro <i>Negative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
B7XV2Q	Phen <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	Radial Diffusion <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
BLQQBW	Visual <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID-Saliva <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
BNVB33	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative

Table 1 - Item 2  
Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
BQJU4Z	Visual LMG HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	AP	Negative	RSID	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
BTQ6Z4	Phenolphthalein HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
BU9Y2D	PEROXIDASES HEXAGON OBTI	Positive/Indicative Positive/Indicative	PSA Micro	Negative Negative		TEST IS NOT PERFORMED
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
C84ZT6	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
CDC42R	RSID Blood	Positive/Indicative	RSID Semen	Negative	RSID Saliva	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
CGD64B	HEXAGON OBTI	Positive/Indicative	PSA Micro	Negative Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
CKDD9K	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response	
CNWLH7	Appearance Kastle Meyer test HemaTrace test	Positive/Indicative Positive/Indicative Positive/Indicative	Acid phosphatase test	Negative	Phadebas paper press test	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
CQFMDR	KM SERATEC Hb	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
CYTE6B	LMG OBTI	Positive/Indicative Positive/Indicative	AP Micro	Inconclusive Negative		
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested	
DH7FGP	KM	Positive/Indicative	AP	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	

Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
DKQAC4	FOB	Positive/Indicative	RSID	Negative	RSID	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
DQFZL7	Visual appearance of Blood	Positive/Indicative	AP	Negative	Phadebas (press test)	Positive/Indicative
	KM	Positive/Indicative				
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
DWKLQ3	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	HemaTrace	Positive/Indicative	PSA	Negative		
			Sperm cells - Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
	<u>Other Tests Performed</u>					
	Nucleated epithelial cells - Micro Negative					
DZQ388	KM	Positive/Indicative	AP	Negative	Phadebas paper press test	Positive/Indicative
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
EBRRPL	LMG	Positive/Indicative	AP	Negative	RSID-Saliva	Positive/Indicative
	RSID-Blood	Positive/Indicative	RSID-Semen	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
ED299W	Visual	Positive/Indicative	AP	Negative	RSID-Saliva	Positive/Indicative
	LMG	Positive/Indicative				
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
EGQD7A	HEXAGON OBTI IMMUNOCHROMATOGRAPHIC RAPID TEST	Positive/Indicative	PSA	Negative	TEST NOT PERFORMED	
	PEROXIDASE TEST	Positive/Indicative	Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
EQF6UM	KM	Positive/Indicative	AP	Negative	RSID	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
FBHAXM	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	Modified Takayama	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
	<u>Other Tests Performed</u>					
	Species - HemaTrace Positive/Indicative					
	Urine - Creatinine Negative					

Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
FDLCGW	KM Bluestar OBTI	Positive/Indicative Positive/Indicative	AP RSID-Semen PSA	Negative Negative Negative	RSID-Saliva	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
FGHUX	KM Takayama	Positive/Indicative Positive/Indicative	Micro PSA	Negative Negative	[Laboratory] does not perform Saliva testing.	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
FNJ9HY	PTMB HemaTrace	Positive/Indicative Positive/Indicative	AP PSA	Negative Negative	Phadebas Tablet	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
FV2LGX	KM Bluestar OBTI	Positive/Indicative Positive/Indicative	AP RSID-Semen PSA Sperm Hy-Liter PI	Negative Negative Negative Negative	RSID-Saliva	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
FWF8TK	kastle meyer KM RSID Blood	Positive/Indicative Positive/Indicative	acid phosphatase AP RSID Semen	Negative Negative	RSID Saliva	Positive/Indicative
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response	
FXTH7R	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
	<u>Other Tests Performed</u>					
	Urine - Jaffe Test Negative					
FZA2QL	O-tol and HemaTrace	Positive/Indicative	AP	Negative	Lab does not test for Saliva	
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested	
G6M6E3	LMG	Positive/Indicative	AP	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
GBR2ZN	KM HemDirect	Positive/Indicative Positive/Indicative	AP Micro	Negative Negative	RSID Saliva	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
GCMLNM	KM SERATEC HB	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	

Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
GFJDUY	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> RSID-Semen <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	RSID-Saliva <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
GTTNRM	Kastle Meyer <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Acid Phosphatase <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas Amylase <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
	<u>Other Tests Performed</u> Urine - Creatinine <i>Negative</i> Species - HemaTrace <i>Positive/Indicative</i>		
GVXEUV	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
H2GYEM	KM <i>Positive/Indicative</i> HemDirect <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID Saliva <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
H3CJ3M	KM <i>Positive/Indicative</i> SERATEC HB <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
HFJHYZ	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
HQZC42	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
HVCZ4X	Phenolphthalein <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
J2ETAD	Visual exam <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Sperm elution and Micro <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas tube test <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
	<u>Other Tests Performed</u> OBTI <i>Positive/Indicative</i>		
J6WQ4V	KM <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative



Table 1 - Item 2  
 Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
JQFVNY	O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
K4X8BR	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> RSID-Semen <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	RSID-Saliva <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
KBH3F3	<i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<i>Negative</i> <b>Conclusion:</b> Negative	<i>Inconclusive</i> <b>Conclusion:</b> Inconclusive
KLZMZY	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
KMTD7D	LMG <i>Positive/Indicative</i> ABA HEME TRACE <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
KQ99JP	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	Our laboratory does not currently have a validated test for Saliva. <b>Conclusion:</b> Not Tested
KTFL8J	KM <i>Positive/Indicative</i> HemaTrace Card <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
KXT4PZ	Hemastix <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
	<u>Other Tests Performed</u> ALS <i>Negative</i>		
KZXWGD	Molecular Serology <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Molecular Serology <i>Negative</i> <b>Conclusion:</b> Negative	Molecular Serology <i>Inconclusive</i> <b>Conclusion:</b> Inconclusive
L2CKH6	ALS BLUE (420-470 nm) with yellow glasses <i>Positive/Indicative</i> Hemastix (SIEMENS) <i>Positive/Indicative</i> OBTI (Bluestar) <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS blue LIGHT (420-470 nm) with yellow glasses <i>Negative</i> ALS UV (350-380 nm) with UV glasses <i>Negative</i> PSA <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested

Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
L73VR3	PEROXIDASES HEXAGON OBTI	Positive/Indicative Positive/Indicative	PSA Micro	Negative Negative	TEST IS NOT PERFORMED.	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
LECJ72	KM	Positive/Indicative	AP	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
LJAZHV	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
M4VWUF	LMG	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
MCLHXT	TMB HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative	SALigAE	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
MFNJZB	LMG	Positive/Indicative	AP, Micro	Negative	Phadebas	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
MUMKGH	Kastle Meyer HemaTrace	Positive/Indicative Positive/Indicative	PSA	Negative	RSID	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
N44NZ2	Bluestar OBTI	Positive/Indicative	PSA	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
NAT4TF	KM Modified Takayama	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
	<u>Other Tests Performed</u>					
	Species - HemaTrace Positive/Indicative					
	Urine - Creatinine Negative					
NBNQ6J	LMG HemaTrace Visual	Positive/Indicative Positive/Indicative Positive/Indicative	AP	Negative	RSID	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	

Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
NDTENE	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	Modified Takayama	Positive/Indicative				
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
	<u>Other Tests Performed</u>					
	Urine-Creatinine		Negative			
NDVZKV	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	HemaTrace	Positive/Indicative	Micro	Negative		
			PSA	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
NGE2EG	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	Seratec Hb	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
NNYWKR	KM	Positive/Indicative	AP	Negative		
	HemaTrace	Positive/Indicative	PSA	Negative		
			Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
NVFQJU	Visual examination	Positive/Indicative	AP	Negative	Phadebas Press test	Positive/Indicative
	KM	Positive/Indicative				
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
NXYREG	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	Seratec Hb	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
P7B7XX	KM	Positive/Indicative	PSA	Negative		
	HEM DIRECT HEMOGLOBIN	Positive/Indicative	Micro - CHRISTMAS TREE	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
P9FVYK	KM		Micro and p30			
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
PFY2VR	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
			Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
	<u>Other Tests Performed</u>					
	Urine - creatinine		Inconclusive			

Table 1 - Item 2  
 Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
PGCTPU	LMG	Positive/Indicative	AP	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
PTJJQG	LMG	Positive/Indicative	AP	Negative		Saliva testing is not performed at this laboratory.
	HemaTrace	Positive/Indicative	PSA	Negative		
			Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
Q9A4JA	RSID	Positive/Indicative	RSID	Negative		RSID Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
Q9UKZC			AP	Negative		
			Micro	Negative		
			PSA	Negative		
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
QCDPKT	O-tol	Positive/Indicative	PSA	Negative		RSID Positive/Indicative
	HemaTrace	Positive/Indicative	Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
QCVL2U	KM	Positive/Indicative	AP	Negative		Phadebas Positive/Indicative
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
QD4XJW	LMG	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
QD6VPL	KM	Positive/Indicative	AP	Negative		RSID-Saliva Positive/Indicative
	Bluestar OBTI	Positive/Indicative	PSA	Negative		
			RSID-Semen	Negative		
			Sperm Hy-Liter PI	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
QDPC8A	darkening under ALS/IR	Positive/Indicative	no fluorescence under ALS	Negative		CSI does not test for Saliva
	visible red/brown stain	Positive/Indicative				
	LMG	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
QEXQEF	KM	Positive/Indicative	AP	Negative		
	HemaTrace	Positive/Indicative	PSA	Negative		
			Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	

Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
QH2EXB	KM <i>Positive/Indicative</i> Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas Amylase <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
	<u>Other Tests Performed</u> Species Origin - HemaTrace <i>Positive/Indicative</i> Urine - Creatinine <i>Negative</i>		
QMUJMA	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
QV428T	TMB <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Microscopy <i>Negative</i> p30 <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
QZ2CEL	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	Saliva Tests are not performed at the [Laboratory] <b>Conclusion:</b> Not Tested
T4XQ3M	Hemastix <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> RSID <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
TCMMGA	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	Laboratory no longer performs this testing <b>Conclusion:</b> Not Tested
THCZHX	Bluestar Obti <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Bluestar Identi-PSA <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> No Response
TLCCEW	Visual <i>Positive/Indicative</i> fluorescence 445 nm to 510 nm <i>Positive/Indicative</i> fluorescence 420 nm to 470 nm <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	Visual <i>Negative</i> fluorescence 445 nm to 510 nm <i>Negative</i> fluorescence 420 nm to 470 nm <i>Negative</i> Bluestar PSA <i>Negative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested
TM7X8N	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested

Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
TMQJZE	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	SERATEC Hb Direct	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
TQPWLH	Phenolphthalein	Positive/Indicative	Acid Phosphatase	Negative	[Laboratory] does not conduct Saliva examinations	
	HemaTrace	Positive/Indicative	p30	Negative		
			Microscopic	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
TRY22R	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
TW9FKC	HemaTrace	Positive/Indicative	p30	Negative	RSID	Positive/Indicative
			Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
TXPCBQ	KM	Positive/Indicative	AP	Negative		
	HemaTrace	Positive/Indicative	PSA	Negative		
			Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
U7VQFP	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
U87FTA	LMG	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
UJKFP4	Visual	Positive/Indicative	Acid phosphatase test of supernatant	Negative	Phadebas tube test of supernatant (Amylase)	Positive/Indicative
	Kastle-Mayer spot test	Positive/Indicative	Microscopy	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
UMZ8EB	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	Seratec Hemoglobin	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
VGR42G	KM	Positive/Indicative	AP	Negative	Not competent / authorised to test for Saliva	
	Bluestar OBTI	Positive/Indicative	RSID-Semen	Negative		
			PSA	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	

Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
VHJYQB	TMB <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
<u>Other Tests Performed</u> Blood - White light Visual examination <i>Positive/Indicative</i> Semen - White light Visual examination <i>Negative</i> Saliva - White light Visual examination <i>Negative</i>			
VMDWF3	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
VXWD89	Kastle Meyer test <i>Positive/Indicative</i> ABAcad HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> No Response	<b>Conclusion:</b> No Response
W4ZMZM	O-tol <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	No current test in laboratory for Semen <b>Conclusion:</b> Not Tested	No current test in laboratory for Saliva <b>Conclusion:</b> Not Tested
WEGF6N	O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Inconclusive</i> Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
WL47RJ	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	ALS <i>Negative</i> Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
WP2REF	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> RSID-Semen <i>Negative</i> PSA <i>Negative</i> Sperm Hy-Liter PI <i>Negative</i> <b>Conclusion:</b> Negative	RSID-Saliva <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
WWQ6HP	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
WZYRRK	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> p30 <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative

Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
X32Q3D	KM	Positive/Indicative	Micro	Negative	Currently, the lab does not have a validated method for testing for Saliva	
		<b>Conclusion:</b> Positive/Indicative	PSA	Negative	<b>Conclusion:</b> Not Tested	
X7HDGH	KM	Positive/Indicative	AP	Negative	Phadebas Positive/Indicative	
	HemaTrace	Positive/Indicative	Micro	Negative		
		<b>Conclusion:</b> Positive/Indicative	PSA	Negative	<b>Conclusion:</b> Positive/Indicative	
XBE2V4	Leucomalachite green Positive/Indicative					
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
XEYPW4	KM	Positive/Indicative	AP	Negative	RSID Positive/Indicative	
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
XPVZZ6	LMG	Positive/Indicative	AP	Negative	No longer performed	
	HemaTrace	Positive/Indicative	PSA	Negative		
		<b>Conclusion:</b> Positive/Indicative	Micro	Negative	<b>Conclusion:</b> Not Tested	
XTUBW6	LMG	Positive/Indicative	AP	Negative	No longer performed	
	HemaTrace	Positive/Indicative	PSA	Negative		
		<b>Conclusion:</b> Positive/Indicative	Micro	Negative	<b>Conclusion:</b> Not Tested	
Y48XHB	Presumptive test- LMG					
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
Y7YU86	RSID Blood Positive/Indicative		RSID Semen Negative		RSID Saliva Positive/Indicative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
Y94GNF	PTMB	Positive/Indicative				
	Bluestar	Positive/Indicative				
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
YDXBLG	KM	Positive/Indicative	AP	Negative	Phadebas Positive/Indicative	
	HemaTrace	Positive/Indicative	Micro	Negative	ALS Negative	
		<b>Conclusion:</b> Positive/Indicative	Prostate specific antigen	Negative		
			ALS	Negative	<b>Conclusion:</b> Positive/Indicative	
			<b>Conclusion:</b> Negative			



Table 1 - Item 2  
 Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
YEUW9F	TMB <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	SALigAE <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
YNE396	Bluestar Hexagon OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Seratec PSA Semiquant <i>Negative</i> <b>Conclusion:</b> Negative	Seratec Amylase Test <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
YNVC27	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
YRECCF	O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
	<u>Other Tests Performed</u> Urine- RSID <i>Negative</i>		
YTQ493	KM <i>Positive/Indicative</i> Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas Amylase <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
	<u>Other Tests Performed</u> Species Origin - HemaTrace <i>Positive/Indicative</i> Urine - Creatinine <i>Negative</i>		
Z7EAZC	KM, HEXAGON OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
Z9YG2K	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
ZADNKG	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> RSID-Semen <i>Negative</i> Sperm Hy-Liter PI <i>Negative</i> <b>Conclusion:</b> Negative	RSID-Saliva <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
ZRC2V3	Hexagon OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	PSA-RSID <i>Negative</i> <b>Conclusion:</b> Negative	RSID Saliva <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
ZYAT4W	Phenolphthalein <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested

### Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
ZZ3LYF	KM	Positive/Indicative	AP	Negative	Not competent / authorised to test for Saliva.	
	Bluestar OBTI	Positive/Indicative	RSID-Semen	Negative		
			PSA	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	

#### Response Summary for Table 1 - Item 2: Four pink plastic handle swabs containing Blood & Saliva Mixture

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	165	1	102
Negative	0	150	0
Inconclusive	0	0	2
Not Tested (NT)	2	17	68
No Response	11	10	6
<b>Total</b>	<b>178</b>	<b>178</b>	<b>178</b>

# Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
2PK689	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
2PZNR7	Visual Examination <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID Semen <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID Saliva <i>Negative</i> <b>Conclusion:</b> Negative
2XF3XH	<b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
38Y4RX	LMG <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
3996DV	KM <i>Negative</i> <b>Conclusion:</b> No Response	AP <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	Phadebas <i>Negative</i> <b>Conclusion:</b> No Response
3H2GPC	PTMB <i>Negative</i> Bluestar <i>Negative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
3RNEGH	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
3URBZK	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Spermatozoa Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
3XDWT9	Visual examination did not detect red or brown stains therefore Blood testing was not performed <b>Conclusion:</b> Not Tested	Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	This laboratory does not have a current validated method for Saliva testing <b>Conclusion:</b> Not Tested
3YL8PC	Not indicated <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not indicated <b>Conclusion:</b> Not Tested

**Table 1 - Item 3**  
**Evidence - Four blue plastic swabs containing Semen**

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
4DGTNC	TMB HemaTrace	Negative Negative	AP Micro PSA	Positive/Indicative Positive/Indicative Positive/Indicative	SALigAE	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
4HDDVG	KM	Negative	AP Micro	Positive/Indicative Positive/Indicative	Phadebas	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
4JLPBX	Phenolphthalein	Negative	Microscopic exam p30	Positive/Indicative Positive/Indicative	Phadebas	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
4N3XXW	KM	Negative	Brentamine (AP) Human seminal protein (p30) Microscopic	Positive/Indicative Positive/Indicative Positive/Indicative	It does not apply according to the scenario	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
4R39RC	TMB	Negative	AP p30 Micro	Positive/Indicative Positive/Indicative Positive/Indicative	SALigAE	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
63V48H	KM SERATEC HemDirect	Negative Negative	AP SERATEC PSA	Positive/Indicative Positive/Indicative	SERATEC AmylaseTest	Negative
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> Negative	
69JRZ3	Visual	Negative	AP RSID Micro	Positive/Indicative Positive/Indicative Positive/Indicative	RSID	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
6BBMPW	benzidine	Negative	acid phosphatase solution PSA kit Microscope	Positive/Indicative Positive/Indicative Positive/Indicative	It does not apply according to the scenario	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
6HUTL4	Visual Examination	Negative	AP RSID Micro	Positive/Indicative Positive/Indicative Positive/Indicative	RSID	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	

Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
6WTV3B	LMG	Negative	AP Micro PSA	Positive/Indicative Positive/Indicative Positive/Indicative	Amylase Diffusion	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
6YVZDB	Hemastix	Negative	AP ABAcad p30 Sperm Search - Micro	Positive/Indicative Positive/Indicative Positive/Indicative	SALlgAE	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
79FQD8	KM	Negative				
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
7C4X2B	KM	Negative	AP PSA Micro	Positive/Indicative Positive/Indicative Positive/Indicative		
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
7FYPA8	KM Micro	Negative Negative	AP Micro	Positive/Indicative Positive/Indicative	RSID	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
7G9YXY	KM	Negative	AP Micro	Positive/Indicative Positive/Indicative	Phadebas	Negative
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
7JZXCW	No brownish traces		PSA	Positive/Indicative	Amylase - Method CNPG3	Negative
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
7LBQEC			Microscopy	Positive/Indicative		
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
7UGHWT			AP Micro	Positive/Indicative Positive/Indicative	RADIAL Diffusion	Negative
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
83AAGT	KM	Negative	AP Micro PSA	Positive/Indicative Positive/Indicative Positive/Indicative		
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
8AV3XY	No visible staining on swabs		PSA Micro	Positive/Indicative Positive/Indicative	RSID	Negative
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response	

Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
8CE398	KM <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
8GNHHW	TMB <i>Negative</i> <b>Conclusion:</b> Negative	PSA <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Seratec Amylase <i>Negative</i> <b>Conclusion:</b> Negative
8Z7GPE	No Blood-like stains observed <b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
9BFE63	Visual <i>Negative</i> LMG <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID Semen <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID Saliva <i>Negative</i> <b>Conclusion:</b> Negative
9EJ86T	No apparent Bloodstains present. <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> Seratec <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Case scenario does not indicate the presence of Saliva. <b>Conclusion:</b> Not Tested
9K9PJ4	Micro <i>Negative</i> KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
9QA7Q4	Not tested based on case scenario and/or no stain observed. <b>Conclusion:</b> Not Tested	RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Generally not tested unless specified. <b>Conclusion:</b> Not Tested
9QVF4E	<b>Conclusion:</b> Negative <u>Other Tests Performed</u> ALS <i>Positive/Indicative</i>	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
9UVNAA	<b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
9VPA2F	<b>Conclusion:</b> Not Tested	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
A7LW93	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative

Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
ACUZGP	<b>Conclusion:</b> Not Tested	PSA Semiquant <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
AEDFMB	LMG <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
AF8ZDG	Visual appearance <i>Negative</i> Bluestar forensic tablets <i>Negative</i> <b>Conclusion:</b> Negative	ALS <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Our lab does not perform screening tests for Saliva <b>Conclusion:</b> Not Tested
ALCYQR	No Visual indications of Blood to test. <b>Conclusion:</b> Not Tested	ALS, AP, Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS, RSID Saliva <i>Negative</i> <b>Conclusion:</b> Negative
AQ64HB	No Blood-like discoloration noted. <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas Amylase Test <i>Negative</i> <b>Conclusion:</b> Negative
ARYUNP	<b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
B7XV2Q	No Red Brown Staining <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Radial Diffusion <i>Negative</i> <b>Conclusion:</b> Negative
BLQQBW	Visual <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID-Saliva <i>Negative</i> <b>Conclusion:</b> Negative
BNVB33	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
BQJU4Z	Visual <i>Negative</i> LMG <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
BTQ6Z4	Phenolphthalein <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested

Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
BU9Y2D			PSA Positive/Indicative Micro Positive/Indicative			
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
C84ZT6			AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative			
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
CDC42R	RSID Blood Negative		RSID Semen Positive/Indicative		RSID Saliva Negative	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
CGD64B	HEXAGON OBTI Negative		PSA Positive/Indicative Micro Positive/Indicative			
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
CKDD9K	KM Negative		AP Positive/Indicative		Phadebas Negative	
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response	
CNWLH7			Acid phosphatase test Positive/Indicative Microscopy Positive/Indicative		Phadebas paper press test Negative	
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
CQFMDR	KM Negative		AP Positive/Indicative Micro Positive/Indicative		Phadebas Negative	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
CYTE6B	LMG Negative OBTI Negative		AP Positive/Indicative Micro Positive/Indicative			
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested	
DH7FGP	No apparent Bloodstains present.		AP Positive/Indicative			
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
DKQAC4	FOB Negative		RSID Positive/Indicative		RSID Negative	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
DQFZL7	Visual appearance of Blood Negative KM Negative		AP Positive/Indicative Micro (Spermatozoa) Positive/Indicative		Phadebas (press test) Negative	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	



Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
DWKLQ3	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Sperm cells - Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
<u>Other Tests Performed</u> Nucleated epithelial cells - Micro <i>Negative</i>			
DZQ388	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Phadebas paper press test <i>Negative</i>  <b>Conclusion:</b> Negative
EBRRPL	LMG <i>Negative</i> RSID-Blood <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	RSID-Saliva <i>Negative</i>  <b>Conclusion:</b> Negative
ED299W	Visual <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID - Semen <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	RSID-Saliva <i>Negative</i>  <b>Conclusion:</b> Negative
EGQD7A	HEXAGON OBTI IMMUNOCHROMATOGRAPHIC RAPID TEST <i>Negative</i> PEROXIDASE TEST <i>Negative</i>  <b>Conclusion:</b> Negative	PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	TEST NOT PERFORMED  <b>Conclusion:</b> Not Tested
EQF6UM	No apparent Bloodstains present.  <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> Seratec <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Not tested due to case scenario.  <b>Conclusion:</b> Not Tested
FBHAXM	KM <i>Negative</i>  <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i>  <b>Conclusion:</b> Negative
FDLCGW	Not tested based on case scenario and/or no stain observed  <b>Conclusion:</b> Not Tested	RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	Generally not tested unless specified  <b>Conclusion:</b> Not Tested
FGHUXX	KM <i>Negative</i>  <b>Conclusion:</b> Negative	Micro <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative	[Laboratory] does not perform Saliva testing.  <b>Conclusion:</b> Not Tested

Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
FNJ9HY	not indicated <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	not indicated <b>Conclusion:</b> Not Tested
FV2LGX	No stain observed <b>Conclusion:</b> Not Tested	RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Generally not tested unless specified <b>Conclusion:</b> Not Tested
FWF8TK	kastle meyer KM <i>Negative</i> <b>Conclusion:</b> No Response	acid phosphatase AP <i>Positive/Indicative</i> RSID Semen <i>Positive/Indicative</i> miroscopy - christmas tree staining <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	RSID Saliva <i>Negative</i> <b>Conclusion:</b> No Response
FXTH7R	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
FZA2QL	<b>Conclusion:</b> Not Tested	AP and p30 <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	Lab does not test for Saliva <b>Conclusion:</b> Not Tested
G6M6E3	LMG <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
GBR2ZN	No red/brown stain present <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	NT based on case scenario <b>Conclusion:</b> Not Tested
GCMLNM	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
GFJDUY	No stain observed <b>Conclusion:</b> Not Tested	RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Generally not tested unless specified <b>Conclusion:</b> Not Tested
GTTNRM	Kastle Meyer <i>Negative</i> <b>Conclusion:</b> Negative	Acid Phosphatase <i>Positive/Indicative</i> Microscopic <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas Amylase <i>Negative</i> <b>Conclusion:</b> Negative
GVXEUV	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative

Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
H2GYEM	no reddish brown stains observed on swabs <b>Conclusion:</b> Not Tested	AP Positive/Indicative Micro Positive/Indicative <b>Conclusion:</b> Positive/Indicative	RSID Saliva Negative <b>Conclusion:</b> Negative
H3CJ3M	KM Negative <b>Conclusion:</b> Negative	AP Positive/Indicative Micro Positive/Indicative <b>Conclusion:</b> Positive/Indicative	Phadebas Negative <b>Conclusion:</b> Negative
HFJHYZ	KM Negative <b>Conclusion:</b> Negative	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative <b>Conclusion:</b> Positive/Indicative	Phadebas Negative <b>Conclusion:</b> Negative
HQZC42	LMG Negative <b>Conclusion:</b> Negative	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
HVCZ4X	Phenolphthalein Negative <b>Conclusion:</b> Negative	AP Positive/Indicative PSA Positive/Indicative <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
J2ETAD	Visual exam Negative KM Negative <b>Conclusion:</b> Negative	AP Positive/Indicative Sperm elution and Micro Positive/Indicative <b>Conclusion:</b> Positive/Indicative	Phadebas tube test Negative <b>Conclusion:</b> Negative
J6WQ4V	KM Negative Micro Negative <b>Conclusion:</b> Negative	AP Positive/Indicative Micro Positive/Indicative <b>Conclusion:</b> Positive/Indicative	RSID Negative <b>Conclusion:</b> Negative
JQFVNY	O-tol Negative <b>Conclusion:</b> Negative	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative <b>Conclusion:</b> Positive/Indicative	RSID Negative <b>Conclusion:</b> Negative
K4X8BR	Not tested based on case scenario and/or no stain observed. <b>Conclusion:</b> Not Tested	RSID-Semen Positive/Indicative PSA Positive/Indicative <b>Conclusion:</b> Positive/Indicative	Generally not tested unless specified. <b>Conclusion:</b> Not Tested
KBH3F3	Negative <b>Conclusion:</b> Negative	Positive/Indicative <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested

Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
KLZMZY			AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative			
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
KMTD7D	LMG Negative		AP Positive/Indicative p30 Positive/Indicative			
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
KQ99JP	KM Negative		Micro Positive/Indicative		Our laboratory does not currently have a validated test for Saliva.	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
KTFL8J	No testing for Blood performed. The scenario as described and the appearance of the swabs (no red-brown staining) contraindicated testing for Blood.		PSA Positive/Indicative Micro Positive/Indicative		No testing for Saliva performed. The scenario as described and the appearance of the swabs contraindicated testing for Saliva.	
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
KXT4PZ			AP Positive/Indicative			
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
	<u>Other Tests Performed</u>					
	ALS Positive/Indicative					
KZXWGD	Molecular Serology Negative		Molecular Serology Positive/Indicative		Molecular Serology Negative	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
L2CKH6	ALS BLUE (420-470 nm) with yellow glasses Negative Hemastix (SIEMENS) Negative		ALS BLUE (420-470 nm) with yellow glasses Positive/Indicative ALS UV (350-380 nm) with UV glasses Positive/Indicative PSA Positive/Indicative AP Positive/Indicative			
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
L73VR3			PSA Positive/Indicative Micro Positive/Indicative			
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
LECJ72	KM Negative		AP Negative			
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
LJAZHV			AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative			
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	

Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
M4VWUF	LMG (no Blood-like stains observed) Negative					
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
MCLHXT	no red/brown staining		AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative		SALigAE Negative	
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
MFNJZB	No Visible Bloodlike staining		AP, Micro Positive/Indicative		Not examined given case circumstances	
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
MUMKGH	reddish-brown staining not observed, so no KM test performed		PSA Positive/Indicative Micro Positive/Indicative		RSID Negative	
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
N44NZ2	Bluestar OBTI Negative		PSA Positive/Indicative			
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
NAT4TF	KM Negative		AP Positive/Indicative Micro Positive/Indicative		Phadebas Negative	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
NBNQ6J	Visual Negative		AP Positive/Indicative RSID Positive/Indicative Micro Positive/Indicative		RSID Negative	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
NDTENE	KM Negative		AP Positive/Indicative Xmas/Micro Positive/Indicative		Phadebas Negative	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
NDVZKV	KM Negative		AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative		Phadebas Negative	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
NGE2EG	KM Negative		AP Positive/Indicative Spermatozoa Identification (Microscopic) Positive/Indicative		Phadebas Negative	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
NNYWKR	KM Negative		AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative			
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	

Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
NVFQJU	Visual examination	Negative	AP	Positive/Indicative	Phadebas Press Test	Negative
		<b>Conclusion:</b> Negative	Micro	Positive/Indicative		<b>Conclusion:</b> Negative
NXYREG	KM	Negative	AP	Positive/Indicative	Phadebas	Negative
		<b>Conclusion:</b> Negative	Micro	Positive/Indicative		<b>Conclusion:</b> Negative
P7B7XX	KM	Negative	PSA	Positive/Indicative		
	HEM DIRECT HEMOGLOBIN	Negative	Micro - CHRISTMAS TREE	Positive/Indicative		<b>Conclusion:</b> Not Tested
P9FVYK			Micro			
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
PFY2VR	KM	Negative	AP	Positive/Indicative		
		<b>Conclusion:</b> Negative	Micro	Positive/Indicative		<b>Conclusion:</b> Not Tested
PGCTPU	LMG	Negative	AP	Positive/Indicative		
		<b>Conclusion:</b> Negative	Micro	Positive/Indicative	PSA	Positive/Indicative
				<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
PTJJQG			AP	Positive/Indicative	Saliva testing is not performed at this laboratory.	
		<b>Conclusion:</b> Not Tested	PSA	Positive/Indicative		<b>Conclusion:</b> Not Tested
			Micro	Positive/Indicative		
Q9A4JA	No RBS present		RSID	Positive/Indicative	RSID	Negative
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
Q9UKZC			AP	Positive/Indicative		
		<b>Conclusion:</b> Not Tested	Micro	Positive/Indicative		<b>Conclusion:</b> Not Tested
			PSA	Positive/Indicative		
QCDPKT	No visible staining		PSA	Positive/Indicative	Not relevant to scenario	
		<b>Conclusion:</b> Not Tested	Micro	Positive/Indicative		<b>Conclusion:</b> Not Tested

Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
QCVL2U			AP Positive/Indicative Micro Exam Positive/Indicative p30 Positive/Indicative			
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
QD4XJW						
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
QD6VPL	Not tested based on case scenario and/or no stain observed.		PSA Positive/Indicative RSID-Semen Positive/Indicative		Generally not tested unless specified.	
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
QDPC8A	No visible appearance of Blood or under ALS/IR Negative		ALS Positive/Indicative AP Positive/Indicative		CSI does not test for Saliva	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
QEXQEF	No Visual characteristics of Blood		AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative			
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
QH2EXB	KM Negative		AP Positive/Indicative Micro Positive/Indicative		Phadebas Amylase Negative	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	
QMUJMA	KM Negative		AP Positive/Indicative			
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
QV428T			Microscopy Positive/Indicative			
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
QZ2CEL	Visual, based on color of swabs testing was not warranted.		AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative		Saliva Tests are not performed at the [Laboratory]	
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
T4XQ3M			AP Positive/Indicative Micro Positive/Indicative RSID Positive/Indicative			
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
TCMMGA	LMG Negative		AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative		Laboratory no longer performs this testing	
	<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	

Table 1 - Item 3  
 Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
THCZHx	Bluestar Obti <i>Negative</i> <b>Conclusion:</b> Negative	Bluestar Identi-PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> No Response
TLCCEW	Visual <i>Negative</i> fluorescence 445 nm to 510 nm <i>Negative</i> fluorescence 420 nm to 470 nm <i>Negative</i> <b>Conclusion:</b> No Response	Visual <i>Negative</i> fluorescence 445 nm to 510 nm <i>Positive/Indicative</i> fluorescence 420 nm to 470 nm <i>Positive/Indicative</i> Bluestar PSA <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested
TM7X8N	<b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
TMQJZE	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not required based on analysis of sample from the scenario <b>Conclusion:</b> Not Tested
TQPWLH	Phenolphthalein <i>Negative</i> HemaTrace <i>Negative</i> <b>Conclusion:</b> Negative	Acid phosphatase <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> Microscopic <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	[Laboratory] does not conduct Saliva examinations <b>Conclusion:</b> Not Tested
TRY22R	no red/brown staining observed <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> Microscopic <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	No indication in the narrative that Saliva was suspected in the stain <b>Conclusion:</b> Not Tested
TW9FKC	No Reddish/Brown Staining Observed <b>Conclusion:</b> Not Tested	p30 <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
TXPCBQ	NT not applicable to case scenario <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested
U7VQFP	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
U87FTA	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested



Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
UJKFP4	Visual <i>Negative</i> Kastle-Mayer test of supernatant <i>Negative</i> <b>Conclusion:</b> Negative	Acid phosphatase test of supernatant <i>Positive/Indicative</i> Microscopy <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas test of supernatant (Amylase) <i>Negative</i> <b>Conclusion:</b> Negative
UMZ8EB	KM <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> Microscopic <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
VGR42G	Not tested based on case scenario and/or no stain observed <b>Conclusion:</b> Not Tested	RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Generally not tested unless specified <b>Conclusion:</b> Not Tested
VHJYQB	TMB <i>Negative</i> <b>Conclusion:</b> Negative	AP <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Blood - White light Visual examination <i>Negative</i> Semen - White light Visual examination <i>Inconclusive</i> Saliva - White light Visual examination <i>Inconclusive</i>		
VMDWF3	no red-brown staining observed <b>Conclusion:</b> Not Tested	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
VXWD89	<b>Conclusion:</b> No Response	brentamine spot test <i>Positive/Indicative</i> ABAcad p30 <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> No Response
W4ZMZM	O-tol <i>Negative</i> <b>Conclusion:</b> Negative	No current test in laboratory for Semen <b>Conclusion:</b> Not Tested	No current test in laboratory for Saliva <b>Conclusion:</b> Not Tested
WEGF6N	Synopsis does not suggest the presence of Blood. <b>Conclusion:</b> Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Synopsis does not suggest the presence of Saliva. <b>Conclusion:</b> Not Tested
WL47RJ	KM <i>Negative</i> <b>Conclusion:</b> Negative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Positive/Indicative</i> Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
WP2REF	No stain observed <b>Conclusion:</b> Not Tested	RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Generally not tested unless specified <b>Conclusion:</b> Not Tested

Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
WVQ6HP	LMG	Negative	AP	Positive/Indicative		
		<b>Conclusion:</b> Negative	Micro	Positive/Indicative		
			PSA	Positive/Indicative		
				<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
WZYRRK			AP	Positive/Indicative		
		<b>Conclusion:</b> Not Tested	Micro	Positive/Indicative		
				<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
X32Q3D	Visual examination did not detect brown stains		Micro	Positive/Indicative	Currently, the lab does not have a validated method for testing for Saliva	
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
X7HDGH	KM	Negative	AP	Positive/Indicative	Phadebas	Negative
		<b>Conclusion:</b> Negative	Micro	Positive/Indicative		
			PSA	Positive/Indicative		
				<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
XBE2V4						
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested
XEYPW4	No apparent Blood stains present		AP	Positive/Indicative	Not tested based on case scenario	
		<b>Conclusion:</b> Not Tested	Seratec	Positive/Indicative		
			Micro	Positive/Indicative		
				<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
XPVZZ6	Scenario does not indicate need for test		AP	Positive/Indicative	No longer performed	
		<b>Conclusion:</b> Not Tested	PSA	Positive/Indicative		
			Micro	Positive/Indicative		
				<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
XTUBW6	Scenario does not indicate need for test		AP	Positive/Indicative	No longer performed	
		<b>Conclusion:</b> Not Tested	PSA	Positive/Indicative		
			Micro	Positive/Indicative		
				<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
Y48XHB	No visible red/brown stains.					
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested
Y7YU86	RSID Blood	Negative	RSID Semen	Positive/Indicative	RSID Saliva	Negative
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
Y94GNF	PTMB	Negative				
	Bluestar	Negative				
		<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested

Table 1 - Item 3  
Evidence - Four blue plastic swabs containing Semen

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
YDXBLG	KM	Negative	AP	Positive/Indicative	Phadebas	Negative
			Micro	Positive/Indicative	ALS	Positive/Indicative
			Prostate specific antigen	Positive/Indicative		
			ALS	Positive/Indicative		
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
YEUW9F	TMB	Negative	AP	Positive/Indicative	SALigAE	Negative
			PSA	Positive/Indicative		
			Micro	Positive/Indicative		
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
YNE396	Bluestar Hexagon OBTI	Negative	Seratec PSA Semiquant	Positive/Indicative	Seratec Amylase Test	Negative
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
YNVC27	KM	Negative	AP	Positive/Indicative	RSID	Negative
			PSA	Positive/Indicative		
			Micro	Positive/Indicative		
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
YRECCF			AP	Positive/Indicative		
			PSA	Positive/Indicative		
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
YTQ493	KM	Negative	AP	Positive/Indicative	Phadebas Amylase	Negative
			Micro	Positive/Indicative		
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
Z7EAZC	LAB Policy		AP	Positive/Indicative	Phadebas	Negative
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
Z9YG2K	LMG	Negative	AP	Positive/Indicative		
			Micro	Positive/Indicative		
			PSA	Positive/Indicative		
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested
ZADNKG	Not tested because no stain observed.		PSA	Positive/Indicative	Generally not tested unless specified.	
			RSID-Semen	Positive/Indicative		
	<b>Conclusion:</b> Not Tested			<b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	
ZRC2V3	Hexagon	Negative	PSA RSID	Positive/Indicative	Amilase	Negative
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative
ZYAT4W	Phenolphthalein	Negative				
		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested

**Table 1 - Item 3**  
*Evidence - Four blue plastic swabs containing Semen*

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
ZZ3LYF	Not tested based on case scenario and/or no stain observed.		RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i>		Generally not tested unless specified.	
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	

**Response Summary for Table 1 - Item 3: Four blue plastic swabs containing Semen**

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	0	154	1
Negative	98	1	74
Inconclusive	0	0	0
Not Tested (NT)	71	15	97
No Response	9	8	6
<b>Total</b>	178	178	178

# Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 4  
Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
2PK689	Img <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
2PZNR7	Visual Examination <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID Saliva <i>Negative</i> <b>Conclusion:</b> Negative
2XF3XH	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
38Y4RX	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
3996DV	KM <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	AP <i>Negative</i> <b>Conclusion:</b> No Response	Phadebas <i>Negative</i> <b>Conclusion:</b> No Response
3H2GPC	PTMB <i>Positive/Indicative</i> Bluestar <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
3RNEGH	KM <i>Positive/Indicative</i> Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
3URBZK	KM <i>Positive/Indicative</i> Seratec Hb <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
3XDWT9	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	This laboratory does not have a current validated method for Saliva testing <b>Conclusion:</b> Not Tested
3YL8PC	PTMB <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not indicated <b>Conclusion:</b> Not Tested	Not indicated <b>Conclusion:</b> Not Tested
4DGTNC	TMB <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	SALIgAE <i>Negative</i> <b>Conclusion:</b> Negative

Table 1 - Item 4  
Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
4HDDVG	KM Modified Takayama <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative	AP <b>Conclusion:</b> Negative	Negative	Phadebas <b>Conclusion:</b> Negative	Negative
4JLPBX	Phenolphthalein <b>Conclusion:</b> Positive/Indicative	Positive/Indicative	Microscopic exam p30 <b>Conclusion:</b> Negative	Negative Negative	Phadebas <b>Conclusion:</b> Negative	Negative
4N3XXW	KM HemaTrace <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative	Brentamine (AP) Human seminal protein (p30) Microscopic <b>Conclusion:</b> Negative	Negative Negative Negative	Human Amylase (alpha-amilosa) Negative <b>Conclusion:</b> Negative	Negative
4R39RC	TMB HemaTrace <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative	AP <b>Conclusion:</b> Negative	Negative	SALigAE <b>Conclusion:</b> Negative	Negative
63V48H	KM SERATEC HemDirect <b>Conclusion:</b> No Response	Positive/Indicative Positive/Indicative	AP SERATEC PSA <b>Conclusion:</b> No Response	Negative Negative	SERATEC AmylaseTest <b>Conclusion:</b> Negative	Negative
69JRZ3	Visual LMG HemaTrace <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative Positive/Indicative	AP <b>Conclusion:</b> Negative	Negative	RSID <b>Conclusion:</b> Negative	Negative
6BBMPW	benzidine Hem Direct kit <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative	acid phosphatase solution <b>Conclusion:</b> Negative	Negative	<b>Conclusion:</b> Not Tested	
6HUTL4	Visual Examination LMG HemaTrace <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative Positive/Indicative	AP <b>Conclusion:</b> Negative	Negative	RSID <b>Conclusion:</b> Negative	Negative
6WTV3B	LMG <b>Conclusion:</b> Positive/Indicative	Positive/Indicative	AP Micro PSA <b>Conclusion:</b> Negative	Negative Negative Negative	Amylase Diffusion <b>Conclusion:</b> Negative	Negative
6YVZDB	Hemastix ABAcad HemaTrace <b>Conclusion:</b> Positive/Indicative	Positive/Indicative Positive/Indicative	AP <b>Conclusion:</b> Negative	Negative	SALigAE <b>Conclusion:</b> Negative	Negative

Table 1 - Item 4  
 Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
79FQD8	KM HemaTrace	Positive/Indicative Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
7C4X2B	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
7FYPA8	KM Micro	Positive/Indicative Positive/Indicative	AP Micro	Negative Negative	RSID	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
7G9YXY	KM Seratec HB	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
7JZXCW	TMB Hexagon OBTI	Positive/Indicative Positive/Indicative	PSA	Negative	Amylase - Method CNPG3	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
7LBQEC	Hemastix HemaTrace	Positive/Indicative Positive/Indicative	Microscopy AP	Negative Negative	RSID-Saliva	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
7UGHWT	KM	Positive/Indicative	AP Micro	Negative Negative	RADIAL Diffusion	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
83AAGT	KM HemaTrace	Positive/Indicative Positive/Indicative	AP Micro PSA	Negative Negative Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
8AV3XY	KM HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative	RSID	Negative
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response	
8CE398	KM HemaTrace	Positive/Indicative Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
8GNHHW	OBTI	Positive/Indicative	PSA	Negative	Seratec Amylase	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	

Table 1 - Item 4  
Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
8Z7GPE	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
9BFE63	Visual <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID Saliva <i>Negative</i> <b>Conclusion:</b> Negative
9EJ86T	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Case scenario does not indicate the presence of Saliva. <b>Conclusion:</b> Not Tested
9K9PJ4	Micro <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
9QA7Q4	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not tested based on case scenario. <b>Conclusion:</b> Not Tested	Generally not tested unless specified. <b>Conclusion:</b> Not Tested
9QVF4E	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
9UVNAA	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
9VPA2F	<b>Conclusion:</b> Not Tested	RSID <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
A7LW93	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
ACUZGP	Phenolphthalein <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
	<u>Other Tests Performed</u> Hemedirect <i>Positive/Indicative</i>		
AEDFMB	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested



Table 1 - Item 4  
Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
AF8ZDG	Visual appearance Bluestar forensic tablets	Positive/Indicative Positive/Indicative	ALS Micro PSA	Negative Negative Negative		Our lab does not perform screening tests for Saliva
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative		<b>Conclusion:</b> Not Tested
ALCYQR	Visual, KM, HemDirect	Positive/Indicative	ALS, AP	Negative		ALS, RSID Saliva
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative		<b>Conclusion:</b> Negative
AQ64HB	KM	Positive/Indicative	AP	Negative		Phadebas Amylase Test
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative		<b>Conclusion:</b> Negative
ARYUNP	KM HemaTrace	Positive/Indicative Positive/Indicative	AP Micro PSA	Negative Negative Negative		
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative		<b>Conclusion:</b> Not Tested
B7XV2Q	Phen	Positive/Indicative	AP Micro	Negative Negative		Radial Diffusion
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative		<b>Conclusion:</b> Negative
BLQQBW	Visual LMG HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	AP	Negative		RSID-Saliva
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative		<b>Conclusion:</b> Negative
BNVB33	KM	Positive/Indicative	AP Micro	Negative Negative		RSID
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative		<b>Conclusion:</b> Negative
BQJU4Z	Visual LMG HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	AP	Negative		RSID
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative		<b>Conclusion:</b> Negative
BTQ6Z4	Phenolphthalein HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative		
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative		<b>Conclusion:</b> Not Tested
BU9Y2D	PEROXIDASES HEXAGON OBTI	Positive/Indicative Positive/Indicative	PSA Micro	Negative Negative		
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative		<b>Conclusion:</b> Not Tested

Table 1 - Item 4  
Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
C84ZT6	KM	Positive/Indicative	AP	Negative		
	HemaTrace	Positive/Indicative	PSA	Negative		
			Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
CDC42R	RSID Blood	Positive/Indicative	RSID Semen	Negative	RSID Saliva	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
CGD64B	HEXAGON OBTI	Positive/Indicative	PSA	Negative		
			Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
CKDD9K	KM	Positive/Indicative	AP	Negative	Phadebas	Negative
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response	
CNWLH7	Appearance	Positive/Indicative	Acid phosphatase test	Negative	Phadebas paper press test	Negative
	Kastle Meyer test	Positive/Indicative				
	HemaTrace test	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
CQFMDR	KM	Positive/Indicative	AP	Negative	Phadebas	Negative
	SERATEC Hb	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
CYTE6B	LMG	Positive/Indicative	AP	Inconclusive		
	OBTI	Positive/Indicative	Micro	Negative		
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested	
DH7FGP	KM	Positive/Indicative	AP	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
DKQAC4	FOB	Positive/Indicative	RSID	Negative	RSID	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
DQFZL7	Visual appearance of Blood	Positive/Indicative	AP	Negative	Phadebas (press test)	Negative
	KM	Positive/Indicative				
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
DWKLQ3	KM	Positive/Indicative	AP	Negative	Phadebas	Inconclusive
	HemaTrace	Positive/Indicative	Sperm cells - Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Inconclusive	
	<u>Other Tests Performed</u>					
	Nucleated epithelial cells - Micro Positive/Indicative					

Table 1 - Item 4  
Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
DZQ388	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas paper press test <i>Negative</i> <b>Conclusion:</b> Negative
EBRRPL	LMG <i>Positive/Indicative</i> RSID-Blood <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> RSID-Semen <i>Negative</i> <b>Conclusion:</b> Negative	RSID-Saliva <i>Negative</i> <b>Conclusion:</b> Negative
ED299W	Visual <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID-Saliva <i>Negative</i> <b>Conclusion:</b> Negative
EGQD7A	HEXAGON OBTI IMMUNOCHROMATOGRAPHIC RAPID TEST <i>Positive/Indicative</i> PEROXIDASE TEST <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	PSA <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	TEST NOT PERFORMED <b>Conclusion:</b> Not Tested
EQF6UM	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Not tested due to case scenario. <b>Conclusion:</b> Not Tested
FBHAXM	KM <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
FDLCGW	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not tested based on case scenario <b>Conclusion:</b> Not Tested	Generally not tested unless specified <b>Conclusion:</b> Not Tested
FGHUXX	KM <i>Positive/Indicative</i> Takayma <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	[Laboratory] does not perform Saliva testing. <b>Conclusion:</b> Not Tested
FNJ9HY	PTMB <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	not indicated <b>Conclusion:</b> Not Tested	not indicated <b>Conclusion:</b> Not Tested
FV2LGX	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not tested based on case scenario <b>Conclusion:</b> Not Tested	Generally not tested unless specified <b>Conclusion:</b> Not Tested
FWF8TK	kastle meyer KM <i>Positive/Indicative</i> RSID Blood <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	acid phosphatase AP <i>Negative</i> RSID Semen <i>Negative</i> <b>Conclusion:</b> No Response	RSID Saliva <i>Negative</i> <b>Conclusion:</b> No Response

Table 1 - Item 4  
 Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
FXTH7R	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
FZA2QL	O-tol and HemaTrace <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested
G6M6E3	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
GBR2ZN	KM <i>Positive/Indicative</i> HemDirect <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	NT based on case scenario <b>Conclusion:</b> Not Tested	RSID Saliva <i>Negative</i> <b>Conclusion:</b> Negative
GCMLNM	KM <i>Positive/Indicative</i> SERATEC HB <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
GFJDUY	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not tested based on case scenario <b>Conclusion:</b> Not Tested	Generally not tested unless specified <b>Conclusion:</b> Not Tested
GTTNRM	Kastle Meyer <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Acid Phosphatase <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas Amylase <i>Negative</i> <b>Conclusion:</b> Negative
GVXEUV	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
H2GYEM	KM <i>Positive/Indicative</i> HemDirect <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	based on scenario <b>Conclusion:</b> Not Tested	RSID Saliva <i>Negative</i> <b>Conclusion:</b> Negative
H3CJ3M	KM <i>Positive/Indicative</i> SERATEC HB <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
HFJHYZ	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
HQZC42	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested

Table 1 - Item 4  
Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
HVCZ4X	Phenolphthalein	Positive/Indicative	AP	Negative		
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
J2ETAD	Visual	Positive/Indicative	AP	Negative	Phadebas tube test-	Inconclusive
	KM	Positive/Indicative	Sperm elution and Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Inconclusive	
	<u>Other Tests Performed</u>					
	OBTI	Positive/Indicative				
J6WQ4V	KM	Positive/Indicative	AP	Negative	RSID	Negative
	Micro	Positive/Indicative	Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
JQFVNY	O-tol	Positive/Indicative	ALS	Negative	RSID	Negative
	HemaTrace	Positive/Indicative	AP	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
K4X8BR	KM	Positive/Indicative	Not tested based on case scenario.		Generally not tested unless specified.	
	Bluestar OBTI	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
KBH3F3		Positive/Indicative		Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
KLZMZY	KM	Positive/Indicative				
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
KMTD7D	LMG	Positive/Indicative	AP	Negative		
	ABA HEME TRACE	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
KQ99JP	KM	Positive/Indicative	Micro	Negative	Our laboratory does not currently have a validated test for Saliva.	
			PSA	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
KTFL8J	KM	Positive/Indicative	No testing for Semen performed. The scenario as described contraindicated testing for Semen.		RSID	Negative
	HemaTrace Card	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Negative	
KXT4PZ	Hemastix	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
	<u>Other Tests Performed</u>					
	ALS	Negative				

Table 1 - Item 4  
Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
KZXWGD	Molecular Serology	Positive/Indicative	Molecular Serology	Negative	Molecular Serology	Negative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
L2CKH6	ALS BLUE (420-470 nm) with yellow glasses	Positive/Indicative	ALS BLUE (420-470 nm) with yellow glasses	Negative		
	Hemastix (SIEMENS)	Positive/Indicative	ALS UV (350-380 nm) with UV glasses	Negative		
	OBTI (Bluestar)	Positive/Indicative	PSA	Negative		
	<b>Conclusion:</b>	Positive/Indicative	AP	Negative	<b>Conclusion:</b>	Not Tested
L73VR3	PEROXIDASES	Positive/Indicative	PSA	Negative		
	HEXAGON OBTI	Positive/Indicative	Micro	Negative		
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Not Tested
LECJ72	KM	Positive/Indicative				
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Not Tested	<b>Conclusion:</b>	Not Tested
LJAZHV	KM	Positive/Indicative				
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Not Tested	<b>Conclusion:</b>	Not Tested
M4VWUF	LMG	Positive/Indicative				
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Not Tested	<b>Conclusion:</b>	Not Tested
MCLHXT	TMB	Positive/Indicative	AP	Negative	SALlgAE	Negative
	HemaTrace	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
MFNJZB	LMG	Positive/Indicative	AP, Micro	Negative	Phadebas	Negative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
MUMKGH	Kastle Meyer	Positive/Indicative	PSA	Negative	RSID	Negative
	HemaTrace	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
N44NZ2	Bluestar OBTI	Positive/Indicative	PSA	Negative		
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Not Tested
NAT4TF	KM	Positive/Indicative	AP	Negative	Phadebas	Negative
	Modified Takayama	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
	<b>Conclusion:</b>	Positive/Indicative				

Table 1 - Item 4  
Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
NBNQ6J	LMG HemaTrace Visual	Positive/Indicative Positive/Indicative Positive/Indicative	AP	Negative	RSID	Negative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
NDTENE	KM Modified Takayama	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Negative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
NDVZKV	KM HemaTrace	Positive/Indicative Positive/Indicative	AP Micro PSA	Negative Negative Negative	Phadebas	Negative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
NGE2EG	KM Seratec Hb	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Negative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
NNYWKR	KM HemaTrace	Positive/Indicative Positive/Indicative	Scenario			
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Not Tested	<b>Conclusion:</b>	Not Tested
NVFQJU	Visual examination KM	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas press test	Negative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
NXYREG	KM Seratec Hb	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Negative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
P7B7XX	KM HEM DIRECT HEMOGLOBIN	Positive/Indicative Positive/Indicative	PSA Micro - CHRISTMAS TREE	Negative Negative		
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Not Tested
P9FVYK	KM		Micro and p30			
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Not Tested
PFY2VR	KM	Positive/Indicative	AP Micro	Negative Negative	Phadebas	Negative
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
PGCTPU	LMG	Positive/Indicative	AP	Negative		
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Not Tested

Table 1 - Item 4  
 Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
PTJJQG	LMG HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative		Saliva testing is not performed at this laboratory.  <b>Conclusion:</b> Not Tested
Q9A4JA	RSID	Positive/Indicative	RSID	Negative		RSID Negative <b>Conclusion:</b> Negative
Q9UKZC			AP Micro PSA	Negative Negative Negative		<b>Conclusion:</b> Not Tested
QCDPKT	O-tol	Positive/Indicative	Low evidentiary value finding Semen on his own hands			Not relevant to the scenario <b>Conclusion:</b> Not Tested
QCVL2U	KM	Positive/Indicative				Phadebas Negative <b>Conclusion:</b> Negative
QD4XJW	LMG	Positive/Indicative				<b>Conclusion:</b> Not Tested
QD6VPL	KM Bluestar OBTI	Positive/Indicative Positive/Indicative	Not tested based on case scenario.			Generally not tested unless specified. <b>Conclusion:</b> Not Tested
QDPC8A	visible red/brown stain stain darkened under ALS/IR LMG	Positive/Indicative Positive/Indicative Positive/Indicative	no fluorescence under ALS	Negative		CSI does not test for Saliva <b>Conclusion:</b> Not Tested
QEXQEF	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative		<b>Conclusion:</b> Not Tested
QH2EXB	KM Takayama	Positive/Indicative Positive/Indicative	AP	Negative		Phadebas Amylase Negative <b>Conclusion:</b> Negative
QMUJMA	KM	Positive/Indicative	AP	Negative		<b>Conclusion:</b> Not Tested



Table 1 - Item 4  
Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
QV428T	TMB	Positive/Indicative	Microscopy	Negative		
		<b>Conclusion:</b> Positive/Indicative	p30	Negative		<b>Conclusion:</b> Not Tested
QZ2CEL	KM	Positive/Indicative	AP	Negative	Saliva Tests are not performed at the [Laboratory]	
	HemaTrace	Positive/Indicative	PSA	Negative		<b>Conclusion:</b> Not Tested
		<b>Conclusion:</b> Positive/Indicative	Micro	Negative		
T4XQ3M	Hemastix	Positive/Indicative			RSID	Negative
	RSID	Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Negative
TCMMGA	LMG	Positive/Indicative	AP	Negative	Laboratory no longer performs this testing	
	HemaTrace	Positive/Indicative	PSA	Negative		<b>Conclusion:</b> Not Tested
		<b>Conclusion:</b> Positive/Indicative	Micro	Negative		
THCZHX	Bluestar Obti	Positive/Indicative	Bluestar Identi-PSA	Negative		
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> No Response
TLCCEW	Visual	Positive/Indicative	Visual	Negative		
	fluorescence 445 nm to 510 nm	Positive/Indicative	fluorescence 445 nm to 510 nm	Negative		
	fluorescence 420 nm to 470 nm	Positive/Indicative	fluorescence 420 nm to 470 nm	Negative		
	Bluestar OBTI	Positive/Indicative	Bluestar PSA	Negative		<b>Conclusion:</b> Not Tested
		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		
TM7X8N	KM	Positive/Indicative	AP	Negative		
	HemaTrace	Positive/Indicative	PSA	Negative		<b>Conclusion:</b> Not Tested
		<b>Conclusion:</b> Positive/Indicative	Micro	Negative		
TMQJZE	KM	Positive/Indicative	AP	Negative	not required based on analysis of sample from the scenario	
	SERATEC Hb Direct	Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
		<b>Conclusion:</b> Positive/Indicative				
TQPWLH	Phenolphthalein	Positive/Indicative	Acid Phosphatase	Negative	[Laboratory] does not conduct Saliva examinations	
	HemaTrace	Positive/Indicative	p30	Negative		<b>Conclusion:</b> Not Tested
		<b>Conclusion:</b> Positive/Indicative	Microscopic	Negative		
TRY22R	KM	Positive/Indicative	ejaculate on the floor and the assault was due to being punched in the face		No indication from the narrative that Saliva should be present	
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested

Table 1 - Item 4  
Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
TW9FKC	HemaTrace	Positive/Indicative	p30	Negative	RSID	Negative
		<b>Conclusion:</b> Positive/Indicative	Micro	Negative		<b>Conclusion:</b> Negative
TXPCBQ	KM	Positive/Indicative	NT not applicable to case scenario			
	HemaTrace	Positive/Indicative				
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested
U7VQFP	KM	Positive/Indicative	AP	Negative	Phadebas	Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
U87FTA	LMG	Positive/Indicative				
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested
UJKFP4	Visual	Positive/Indicative	Acid phosphatase test of supernatant		Phadebas tube test of supernatant	
	Kastle-Mayer spot test	Positive/Indicative		Negative	(Amylase) Negative	
		<b>Conclusion:</b> Positive/Indicative	Microscopy	Negative		<b>Conclusion:</b> Negative
UMZ8EB	KM	Positive/Indicative	Ap	Negative	Phadebas	Negative
	Seratec Hemoglobin	Positive/Indicative				
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
VGR42G	KM	Positive/Indicative	Not tested based on case scenario		Generally not tested unless specified	
	Bluestar OBTI	Positive/Indicative				
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested
VHJYQB	TMB	Positive/Indicative	AP	Negative	RSID	Negative
	RSID	Positive/Indicative				
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
	<u>Other Tests Performed</u>					
	Blood - White light Visual examination Positive/Indicative					
	Semen - White light Visual examination Negative					
	Saliva - White light Visual examination Negative					
VMDWF3	RSID	Positive/Indicative	RSID	Negative	RSID	Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
VXWD89	Kastle Meyer test	Positive/Indicative				
	ABACard HemaTrace	Positive/Indicative				
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response
W4ZMZM	O-tol	Positive/Indicative	No current test in laboratory for Semen		No current test in laboratory for Saliva	
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested

Table 1 - Item 4  
 Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
WEGF6N	O-tol HemaTrace	Positive/Indicative Positive/Indicative	AP Micro PSA	Inconclusive Negative Negative	RSID	Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
WL47RJ	KM HemaTrace	Positive/Indicative Positive/Indicative	ALS AP	Negative Negative	ALS Phadebas	Negative Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
WP2REF	KM Bluestar OBTI	Positive/Indicative Positive/Indicative	Not tested based on case scenario		Generally not tested unless specified	
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested
WVQ6HP	LMG	Positive/Indicative	AP	Negative		
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
WZYRRK	KM	Positive/Indicative				
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested
X32Q3D	KM	Positive/Indicative	Micro PSA	Negative Negative	Currently, the lab does not have a validated method for testing for Saliva	
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
X7HDGH	KM HemaTrace	Positive/Indicative Positive/Indicative	AP Micro	Negative Negative	Phadebas	Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
XBE2V4	Leucomalachite green	Positive/Indicative				
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested
XEYPW4	KM HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative	Not tested based on case scenario	
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
XPVZZ6	LMG	Positive/Indicative	AP	Negative	No longer performed	
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
XTUBW6	LMG HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative	No longer performed	
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
Y48XHB	Presumptive test- LMG					
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested

Table 1 - Item 4  
Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
Y7YU86	RSID Blood	Positive/Indicative	RSID Semen	Negative	RSID Saliva	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
Y94GNF	PTMB	Positive/Indicative				
	Bluestar	Positive/Indicative				
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
YDXBLG	KM	Positive/Indicative	AP	Negative	Phadebas	Inconclusive
	HemaTrace	Positive/Indicative	ALS	Negative	ALS	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Inconclusive	
YEUW9F	TMB	Positive/Indicative	AP	Negative	SALigAE	Negative
	HemaTrace	Positive/Indicative	PSA	Negative		
			Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
YNE396	Bluestar Hexagon OBTI	Positive/Indicative	Seratec PSA semiquant	Inconclusive	Seratec Amylase Test	Negative
			Seratec PSA semiquant	Negative		
			Seratec PSA semiquant	Negative		
			RSID Semen	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
YNVC27	KM	Positive/Indicative	AP	Negative	RSID	Negative
	HemaTrace	Positive/Indicative	PSA	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
YRECCF	O-tol	Positive/Indicative			RSID	Negative
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Negative	
YTQ493	KM	Positive/Indicative	AP	Negative	Phadebas Amylase	Negative
	Takayama	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
Z7EAZC	KM HEXAGON OBTI	Positive/Indicative	case scenario		Phadebas	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Negative	
Z9YG2K	LMG	Positive/Indicative	AP	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
ZADNKG	KM	Positive/Indicative	Not tested based on case scenario.		Generally not tested unless specified.	
	Bluestar OBTI	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	

Table 1 - Item 4  
 Evidence - Four yellow plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
ZRC2V3	Hexagon OBTI	Positive/Indicative	PSA	Negative	RSID Saliva	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
ZYAT4W	Phenolphthalein	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
ZZ3LYF	KM	Positive/Indicative	Not tested based on case scenario.		Generally not tested unless specified.	
	Bluestar OBTI	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	

**Response Summary for Table 1 -  
 Item 4: Four yellow plastic swabs containing Blood**

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	166	0	2
Negative	0	124	80
Inconclusive	0	0	3
Not Tested (NT)	2	45	87
No Response	10	9	6
<b>Total</b>	<b>178</b>	<b>178</b>	<b>178</b>

# Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
2PK689	Img <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
2PZNR7	Visual Examination <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID Saliva <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
2XF3XH	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
38Y4RX	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
3996DV	KM <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	AP <i>Negative</i> <b>Conclusion:</b> No Response	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> No Response
3H2GPC	PTMB <i>Positive/Indicative</i> Blueastar <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
3RNEGH	KM <i>Positive/Indicative</i> Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
3URBZK	KM <i>Positive/Indicative</i> Seratec Hb <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
3XDWT9	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	This laboratory does not have a current validated method for Saliva testing <b>Conclusion:</b> Not Tested
3YL8PC	PTMB <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not indicated <b>Conclusion:</b> Not Tested	Phadebus Tablet <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
4DGTNC	TMB <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	SALigAE <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
	<u>Other Tests Performed</u> Body Fluid - Test Performed: ALS <i>Negative</i>		
4HDDVG	KM <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
	4JLPBX	Phenolphthalein <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Microscopic exam <i>Negative</i> p30 <i>Negative</i> <b>Conclusion:</b> Negative
4N3XXW	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Does not apply according to the scenario <b>Conclusion:</b> Not Tested	Human Amylase (alpha-amilosa) <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
4R39RC	TMB <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	SALigAE <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
	<u>Other Tests Performed</u> ALS <i>Negative</i>		
63V48H	KM <i>Positive/Indicative</i> SERATEC HemDirect <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	AP <i>Negative</i> SERATEC PSA <i>Negative</i> <b>Conclusion:</b> No Response	SERATEC AmylaseTest <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
69JRZ3	Visual <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Positive/Indicative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	ALS <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
6BBMPW	benzidine <i>Positive/Indicative</i> Hem Direct kit <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	acid phosphatase solution <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
6HUTL4	Visual Examination <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Positive/Indicative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	ALS <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
6WTV3B	LMG Positive/Indicative ALS Positive/Indicative		AP Positive/Indicative Micro Negative AP Mapping Negative ALS Negative		Amylase Diffusion Positive/Indicative Amylase Mapping Positive/Indicative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
6YVZDB	Hemastix Positive/Indicative ABAcard HemaTrace Positive/Indicative		AP Negative		SALigAE Positive/Indicative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
	<u>Other Tests Performed</u> Body fluid - ALS Negative					
79FQD8	KM Positive/Indicative HemaTrace Positive/Indicative					
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
7C4X2B	KM Positive/Indicative HemaTrace Positive/Indicative		AP Negative			
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
7FYPA8	KM Positive/Indicative Micro Positive/Indicative		AP Negative Micro Negative ALS Negative		RSID Positive/Indicative ALS Negative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
7G9YXY	KM Positive/Indicative Seratec HB Positive/Indicative		AP Negative		Phadebas Positive/Indicative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
7JZXCW	TMB Positive/Indicative Hexagon OBTI Positive/Indicative		PSA Negative		Amylase - Method CNPG3 Positive/Indicative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
	<u>Other Tests Performed</u> ALS Positive/Indicative					
7LBQEC	Hemastix Positive/Indicative				Phadebas Positive/Indicative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative	
7UGHWT	KM Positive/Indicative		AP Negative Micro Negative		RADIAL Diffusion Positive/Indicative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
83AAGT	KM Positive/Indicative HemaTrace Positive/Indicative		ALS Inconclusive AP Negative			
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	



Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
8AV3XY	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	AP <i>Negative</i> <b>Conclusion:</b> No Response	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> No Response
	<u>Other Tests Performed</u> ALS <i>Negative</i>		
8CE398	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
8GNHHW	OBTI <i>Positive/Indicative</i> TMB <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	PSA <i>Negative</i> <b>Conclusion:</b> Negative	Seratec Amylase <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
8Z7GPE	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
9BFE63	Visual <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID Saliva <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
9EJ86T	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	No areas of biological fluorescence observed. <b>Conclusion:</b> Not Tested	No areas of biological fluorescence observed. <b>Conclusion:</b> Not Tested
9K9PJ4	Micro <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
9QA7Q4	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not tested based on case scenario. <b>Conclusion:</b> Not Tested	Generally not tested unless specified. <b>Conclusion:</b> Not Tested
9QVF4E	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
	<u>Other Tests Performed</u> ALS <i>Negative</i>		
9UVNAA	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
9VPA2F	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
A7LW93	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
ACUZGP	Phenolphthalein <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
	<u>Other Tests Performed</u> Hemdirect <i>Positive/Indicative</i>		
AEDFMB	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
AF8ZDG	Visual appearance <i>Positive/Indicative</i> Bluestar forensic tablets <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	Our lab does not perform screening tests for Saliva <b>Conclusion:</b> Not Tested
ALCYQR	Visual, KM, HemDirect <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS, AP <i>Negative</i> <b>Conclusion:</b> Negative	ALS, RSID Saliva <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
AQ64HB	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas Amylase <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
ARYUNP	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
B7XV2Q	Phen <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	Radial Diffusion <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
BLQQBW	Visual <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	ALS <i>Negative</i> RSID-Saliva <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
BNVB33	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
BQJU4Z	Visual LMG HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	ALS AP	Negative Negative	ALS RSID	Positive/Indicative Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
BTQ6Z4	Phenolphthalein HemaTrace	Positive/Indicative Positive/Indicative	ALS	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
BU9Y2D	PEROXIDASES HEXAGON OBTI	Positive/Indicative Positive/Indicative	ALS			
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
C84ZT6	KM HemaTrace	Positive/Indicative Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
CDC42R	RSID Blood	Positive/Indicative	RSID Semen	Negative	RSID Saliva	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
CGD64B	Hydrogen peroxide HEXAGON OBTI	Positive/Indicative Positive/Indicative	ALS	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
CKDD9K	KM	Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response	
CNWLH7	Appearance Kastle Meyer test HemaTrace test	Positive/Indicative Positive/Indicative Positive/Indicative	Acid phosphatase	Negative	Phadebas paper press test	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
CQFMDR	KM SERATEC Hb	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
CYTE6B	LMG OBTI	Positive/Indicative Positive/Indicative		not relevant, based on scenario description		
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
DH7FGP	KM	Positive/Indicative		No areas of biological fluorescence observed.		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
DKQAC4	FOB	Positive/Indicative	RSID	Negative	RSID	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
DQFZL7	Visual appearance of Blood <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i>  <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative
DWKLQ3	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Sperm cells - Micro <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative
	<u>Other Tests Performed</u> Nucleated epithelial cells - Micro <i>Negative</i> ALS <i>Negative</i>		
DZQ388	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i>  <b>Conclusion:</b> Negative	Phadebas paper press test <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative
EBRRPL	LMG <i>Positive/Indicative</i> RSID-Blood <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> RSID-Semen <i>Negative</i> <b>Conclusion:</b> Negative	RSID-Semen <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative
ED299W	Visual <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i>  <b>Conclusion:</b> Negative	RSID-Saliva <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative
EGQD7A	HEXAGON OBTI IMMUNOCHROMATOGRAPHIC RAPID TEST <i>Positive/Indicative</i> PEROXIDASE TEST <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Positive/Indicative</i> PSA <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	TEST NOT PERFORMED  <b>Conclusion:</b> Not Tested
EQF6UM	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	No areas of biological fluorescence observed.  <b>Conclusion:</b> Not Tested	No areas of biological fluorescence observed.  <b>Conclusion:</b> Not Tested
FBHAXM	KM <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i>  <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i>  <b>Conclusion:</b> Positive/Indicative
FDLCGW	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not tested based on case scenario  <b>Conclusion:</b> Not Tested	Generally not tested unless specified  <b>Conclusion:</b> Not Tested

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
FGHUXX	KM Positive/Indicative Takayama Negative HemaTrace Positive/Indicative		ALS Negative AP Negative		[Laboratory] does not conduct Saliva testing	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
FNJ9HY	PTMB Positive/Indicative		not indicated		Phadebas Tablet Positive/Indicative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative	
FV2LGX	KM Positive/Indicative Bluestar OBTI Positive/Indicative		Not tested based on case scenario		Generally not tested unless specified	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
FWF8TK	kastle meyer KM Positive/Indicative RSID Blood Positive/Indicative		acid phosphatase AP Negative RSID Semen Negative		RSID Saliva Positive/Indicative	
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response	
FXTH7R	KM Positive/Indicative		AP Negative		Phadebas Positive/Indicative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
FZA2QL	O-tol and HemaTrace Positive/Indicative				Lab does not test for Saliva	
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
G6M6E3	LMG Positive/Indicative		AP Negative			
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
GBR2ZN	KM Positive/Indicative HemDirect Positive/Indicative		ALS Negative		RSID Saliva Positive/Indicative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
GCMLNM	KM Positive/Indicative SERATEC HB Positive/Indicative		AP Negative		Phadebas Positive/Indicative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
GFJDUY	KM Positive/Indicative Bluestar OBTI Positive/Indicative		Not tested based on case scenario		Generally not tested unless specified	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
GTTNRM	Kastle Meyer Positive/Indicative Modified Takayama Positive/Indicative		Acid Phosphatase Negative		Phadebas Amylase Positive/Indicative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
GVXEUV	KM Positive/Indicative		AP Negative Micro Negative		RSID Positive/Indicative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
H2GYEM	KM <i>Positive/Indicative</i> HemDirect <i>Positive/Indicative</i>		no fluorescence under ALS and case scenario		no fluorescence under ALS and case scenario	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
H3CJ3M	KM <i>Positive/Indicative</i> SERATEC HB <i>Positive/Indicative</i>		AP <i>Negative</i>		Phadebas <i>Positive/Indicative</i>	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
HFJHYZ	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i>		AP <i>Negative</i>		Phadebas <i>Positive/Indicative</i>	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
HQZC42	LMG <i>Positive/Indicative</i>		AP <i>Negative</i>			
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
HVCZ4X	Phenolphthalein <i>Positive/Indicative</i>		ALS <i>Negative</i> AP <i>Negative</i>			
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
J2ETAD	Visual exam <i>Positive/Indicative</i> KM <i>Positive/Indicative</i>		AP <i>Negative</i>		Phadebas paper test <i>Positive/Indicative</i>	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
	<u>Other Tests Performed</u>					
	OBTI <i>Positive/Indicative</i>					
J6WQ4V	KM <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>		AP <i>Negative</i> Micro <i>Negative</i>		RSID <i>Positive/Indicative</i>	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
JQFVNY	O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i>		ALS <i>Negative</i> AP <i>Negative</i>		RSID <i>Positive/Indicative</i>	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
K4X8BR	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i>		Not tested based on case scenario.		Generally not tested unless specified.	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
KBH3F3	<i>Positive/Indicative</i>		<i>Negative</i>			
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
KLZMZY	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i>					
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
KMTD7D	LMG ABA HEME TRACE	Positive/Indicative Positive/Indicative	AP	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
KQ99JP	KM	Positive/Indicative	ALS	Negative		Our laboratory does not currently have a validated test for Saliva.
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
KTFL8J	KM HemaTrace card	Positive/Indicative Positive/Indicative	No testing for Semen performed. The scenario as described contraindicated testing for Semen.		Phadebas Amylase Mapping	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		RSID	Positive/Indicative
<b>Other Tests Performed</b>						
ALS was used to search for staining areas. Nothing further detected.						
KXT4PZ	Hemastix	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
<b>Other Tests Performed</b>						
ALS Negative						
KZXWGD	Molecular Serology	Positive/Indicative	Molecular Serology	Negative	Molecular Serology	Inconclusive
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Inconclusive	
L2CKH6	ALS BLUE (420-470 nm) with yellow glasses Hemastix (SIEMENS) OBTI (Bluestar)	Positive/Indicative Positive/Indicative Positive/Indicative	ALS BLUE (420-470 nm) with yellow glasses ALS UV (350-380 nm) with UV glasses	Negative Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
L73VR3	PEROXIDASES HEXAGON OBTI	Positive/Indicative Positive/Indicative	ALS	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested	
LECJ72	KM	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
LJAZHV	KM HemaTrace	Positive/Indicative Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
M4VWUF	LMG	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
MCLHXT	TMB <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	SALigAE <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
	<u>Other Tests Performed</u> Body Fluid - ALS <i>Negative</i>		
MFNJZB	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not examined given case circumstances <b>Conclusion:</b> Not Tested	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
MUMKGH	Kastle Meyer <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	scenario does not associate lanyard with possible Semen, so no Semen testing performed <b>Conclusion:</b> Not Tested	Phadebas Paper <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
	<u>Other Tests Performed</u> ALS @ 415nm, 450nm, 470nm & 490nm <i>Negative</i>		
N44NZ2	Bluestar Obti <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
NAT4TF	KM <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
NBNQ6J	Visual <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	ALS <i>Negative</i> RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
NDTENE	KM <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
NDVZKV	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
NGE2EG	KM <i>Positive/Indicative</i> Seratec Hb <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
NNYWKR	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Scenario <b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested



Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
NVFQJU	Visual examination <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas Press test <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
NXYREG	KM <i>Positive/Indicative</i> Seratec Hb <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
P7B7XX	KM <i>Positive/Indicative</i> HEM DIRECT HEMOGLOBIN <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	PSA <i>Negative</i> Micro - CHRISTMAS TREE <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
P9FVYK	KM <b>Conclusion:</b> Positive/Indicative	ALS, Micro and p30 <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
PFY2VR	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
PGCTPU	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
PTJJQG	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	Saliva testing is not performed at this laboratory. <b>Conclusion:</b> Not Tested
Q9A4JA	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
Q9UKZC	<b>Conclusion:</b> Not Tested	AP <i>Negative</i> Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
QCDPKT	O-tol <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not relevant to scenario <b>Conclusion:</b> Not Tested	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
QCVL2U	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
QD4XJW	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
QD6VPL	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not tested based on case scenario. <b>Conclusion:</b> Not Tested	Generally not tested unless specified. <b>Conclusion:</b> Not Tested
QDPC8A	visible red/brown stain <i>Positive/Indicative</i> stain darkened under ALS/IR <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	no fluorescence under ALS <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
QEXQEF	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> ALS <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> ALS <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
QH2EXB	KM <i>Positive/Indicative</i> Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas Amylase <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
QMUJMA	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
QV428T	TMB <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
QZ2CEL	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Based on scenario testing is not warranted <b>Conclusion:</b> Not Tested	Saliva Tests are not performed at the [Laboratory] <b>Conclusion:</b> Not Tested
T4XQ3M	Hemastix <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	ALS <i>Negative</i> RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
TCMMGA	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Scenario did not indicate need for test <b>Conclusion:</b> Not Tested	Laboratory no longer performs this testing <b>Conclusion:</b> Not Tested
THCZHx	Bluestar Obti <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Bluestar Identi-PSA <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> No Response

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
TLCCEW	Visual <i>Positive/Indicative</i> fluorescence 445 nm to 510 nm <i>Positive/Indicative</i> fluorescence 420 nm to 470 nm <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	Visual <i>Negative</i> fluorescence 445 nm to 510 nm <i>Negative</i> fluorescence 420 nm to 470 nm <i>Negative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested
TM7X8N	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
TMQJZE	KM <i>Positive/Indicative</i> SERATEC Hb Direct <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
TQPWLH	Phenolphthalein <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Alternate Light Source <i>Negative</i> Acid phosphatase <i>Negative</i> p30 <i>Negative</i> Microscopic <i>Negative</i> <b>Conclusion:</b> Negative	[Laboratory] does not conduct Saliva examinations <b>Conclusion:</b> Not Tested
TRY22R	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	no indication of possible Semen to the face or lanyard <b>Conclusion:</b> Not Tested	Phadebas on beige/tan stain <i>Positive/Indicative</i> Phadebas on red/brown stain <i>Negative</i> <b>Conclusion:</b> Positive/Indicative
TW9FKC	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative  <u>Other Tests Performed</u> ALS <i>Negative</i>	AP <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
TXPCBQ	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	NT not applicable to case scenario <b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
U7VQFP	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Based on scenario given- lanyard used to tie up victim so not tested for Semen <b>Conclusion:</b> Not Tested	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
U87FTA	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
UJKFP4	Visual <i>Positive/Indicative</i> Kastle-Mayer spot-test <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Acid phosphatase screen of item <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas screen of item (Amylase) <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
<u>Other Tests Performed</u> Crime-lite (no significant fluorescence observed) <i>Inconclusive</i>			
UMZ8EB	KM <i>Positive/Indicative</i> Seratec Hemoglobin <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
VGR42G	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not tested based on case scenario <b>Conclusion:</b> Not Tested	Generally not tested unless specified <b>Conclusion:</b> Not Tested
VHJYQB	ALS <i>Positive/Indicative</i> TMB <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	ALS <i>Negative</i> RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
<u>Other Tests Performed</u> Blood - White light Visual examination <i>Positive/Indicative</i> Semen - White light Visual examination <i>Negative</i> Saliva - White light Visual examination <i>Negative</i>			
VMDWF3	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
VXWD89	Kastle Meyer <i>Positive/Indicative</i> ABACard HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> No Response	<b>Conclusion:</b> No Response
W4ZMZM	O-tol <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	No current test in laboratory for Semen <b>Conclusion:</b> Not Tested	No current test in laboratory for Saliva <b>Conclusion:</b> Not Tested
WEGF6N	O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Synopsis does not suggest the presence of Semen. <b>Conclusion:</b> Not Tested	ALS <i>Negative</i> RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
WL47RJ	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	ALS <i>Negative</i> Phadebas <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative
WP2REF	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not tested based on case scenario <b>Conclusion:</b> Not Tested	Generally not tested unless specified <b>Conclusion:</b> Not Tested

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
WVQ6HP	LMG	Positive/Indicative	AP	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
WZYRRK	KM	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
X32Q3D	KM	Positive/Indicative	Micro	Negative		Currently, the lab does not have a validated method for testing for Saliva
			PSA	Negative		
			ALS	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
X7HDGH	KM	Positive/Indicative	ALS	Negative		ALS Negative
	HemaTrace	Positive/Indicative	AP	Negative		Phadebas Positive/Indicative
			Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
XBE2V4	Leucomalachite green	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
XEYPW4						
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
XPVZZ6	LMG	Positive/Indicative		Scenario does not indicate need for test		No longer performed
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
XTUBW6	LMG	Positive/Indicative		Scenario does not indicate need for test		No longer performed
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
Y48XHB	Presumptive test- LMG					
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
Y7YU86	RSID Blood	Positive/Indicative	RSID Semen	Negative		RSID Saliva Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
Y94GNF	PTMB	Positive/Indicative				
	Bluestar	Positive/Indicative				
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
YDXBLG	KM	Positive/Indicative	AP	Negative		Phadebas Positive/Indicative
	HemaTrace	Positive/Indicative	ALS	Positive/Indicative		ALS Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	

Table 1 - Item 5  
Evidence - Lanyard containing Blood & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
YEUW9F	TMB HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro ALS	Negative Negative Negative Negative	SALigAE	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
YNE396	Bluestar Hexagon OBTI	Positive/Indicative	Seratec PSA Semiquant	Negative	Seratec Amylase Test	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
YNVC27	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA	Negative Negative	Phadebas Press RSID	Positive/Indicative Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
	<u>Other Tests Performed</u> Non-specific body fluid - ALS Negative					
YRECCF	O-tol HemaTrace	Positive/Indicative Positive/Indicative			RSID	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative	
YTQ493	KM Takayama	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas Amylase	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
Z7EAZC	KM HEXAGON OBTI	Positive/Indicative	case scenario		Phadebas	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Positive/Indicative	
Z9YG2K	LMG	Positive/Indicative	AP	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
ZADNKG	KM Bluestar OBTI	Positive/Indicative Positive/Indicative	Not tested based on case scenario.		Generally not tested unless specified.	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
ZRC2V3	Hexagon OBTI	Positive/Indicative	PSA RSID	Negative	RSID Saliva	Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
ZYAT4W	Phenolphthalein	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
ZZ3LYF	KM Bluestar OBTI	Positive/Indicative Positive/Indicative	Not tested based on case scenario.		Generally not tested unless specified.	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	

**Response Summary for Table 1 -  
Item 5: Lanyard containing Blood & Saliva Mixture**

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	164	0	88
Negative	0	102	0
Inconclusive	0	0	1
Not Tested (NT)	3	68	83
No Response	11	8	6
<hr/>			
Total	178	178	178

# Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
2PK689	Img <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
	<u>Other Tests Performed</u> ALS <i>Positive/Indicative</i>		
2PZNR7	Visual Examination <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID Saliva <i>Negative</i> <b>Conclusion:</b> Negative
2XF3XH	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
38Y4RX	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
3996DV	KM <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	AP <i>Negative</i> <b>Conclusion:</b> No Response	Phadebas <i>Negative</i> <b>Conclusion:</b> No Response
3H2GPC	PTMB <i>Positive/Indicative</i> Bluestar <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
3RNEGH	KM <i>Positive/Indicative</i> Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
3URBZK	KM <i>Positive/Indicative</i> Seratec Hb <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
3XDWT9	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	This laboratory does not have a current validated method for Saliva testing <b>Conclusion:</b> Not Tested
3YL8PC	PTMB <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not indicated <b>Conclusion:</b> Not Tested	Not indicated <b>Conclusion:</b> Not Tested



Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
4DGTNC	TMB	Positive/Indicative	AP	Negative	SALigAE	Negative
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
	<u>Other Tests Performed</u>					
	Body Fluid - Test Performed: ALS Negative					
4HDDVG	KM	Positive/Indicative	AP	Negative	Phadebas	Negative
	Modified Takayama	Positive/Indicative	Micro	Negative		
			PSA	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
4JLPBX	Phenolphthalein	Positive/Indicative	Microscopic exam	Negative	Phadebas	Negative
	HemaTrace	Positive/Indicative	p30	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
4N3XXW	KM	Positive/Indicative	Brentamine (AP)	Negative	Human Amylase (alpha-amilosa)	Positive/Indicative
	HemaTrace	Positive/Indicative	Human seminal protein (p30)	Negative		
			Microscopic	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
4R39RC	TMB	Positive/Indicative	AP	Negative	SALigAE	Negative
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
	<u>Other Tests Performed</u>					
	ALS Negative					
63V48H	KM	Positive/Indicative	AP	Negative	SERATEC AmylaseTest	Negative
	SERATEC HemDirect	Positive/Indicative	SERATEC PSA	Negative		
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> Negative	
69JRZ3	Visual	Positive/Indicative	ALS	Negative	ALS	Negative
	LMG	Positive/Indicative	AP	Negative	RSID	Negative
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
6BBMPW	benzidine	Positive/Indicative	acid phosphatase solution	Negative		
	Hem Direct kit	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
6HUTL4	Visual Examination	Positive/Indicative	ALS	Negative	ALS	Negative
	LMG	Positive/Indicative	AP	Negative	RSID	Negative
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
6WTV3B	LMG ALS	Positive/Indicative Positive/Indicative	AP Micro ALS	Negative Negative Negative	Amylase Diffusion	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
6YVZDB	Hemastix ABAcad HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative	SALlgAE	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
	<u>Other Tests Performed</u>					
	Body Fluid - ALS Negative					
79FQD8	KM HemaTrace	Positive/Indicative Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
7C4X2B	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Inconclusive Negative Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
7FYPA8	KM Micro	Positive/Indicative Positive/Indicative	AP Micro ALS	Negative Negative Positive/Indicative	RSID ALS	Negative Positive/Indicative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
7G9YXY	KM Seratec HB	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
7JZXCW	TMB Hexagon OBTI	Positive/Indicative Positive/Indicative	PSA	Negative	Amylase - Method CNPG3	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
	<u>Other Tests Performed</u>					
	ALS Positive/Indicative					
7LBQEC	Hemastix	Positive/Indicative			Phadebas	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Negative	
7UGHWT	KM	Positive/Indicative	AP Micro	Negative Negative	RADIAL Diffusion	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
83AAGT	KM HemaTrace	Positive/Indicative Positive/Indicative	ALS AP	Negative Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
8AV3XY	KM	Positive/Indicative	AP	Negative	RSID	Negative
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response	
	<u>Other Tests Performed</u>					
	ALS Negative					
8CE398	KM	Positive/Indicative				
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
8GNHHW	OBTI	Positive/Indicative	PSA	Negative	Seratec Amylase	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
8Z7GPE	LMG	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
9BFE63	Visual	Positive/Indicative	ALS	Negative	RSID Saliva	Negative
	LMG	Positive/Indicative	AP	Negative		
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
9EJ86T	KM	Positive/Indicative	No areas of biological fluorescence observed.		No areas of biological fluorescence observed.	
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
9K9PJ4	Micro	Positive/Indicative	ALS	Positive/Indicative	ALS	Positive/Indicative
	KM	Positive/Indicative	AP	Negative	RSID	Negative
			Micro	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
9QA7Q4	KM	Positive/Indicative	Not tested based on case scenario.		Generally not tested unless specified.	
	Bluestar	OBTI Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
9QVF4E	LMG	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
	<u>Other Tests Performed</u>					
	ALS Positive/Indicative					
9UVNAA	KM	Positive/Indicative				
	HemaTrace	Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
9VPA2F						
	<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
A7LW93	KM	Positive/Indicative	AP	Negative	RSID	Negative
		<b>Conclusion:</b> Positive/Indicative	Micro	Negative		<b>Conclusion:</b> Negative
ACUZGP	Phenolphthalein	Positive/Indicative	Alternative light source	Negative		
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
	<u>Other Tests Performed</u>					
	Hemedirect	Positive/Indicative				
AEDFMB	LMG	Positive/Indicative	AP	Negative		
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
AF8ZDG	Visual appearance	Positive/Indicative	ALS	Negative		Our lab does not perform screening tests for Saliva
	Bluestar forensic tablets	Positive/Indicative	Micro	Negative		
		<b>Conclusion:</b> Positive/Indicative	PSA	Negative		
ALCYQR	Visual, KM, HemDirect	Positive/Indicative	ALS, AP	Negative	ALS, RSID Saliva	Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
AQ64HB	KM	Positive/Indicative	AP	Negative	Phadebas Amylase Test - Blood-like stain	Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative	Phadebas Amylase Test - Press Test of entire exhibit	Negative
						<b>Conclusion:</b> Negative
ARYUNP	KM	Positive/Indicative	ALS	Positive/Indicative		
	HemaTrace	Positive/Indicative	AP	Negative		
		<b>Conclusion:</b> Positive/Indicative	Micro	Negative		
			PSA	Negative		
				<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
B7XV2Q	Phen	Positive/Indicative	AP	Negative	Radial Diffusion	Negative
		<b>Conclusion:</b> Positive/Indicative	Micro	Negative		<b>Conclusion:</b> Negative
BLQQBW	Visual	Positive/Indicative	ALS	Negative	ALS	Negative
	LMG	Positive/Indicative	AP	Negative	RSID-Semen	Negative
	HemaTrace	Positive/Indicative				
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
BNVB33	KM	Positive/Indicative	AP	Negative	RSID	Negative
		<b>Conclusion:</b> Positive/Indicative	Micro	Negative		<b>Conclusion:</b> Negative
				<b>Conclusion:</b> Negative		

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
BQJU4Z	Visual LMG HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	ALS AP	Negative Negative	ALS RSID	Negative Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
BTQ6Z4	Phenolphthalein HemaTrace	Positive/Indicative Positive/Indicative	ALS	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
BU9Y2D	PEROXIDASES HEXAGON OBTI	Positive/Indicative Positive/Indicative	ALS	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested	
C84ZT6	KM HemaTrace	Positive/Indicative Positive/Indicative				
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
CDC42R	RSID Blood	Positive/Indicative	RSID Semen	Negative	RSID Saliva	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
CGD64B	Hydrogen peroxide HEXAGON OBTI	Positive/Indicative Positive/Indicative	ALS	Negative		
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
CKDD9K	KM	Positive/Indicative	AP	Negative	Phadebas	Negative
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response	
CNWLH7	Appearance Kastle Meyer test HemaTrace test	Positive/Indicative Positive/Indicative Positive/Indicative	Acid phosphatase tset	Negative	Phadebas paper press test	Inconclusive
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Inconclusive	
CQFMDR	KM SERATEC Hb	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Negative
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
CYTE6B	LMG OBTI	Positive/Indicative Positive/Indicative	not relevant, based on scenario description			
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
	<u>Other Tests Performed</u>					
	Forensic Light source Crime Lite Negative					
DH7FGP	KM	Positive/Indicative	Not tested due to case scenario.			
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
DKQAC4	FOB <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
DQFZL7	Visual appearance of Blood <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
DWKLQ3	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> Sperm cell - Micro <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Nucleated epithelial cells - Micro <i>Negative</i> ALS <i>Negative</i>		
DZQ388	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas paper press test <i>Negative</i> <b>Conclusion:</b> Negative
EBRRPL	LMG <i>Positive/Indicative</i> RSID-Blood <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> RSID-Semen <i>Negative</i> <b>Conclusion:</b> Negative	RSID-Saliva <i>Negative</i> <b>Conclusion:</b> Negative
ED299W	Visual <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	RSID-Saliva <i>Negative</i> <b>Conclusion:</b> Negative
EGQD7A	HEXAGON OBTI IMMUNOCHROMATOGRAPHIC RAPID TEST <i>Positive/Indicative</i> PEROXIDASE TEST <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Positive/Indicative</i> PSA <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	TEST NOT PERFORMED <b>Conclusion:</b> Not Tested
EQF6UM	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	No areas of biological fluorescence observed. <b>Conclusion:</b> Not Tested	No areas of biological fluorescence observed. <b>Conclusion:</b> Not Tested
FBHAXM	KM <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
FDLCGW	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not tested based on case scenario <b>Conclusion:</b> Not Tested	Generally not tested unless specified <b>Conclusion:</b> Not Tested

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
FGHUXX	KM <i>Positive/Indicative</i> Takayama <i>Positive/Indicative</i>		ALS <i>Negative</i> AP <i>Positive/Indicative</i> Micro <i>Negative</i> PSA <i>Negative</i>		[Laboratory] does not conduct Saliva testing	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
FNJ9HY	PTMB <i>Positive/Indicative</i>		not indicated		Phadebas Tablet <i>Negative</i>	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Negative	
FV2LGX	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i>		Not tested based on case scenario		Generally not tested unless specified	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
FWF8TK	kastle meyer KM <i>Positive/Indicative</i> RSID Blood <i>Positive/Indicative</i>		acid phosphatase AP <i>Negative</i> RSID Semen <i>Negative</i>		RSID Saliva <i>Negative</i>	
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response		<b>Conclusion:</b> No Response	
FXTH7R	KM <i>Positive/Indicative</i>		AP <i>Negative</i>		Phadebas <i>Positive/Indicative</i>	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Positive/Indicative	
FZA2QL	O-tol and HemaTrace <i>Positive/Indicative</i>				Lab does not test for Saliva	
	<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
G6M6E3	LMG <i>Positive/Indicative</i>		AP <i>Negative</i>			
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
GBR2ZN	KM <i>Positive/Indicative</i> HemDirect <i>Positive/Indicative</i>		ALS <i>Negative</i>		RSID Saliva <i>Negative</i>	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
GCMLNM	KM <i>Positive/Indicative</i> SERATEC HB <i>Positive/Indicative</i>		AP <i>Negative</i>		Phadebas <i>Negative</i>	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
GFJDUY	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i>		Not tested based on case scenario		Generally not tested unless specified	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
GTTNRM	Kastle Meyer <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i>		Acid Phosphatase <i>Negative</i>		Phadebas Amylase <i>Negative</i>	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
GVXEUV	KM <i>Positive/Indicative</i>		AP <i>Negative</i> Micro <i>Negative</i>		RSID <i>Negative</i>	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
H2GYEM	KM Positive/Indicative HemDirect Positive/Indicative		case scenario		RSID Saliva Negative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Negative	
H3CJ3M	KM Positive/Indicative SERATEC HB Positive/Indicative		AP Negative		Phadebas Negative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
HFJHYZ	KM Positive/Indicative HemaTrace Positive/Indicative		AP Negative		Phadebas Negative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
HQZC42	LMG Positive/Indicative		AP Negative			
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
HVCZ4X	Phenolphthalein Positive/Indicative		ALS Negative AP Negative			
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested	
J2ETAD	Visual exam Positive/Indicative KM Positive/Indicative		AP Negative		Phadebas paper test Negative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
	<u>Other Tests Performed</u>					
	OBTI Positive/Indicative					
J6WQ4V	KM Positive/Indicative Micro Positive/Indicative		AP Negative Micro Negative		RSID Negative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
JQFVNY	O-tol Positive/Indicative HemaTrace Positive/Indicative		ALS Negative AP Negative		RSID Negative	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative	
K4X8BR	KM Positive/Indicative Bluestar OBTI Positive/Indicative		Not tested based on case scenario.		Generally not tested unless specified.	
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	
KBH3F3	Positive/Indicative		Positive/Indicative			
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested	
KLZMZY	KM Positive/Indicative HemaTrace Positive/Indicative					
	<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested	



Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
KMTD7D	LMG <i>Positive/Indicative</i> ABA HEME TRACE <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
KQ99JP	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Positive/Indicative</i> Micro <i>Negative</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	Our laboratory does not currently have a validated test for Saliva. <b>Conclusion:</b> Not Tested
KTFL8J	KM on Stain 1 <i>Positive/Indicative</i> HemaTrace Card on Stain 1 <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	No testing for Semen performed. The scenario as described contraindicated testing for Semen. Also, areas for Semen testing not identified using ALS. <b>Conclusion:</b> Not Tested	Phadebas Amylase Mapping <i>Negative</i> RSID <i>Negative</i> <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> ALS was used to search for staining areas. Nothing further detected.		
KXT4PZ	Hemastix <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
	<u>Other Tests Performed</u> ALS <i>Positive/Indicative</i>		
KZXWGD	Molecular Serology <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Molecular Serology <i>Negative</i> <b>Conclusion:</b> Negative	Molecular Serology <i>Negative</i> <b>Conclusion:</b> Negative
L2CKH6	ALS BLUE (420-470 nm) with yellow glasses <i>Positive/Indicative</i> Hemastix (SIEMENS) <i>Positive/Indicative</i> OBTI (Bluestar) <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS BLUE (420-470 nm) with yellow glasses <i>Negative</i> ALS UV (350-380 nm) with UV glasses <i>Negative</i> PSA <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
L73VR3	PEROXIDASES <i>Positive/Indicative</i> HEXAGON OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Positive/Indicative</i> PSA <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
LECJ72	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
LJAZHV	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
M4VWUF	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
MCLHXT	TMB <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	SALigAE <i>Negative</i> <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> BodyFluid - ALS <i>Negative</i>		
MFNJZB	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not examined given case circumstances <b>Conclusion:</b> Not Tested	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
MUMKGH	Kastle Meyer <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Acid Phosphatase <i>Inconclusive</i> PSA <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas Paper <i>Inconclusive</i> RSID <i>Negative</i> <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> ALS @ 415nm, 450nm, 470nm & 490nm <i>Positive/Indicative</i>		
N44NZ2	Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
NAT4TF	KM <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
NBNQ6J	Visual <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	ALS <i>Negative</i> RSID <i>Negative</i> <b>Conclusion:</b> Negative
NDTENE	KM <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
NDVZKV	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
NGE2EG	KM <i>Positive/Indicative</i> Seratec Hb <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
NNYWKR	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Scenario <b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
NVFQJU	Visual examination KM	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas press test	Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
NXYREG	KM Seratec Hb	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
P7B7XX	KM HEM DIRECT HEMOGLOBINE	Positive/Indicative Positive/Indicative	PSA Micro - CHRISTMAS TREE	Negative Negative		
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
P9FVYK	KM		ALS, Micro and p30			
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
PFY2VR	KM	Positive/Indicative			Phadebas	Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Negative
PGCTPU	LMG	Positive/Indicative	AP	Negative		
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
PTJJQG	LMG HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative	Saliva testing is not performed at this laboratory.	
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
Q9A4JA	RSID	Positive/Indicative	RSID	Negative	RSID	Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
Q9UKZC			AP Micro PSA	Negative Negative Negative		
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
QCDPKT	O-tol	Positive/Indicative	Not relevant to scenario		RSID	Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Negative
QCVL2U	KM	Positive/Indicative			Phadebas	Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Negative
QD4XJW	LMG	Positive/Indicative				
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
QD6VPL	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not tested based on case scenario. <b>Conclusion:</b> Not Tested	Generally not tested unless specified. <b>Conclusion:</b> Not Tested
QDPC8A	visible red/brown stain <i>Positive/Indicative</i> darkened under ALS/IR <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS (faint lightening around perimeter) <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
QEXQEF	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> ALS <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> ALS <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
QH2EXB	KM <i>Positive/Indicative</i> Takayama <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas Amylase <i>Negative</i> <b>Conclusion:</b> Negative
QMUJMA	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
QV428T	TMB <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
QZ2CEL	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Positive/Indicative</i> AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	Saliva Tests are not performed at the [Laboratory] <b>Conclusion:</b> Not Tested
T4XQ3M	Hemastix <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	ALS <i>Negative</i> RSID <i>Negative</i> <b>Conclusion:</b> Negative
TCMMGA	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Scenario did not indicate need for test <b>Conclusion:</b> Not Tested	Laboratory no longer performs this testing <b>Conclusion:</b> Not Tested
THCZHX	Bluestar Obti <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Bluestar Identi-PSA <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> No Response

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
TLCCEW	Visual <i>Positive/Indicative</i> fluorescence 445 nm to 510 nm <i>Positive/Indicative</i> fluorescence 420 nm to 470 nm <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> No Response	Visual <i>Negative</i> fluorescence 445 nm to 510 nm <i>Negative</i> fluorescence 420 nm to 470 nm <i>Negative</i> Bluestar PSA <i>Negative</i> <b>Conclusion:</b> No Response	<b>Conclusion:</b> Not Tested
TM7X8N	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
TMQJZE	KM <i>Positive/Indicative</i> SERATEC Hb Direct <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phedebas <i>Negative</i> <b>Conclusion:</b> Negative
TQPWLH	Phenolphthalein <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Alternate Light Source <i>Negative</i> Acid phosphatase <i>Negative</i> p30 <i>Negative</i> Microscopic <i>Negative</i> <b>Conclusion:</b> Negative	[Laboratory] does not conduct Saliva examinations <b>Conclusion:</b> Not Tested
TRY22R	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	no stains observed other than the visible red/brown stain <b>Conclusion:</b> Not Tested	no stains observed other than the visible red/brown stain <b>Conclusion:</b> Not Tested
TW9FKC	HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative <u>Other Tests Performed</u> ALS <i>Negative</i>	p30 <i>Negative</i> Micro <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas Press Test <i>Negative</i> <b>Conclusion:</b> Negative
TXPCBQ	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	NT not applicable to case scenario <b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
U7VQFP	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
U87FTA	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	<b>Conclusion:</b> Not Tested	<b>Conclusion:</b> Not Tested
UJKFP4	Visual <i>Positive/Indicative</i> Kastle-Mayer spot test <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative <u>Other Tests Performed</u> Crime-lite (no significant fluorescence observed) <i>Inconclusive</i>	Acid phosphatase screen of item <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas screen of item (Amylase) <i>Inconclusive</i> <b>Conclusion:</b> Inconclusive

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
UMZ8EB	KM <i>Positive/Indicative</i> Seratec Hemoglobin <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
VGR42G	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not tested based on case scenario <b>Conclusion:</b> Not Tested	Generally not tested unless specified <b>Conclusion:</b> Not Tested
VHJYQB	ALS <i>Positive/Indicative</i> TMB <i>Positive/Indicative</i> RSID <i>Negative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Positive/Indicative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	ALS <i>Positive/Indicative</i> RSID <i>Negative</i> <b>Conclusion:</b> Negative
	<u>Other Tests Performed</u> Blood - White light Visual examination <i>Positive/Indicative</i> Semen - White light Visual examination <i>Positive/Indicative</i> Saliva - White light Visual examination <i>Positive/Indicative</i>		
VMDWF3	RSID <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	RSID <i>Negative</i> <b>Conclusion:</b> Negative	RSID <i>Negative</i> <b>Conclusion:</b> Negative
VXWD89	Kastle Meyer <i>Positive/Indicative</i> ABACard HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	brentamine spot test <i>Negative</i> ABACard p30 <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> No Response
W4ZMZM	O-tol <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	No current test in laboratory for Semen <b>Conclusion:</b> Not Tested	No current test in laboratory for Saliva <b>Conclusion:</b> Not Tested
WEGF6N	O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Synopsis does not suggest the presence of Semen. <b>Conclusion:</b> Not Tested	ALS <i>Negative</i> RSID <i>Negative</i> <b>Conclusion:</b> Negative
WL47RJ	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> <b>Conclusion:</b> Negative	ALS <i>Negative</i> Phadebas <i>Negative</i> <b>Conclusion:</b> Negative
WP2REF	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	Not tested based on case scenario <b>Conclusion:</b> Not Tested	Generally not tested unless specified <b>Conclusion:</b> Not Tested
WVQ6HP	LMG <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested
WZYRRK	KM <i>Positive/Indicative</i> <b>Conclusion:</b> Positive/Indicative	AP <i>Negative</i> <b>Conclusion:</b> Negative	<b>Conclusion:</b> Not Tested

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
X32Q3D	KM	Positive/Indicative	Micro	Negative		Currently, the lab does not have a validated method for testing for Saliva
			PSA	Negative		
			ALS	Negative		
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
X7HDGH	KM	Positive/Indicative	ALS	Negative		ALS Negative
	HemaTrace	Positive/Indicative	AP	Negative		Phadebas Negative
			Micro	Negative		
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
XBE2V4	Leucomalachite green	Positive/Indicative				
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested
XEYPW4						
		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested
XPVZZ6	LMG	Positive/Indicative	AP	Negative		No longer performed
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Not Tested
XTUBW6	LMG	Positive/Indicative		Scenario does not indicate need for test		No longer performed
	HemaTrace	Positive/Indicative				
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested
Y48XHB	Presumptive test-	LMG				
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested
Y7YU86	RSID Blood	Positive/Indicative	RSID Semen	Negative		RSID Saliva Negative
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative
Y94GNF	PTMB	Positive/Indicative				
	Bluestar	Positive/Indicative				
		<b>Conclusion:</b> No Response		<b>Conclusion:</b> Not Tested		<b>Conclusion:</b> Not Tested
YDXBLG	KM	Positive/Indicative	AP	Negative		Phadebas Inconclusive
	HemaTrace	Positive/Indicative	ALS	Inconclusive		ALS Inconclusive
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Inconclusive
YE UW9F	TMB	Positive/Indicative	AP	Negative		SALigAE Negative
	HemaTrace	Positive/Indicative	PSA	Negative		
			Micro	Negative		
			ALS	Negative		
		<b>Conclusion:</b> Positive/Indicative		<b>Conclusion:</b> Negative		<b>Conclusion:</b> Negative

Table 1 - Item 6  
Evidence - Washcloth containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
YNE396	Bluestar Hexagon OBTI <i>Positive/Indicative</i>		Seratec PSA Semiquant RSID Semen	<i>Negative</i> <i>Negative</i>	Seratec Amylase Test	<i>Negative</i>
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
YNVC27	KM HemaTrace	<i>Positive/Indicative</i> <i>Positive/Indicative</i>	AP PSA	<i>Negative</i> <i>Negative</i>	Phadebas Press* RSID	<i>Positive/Indicative</i> <i>Negative</i>
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
	<u>Other Tests Performed</u>					
	Non-specific body fluid - ALS <i>Negative</i>					
YRECCF	O-tol HemaTrace	<i>Positive/Indicative</i> <i>Positive/Indicative</i>	AP	<i>Negative</i>	RSID	<i>Negative</i>
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
YTQ493	KM Takayama	<i>Positive/Indicative</i> <i>Positive/Indicative</i>	AP	<i>Negative</i>	Phadebas Amylase	<i>Negative</i>
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
Z7EAZC	KM HEXAGON OBTI	<i>Positive/Indicative</i>	AP	<i>Negative</i>	Phadebas	<i>Negative</i>
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
Z9YG2K	LMG	<i>Positive/Indicative</i>	AP	<i>Negative</i>		
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Not Tested
ZADNKG	KM Bluestar OBTI	<i>Positive/Indicative</i> <i>Positive/Indicative</i>	Not tested based on case scenario.		Generally not tested unless specified.	
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Not Tested	<b>Conclusion:</b>	Not Tested
ZRC2V3	Hexagon OBTI	<i>Positive/Indicative</i>	PSA RSID	<i>Negative</i>	RSID Saliva	<i>Negative</i>
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Negative	<b>Conclusion:</b>	Negative
ZYAT4W	Phenolphthalein	<i>Positive/Indicative</i>				
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Not Tested	<b>Conclusion:</b>	Not Tested
ZZ3LYF	KM Bluestar OBTI	<i>Positive/Indicative</i> <i>Positive/Indicative</i>	Not tested based on case scenario.		Generally not tested unless specified.	
	<b>Conclusion:</b>	Positive/Indicative	<b>Conclusion:</b>	Not Tested	<b>Conclusion:</b>	Not Tested



**Response Summary for Table 1 -  
Item 6: Washcloth containing Blood**

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	164	1	2
Negative	0	115	83
Inconclusive	0	0	3
Not Tested (NT)	3	55	84
No Response	11	7	6
<hr/>			
Total	178	178	178

# Additional Comments

## TABLE 2

WebCode	Additional Comments
2PK689	Items recovered at a crime scene with no visual staining would be collected as a whole item and forwarded to the biological sciences section for further analysis. Presumptive testing other than for possible blood is not performed at a crime scene. ALS may be used on a case-by-case basis, but typically the whole item would be collected for further analysis by the biological sciences section. Presumptive testing of items 01.02.01, 01.04.01, 01.05.01, and 01.06.01 for the presence of blood was positive. A positive presumptive test is not a confirmation of the identity of blood.
2XF3XH	RESULTS AND CONCLUSIONS: Items 2, 4, 5 and 6 contain blood. Presumptive* testing for human origin was positive. No blood was found Item 1. Items 1 and 3 contain semen. No semen was found Items 2 or 4. No analysis was performed on Item 7. *Presumptive tests do not confirm the presence of bodily fluids or human origin.
3H2GPC	Crime scene personnel not trained on semen or saliva screening.
3XDWT9	Presumptive testing detected the possible presence of blood on Items 2, 4, 5 and 6. Spermatozoa were observed on Items 1 and 3, which confirms the presence of semen. Semen was not detected on Items 2, 4, 5, and 6.
3YL8PC	No testing was performed on the control swabs
4JLPBX	Viewed Item 5 and 6 with ALS - blue light
6HUTL4	All microscopic examinations were performed after the Kernechtrot-Picroindigocarmine staining procedure. No hairs were observed on the item 5 lanyard or the item 6 washcloth.
79FQD8	Any evidence items that need testing for latent stains such as semen or saliva would need to be forwarded for DNA testing. Our lab currently only performs blood identification.
7C4X2B	No blood was found on Items 1 or 3. Items 2, 4, 5, and 6 contain blood. Presumptive* testing for human origin was positive. Items 1 and 3 contain semen. Due to the color of the item, presumptive* testing for semen was inconclusive on Item 6. Other testing for semen was negative. No semen was found on Items 2, 4 or 5. No analysis was conducted on Item 7. *presumptive tests do not confirm the presence of bodily fluids or human origin
7JZXCW	ITEM 5: - Slight luminescence was detected. - Amylase positive on blood trace. ITEM 6: -Very small luminescence was detected.
8CE398	This laboratory does not test for semen. If a sexual assault is suspected the evidence would be sent to another laboratory for screening and DNA testing. This laboratory does not test for saliva. If saliva is suspected on an evidence item, it would be sent to another laboratory for DNA testing.
8Z7GPE	Only LMG testing is performed at this laboratory.
9BFE63	No hairs observed on Items 1-7.
9K9PJ4	RSID positives for Item 3 and Item 4 were weak positives.
9QVF4E	We do not possess any tests for semen and saliva. These tests are only performed at the lab in [Country]. The only preliminary-test we perform in [Country] is LMG for blood.
9UVNAA	No blood was found on Item 1. Items 2, 4, 5 and 6 contain blood. Presumptive* testing for human origin was positive. No semen was found on Item 2. Items 1 and 3 contain semen. No analysis was performed on Item 7. *presumptive tests do not confirm the presence of bodily fluids or human origin
AEDFMB	My laboratory does not perform screening tests for saliva.
AF8ZDG	As a confirmatory test for human blood our lab uses Bluestar Identi-hem. Unfortunately, it is currently unavailable at our lab, so our lab proceeds directly to DNA profiling. Items submitted to our lab for saliva testing are sampled and forwarded to our DNA department for direct DNA profiling.
ALCYQR	Items 1 and 3: No visual indications of blood were detected, so no testing performed. Items 2, 4, 5, and 6: These items would be reported as having "No indications of seminal fluid", as the analyses stopped after AP testing.
ARYUNP	Serology only

TABLE 2

WebCode	Additional Comments
BLQQBW	No hairs were observed on the submitted items.
BQJU4Z	No hairs were observed on items 5 and 6.
BTQ6Z4	Our lab section does not perform body fluid ID for saliva.
BU9Y2D	ITEM 1: SEMEN IS DETECTED. ITEM 2: HUMAN BLOOD IS DETECTED. ITEM 3: SEMEN IS DETECTED. ITEM 4: HUMAN BLOOD IS DETECTED. ITEM 5: HUMAN BLOOD IS DETECTED. ITEM 6: HUMAN BLOOD IS DETECTED.
C84ZT6	No blood was found on Item 1. Item 1 and 3 contains semen. Items 2, 4, 5 and 6 contain blood. Presumptive* testing for human origin was positive. No semen was found on Items 2 or 4. No analysis was performed on Item 7. *presumptive tests do not confirm the presence of bodily fluids or human origin
CNWLH7	The results and conclusions provided below form my expert opinion, which is based on my scientific knowledge, experience and training. The reporting of blood is based on its appearance and a positive result for a presumptive chemical test for blood. A presumptive test is used to screen for a substance, but cannot confirm its presence. Item 1, Swabs: This item comprised an envelope containing four duplicate swabs. One swab was selected at random for testing. No staining was observed on the swab head. It was tested for the presence of semen and saliva. Semen, containing sperm, was detected on the swab head. No saliva was detected. Item 2, Swabs: This item comprised an envelope containing four duplicate pink-handled swabs. One swab was selected at random for testing. Red-brown staining, with the appearance of blood, was observed on the swab head. It was tested for the presence of blood, semen and saliva. Blood was detected on the swab head. Further testing was carried out on a sample of this bloodstaining to determine if the blood was human in origin. This sample gave a positive result. This test gives positive results for human blood, but can also give positive results for higher primate blood and for blood from members of the Mustelidae family, e.g. ferrets, stoats and weasels. No semen was detected on the swab head. A positive result was obtained in a test for the presence of alpha-amylase, which is a component of saliva. It should be noted that alpha-amylase is present in smaller quantities in some other bodily fluids. Therefore a positive result with this test is referred to as indicating the probable presence of saliva. In my opinion the swab head contained probable saliva. Item 3, Swabs: This item comprised an envelope containing four duplicate blue-handled swabs. One swab was selected at random for testing. No staining was observed on the swab head. It was tested for the presence of semen and saliva. Semen, containing sperm, was detected on the swab head. No saliva was detected. Item 4, Swabs: This item comprised an envelope containing four duplicate yellow-handled swabs. One swab was selected at random for testing. Red-brown staining, with the appearance of blood was observed on the swab. Blood was detected in this stained area. Further testing was carried out on a sample of this bloodstaining to determine if the blood was human in origin. This sample gave a positive result. Therefore the blood was either human blood, higher primate blood, or blood from members of the Mustelidae family. Neither semen nor saliva was detected on the swab. Item 5, Lanyard: This item was an orange lanyard. Red-brown staining, with the appearance of blood and measuring approximately 10 by 52 millimetres, was observed on one side of the lanyard. Blood was detected in this stained area. Further testing was carried out on a sample of this bloodstaining to determine if the blood was human in origin. This sample gave a positive result. Therefore the blood was either human blood, higher primate blood, or blood from members of the Mustelidae family. Due to the circumstances of this case, in my opinion the blood was most likely human blood. No semen was detected on the lanyard, when it was tested. A positive result was obtained for the presence of alpha-amylase in the same area as the bloodstaining. Therefore this area contained probable saliva. In my opinion, given the close correspondence in location of the bloodstaining and probable saliva staining, it is most likely that these two fluids were applied as a mixture. Item 6, Washcloth: This item comprised an envelope containing a blue fabric washcloth, measuring approximately 340 by 300 millimetres at its widest dimensions. A circular red-brown stain, measuring approximately 20 millimetres in diameter, was present in one corner of the washcloth. Testing showed that the red-brown staining was blood. Both sides of the washcloth were tested for the presence of semen; no semen was detected. Both sides of the washcloth were also tested for the presence of saliva. A very weak positive result was obtained with the test used, on the corner diagonally opposite the bloodstaining. Due to the very low strength of this reaction, in my opinion it is inconclusive as to whether or not saliva is present on the washcloth.

TABLE 2

WebCode	Additional Comments
CYTE6B	Saliva tests are not included in our range of tests carried out for routine cases.
DH7FGP	Saliva testing is outside the scope of participant's job function and was not performed on items 1-6.
DWKLQ3	Item 1 - is a swab with light yellow stain and per lab policy, visible blood-like stain (reddish in color) must be present on evidence to infer that blood is indicated. Item 4 - Phadebas result was deemed inconclusive because the detection of amylase was obscured by the red-brown stain color of the extract.
EGQD7A	THE LABORATORY DOES NOT PERFORM SALIVA ANALYSIS. ITEM1 CONCLUSION: NEGATIVE TO HUMAN BLOOD, POSITIVE TO SEMEN. ITEM2: CONCLUSION: POSITIVE TO HUMAN BLOOD, NEGATIVE TO SEMEN. ITEM3: CONCLUSION: NEGATIVE TO HUMAN BLOOD, POSITIVE SEMEN. ITEM 4: CONCLUSION: POSITIVE TO HUMAN BLOOD, NEGATIVE TO SEMEN. ITEM 5: CONCLUSION: POSITIVE TO HUMAN BLOOD, NEGATIVE ALTERNATIVE LIGHT SOURCE, NEGATIVE TO SEMEN. ITEM 6: CONCLUSION: POSITIVE TO HUMAN BLOOD, POSITIVE ALTERNATIVE LIGHT SOURCE, NEGATIVE TO SEMEN.
FGHUXX	Sperm Cell(s) were identified on items Q01 and Q03. Blood was identified on items Q02, Q04, Q05, and Q06. No semen was identified on items Q02, Q04, Q05, and Q06. No blood was identified on items Q01 and Q03. Appropriate samples are being submitted for DNA testing. If a DNA profile is developed, known oral swabs are needed from the victim and any developed suspect(s). RETAINED SAMPLES: Item(s) Q01, Q02, Q03, Q04, Cutting(s) from Q05 and Q06, Tape lift(s) from Q05 and Q06. Also submitted, but not examined: Control Swabs
G6M6E3	My laboratory does not perform screening tests for saliva.
H2GYEM	The control swabs were not opened, not tested. The result for Item 2 should be negative for indications of seminal fluid based on [Laboratory] reporting of results.
HQZC42	My laboratory does not perform screening tests for saliva.
HVCZ4X	the crime scene section of this lab does not perform saliva testing. Items not tested.
J2ETAD	Item 4- possible masking of reaction from Phadebas tube test by heavy blood staining therefore inconclusive. In casework situations identification of saliva would not rely on Phadebas reactions alone (given other substances/ body fluids can contain amylase). The case circumstances , location of Phadebas reactions and results of DNA analysis would also be considered.
JQFVNY	On casework, I would not normally test for PSA if sperm cells were identified.
KLZMZY	Items 2, 4, 5, and 6 contain blood. Presumptive* testing for human origin was positive. No blood was found on Item 1. Items 1 and 3 contain semen. No semen was found on Item 2. No analysis was performed on Item 7. *Presumptive tests do not confirm the presence of bodily fluids or human origin. Results relate only to the items or sampling of items tested.
KMTD7D	No saliva test was performed due to reagents not being available at the lab. Please also note that saliva testing is not a routine test conducted. It is done based on the case scenario.
KQ99JP	Visual examination with an Alternate Light Source was performed on Items 5 and 6. Presumptive testing detected the possible presence of blood on Items 2, 4, 5, and 6. Blood was not detected on Item 1 or 3. Spermatozoa were observed on Items 1 and 3, which confirms the presence of semen. Semen was not detected on Item 2, 4, 5, or 6.
KTFL8J	Item E1: Sample remained. Item E2: Sample remained. Item E3: Sample remained. Item E4: Sample remained. Item E5: ALS was used to search for staining areas. Nothing further detected using ALS. Sample remained. Item E6: ALS was used to search for staining areas. Nothing further detected using ALS. Sample remained.
KXT4PZ	I tested a red-brown stained area and a separate fluorescent area on Item 6 that wasn't visible without the use of the ALS.
L73VR3	ITEM 1.- SEMEN IS DETECTED. ITEM 2.- HUMAN BLOOD IS DETECTED. ITEM 3.- SEMEN IS DETECTED. ITEM 4.- HUMAN BLOOD IS DETECTED. ITEM 5.- HUMAN BLOOD IS DETECTED. ITEM 6.- HUMAN BLOOD IS DETECTED.

TABLE 2

WebCode	Additional Comments
LJAZHV	Items 2, 4, 5, and 6 contain blood. Presumptive* testing for human origin was positive. No blood was found on Item 1. Items 1 and 3 contain semen. No semen was found on Items 2 or 6. No analysis was performed on Item 7. Results relate only to the items or sampling of items tested. *presumptive tests do not confirm the presence of bodily fluids or human origin
MUMKGH	Phadebas paper was used on item 6, the blue washcloth. After allowing the phadebas paper to thoroughly dry overnight, results were inconclusive. I arrived at this conclusion because the washcloth is also blue, and throughout the dry phadebas paper were small specks of blue sporadically throughout the phadebas paper~1-2mm in size. At the site of the RBS and light fluorescence observed under ambient light and ALS respectively, there was little blue coloration present. RSID testing was performed on a small cutting from the center of the stain previously observed. Results from the RSID test were Negative.
N44NZ2	We don't have tests available for Saliva in our lab, so we did not test this. Item 5 and 6 are checked with blue light ~420- 470 nm. There were no luminescent stains on the items, wich indicates the abscent of semen.
NBNQ6J	No hairs were observed on the lanyard or washcloth.
NGE2EG	The control swabs were tested for blood, semen and saliva. Each test came back negative.
NNYWKR	Items 2, 4, 5 and 6 contain blood. Presumptive* testing for human origin was positive. No blood was found on Items 1 or 3. Items 1 and 3 contain semen. No semen was found on Item 2. *Presumptive tests do not confirm the presence of bodily fluids or human origin.
P9FVYK	The laboratory does not have a validated saliva test. Items number 1 and 3 were not red stained, therefore no blood test were performed. Presumptive testing detected the possible presence of blood on Items 1B and 1D and Samples 1E1 and 1F1. Spermatozoa were observed on Items 1A and 1C, which confirms the presence of semen. Semen was not detected on Items 1B and 1D and Samples 1E1 and 1F1. Visual examination with an alternate light source was used on Items 1E and 1F.
PGCTPU	My laboratory does not perform screening tests for saliva.
QD4XJW	The only pre-test we do at our laboratory is LMG for blood.
QEXQEF	Saliva testing not performed in this lab.
QZ2CEL	Item 6 had one area that fluoresced. Due to the victim wiping his mouth it was possibly saliva. The [Laboratory] does not perform saliva testing.
TM7X8N	Items 2, 4, 5, and 6 contain blood. Presumptive* testing for human origin was positive. No blood was found on Item 1. Items 1 and 3 contain semen. No semen was found on Items 2 or 4. No analysis was conducted on Item 7. Results relate only to the items or sampling of items tested. *presumptive tests do not confirm the presence of bodily fluids or human origin
TQPWLH	Due to the scenario given in the case, I took a general swabbing of the lanyard (sub# 5) and preserved that swabbing (5-T2) for possible skin cells from the suspect. Nothing was done with the "Control Swab". It was placed back inside the packaging with the other submissions.
TRY22R	ALS examination of Item 5 showed a faint beige/tan stain between two crease marks in the middle of the lanyard cord. Phadebas testing was performed on both the visible Red/brown stain as well as the stain revealed by ALS due to the narrative of the use of this item.
TXPCBQ	Items 2, 4, 5 and 6 contain blood. Presumptive* testing for human origin was positive. Items 1 and 3 contains semen. No blood was found on Item 1. No semen was found on Item 2. No analysis was conducted on Item 7. *Presumptive tests do not confirm the presence of bodily fluids or human origin.

TABLE 2

WebCode	Additional Comments
UJKFP4	Item 4 swabs were heavily bloodstained and this may have acted to mask any low-level indications of amylase in the Phadebas tube test used. We could potentially have reported this as inconclusive for saliva but in the context of the trial being for 'Body Fluid Identification' purposes, felt it would be more appropriate to report this as negative given there was no clear indication of amylase (and therefore saliva) from this sample. We examined items 5 and 6 using specialized lighting (i.e. a crime-lite) but no significant fluorescence was observed to indicate any body fluids away from the areas of visible staining. We do however acknowledge that body fluids do not consistently fluoresce and so deemed this test inconclusive. We applied mini-tapes to item 5 in keeping with an approach we would take in casework. Item 6's amylase test was considered inconclusive as although a weak reaction was observed, this was closely associated with the blood and it was considered potentially to be a false positive. NB. The items were not subjected to DNA testing which may have indicated that mixtures of body fluids were present. This would be something we would consider in real casework, e.g. in relation to the item 4 swabs where it is arguable the blood could mask visual indications of amylase, but was not felt appropriate for the purposes of the trial. Item 6 could similarly have been subjected to DNA testing to see if we might have indications of body fluids from two different individuals (which may have tended to support that the amylase weak reactions observed were of a body fluid). Again however, in the context of this being a trial, and with uncertainty as to how the trial may have been set up, this was not undertaken and reporting the test as inconclusive seemed more appropriate.
WL47RJ	In regards to the Blood Test for Item 1, the result was reported as Inconclusive because the swabs were documented to have an off-white appearance and a negative KM test was obtained. Even though the HemaTrace result was positive, per our Technical Manual: "The HemaTrace test can produce positive results from other human body fluids that may contain hemoglobin at low levels. Analysts should be cautious of inferring that blood is indicated from samples that do not exhibit some other indication (i.e. blood-like color) of blood being present."
WVQ6HP	My laboratory does not perform screening tests for saliva.
X32Q3D	Kastle-Meyer Test for Blood: Presumptive testing detected the possible presence of blood on Items 2, 4, 5, and 6. Visual examination did not detect the presence of brown stains on items 1 and 3. Microscopic Examination for Spermatozoa and p30 Testing: Spermatozoa were observed on Items 1 and 3, which confirms the presence of semen. Semen was not detected on Items 2, 4, 5, and 6. Notes: Visual examination with an alternate light source was performed on items 5 and 6.
X7HDGH	Reporting item 1 as inconclusive for blood because it tested positive for Hematrace (twice) while testing negative with KM and not having any blood-like discoloration.
XEYPW4	Examination of Items 5 and 6 are outside the scope of participant's job function.
Y94GNF	Crime scene personnel not trained on saliva or semen screening.
YDXBLG	Regarding Item 5 and ALS testing, Item 5 had two samples taken from it: a swabbing of the plastic piece and a cutting of the red-brown-stained fabric. Although the results are positive for ALS, it was the plastic piece that fluoresced, and the plastic piece was negative for both semen and saliva; the red-brown stain did not fluoresce, despite it being positive for saliva.....Regarding Item 6 and ALS testing, Item 6 had three samples taken from it: two different cuttings of non-red-brown-stained areas, whose fluorescence was difficult to discern with certainty, and a cutting of the red-brown-stained area, which seemed to be possibly fluorescent around its outer edge. These areas were AP-, and therefore, not semen. These areas were inconclusive for saliva.....For phadebas (saliva) testing, the blood from Items 4 and 6 (the red-brown-stained area cutting) created supernatants that were brownish-colored, which could mask any blue color that could be there from a possibly positive stain. Additionally, the two cuttings of the non-red-brown-stained areas on Item 6 created supernatants that were blue, likely from the leached dye in the substrate, which could mask any blue color that could be there from a possibly positive stain.
YNE396	Regarding Item 6 (washcloth), we had detected 3 different traces. Two of them are latent ones and were observed through Crimescope CS -16 (420-470 nm wavelength), while the third is an evident blood stain (tested positive with Bluestar Hexagon OBTI, while tested negative with PSA and Amylase as described). Both of latent traces had been negative to all of the 3 specific tests (OBTI, PSA and Amylase as described).

TABLE 2

WebCode	Additional Comments
YNVC27	*Item 6, Phadebas Press Test. Atypical results: very weak, with no visible feathering or bleeding of the dye around the reaction sites.
Z9YG2K	My laboratory does not perform screening tests for saliva.
ZYAT4W	Our laboratory does not test for semen or saliva. The only body fluid we test for is blood, therefore the only test I performed was a blood test (Phenolphthalein).

-End of Report-  
(Appendix may follow)

## Test No. 23-5782: Body Fluid Identification

DATA MUST BE SUBMITTED BY **Nov. 20, 2023, 11:59 p.m. EST** TO BE INCLUDED IN THE REPORT

Participant Code: U1234A

WebCode: U6MEYN

The Accreditation Release section can be accessed by using the "Continue to Final Submission" button above. This information can be entered at any time prior to submitting to CTS.

**Scenario (this scenario applies only to Items 3-6):**

Police are investigating a sexual assault involving a male victim and male suspect. The victim reports they were together in the suspect's office. The victim was punched in the face by the suspect and the suspect tied his lanyard around the victim's mouth. The suspect then masturbated and ejaculated on the floor. The victim was able to remove the lanyard, wipe his face with a washcloth from the suspect's office, and escape. The following items are being submitted to your laboratory for examination: One envelope containing four duplicate swabs (Item 3) taken from a stain on the office floor. One envelope containing four duplicate swabs (Item 4) taken from a stain on the suspect's hand. One stained lanyard (Item 5) used by the suspect to tie the victim. One washcloth (Item 6) used by the victim directly following the assault.

**Items Submitted (Sample Pack BF12):**

Item 1: Swab with white handle (provided as 4 replicates) - please examine this item for body fluids tested in normal laboratory casework; scenario not applicable

Item 2: Swab with pink handle (provided as 4 replicates) - please examine this item for body fluids tested in normal laboratory casework; scenario not applicable

Item 3: Swab with blue handle (provided as 4 replicates) - taken from a stain on the office floor

Item 4: Swab with yellow handle (provided as 4 replicates) - taken from a stain on the suspect's hand

Item 5: One lanyard used by suspect to tie victim during the assault

Item 6: One washcloth used by victim after the assault

Control swabs (unused): One package to be used for substrate control testing, as necessary.



**Instructions:**

Please indicate the overall conclusion for each body fluid using the check boxes provided. List the tests performed and the results of each test. Item descriptions (as applicable) are provided at the top of each Item page.

For the purposes of this test, all swabs within an Item have been prepared as duplicates in order to provide sufficient sample quantity. A package of unused swabs is also provided for control purposes.

Please use the abbreviations listed in this response key to fill in the Screening Test tables. This is not an all inclusive list of tests, and should not be used to determine what tests should be performed. TESTS NOT ON THIS LIST MAY BE USED FOR SCREENING.

Test	Abbreviation	Test	Abbreviation
Acid Phosphatase	AP	Ortho-tolidine	O-tol
Alternate Light Source	ALS	Phenolphthalein-Tetramethyl benzidine	PTMB
Kastle Meyer	KM	Prostate Specific Antigen	PSA
Leucomalachite Green	LMG	Rapid Stain Identification	RSID
Microscopic	Micro	Tetramethyl benzidine	TMB

**Example**

**Blood Tests**

		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for blood?	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	HemaTrace	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	
	KM	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Pos./Ind. - Positive/Indicative; Neg. - Negative; Inc. - Inconclusive; NT - Not Tested

\*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

### Item 1: Cotton tipped swabs with white plastic handle

Four replicate white handled swabs. Scenario is not applicable for this Item. Utilizing your laboratory policies, please examine this item for body fluids tested in normal laboratory casework.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for blood?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for semen?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for saliva?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

\*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Pos./Ind.	Neg.	Inc.
Please report Other Tests as 'Body Fluid - Test Performed'				
	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

### Item 2: Cotton tipped swabs with pink plastic handle

Four replicate pink handled swabs. Scenario is not applicable for this Item. Utilizing your laboratory policies, please examine this item for body fluids tested in normal laboratory casework.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for blood?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for semen?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for saliva?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

\*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Pos./Ind.	Neg.	Inc.
Please report Other Tests as 'Body Fluid - Test Performed'		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	<input type="text"/>			

### Item 3: Cotton tipped swabs with blue plastic handle

Four replicate blue handled swabs taken from a stain on the office floor.

<b>Blood Tests</b>		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for blood?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

<b>Semen Tests</b>		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for semen?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

<b>Saliva Tests</b>		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for saliva?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

\*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

<b>Other Tests</b>		Pos./Ind.	Neg.	Inc.
	Please report Other Tests as 'Body Fluid - Test Performed'			
	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

### Item 4: Cotton tipped swabs with yellow plastic handle

Four replicate yellow handled swabs taken from a stain on the suspect's hand.

<b>Blood Tests</b>		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for blood?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

<b>Semen Tests</b>		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for semen?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

<b>Saliva Tests</b>		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for saliva?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

\*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

<b>Other Tests</b>		Pos./Ind.	Neg.	Inc.
Please report Other Tests as 'Body Fluid - Test Performed'				
	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

**Item 5: One lanyard used by the suspect to tie the victim during the assault**

The item of evidence is a lanyard used by the suspect to tie the victim during the assault.

<b>Blood Tests</b>		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for blood?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

<b>Semen Tests</b>		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for semen?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

<b>Saliva Tests</b>		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for saliva?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

\*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

<b>Other Tests</b>		Pos./Ind.	Neg.	Inc.
Please report Other Tests as 'Body Fluid - Test Performed'				
	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

**Item 6: One washcloth used by the victim after the assault**

The item of evidence is a washcloth used by the victim after the assault.

<b>Blood Tests</b>		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for blood?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

<b>Semen Tests</b>		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for semen?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

<b>Saliva Tests</b>		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for saliva?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

\*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

<b>Other Tests</b>		Pos./Ind.	Neg.	Inc.
	Please report Other Tests as 'Body Fluid - Test Performed'			
	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

### Additional Comments

*Please note: Any additional formatting applied in the free form space below will not transfer to the Summary Report and may cause your information to be illegible. This includes additional spacing and returns that present your responses in lists and tabular formats.*



## RELEASE OF DATA TO ACCREDITATION BODIES

The Accreditation Release is accessed by pressing the "Continue to Final Submission" button online and can be completed at any time prior to submission to CTS.

CTS submits external proficiency test data directly to ASCLD/LAB, ANAB, and/or A2LA. Please select one of the following statements to ensure your data is handled appropriately.

- This participant's data is intended for submission to ASCLD/LAB, ANAB, and/or A2LA. (Accreditation Release section below must be completed.)
- This participant's data is **not** intended for submission to ASCLD/LAB, ANAB, and/or A2LA.

Have the laboratory's designated individual complete the following steps only if your laboratory is accredited in this testing/calibration discipline by one or more of the following Accreditation Bodies.

### Step 1: Provide the applicable Accreditation Certificate Number(s) for your laboratory.

ANAB Certificate No.   
(Include ASCLD/LAB Certificate here)

A2LA Certificate No.

### Step 2: Complete the Laboratory Identifying Information in its entirety.

Authorized Contact Person and Title

Laboratory Name

Location (City/State)