



Body Fluid Identification Test No. 14-579 Summary Report

This proficiency test was sent to 219 participants. Each participant received a sample pack containing six items for examination. Participants were asked to examine each item and report the body fluid(s) identified. Data were returned from 151 participants (69% response rate) and are compiled into the following tables:

	<u>Page</u>
<u>Manufacturer's Information</u>	2
<u>Summary Comments</u>	3
<u>Key to Screening Results</u>	4
<u>Table 1: Screening Conclusions & Test(s) Performed</u>	5
<u>Table 2: Serology Screening Results</u>	101
<u>Table 3: Additional Comments</u>	189
<u>Appendix: Data Sheet</u>	195

This report contains the data received from the participants in this test. Since these participants are located in many countries around the world, and it is their option how the samples are to be used (e.g., training exercise, known or blind proficiency testing, research and development of new techniques, etc.), the results compiled in the Summary Report are not intended to be an overview of the quality of work performed in the profession and cannot be interpreted as such. The Summary Comments are included for the benefit of participants to assist with maintaining or enhancing the quality of their results. These comments are not intended to reflect the general state of the art within the profession.

Participant results are reported using a randomly assigned "WebCode". This code maintains participant's anonymity, provides linking of the various report sections, and will change with every report.

Manufacturer's Information

Each sample pack consisted of six items for testing, four of which were swabs and two of which were evidence items. Participants were asked to examine Items 1 and 2 for the body fluids that they typically encounter in casework. Items 3-6 were to be examined for scenario specific body fluids.

SAMPLE PREPARATION: The stains on Item 1 (white swabs) and Item 4 (yellow swabs) were prepared using human semen from a single donor. The stain on Item 2 (pink swabs) was prepared using human saliva. The stains on Item 3 (blue swabs) and Item 5 (black headpiece) were prepared using human whole blood which was drawn into citric acid preservative bloodbank bags. The stain on Item 6 (red women's underwear) was a mixture prepared from human whole blood and semen. The raw semen was from a single donor, stored frozen at a temperature below -20°C. The semen was thawed and mixed 1:1 with TAE buffer, and then mixed 1:1 with the blood. This mixture was used to spot the Item 6 substrate. Stains from different sources were prepared at separate times and were packaged once they were thoroughly dried.

SAMPLE SET ASSEMBLY: Each sample set contained a package of unused swabs, questioned swabs (Items 1, 2, 3 and 4) and two questioned items of evidence, a black headpiece (Item 5) and a pair of red women's underwear (Item 6) which were placed in a pre-labeled sample pack envelope. The sealed sample pack envelopes were then packaged in pre-labeled Heat Seal envelopes and sealed. This process was repeated until all of the sample sets were prepared.

VERIFICATION: Laboratories that conducted predistribution analysis of the samples reported consistent body fluid conclusions.

Sample Preparation Schedule

<u>Item</u>	<u>Substrate</u>	<u>Fluid</u>	<u>Volume</u>
1	4 White swabs	Semen	100 µl/swab
2	4 Pink swabs	Saliva	100 µl/swab
3	4 Blue swabs	Blood	100 µl/swab
4	4 Yellow swabs	Semen	100 µl/swab
5	Black Du-rag	Blood	100 µl
6	Red women's underwear	Blood and Semen	100 µl

Completed sample sets were stored at -20°C until shipment on September 09, 2014, following receipt of predistribution results.

Summary Comments

The 14-579 Body Fluid Identification test was designed to allow participants to assess their proficiency in the identification of stains that may contain dried body fluids. Each participant received six items for testing, four of which were swabs and two were items of clothing. No scenario was provided for Items 1 and 2; participants were requested to test for body fluids that they typically encounter in casework. Scenarios were provided for the remaining items; participants were asked to treat these as they would casework.

Item 1 consisted of four swabs containing semen. For blood, the majority (91%) of participants reporting a conclusion were negative, 10 were positive/indicative and one was inconclusive. For semen, the majority (98%) of participants reporting a conclusion for semen were positive/indicative and three were negative. For saliva, a consensus conclusion could not be established, there were 85 participants that were negative, 20 that were positive/indicative and 10 that were inconclusive.

Item 2 consisted of four swabs containing saliva. For blood, the majority (98%) of participants reporting a conclusion were negative, one was positive/indicative and one was inconclusive. For semen, the majority (99%) of participants reporting a conclusion were negative and one was positive/indicative. For saliva, all participants reporting a conclusion were positive/indicative.

Item 3 consisted of four swabs containing blood. For blood, the majority (98%) of participants reporting a conclusion were positive/indicative and three were negative. For semen, all participants reporting a conclusion were negative. For saliva, the majority (92.6%) of participants reporting a conclusion were negative, five were inconclusive and one was positive/indicative.

Item 4 consisted of four swabs containing semen. For blood, the majority (94%) of participants reporting a conclusion were negative, five were positive/indicative and one was inconclusive. For semen, the majority (98%) of participants reporting a conclusion were positive/indicative and three were negative. For saliva, the majority (80%) of participants reporting a conclusion were negative, 11 were positive/indicative and nine were inconclusive.

Item 5 consisted of one black headpiece containing blood. For blood, all participants reporting a conclusion were positive/indicative. For semen, the majority (97%) of participants reporting a conclusion were negative, two were positive/indicative and one was inconclusive. For saliva, the majority (97%) of participants reporting a conclusion were negative and two were inconclusive.

Item 6 consisted of one red women's underwear containing a mixture of blood and semen. For blood, all participants reporting a conclusion were positive/indicative. For semen, the majority (97%) of participants reporting a conclusion were positive/indicative, three were negative and one was inconclusive. For saliva, the majority (95%) of participants reporting a conclusion were negative, two were positive/indicative and two were inconclusive.

Key for Screening Tests Used

Participants were asked to use, where possible, the following chart of abbreviated screening test names. This was not an all inclusive list and was not designed to determine what tests should be performed. Participants were advised that tests not on this list may be used for screening.

<u>Test</u>	<u>Abbreviation</u>
Acid Phosphatase	AP
Kastle Meyer	KM
Leucomalachite Green	LMG
Microscopic	Micro
Ortho-tolidine	O-tol
Phenolphthalein-Tetramethyl benzidin	PTMB
Prostate Specific Antigen	PSA
Tetramethyl benzidine	TMB

Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
27CG2C	KM	AP MICRO	PHADEBAS
28BLGK	Negative LMG	Positive/ Indicative AP Micro P30	Not Tested
2EMTZU	Negative KM	Positive/ Indicative AP Micro	Positive/ Indicative RSID
2FHERZ	Negative O-tol	Positive/ Indicative micro	Negative RSID Saliva
2PQVGK	Negative KM <u>Other Tests Performed</u> Sperm Search	Positive/ Indicative PSA	Negative RSID
2VEVGY	Negative O-tol	Positive/ Indicative AP Micro	Positive/ Indicative RSID
3FF4BN	Negative KM	Positive/ Indicative AP Micro	Negative Phadebas
3MXG82	Negative KM <u>Other Tests Performed</u> Urine - Creatinine Urine - DMAC Faeces - Edelmanns	Positive/ Indicative AP Micro	Inconclusive Phadebas tube test
3P4YCN	Negative KM	Positive/ Indicative AP Micro	Negative RSID-Saliva

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
4APZYB	Negative KM	Positive/ Indicative AP Micro	Negative Amylase (radial diffusion)
4BJLQG	Negative LMG	Positive/ Indicative AP PSA Micro	Not Tested
4HXUT4	Negative LMG	Positive/ Indicative AP	Not Tested
4JVXRK	Negative KM	Positive/ Indicative AP Micro	Positive/ Indicative RSID
68YCXV	Negative OT	Positive/ Indicative Micro	Negative RSID
6HVCDK	Negative Visual	Positive/ Indicative AP micro	Negative RSID-Saliva
6PZ4XB	Negative KM	Positive/ Indicative AP Micro	Positive/ Indicative Phadebas
6UWLGW	Negative O-tol	Positive/ Indicative Micro	Inconclusive RSID
6ZGMFL	Not Tested Swab was not red stained	Positive/ Indicative Micro	Not Tested
7M6VAR	Not Tested No visible stain observed.	Positive/ Indicative AP, Micro	Negative Amylase
7MNLCK	Not Tested Swabs have no red stains	Positive/ Indicative Microscopic	Not Tested Per Agency's Policy

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
7XNELB	Negative KM	Positive/ Indicative AP P30 KPIC	Negative Phadebas
8C6JGB	Negative KM <u>Other Tests Performed</u> Urine - Creatinine Feces - Urobilinogen	Positive/ Indicative AP Micro	Negative Phadebas
8NUMDA	Positive/ Indicative LMG HemaTrace	Positive/ Indicative AP Micro PSA ALS	Not Tested This laboratory does not perform testing for saliva.
8PP7HW	Negative LMG test	Not Tested	Inconclusive Phadebas tube test
8XLW8D	Negative LMG Visual	Positive/ Indicative AP P30 Micro	Negative RSID - Saliva
8YD3JW	Not Tested Yellowish staining observed, but no visible blood-like staining observed, therefore no presumptive test for blood performed	Positive/ Indicative AP: x4 AP positive Micro: 2+ heads & tails, negative for NECs	Negative Phadebas: x4 Phadebas negative
9B677N	Negative KM	Positive/ Indicative AP Micro	Positive/ Indicative RSID
9GQWTH	Negative KM	Positive/ Indicative AP Micro	Negative RSID
9HLHJN	Negative KM	Positive/ Indicative AP Micro	Positive/ Indicative RSID

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
9JH363	Negative KM	Positive/ Indicative AP PSA	Negative Phadebas Press
9Q3XCX	Negative KM	Positive/ Indicative AP MICRO	Positive/ Indicative AMYLASE DIFFUSION
9UE4NY	Not Tested Visual exam - no staining observed	Positive/ Indicative AP Micro P30	Negative Amylase Diffusion
9W3CXZ	Negative KM	Positive/ Indicative AP PSA Semiquant	Negative Phadebas Press
AWHZRY	Negative KM	Positive/ Indicative AP PSA Micro	Negative Phadebas
B2DAF7	Positive/ Indicative Hemastix - wk pos 25 secs Hematrace - pos 3 mins	Positive/ Indicative AP - pos 8 secs Micro - pos p30 - pos 1 min	Negative RSID - neg 8.5 min
B33FWW	Not Tested Not tested: No red-brown staining observed on the Items.	Positive/ Indicative AP RSID Micro	Not Tested
BA8DNA	Negative Visual LMG	Positive/ Indicative AP P30 Spermatozoa (Micro)	Negative RSID
BHXXRM	Not Tested Not tested due to lack of reddish-brown color on swabs.	Positive/ Indicative Item 1 AnB: AP, Micro Item 1 CnD: AP, Micro	Negative Item 1 AnB: RSID Item 1 CnD: RSID
BMUKKB	Negative KM	Positive/ Indicative AP Seratec Micro	Negative RSID

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
BNKQXU	Positive/ Indicative Hematrace	Positive/ Indicative PSA Micro	Negative RSID
BU8H89	Negative LMG	Positive/ Indicative AP PSA Micro	Not Tested
BYM4XV	Not Tested no apparent bloodstain observed	Positive/ Indicative AP Micro	Negative Amylase diffusion
C4UYB4	Positive/ Indicative LMG HemaTrace	Positive/ Indicative ALS AP Micro P30	Not Tested Saliva testing is not offered as a service in this laboratory.
CGLJZ8	Negative KM <u>Other Tests Performed</u> Urine - Creatinine Urine - DMAC Fecal material - Urobilinogen	Positive/ Indicative AP Micro	Negative Phadebas Amylase
DL73AZ	Negative KM	Positive/ Indicative AP PSA Micro	Negative Phadebas Press
EB38QV	Negative KM	Positive/ Indicative PSA MICRO	Negative PHADEBAS
EB9GR9	Negative TMB	Positive/ Indicative AP	Negative RSID Saliva Field Kit
EFF3F6	Negative KM <u>Other Tests Performed</u> Urine - Creatinine, DMAC Fecal Material - Urobilinogen	Positive/ Indicative AP (Seri AP) Micro (Xmas Stain)	Negative Amylase (Phadebas)

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
ERDV3B	Negative KM	Positive/ Indicative AP PSA Micro	Negative RSID
EVPZDB	Not Tested White scene or kit swabs not tested for blood in this lab	Positive/ Indicative Micro	Not Tested Saliva testing not currently being done in this lab
EXGQ74	Negative KM <u>Other Tests Performed</u> Urine - Creatinine Fecal Material - Urobilinogen	Positive/ Indicative AP (SAP) PSA (ABA) Micro (Xmas)	Negative Phadebas Amylase
F42C9M	Negative Phenolphthalein	Negative RSID	Negative SaligaE
FVMNEL	Negative O-Tol	Positive/ Indicative Micro	Negative RSID
FXB6NT	Negative KM	Positive/ Indicative AP P30 Micro	Negative Phadebas
G74RZB	Negative O-tol	Positive/ Indicative AP Micro	Positive/ Indicative RSID
G8XFEZ	Positive/ Indicative TMB HemaTrace Haemochromogen <u>Other Tests Performed</u> Ouchterlony (Human origin)	Positive/ Indicative AP PSA Microscopic	Not Tested
G9NLVQ	Negative LMG	Positive/ Indicative AP	Negative PHADEBAS PAPER TEST
GFC6FU	Negative KM	Positive/ Indicative AP PSA	Negative Phadebas

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
GPVJEP	Negative KM Visual exam	Positive/ Indicative AP H&E Slide	Negative Phadebas paper
GV43GN	Not Tested no staining on swab	Positive/ Indicative AP Micro	Negative Amylase diffusion
GZXZM8	Not Tested *swabs were not red stained	Positive/ Indicative Micro	Not Tested
H7J7Z	Negative KM <u>Other Tests Performed</u> Urine - Creatinine Fecal Material - Urobilinogen	Positive/ Indicative AP (Map) Micro (Xmas)	Negative a-Amylase
H9ADNB	Negative O-tol	Positive/ Indicative ALS AP Micro PSA	Positive/ Indicative RSID
HDYUYF	Positive/ Indicative KM Abacard Hematrace	Positive/ Indicative AP Micro Abacard p30	Not Tested
HEDJ2T	Negative KM ABAcad HemaTrace	Positive/ Indicative AP ABAcad P30 Micro	Negative Phadebas test
HHBZEM	Not Tested A red stain was not observed.	Positive/ Indicative AP Micro PSA ALS	Negative Amylase diffusion
HQ63HJ	Negative O-tol	Positive/ Indicative Micro	Negative RSID
HWRXDW	Negative KM	Positive/ Indicative AP Micro	Negative Radial Diffusion

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
J23ATA	Negative KM	Positive/ Indicative RSID-Semen PSA	Not Tested Saliva test is not done if item is tested positive for blood or semen. No saliva test was done for Item 1 as it was tested positive for semen
J94HFP	Not Tested	Positive/ Indicative AP Micro	Negative Phadebas
JHVYAZ	Negative KM <u>Other Tests Performed</u> Urine - Creatinine Feces - Urobilinogen	Positive/ Indicative AP Micro	Negative Phadebas
JJM2UU	Negative KM <u>Other Tests Performed</u> Urine - DMAC	Positive/ Indicative AP Micro	Inconclusive Phadebas tube test Micro
JNFWZ9	Negative K-M <u>Other Tests Performed</u> Urine - Visual & DMAC Feces - Visual & Microscopy	Positive/ Indicative AP Sperm search w/ microscope	Negative RSID-Saliva
JPQ3UN	Negative LMG	Positive/ Indicative AP	Negative PHADEBAS PAPER TEST
JRWLWQ	Negative Visual Exam	Positive/ Indicative Fast blue (+ve) (AP) Micro (+ve)	Not Tested Current protocol does not allow for Saliva testing on swabs.
JTMRDG	Negative KM	Positive/ Indicative ALS AP RSID-Semen Microscopy	Negative ALS RSID-Saliva

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
JYWP6A	Negative KM	Positive/ Indicative RSID PSA	Not Tested Saliva test is not done if item is tested positive for either blood or semen. Item 1 is tested positive for semen, thus no saliva test was performed.
K6XJEV	Negative KM	Positive/ Indicative AP Micro	Positive/ Indicative Phadebas
K7UYQR	Negative Visual KM	Positive/ Indicative AP Micro	Positive/ Indicative Phadebas
KQ3XGF	Negative KM	Positive/ Indicative AP Micro ABAcad p30	Inconclusive radial diffusion
KUE3RF	Negative O-TOL	Positive/ Indicative AP MICRO	Negative RSID
KXYWNG	Negative HemaTrace	Positive/ Indicative ABAcad P30	Negative Saligae
KZ3R7Y	Negative LMG	Positive/ Indicative AP PSA Microscopic	Not Tested
LXURT6	Negative Visual exam	Positive/ Indicative AP Micro	Negative RSID
M9CTN6	Negative Visual KM	Positive/ Indicative AP Micro	Negative RSID
MED8A7	Negative O-Tol	Positive/ Indicative ALS AP Micro	Negative ALS RSID

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
MEWRED	Not Tested Swabs were not tested for blood - does not display the characteristic reddish-brown stain associated with blood.	Positive/ Indicative AP, Micro AP, Micro AP, Micro AP, Micro	Negative RSID Saliva RSID Saliva RSID Saliva RSID Saliva
MQAX88	Negative KM	Positive/ Indicative RSID PSA	Not Tested
MY39JM	Negative KM	Positive/ Indicative AP	Not Tested *saliva testing not performed in this lab
MZHQHM	Negative Hemastix	Positive/ Indicative Fast Blue (AP) Micro	Not Tested Current protocol does not allow for Saliva testing on swabs.
N3ZCLU	Negative KM	Positive/ Indicative AP p30 KPIC	Inconclusive Phadebas
N9PM3J	Not Tested No red-brown discoloration	Positive/ Indicative AP Micro PSA	Negative Amylase Diffusion
N9PRB2	Negative KM	Positive/ Indicative AP Micro (sperm search) PSA	Negative RSID
NQ6YNH	Not Tested	Positive/ Indicative AP Micro	Negative RSID-Saliva
P8P4F9	Not Tested No reddish/brown stain observed during visual inspection.	Positive/ Indicative AP Microscopy	Negative RSID Saliva
PCGF7Y	Negative KM	Positive/ Indicative AP Micro	Negative Amylase testing by Radial Diffusion

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
PDCWKG	Negative KM O-tol ABA Card	Positive/ Indicative AP	Negative SALlgAE
PJ2EYR	Positive/ Indicative LMG Hematrace	Positive/ Indicative AP Micro P30 ALS	Not Tested Testing not offered.
PJ7RQU	Negative KM	Positive/ Indicative AP Seratec Micro	Negative RSID
PKWW3E	Negative KM Visual	Positive/ Indicative AP Slide made (Micro)	Inconclusive Phadebas
PMLJF3	Not Tested *NT because no red brown stain observed.	Positive/ Indicative AP Micro	Positive/ Indicative RSID
PNF7UR	Inconclusive TMB HemaTrace <u>Other Tests Performed</u> Human - Ouchterlony Faeces - Price (urobilin)	Positive/ Indicative AP PSA Micro	Negative Phadebas (amylase)
QADXPB	Negative Hemastix (TMB)	Positive/ Indicative Microscopy	Negative RSID a-amylase
QP6VRD	Negative KM	Positive/ Indicative BCIP Microscopic Sperm Search	Positive/ Indicative Amylase Diffusion
R9VBPE	Negative LMG	Positive/ Indicative Microscopy AP	Not Tested

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
RBKG8Q	Negative KM	Positive/ Indicative AP P-30 KPIC	Negative Phadebas
RDTPGB	Negative O-tol	Positive/ Indicative Micro	Negative RSID
RNACNQ	Positive/ Indicative TMB HemaTrace Haemachromogen[sic] <u>Other Tests Performed</u> Faeces - Prices	Positive/ Indicative AP PSA Micro	Negative Amylase
RNQTHJ	Negative Visually	Positive/ Indicative AP Micro	Not Tested
RVQ877	Not Tested No visual apparent blood staining	Positive/ Indicative AP Micro PSA	Negative Amylase diffusion
T2V4RR	Negative KM	Positive/ Indicative AP Seratec Micro	Negative RSID
T7LPCY	Negative KM	Positive/ Indicative Microscopic	Not Tested No currently validated procedures for saliva test.
TJBFHW	Negative KM	Positive/ Indicative AP Micro	Positive/ Indicative RSID
TKLJQE	Negative LMG	Positive/ Indicative AP Micro	Negative Amylase Diffusion
TX4EYX	Negative KM	Positive/ Indicative Micro	Not Tested

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
U3J3HU	Negative Kastle-Meyer Test	Positive/ Indicative Acid Phosphatase Microscopic	Positive/ Indicative RSID Saliva
UAYML8	Negative KM test	Positive/ Indicative ALS AP Test RSID-Semen Microscopy	Negative ALS RSID-Saliva
UAZKRV	Negative KM	Positive/ Indicative Micro	Not Tested
UF9YH8	Negative o-tol	Positive/ Indicative ALS AP Micro	Negative ALS RSID
UHA7HB	Negative KM	Positive/ Indicative AP	Negative Phadebas
UPT6GW	Negative Visual examination	Positive/ Indicative AP PSA	Positive/ Indicative Phadebas amylase
UPU6B6	Not Tested	Positive/ Indicative AP MICRO	Inconclusive AMYLASE PLATE
V4BH9W	Negative KM	Positive/ Indicative ALS AP Micro	Negative RSID ALS
V7F8B6	Negative O-TOL	Positive/ Indicative AP Micro	Positive/ Indicative RSID
V9Z287	Negative KM	Positive/ Indicative AP	Negative Phadebas
VKV7HR	Negative Tetrabase	Positive/ Indicative Micro	Negative RSID

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
VM2TTB	Negative LMG	Positive/ Indicative AP Micro PSA	Negative Amylase Diffusion
W4LNZR	Negative KM	Positive/ Indicative AP Micro exam	Negative Amylase with radial diffusion
WBNQ8D	Negative KM	Positive/ Indicative AP Micro P30	Negative Phadebas
WF2Q7P	Negative KM TEST HemaTrace <u>Other Tests Performed</u> Urine RSID	Positive/ Indicative AP Micro H+E stain	Negative Rsid saliva
WM2X2R	Negative PTMB Slight yellowish stain present. <u>Other Tests Performed</u> Urine - Creatinine	Positive/ Indicative Acid Phosphatase (AP) microscopic (micro)	Negative Phadebas
WXGZ8Q	Negative K-M	Positive/ Indicative AP Micro	Positive/ Indicative RSID Saliva
XD83L9	Negative KM	Positive/ Indicative AP	Not Tested Not competent in saliva testing
XHZBMU	Not Tested	Positive/ Indicative Micro	Not Tested
XJC9UK	Positive/ Indicative KM	Positive/ Indicative AP Micro	
XXUUMQ	Negative KM	Positive/ Indicative AP Micro	Positive/ Indicative RSID

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
Y28ZFH	Positive/ Indicative LMG HemaTrace	Positive/ Indicative AP Micro PSA ALS	Not Tested This laboratory does not currently offer saliva testing
Y8GAQK	Negative KM <u>Other Tests Performed</u> Urine - Creatinine Urine - Urease Fecal Material - Urobilinogen	Positive/ Indicative AP Micro - Xmas tree stain	Inconclusive Phadebas Amylase
Y8WHUU	Negative LMG	Positive/ Indicative AP PSA Micro	Not Tested
Y9VLM6	Negative Phenolphthalein	Negative RSID	Negative SaligaE
YBDZQ6	Negative Visual KM	Positive/ Indicative AP Microscopic exam - spermatazoa[sic] identified	Not Tested Phadebas tube test not accredited
YLWGGQ	Not Tested Not a brown stains[sic]	Positive/ Indicative Microscopic	Not Tested
YMCGYU	Negative KM	Positive/ Indicative RSID-Semen PSA	Not Tested Not tested for saliva since item tested positive for semen.
YN3XWA	Negative Hemastix -ve @ 60 secs	Positive/ Indicative Micro	Negative RSID -ve @ 8½ mins
YQT4CZ	Not Tested No red brown staining	Positive/ Indicative AP Micro	Inconclusive Radial Diffusion
YTUAPK	Negative LMG	Positive/ Indicative AP PSA Micro	Not Tested

Table 1 - Item 1

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
YZDAA4	Negative Phenolphthalein	Negative RSID	Negative Saligae
ZCTE3K	Negative Visual	Positive/ Indicative AP P30 Microscopic	Negative RSID
ZKJZ9J	Negative KM	Positive/ Indicative AP Seratec Micro	Negative RSID
ZRYNXE	Negative KM <u>Other Tests Performed</u> URINE, TEST PERFORMED CREATININE - FAINT	Positive/ Indicative AP MICRO	Positive/ Indicative PHADEBAS (TUBE METHOD)
ZYFE47	Not Tested No red discoloration	Positive/ Indicative AP Micro PSA	Negative Amylase Diffusion

Response Summary for Table 1 - Item 1

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	10	146	20
Negative	115	3	85
Inconclusive	1	0	10
Not Tested (NT)	24	1	34

Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
27CG2C	KM	AP MICRO	PHADEBAS
28BLGK	Negative LMG	Negative AP	Not Tested
2EMTZU	Negative KM	Negative AP Micro	Positive/ Indicative RSID
2FHERZ	Negative O-tol	Negative RSID Semen Micro	Positive/ Indicative RSID Saliva
2PQVGK	Negative KM	Negative PSA	Positive/ Indicative RSID
2VEVGY	Negative O-tol	Negative AP Micro RSID	Positive/ Indicative RSID
3FF4BN	Negative KM	Negative AP Micro	Positive/ Indicative Phadebas
3MXG82	Negative KM	Negative AP Micro	Positive/ Indicative Phadebas
	<u>Other Tests Performed</u>		
	Urine - Creatinine		
	Urine - DMAC		
	Faeces - Edelmanns		
3P4YCN	Negative KM	Negative AP Micro	Positive/ Indicative RSID-Saliva
4APZYB	Negative KM	Negative AP Micro	Positive/ Indicative Amylase (radial diffusion)

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
4BJLQG	Negative LMG	Negative AP	Not Tested
4HXUT4	Negative LMG	Negative AP	Not Tested
4JVXRK	Negative KM	Negative AP Micro	Positive/ Indicative RSID
68YCXV	Negative OT	Negative Micro RSID	Positive/ Indicative RSID
6HVCDK	Negative Visual	Negative AP	Positive/ Indicative RSID-Saliva
6PZ4XB	Negative KM	Negative AP Micro	Positive/ Indicative Phadebas
6UWLGW	Negative O-tol	Negative Micro RSID	Positive/ Indicative RSID
6ZGMFL	Not Tested Swab was not red stained	Negative Micro P-30	Not Tested
7M6VAR	Not Tested No visible stain observed	Negative AP	Positive/ Indicative Amylase
7MNLCK	Not Tested Swabs have no red stains	Negative Microscopic P30	Not Tested Per Agency's Policy
7XNELB	Negative KM	Negative AP P30 KPIC	Positive/ Indicative Phadebas

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
8C6JGB	Negative KM <u>Other Tests Performed</u> Urine - Creatinine Feces - Urobilinogen	Negative AP	Positive/ Indicative Phadebas
8NUMDA	Negative LMG HemaTrace	Negative AP Micro PSA ALS	Not Tested This laboratory does not perform testing for saliva.
8PP7HW	Negative LMG test	Not Tested	Positive/ Indicative Phadebas tube test
8XLW8D	Negative Visual	Negative AP	Positive/ Indicative RSID-Saliva
8YD3JW	Not Tested No visible blood-like staining observed, therefore no presumptive test for blood performed	Negative AP: x4 AP negative	Positive/ Indicative Phadebas: x4 Phadebas positive
9B677N	Negative KM	Negative AP	Positive/ Indicative RSID
9GQWTH	Negative KM	Negative AP Micro	Positive/ Indicative RSID
9HLHJN	Negative KM	Negative AP	Positive/ Indicative RSID
9JH363	Negative KM	Negative AP Micro	Positive/ Indicative Phadebas Press
9Q3XCX	Negative KM	Negative AP	Positive/ Indicative AMYLASE DIFFUSION
9UE4NY	Not Tested Visual exam - no staining observed	Negative AP Micro	Positive/ Indicative Amylase Diffusion

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
9W3CXZ	Negative KM	Negative AP Micro	Positive/ Indicative Phadebas press
AWHZRY	Negative KM	Negative AP PSA Micro	Positive/ Indicative Phadebas
B2DAF7	Negative Hemastix - neg 60 secs Hematrace - neg 10 mins	Negative AP - neg 60 secs Micro - neg p30 - neg 10 mins	Positive/ Indicative RSID - pos 1 min
B33FWW	Not Tested Not tested: No red-brown staining observed on the Items.	Negative AP	Positive/ Indicative RSID
BA8DNA	Negative Visual	Negative AP	Positive/ Indicative RSID
BHXXRM	Not Tested Not tested due to lack of reddish-brown color on swabs.	Negative Item 2 AnB: AP, Micro Item 2 CnD: AP, Micro	Positive/ Indicative Item 2 AnB: RSID Item 2 CnD: RSID
BMUKKB	Negative KM	Negative AP Seratec Micro	Positive/ Indicative RSID
BNKQXU	Negative Hematrace	Negative PSA Micro	Positive/ Indicative RSID
BU8H89	Negative LMG	Negative AP	Not Tested
BYM4XV	Not Tested no apparent bloodstain observed	Negative AP Micro	Positive/ Indicative Amylase diffusion

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
C4UYB4	Negative LMG HemaTrace	Negative ALS AP Micro P30	Not Tested Saliva testing is not offered as a service in this laboratory.
CGLJZ8	Negative KM <u>Other Tests Performed</u> Urine - Creatinine Urine - DMAC Fecal Material - Urobilinogen	Negative AP PSA	Positive/ Indicative Phadebas Amylase
DL73AZ	Negative KM	Negative AP Micro	Positive/ Indicative Phadebas Press
EB38QV	Negative KM	Negative PSA MICRO	Positive/ Indicative PHADEBAS
EB9GR9	Negative TMB	Negative AP	Positive/ Indicative RSID Saliva Field Kit
EFF3F6	Negative KM <u>Other Tests Performed</u> Urine - Creatinine, DMAC Fecal Material - Urobilinogen	Negative AP (Seri AP)	Positive/ Indicative Amylase (Phadebas)
ERDV3B	Negative KM	Negative AP PSA	Positive/ Indicative RSID
EVPZDB	Not Tested White scene or kit swabs are not tested for blood in this lab	Positive/ Indicative Micro (one (1) head observed) P30	Not Tested Saliva testing not being done in this lab at this time
EXGQ74	Negative KM <u>Other Tests Performed</u> Urine - Creatinine Fecal Material - Urobilinogen	Negative AP (SAP) PSA (ABA) Micro (Xmas)	Positive/ Indicative Phadebas Amylase

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
F42C9M	Negative Phenolphthalein	Negative RSID	Positive/ Indicative SaligaE
FVMNEL	Negative O-Tol	Negative Micro RSID	Positive/ Indicative RSID
FXB6NT	Negative KM	Negative AP PSA Micro	Positive/ Indicative Phadebas
G74RZB	Negative O-tol	Negative AP	Positive/ Indicative RSID
G8XFEZ	Inconclusive TMB HemaTrace Haemochromogen <u>Other Tests Performed</u> Ouchterlony (Human origin)	Negative AP PSA Microscopic	Not Tested
G9NLVQ	Negative LMG	Negative AP	Positive/ Indicative PHADEBAS PAPER TEST
GFC6FU	Negative KM	Negative AP Micro	Positive/ Indicative Phadebas
GPVJEP	Negative Visual exam KM	Negative AP H&E Slide	Positive/ Indicative Phadebas Paper H&E slide
GV43GN	Not Tested no staining on swab	Negative AP Micro	Positive/ Indicative Amylase diffusion
GZXZM8	Not Tested *swabs were not red stained	Negative Micro P30	Not Tested

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
H7J77Z	Negative KM <u>Other Tests Performed</u> Urine - Creatinine Fecal Material - Urobilinogen	Negative AP (Map)	Positive/ Indicative a-Amylase
H9ADNB	Negative O-tol	Negative ALS AP Micro PSA	Positive/ Indicative RSID
HDYUYF	Negative KM Abacard Hematrace	Negative AP Micro Abacard p30	Not Tested
HEDJ2T	Negative KM ABAcad HemaTrace	Negative AP ABAcad P30 Micro	Positive/ Indicative Phadebas test
HHBZEM	Not Tested A red stain was not observed.	Negative AP PSA Micro ALS	Positive/ Indicative Amylase diffusion
HQ63HJ	Negative O-tol	Negative Micro RSID	Positive/ Indicative RSID
HWRXDW	Negative KM	Negative AP Micro	Positive/ Indicative Radial Diffusion
J23ATA	Negative KM	Negative RSID-Semen PSA	Positive/ Indicative Starch-Iodine Assay
J94HFP	Not Tested	Negative AP	Positive/ Indicative Phadebas

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
JHVYAZ	Negative KM <u>Other Tests Performed</u> Urine - Creatinine Feces - Urobilinogen	Negative AP	Positive/ Indicative Phadebas
JJM2UU	Negative KM <u>Other Tests Performed</u> Urine - DMAC	Negative AP Micro	Positive/ Indicative Phadebas tube test Micro
JNFWZ9	Negative K-M <u>Other Tests Performed</u> Urine - Visual & DMAC Feces - Visual & Microscopy	Negative AP Microscopic sperm search	Positive/ Indicative RSID-Saliva
JPQ3UN	Negative LMG	Negative AP	Positive/ Indicative PHADEBAS PAPER TEST
JRWLWQ	Negative Visual Exam	Negative Fast Blue (Negative) (AP) Micro (Negative) P30 (Negative) (PSA)	Not Tested Current protocol does not allow for Saliva testing on swabs.
JTMRDG	Negative KM	Negative ALS AP	Positive/ Indicative ALS RSID-Saliva
JYWP6A	Negative KM	Negative RSID PSA	Positive/ Indicative Starch-iodine assay
K6XJEV	Negative KM	Negative AP Micro	Positive/ Indicative Phadebas
K7UYQR	Negative Visual KM	Negative AP Micro	Positive/ Indicative Phadebas
KQ3XGF	Negative KM	Negative AP Micro ABAcad p30	Positive/ Indicative radial diffusion Micro (for nucleated epithelial cells)

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
KUE3RF	Negative O-TOL	Negative AP MICRO RSID	Positive/ Indicative RSID
KXYWNG	Negative HemaTrace	Negative ABACard P30	Positive/ Indicative Saligae
KZ3R7Y	Negative LMG	Negative AP	Not Tested
LXURT6	Negative Visual exam	Negative AP	Positive/ Indicative RSID
M9CTN6	Negative Visual	Negative ALS AP Micro	Positive/ Indicative RSID
MED8A7	Negative O-Tol	Negative ALS AP Micro RSID	Positive/ Indicative ALS RSID
MEWRED	Not Tested Swabs were not tested for blood - does not display the characteristic reddish-brown stain associated with blood.	Negative AP, Micro AP, Micro AP, Micro AP, Micro	Positive/ Indicative RSID Saliva RSID Saliva RSID Saliva RSID Saliva
MQAX88	Negative KM	Negative RSID PSA	Positive/ Indicative Starch-iodine
MY39JM	Not Tested	Negative AP	Not Tested *saliva testing not performed in this lab
MZHQHM	Negative Visual Exam	Negative Fast Blue (AP) Micro p30 (PSA)	Not Tested Current protocol does not allow for Saliva testing on swabs.

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
N3ZCLU	Negative KM	Negative AP p30 KPIC	Positive/ Indicative Phadebas
N9PM3J	Not Tested No red-brown discoloration	Negative AP Micro	Positive/ Indicative Amylase Diffusion
N9PRB2	Negative KM	Negative AP PSA	Positive/ Indicative RSID
NQ6YNH	Not Tested	Negative AP RSID-Semen	Positive/ Indicative RSID-Saliva
P8P4F9	Not Tested No reddish/brown stain observed during visual inspection.	Negative AP Microscopy	Positive/ Indicative RSID Saliva
PCGF7Y	Negative KM	Negative AP Micro	Positive/ Indicative Amylase Testing by Radial Diffusion
PDCWKG	Negative KM O-tol ABA Card	Negative AP	Positive/ Indicative SALlgAE
PJ2EYR	Negative LMG Hematrace	Negative ALS AP micro P30	Not Tested Testing not offered.
PJ7RQU	Negative KM	Negative AP Seratec Micro	Positive/ Indicative RSID
PKWW3E	Negative KM Visual	Negative AP Slide made (micro)	Positive/ Indicative Phadebas

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
PMLJF3	Not Tested *NT because no red brown stain observed.	Negative AP	Positive/ Indicative RSID
PNF7UR	Negative TMB HemaTrace <u>Other Tests Performed</u> Human - Ouchterlony Faeces - Price (urobilin)	Negative AP PSA Micro	Positive/ Indicative Phadebas (amylase)
QADXPB	Negative Hemastix (TMB)	Negative AP negative @ 60 secs Microscopy P30	Positive/ Indicative RSID (+)ve @ 36 secs
QP6VRD	Negative KM	Negative BCIP	Positive/ Indicative Amylase Diffusion
R9VBPE	Negative LMG	Negative AP	Not Tested
RBKG8Q	Negative KM	Negative AP P-30 KPIC	Positive/ Indicative Phadebas
RDTPGB	Negative O-tol	Negative Micro RSID	Positive/ Indicative RSID
RNACNQ	Negative TMB HaemaTrace[sic] Haemachromogen[sic] <u>Other Tests Performed</u> Faeces - Prices	Negative AP PSA Micro	Positive/ Indicative Amylase Micro
RNQTHJ	Negative Visually <u>Other Tests Performed</u> DNA-analysis	Negative AP	Not Tested
RVQ877	Not Tested No visual apparent blood staining	Negative AP Micro	Positive/ Indicative Amylase diffusion

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
T2V4RR	Negative KM	Negative AP Seratec	Positive/ Indicative RSID
T7LPCY	Negative KM	Negative Microscopic P30	Not Tested No currently validated procedures for saliva test.
TJBFHW	Negative KM	Negative AP	Positive/ Indicative RSID
TKLJQE	Not Tested No stain consistent with blood	Negative AP Micro PSA	Positive/ Indicative Amylase Diffusion
TX4EYX	Negative KM	Negative Micro P30	Not Tested
U3J3HU	Negative KM	Negative AP	Positive/ Indicative RSID
UAYML8	Negative KM Test	Negative ALS AP Test	Positive/ Indicative ALS RSID-Saliva
UAZKRV	Not Tested Swabs were not stained	Negative Micro P30	Not Tested
UF9YH8	Negative o-tol	Negative ALS AP Micro RSID-semenogelin	Positive/ Indicative ALS RSID
UHA7HB	KM	Negative AP	Positive/ Indicative Phadebas
UPT6GW	Negative Visual examination	Negative AP	Positive/ Indicative Phadebas amylase

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
UPU6B6	Not Tested	Negative AP MICRO	Positive/ Indicative AMYLASE PLATE
V4BH9W	Negative KM	Negative ALS AP Micro	Positive/ Indicative ALS RSID
V7F8B6	Negative O-Tol	Negative AP RSID Micro	Positive/ Indicative RSID
V9Z287	Negative KM	Negative AP	Positive/ Indicative PHADEBAS
VKV7HR	Negative tetrabase	Negative Micro PSA RSID	Positive/ Indicative RSID
VM2TTB	Not Tested No apparent discoloration to swabs	Negative AP Micro	Positive/ Indicative Amylase Diffusion
W4LNZR	Negative KM	Negative AP Micro exam	Positive/ Indicative Amylase with radial diffusion
WBNQ8D	Negative KM	Negative AP Micro	Positive/ Indicative Phadebas
WF2Q7P	Negative KM HemaTrace <u>Other Tests Performed</u> Rsid Urine	Negative AP Micro Rsid semen	Positive/ Indicative Rsid saliva
WM2X2R	Not Tested no visible dried red stain present <u>Other Tests Performed</u> Urine - Creatinine	Negative Acid Phosphatase (AP) P30/PSA Micro	Positive/ Indicative Phadebas (amylase)

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
WXGZ8Q	Negative K-M	Negative AP	Positive/ Indicative RSID Saliva
XD83L9	Negative KM	Negative AP	Not Tested Not competent in saliva testing
XHZBMU	Not Tested no red stain observed	Negative Micro P30	Not Tested
XJC9UK	Positive/ Indicative KM	Negative AP	Positive/ Indicative Amylase
XXUUMQ	Negative KM	Negative AP	Positive/ Indicative RSID
Y28ZFH	Negative LMG HemaTrace	Negative AP Micro PSA ALS	Not Tested This laboratory does not currently offer saliva testing
Y8GAQK	Negative KM <u>Other Tests Performed</u> Urine - Creatinine Fecal Material - Urobilinogen	Negative AP	Positive/ Indicative Phadebas Amylase
Y8WHUU	Negative LMG	Negative AP	Not Tested
Y9VLM6	Negative Phenolphthalein	Negative RSID	Positive/ Indicative SaligaE
YBDZQ6	Negative Visual KM	Negative AP Microscopic exam - no spermatzoa[sic]	Not Tested Phadebas tube test not accredited
YLWGGQ	Not Tested Not brown stains	Negative Microscopic P30	Not Tested

Table 1 - Item 2

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
YMCGYU	Negative KM	Negative RSID-Semen PSA	Positive/ Indicative Starch-iodine test
YN3XWA	Negative Hemastix -ve @ 60 secs	Negative Micro AP -ve @ 60 secs	Positive/ Indicative RSID (+)ve @ 1 min
YQT4CZ	Not Tested No red brown staining	Negative AP Micro	Positive/ Indicative Radial Diffusion
YTUAPK	Negative LMG	Negative AP	Not Tested
YZDAA4	Negative Phenolphthalein	Negative RSID	Positive/ Indicative SaligaE
ZCTE3K	Negative Visual	Negative AP	Positive/ Indicative RSID
ZKJZ9J	Negative KM	Negative AP Seratec Micro	Positive/ Indicative RSID
ZRYNXE	Negative KM	Negative AP MICRO	Positive/ Indicative PHADEBAS (TUBE METHOD)
ZYFE47	Not Tested No red discoloration	Negative AP Micro	Positive/ Indicative Amylase Diffusion

Response Summary for Table 1 - Item 2

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	1	1	121
Negative	118	148	0
Inconclusive	1	0	0
Not Tested (NT)	29	1	29

Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
27CG2C	KM	AP MICRO	PHADEBAS
28BLGK	Positive/ Indicative LMG HemaTrace Ouchterlony	Negative AP	Not Tested
2EMTZU	Positive/ Indicative KM	Negative AP	Not Tested Not required per policy
2FHERZ	Positive/ Indicative O-tol ABACARD Hematrace	Negative micro RSID Semen	Negative RSID Saliva
2PQVGK	Positive/ Indicative KM <u>Other Tests Performed</u> Sperm Search	Negative PSA	Negative RSID
2VEVGY	Positive/ Indicative O-tol HemaTrace	Negative AP Micro RSID	Negative RSID
3FF4BN	Positive/ Indicative KM Hematrace	Negative AP Micro	Negative Phadebas
3MXG82	Positive/ Indicative KM	Negative AP Micro	Inconclusive Phadebas tube test
3P4YCN	Positive/ Indicative KM	Not Tested The swabs were heavily stained reddish-brown and there is no reason to believe semen would also be present on these swabs.	Not Tested The swabs were heavily stained reddish-brown and there is no reason to believe saliva would also be present on these swabs.

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
4APZYB	Positive/ Indicative KM HemaTrace	Negative AP Micro	Negative Amylase (radial diffusion)
4BJLQG	Positive/ Indicative LMG HemaTrace Ouchterlony	Negative AP	Not Tested
4HXUT4	Positive/ Indicative LMG	Negative AP	Not Tested
4JVXRK	Positive/ Indicative KM	Negative AP	Not Tested Not tested per policy *Not Required
68YCXV	Positive/ Indicative OT HemaTrace	Negative Micro RSID	Negative RSID
6HVCDK	Positive/ Indicative Visual KM HemaTrace	Negative AP	Negative RSID-Saliva
6PZ4XB	Positive/ Indicative KM	Negative AP	Not Tested Case scenario did not require testing for saliva
6UWLGW	Positive/ Indicative O-tol Hema Trace	Negative Micro RSID	Negative RSID
6ZGMFL	Positive/ Indicative Kastle-Meyer	Negative Micro P30	Not Tested
7M6VAR	Positive/ Indicative LMG	Negative AP, Micro, P30	Negative Amylase
7MNLCK	Positive/ Indicative Kastle-Meyer Test	Negative Microscopic P30	Not Tested Per Agency's Policy

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
7XNELB	Positive/ Indicative KM	Negative AP P30 KPIC	Negative Phadebas
8C6JGB	Positive/ Indicative KM Modified Takayama	Negative AP	Negative Phadebas
8NUMDA	Positive/ Indicative LMG HemaTrace	Not Tested Testing was not performed due to case scenario.	Not Tested This laboratory does not perform testing for saliva.
8PP7HW	Positive/ Indicative LMG test	Not Tested	Inconclusive Phadebas tube test
8XLW8D	Positive/ Indicative LMG HemaTrace Visual	Negative AP	Negative RSID-Saliva
8YD3JW	Positive/ Indicative LMG: Heavy blood-like staining observed, LMG presumptive test performed, x4 LMG positive	Not Tested As visible blood-like staining that gave an LMG positive result was observed, no presumptive test for semen performed	Not Tested As visible blood-like staining that gave an LMG positive result was observed, no presumptive test for saliva performed
9B677N	Positive/ Indicative KM	Negative AP	Not Tested Not tested per policy.
9GQWTH	Positive/ Indicative KM	Not Tested Not applicable to case scenario	Not Tested Not applicable to case scenario
9HLHJN	Positive/ Indicative KM	Negative AP	Not Tested Per laboratory policy, only intimate samples are saliva tested.
9JH363	Positive/ Indicative KM ABACard HemaTrace	Negative AP Micro	Negative Phadebas Press

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
9Q3XCX	Positive/ Indicative KM ABACARD HEMATRACE	Not Tested NOT TESTED DUE TO SCENARIO, LOCATION OF STAIN, AND POSITIVE KM TEST.	Not Tested NOT TESTED DUE TO SCENARIO, LOCATION OF STAIN, AND POSITIVE KM TEST.
9UE4NY	Positive/ Indicative LMG	Negative AP Micro P30	Negative Amylase Diffusion
9W3CXZ	Positive/ Indicative KM ABAcad Hematrace	Negative AP Micro	Negative Phadebas press
AWHZRY	Positive/ Indicative KM	Negative AP P30 Micro	Negative Phadebas
B2DAF7	Positive/ Indicative Hemastix - pos 5 secs Hematrace - pos 2 mins	AP - neg 60 secs Micro - neg	Not Tested No indication in scenario of possibility of saliva
B33FWW	Positive/ Indicative KM RSID	Negative AP	Not Tested
BA8DNA	Positive/ Indicative LMG HemaTrace Visual	Negative AP P30 Spermatozoa (Micro)	Negative RSID
BHXXRM	Positive/ Indicative Item 3 AnB: O-tol Item 3 CnD: O-tol	Negative Item 3 AnB: AP, Micro Item 3 CnD: AP, Micro	Not Tested Not tested due to swab type.
BMUKKB	Positive/ Indicative KM Hematrace	Negative AP Seratec Micro	Negative RSID
BNKQXU	Positive/ Indicative Hematrace	Negative PSA MICRO	Negative RSID

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
BU8H89	Positive/ Indicative LMG HemaTrace Ouchterlony Species Determination	Negative AP	Not Tested
BYM4XV	Positive/ Indicative LMG HemaTrace	Negative AP Micro	Negative Amylase diffusion
C4UYB4	Positive/ Indicative LMG HemaTrace	Negative ALS Micro P30	Not Tested Saliva testing is not offered as a service in this laboratory.
CGLJZ8	Positive/ Indicative KM Modified Takayama	Not Tested Not tested based on case scenario	Not Tested Not tested based on case scenario
DL73AZ	Positive/ Indicative KM ABACard HemaTrace	Negative AP Micro	Negative Phadebas Press
EB38QV	Positive/ Indicative KM	Negative PSA MICRO	Negative PHADEBAS
EB9GR9	Positive/ Indicative TMB	Negative AP	Negative RSID Saliva Field Kit
EFF3F6	Positive/ Indicative KM Takayama	Negative AP (Seri AP)	Negative Amylase (Phadebas)
ERDV3B	Positive/ Indicative KM	Negative AP	Not Tested Not routinely performed on KM positive stains
EVPZDB	Positive/ Indicative KM	Not Tested Brown or red stained scene swabs are not tested for semen in this lab	Not Tested Saliva testing not being done in this lab at this time

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
EXGQ74	Positive/ Indicative KM Mod Taka	Negative AP (SAP) PSA (ABA) Micro (Xmas)	Negative Phadebas Amylase
F42C9M	Negative Phenolphthalein	Negative RSID	Negative SaligaE
FVMNEL	Positive/ Indicative O-Tol Hematrace	Negative Micro RSID	Negative RSID
FXB6NT	Positive/ Indicative KM	Negative AP PSA Micro	Negative Phadebas
G74RZB	Positive/ Indicative O-tol HemaTrace	Not Tested No information to indicate that the human bloodstain on the wall would contain semen and no specific request for this testing	Not Tested No information to indicate oral contact and no specific request for this testing
G8XFEZ	Positive/ Indicative TMB HemaTrace Haemochromogen	Negative AP	Not Tested
G9NLVQ	Positive/ Indicative LMG	Negative AP	Not Tested
GFC6FU	Positive/ Indicative KM ABA	Negative AP Micro	Negative Phadebas
GPVJEP	Positive/ Indicative Visual exam KM	Negative AP H&E Slide	Negative Phadebas Paper
GV43GN	Positive/ Indicative LMG	Negative AP Micro	Negative Amylase diffusion

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
GZXZM8	Positive/ Indicative KM	Negative Micro P30	Not Tested
H7J77Z	Positive/ Indicative KM Takayama	Not Tested Analyst's discretion (based on scenario provided).	Not Tested Analyst's discretion (based on scenario provided).
H9ADNB	Positive/ Indicative O-tol HemaTrace	Negative ALS AP	Not Tested The swabs were very dark and there was no evidence of dilution
HDYUYF	Positive/ Indicative KM Abacard Hematrace	Negative AP Micro Abacard p30	Not Tested
HEDJ2T	Positive/ Indicative KM ABAcard HemaTrace	Negative AP Micro	Negative Phadebas test
HHBZEM	Positive/ Indicative LMG Hematrace	Negative AP PSA Micro ALS	Negative Amylase diffusion
HQ63HJ	Positive/ Indicative O-tol Hematrace	Negative Micro RSID	Negative RSID
HWRXDW	Positive/ Indicative KM HemaTrace	Negative AP Micro	Negative Radial Diffusion
J23ATA	Positive/ Indicative KM HO	Negative RSID-Semen PSA	Not Tested
J94HFP	Positive/ Indicative TMB	Negative AP	Negative Phadebas

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
JHVYAZ	Positive/ Indicative KM Modified Takayama	Negative AP	Negative Phadebas
JJM2UU	Positive/ Indicative KM <u>Other Tests Performed</u> Urine - DMAC	Negative AP Micro	Negative Phadebas tube test Micro
JNFWZ9	Positive/ Indicative K-M	Not Tested No scene information was provided that made it seem likely this stain would contain more than one body fluid	Not Tested No scene information was provided that made it seem likely this stain would contain more than one body fluid
JPQ3UN	Positive/ Indicative LMG	Negative AP	Not Tested
JRWLWQ	Positive/ Indicative Hemochromogen (Positive)	Negative Fast Blue (-ve) (AP)	Not Tested Current protocols do not allow for Saliva testing on swabs.
JTMRDG	Positive/ Indicative KM RSID-Blood	Negative ALS AP	Negative ALS RSID-Saliva
JYWP6A	Positive/ Indicative KM Hexagon OBTI	Negative RSID PSA	Not Tested
K6XJEV	Positive/ Indicative KM	Negative AP Micro	Inconclusive Phadebas
K7UYQR	Positive/ Indicative Visual KM	Negative AP Micro	Inconclusive Phadebas
KQ3XGF	Positive/ Indicative KM Takayama	Negative AP Micro ABAcad p30	radial diffusion

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
KUE3RF	Positive/ Indicative O-TOL HEMATRACE	Negative AP MICRO RSID	Negative RSID
KXYWNG	Positive/ Indicative HemaTrace	Negative ABAcad P30	Negative Saligae
KZ3R7Y	Positive/ Indicative LMG HemaTrace Ouchterlony	Negative AP PSA Microscopic	Not Tested
LXURT6	Positive/ Indicative Visual exam KM HemaTrace	Negative AP	Negative RSID
M9CTN6	Positive/ Indicative Visual KM HemaTrace	Negative AP Micro	Negative RSID
MED8A7	Positive/ Indicative O-Tol Hematrace	Negative ALS Micro RSID	Negative ALS RSID
MEWRED	Positive/ Indicative O-tol O-tol O-tol O-tol	Negative AP, Micro AP, Micro AP, Micro AP, Micro	Not Tested Did not test for saliva based on case scenario.
MQAX88	Positive/ Indicative KM Hexagon OBTI	Negative RSID PSA	Not Tested
MY39JM	Positive/ Indicative KM	Negative AP	Not Tested *saliva testing not performed in this lab
MZHQHM	Positive/ Indicative Hemochromogen	Negative Fast Blue (AP)	Not Tested Current protocols do not allow for Saliva testing on swabs.

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
N3ZCLU	Positive/ Indicative KM	Negative AP p30 KPIC	Negative Phadebas
N9PM3J	Positive/ Indicative LMG HemaTrace	Negative AP Micro	Negative Amylase Diffusion
N9PRB2	Positive/ Indicative KM	Negative AP	Not Tested KM positive stains are not routinely tested for saliva
NQ6YNH	Positive/ Indicative KM RSID-Blood	Negative AP RSID-Semen	Not Tested
P8P4F9	Positive/ Indicative O-tol	Negative AP Microscopy	Negative RSID Saliva
PCGF7Y	Positive/ Indicative KM HemaTrace	Not Tested Scenario didn't indicate the need for this testing.	Not Tested Scenario didn't indicate the need for this testing.
PDCWKG	Positive/ Indicative KM O-tol ABA Card	Negative AP	Negative SALlgAE
PJ2EYR	Positive/ Indicative LMG HemaTrace	Negative AP ALS	Not Tested Testing not offered
PJ7RQU	Positive/ Indicative KM Hematrace	Negative AP Seratec Micro	Negative RSID
PKWW3E	Positive/ Indicative KM Visual	Not Tested Area recovered - wasn't described as sex swab	Not Tested Area recovered from - was deemed unnecessary (& not dilute looking)

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
PMLJF3	Positive/ Indicative O-tol HemaTrace	Not Tested Item 3 was not tested for semen because of the location of the stain.	Not Tested Item 3 was not tested for saliva because of the location of the stain.
PNF7UR	Positive/ Indicative TMB HemaTrace Haemochromogen	Negative AP PSA Micro	Negative Phadebas (amylase)
QADXPB	Positive/ Indicative Hemastix (TMB) HemaTrace	Negative AP Test Microscopy P30	Negative RSID a-amylase
QP6VRD	Positive/ Indicative KM Hematrace	Negative BCIP	Positive/ Indicative Amylase Diffusion
R9VBPE	Positive/ Indicative HemaTrace LMG	Negative AP	Not Tested Scenario specific - case details does not declare biting nor oral contact
RBKG8Q	Positive/ Indicative KM	Negative AP P-30 KPIC	Negative Phadebas
RDTPGB	Positive/ Indicative O-tol HemaTrace	Negative Micro RSID	Negative RSID
RNACNQ	Positive/ Indicative TMB HemaTrace	Negative AP	Negative Amylase
RNQTHJ	Positive/ Indicative LMG OBTI	Negative AP Micro	Not Tested
RVQ877	Positive/ Indicative LMG HemaTrace	Negative AP Micro	Negative Amylase diffusion

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
T2V4RR	Positive/ Indicative KM Hematrace	Negative AP Seratec	Negative RSID
T7LPCY	Positive/ Indicative KM	Negative Microscopic P30	Not Tested No currently validated procedures for saliva test.
TJBFHW	Positive/ Indicative KM	Negative AP	Not Tested NT (Per Policy) - Item not eligible for saliva testing
TKLJQE	Positive/ Indicative LMG HemaTrace	Negative AP Micro	Negative Amylase Diffusion
TX4EYX	Positive/ Indicative KM	Negative Micro P30	Not Tested
U3J3HU	Positive/ Indicative KM	Negative AP	Not Tested Per Policy; not required to test
UAYML8	Positive/ Indicative KM Test RSID-Blood	Negative ALS AP Test	Negative ALS RSID-Saliva
UAZKRV	Positive/ Indicative KM	Negative Micro P-30	Not Tested
UF9YH8	Positive/ Indicative ALS O-tol Hematrace	Negative ALS Micro RSID-semenogelin	Negative ALS RSID
UHA7HB	Positive/ Indicative KM	Negative AP	
UPT6GW	Positive/ Indicative KM	Negative AP	Negative Phadebas amylase

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
UPU6B6	Positive/ Indicative KM	Not Tested	Not Tested
V4BH9W	Positive/ Indicative KM	Negative ALS AP	Not Tested Case scenario & lab policy/procedures
V7F8B6	Positive/ Indicative O-Tol RSID	Negative AP RSID Micro	Negative RSID
V9Z287	Positive/ Indicative KM	Negative AP	Negative PHADEBAS
VKV7HR	Positive/ Indicative Tetrabase RSID	Negative Micro PSA RSID	Negative RSID
VM2TTB	Positive/ Indicative LMG Hematrace	Negative AP Micro	Negative Amylase Diffusion
W4LNZR	Positive/ Indicative KM HemaTrace	Negative AP Micro exam	Negative Amylase with radial diffusion
WBNQ8D	Positive/ Indicative KM HemaTrace	Negative AP Micro	Negative Phadebas
WF2Q7P	Positive/ Indicative KM Hematrace <u>Other Tests Performed</u> Rsid Urine	Negative AP Micro Rsid Semen	Negative Rsid Saliva
WM2X2R	Positive/ Indicative PTMB HemaTrace <u>Other Tests Performed</u> Urine - Creatinine Fecal - Micro Fecal - Urobilinogen	Negative Acid Phosphatase (AP)	Negative Phadebas (Amylase)

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
WXGZ8Q	Positive/ Indicative K-M	Negative AP	Not Tested Per policy, item does not meet criteria to be tested for this body fluid - therefore not required to test.
XD83L9	Positive/ Indicative KM	Negative AP	Not Tested Not competent in saliva testing
XHZBMU	Positive/ Indicative KM	Negative Micro P30	Not Tested
XJC9UK	Positive/ Indicative KM	Negative AP	
XXUUMQ	Positive/ Indicative KM	Negative AP	Not Tested Saliva test not performed on this item per laboratory policy - not required to test
Y28ZFH	Positive/ Indicative LMG HemaTrace	Not Tested Testing not performed due to the probative value of the sample based on the scenario given.	Not Tested This laboratory does not offer saliva testing.
Y8GAQK	Positive/ Indicative KM Modified Takayama	Negative AP	Negative Phadebas Amylase
Y8WHUU	Positive/ Indicative LMG HemaTrace Ouchterlony - Human	Negative AP	Not Tested
Y9VLM6	Negative Phenolphthalein	Negative RSID	Negative SaligaE
YBDZQ6	Positive/ Indicative Visual KM	Negative AP Microscopic exam - no spermatazoa[sic]	Not Tested Phadebas tube test not accredited

Table 1 - Item 3

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
YLWGGQ	Positive/ Indicative Kastle-Meyer Test	Negative Microscopic P30	Not Tested
YMCGYU	Positive/ Indicative KM Hexagon OBTI	Negative RSID-Semen PSA	Not Tested
YN3XWA	Positive/ Indicative Hemastix (+)ve @ 5 secs Hematrace (+)ve @ 20 secs	Negative Micro AP -ve @ 60 secs	Negative RSID -ve @ 8½ mins
YQT4CZ	Positive/ Indicative KM	Not Tested No senario[sic] to indicate other bodily fluids present	Not Tested No scenario to indicate other bodily fluids present
YTUAPK	Positive/ Indicative LMG HemaTrace Ouchterlony (human)	Negative AP	Not Tested
YZDAA4	Negative Phenolphthalein	Negative RSID	Negative SaligaE
ZCTE3K	Positive/ Indicative LMG HemaTrace Visual	Negative AP P30 Microscopic	Negative RSID
ZKJZ9J	Positive/ Indicative KM Hematrace	Negative AP Seratec Micro	Negative RSID
ZRYNXE	Positive/ Indicative KM	Negative AP MICRO CHOLINE	Inconclusive PHADEBAS (TUBE METHOD)
ZYFE47	Positive/ Indicative LMG Hematrace	Negative AP Micro	Negative Amylase Diffusion

Response Summary for Table 1 - Item 3

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	147	0	1
Negative	3	132	75
Inconclusive	0	0	5
Not Tested (NT)	0	17	66

Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
27CG2C	KM	AP MICRO	PHADEBAS
28BLGK	Negative, Not Tested LMG	Positive/ Indicative AP P30 Micro	Not Tested
2EMTZU	Not Tested Not required per policy	Positive/ Indicative AP Micro	Not Tested Not required per policy
2FHERZ	Not Tested Swabs were not stained red-brown	Positive/ Indicative micro	Negative RSID Saliva
2PQVGK	Negative KM <u>Other Tests Performed</u> Sperm Search	Positive/ Indicative PSA	Negative RSID
2VEVGY	Negative O-tol	Positive/ Indicative AP Micro	Positive/ Indicative RSID
3FF4BN	Negative KM	Positive/ Indicative AP Micro	Negative Phadebas
3MXG82	Negative KM	Positive/ Indicative AP Micro	Inconclusive Phadebas
3P4YCN	Not Tested Light yellow staining was noted on each swab; no reason to test for blood.	Positive/ Indicative AP Micro	Not Tested Semen was detected and there is no reason to believe saliva would also be present on these swabs.

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
4APZYB	Not Tested No characteristic staining observed, would not test based on staining & scenario	Positive/ Indicative AP Micro	Negative Amylase (radial diffusion)
4BJLQG	Negative LMG	Positive/ Indicative AP PSA Micro	Not Tested
4HXUT4	Negative LMG	Positive/ Indicative AP	Not Tested
4JVXRK	Not Tested No reddish-brown stain. Not tested per policy *Not Required	Positive/ Indicative AP Micro	Not Tested Not tested per policy *Not Required
68YCXV	Negative OT	Positive/ Indicative Micro	Negative RSID
6HVCDK	Negative Visual	Positive/ Indicative AP micro	Negative RSID-Saliva
6PZ4XB	Not Tested Case scenario did not require testing for blood, absence of red-brown staining	Positive/ Indicative AP Micro	Not Tested Case scenario did not require testing for saliva
6UWLGW	Negative O-tol	Positive/ Indicative Micro	Positive/ Indicative RSID
6ZGMFL	Not Tested Swab not red stained	Positive/ Indicative Micro	Not Tested
7M6VAR	Not Tested no visible stain observed	Positive/ Indicative ACP, Micro	Negative Amylase
7MNLCK	Swabs have no red stains	Positive/ Indicative Microscopic	Not Tested Per Agency's Policy

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
7XNELB	Negative KM	Positive/ Indicative AP P30 KPIC	Negative Phadebas
8C6JGB	Negative KM	Positive/ Indicative AP Micro	Negative Phadebas
8NUMDA	Not Tested Testing was not performed due to case scenario.	Positive/ Indicative AP Micro ALS	Not Tested This laboratory does not perform testing for saliva.
8PP7HW	Negative LMG test	Not Tested	Inconclusive Phadebas tube test
8XLW8D	Negative LMG Visual	Positive/ Indicative AP P30 Micro	Negative RSID-Saliva
8YD3JW	Not Tested No visible blood-like staining observed, therefore no presumptive test for blood performed	Positive/ Indicative AP: x4 AP positive Micro: 2+ heads & tails, negative for NECs	Negative Phadebas: x4 Phadebas negative
9B677N	Negative KM	Positive/ Indicative AP Micro	Not Tested Not tested per policy.
9GQWTH	Not Tested No staining consistent with blood observed	Positive/ Indicative AP Micro	Not Tested Not applicable to case scenario
9HLHJN	Negative KM	Positive/ Indicative AP Micro	Not Tested Per laboratory policy, only intimate swabs from the vaginal area are saliva tested.
9JH363	Negative KM	Positive/ Indicative AP PSA	Negative Phadebas Press

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
9Q3XCX	Not Tested NOT TESTED DUE TO LOCATION OF SWABBING, NO VISIBLE STAINING, AND POSITIVE SEMEN TEST.	Positive/ Indicative AP MICRO	Not Tested NO TESTING DUE TO POS INDICATION OF SEMEN.
9UE4NY	Not Tested Visual exam - no staining observed	Positive/ Indicative AP Micro	Negative Amylase Diffusion
9W3CXZ	Negative KM	Positive/ Indicative AP PSA Semiquant	Negative Phadebas press
AWHZRY	Negative KM	Positive/ Indicative AP PSA Micro	Negative Phadebas
B2DAF7	Inconclusive Hemastix - wk pos 25 secs. Did not confirm presence of blood as semen positive on victim.	Positive/ Indicative AP - pos 8 secs Micro - pos	Negative RSID - neg 8.5 mins
B33FWW	Not Tested Not tested: No red-brown staining observed on the items.	Positive/ Indicative AP RSID Micro	Negative RSID
BA8DNA	Negative Visual LMG	Positive/ Indicative AP P30 Spermatozoa (Micro)	Negative RSID
BHXXRM	Not Tested Not tested due to lack of reddish-brown color on swabs.	Positive/ Indicative Item 4 AnB: AP, Micro Item 4 CnD: AP, Micro	Negative Item 4 AnB: RSID Item 4 CnD: RSID
BMUKKB	Negative KM	Positive/ Indicative AP Seratec Micro	Negative RSID

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
BNKQXU	Positive/ Indicative HEMATRACE	Positive/ Indicative PSA MICRO	Negative RSID
BU8H89	Negative LMG	Positive/ Indicative AP PSA Micro	Not Tested
BYM4XV	Not Tested no apparent bloodstain observed	Positive/ Indicative AP Micro PSA	Negative Amylase diffusion
C4UYB4	Not Tested No red/brown staining was visualized on the samples; therefore, blood testing was not performed.	Positive/ Indicative ALS AP Micro	Not Tested Saliva testing is not offered as a service in this laboratory.
CGLJZ8	Negative KM	Positive/ Indicative AP Micro	Negative Phadebas Amylase
DL73AZ	Negative KM	Positive/ Indicative AP PSA Micro	Negative Phadebas Press
EB38QV	Negative KM	Positive/ Indicative PSA MICRO	Negative PHADEBAS
EB9GR9	Negative TMB	Positive/ Indicative AP	Negative RSID Saliva Field Kit
EFF3F6	Negative KM	Positive/ Indicative AP (Seri AP) Micro (Xmas Stain)	Negative Amylase (Phadebas)
ERDV3B	Not Tested Not routinely performed on DRS-negative stains	Positive/ Indicative PSA AP Micro	Negative RSID

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
EVPZDB	Not Tested White scene or kit swabs are not tested for blood in this lab	Positive/ Indicative Micro	Not Tested Saliva testing not being done in this lab at this time
EXGQ74	Negative KM	Positive/ Indicative AP (SAP) PSA (ABA) Micro (Xmas)	Inconclusive Phadebas Amylase
F42C9M	Negative Phenolphthalein	Negative RSID	Negative SaligaE
FVMNEL	Not Tested Not applicable to case history	Positive/ Indicative Micro	Negative RSID
FXB6NT	Negative KM	Positive/ Indicative AP PSA Micro	Negative Phadebas
G74RZB	Negative O-tol	Positive/ Indicative AP Micro	Not Tested No information to indicate oral contact by suspect and semen is present in this stain.
G8XFEZ	Positive/ Indicative TMB HemaTrace	Positive/ Indicative AP PSA Microscopic	Not Tested
G9NLVQ	Negative LMG	Positive/ Indicative AP	Negative PHADEBAS PAPER TEST
GFC6FU	Negative KM	Positive/ Indicative AP PSA	Negative Phadebas
GPVJEP	Negative Visual Exam KM	Positive/ Indicative AP H&E Slide	Negative Phadebas Paper
GV43GN	Not Tested no staining on swab	Positive/ Indicative AP Micro	Negative Amylase diffusion

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
GZXZM8	Not Tested *swabs were not red stained	Positive/ Indicative Micro	Not Tested
H7J77Z	Not Tested Analyst's discretion (based on scenario provided).	Positive/ Indicative AP (Map) Micro (Xmas)	Negative a-Amylase
H9ADNB	Negative O-tol	Positive/ Indicative ALS AP Micro	Positive/ Indicative RSID
HDYUYF	Positive/ Indicative KM Abacard Hematrace	Positive/ Indicative AP Micro Abacard p30	Not Tested
HEDJ2T	Negative KM	Positive/ Indicative AP Micro	Negative Phadebas test
HHBZEM	Not Tested A red stain was not observed.	Positive/ Indicative AP Micro ALS	Negative Amylase diffusion
HQ63HJ	Negative O-tol	Positive/ Indicative Micro	Positive/ Indicative RSID
HWRXDW	Negative KM	Positive/ Indicative AP Micro	Negative Radial Diffusion
J23ATA	Negative KM	Positive/ Indicative RSID-Semen PSA	Not Tested
J94HFP	Not Tested	Positive/ Indicative AP Micro	Negative Phadebas
JHVYAZ	Negative KM	Positive/ Indicative AP Micro	Negative Phadebas

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
JJM2UU	Negative KM <u>Other Tests Performed</u> Urine - DMAC	Positive/ Indicative AP Micro	Inconclusive Phadebas tube test Micro
JNFWZ9	Not Tested No visual indication that blood was present	Positive/ Indicative AP Microscopic sperm search	Not Tested No senario[sic] information provided that supports likelihood of a mixed semen/saliva stain
JPQ3UN	Negative LMG	Positive/ Indicative AP	Negative PHADEBAS PAPER TEST
JRWLWQ	Negative Visual Exam	Positive/ Indicative Fast Blue (+ve) (AP) Micro (+ve)	Not Tested Current protocols do not allow for Saliva testing on swabs.
JTMRDG	Negative KM	Positive/ Indicative ALS AP RSID-Semen Microscopy	Negative ALS RSID-Saliva
JYWP6A	Negative KM	Positive/ Indicative RSID PSA	Not Tested
K6XJEV	Negative KM	Positive/ Indicative AP Micro	Positive/ Indicative Phadebas
K7UYQR	Negative Visual KM	Positive/ Indicative AP Micro	Positive/ Indicative Phadebas
KQ3XGF	Negative KM	Positive/ Indicative AP Micro ABAcad p30	Inconclusive radial diffusion
KUE3RF	Negative O-TOL	Positive/ Indicative AP MICRO	Negative RSID

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
KXYWNG	Negative HemaTrace	Positive/ Indicative ABACard P30	Negative Saligae
KZ3R7Y	Negative LMG	Positive/ Indicative AP PSA Microscopic	Not Tested
LXURT6	Negative Visual exam	Positive/ Indicative AP Micro	Negative RSID
M9CTN6	Negative Visual	Positive/ Indicative ALS AP Micro	Negative RSID
MED8A7	Negative O-Tol	Positive/ Indicative ALS Micro	Negative ALS RSID
MEWRED	Not Tested Swabs were not tested for blood - does not display the characteristic reddish-brown stain associated with blood.	Positive/ Indicative AP, Micro AP, Micro AP, Micro AP, Micro	Negative RSID Saliva RSID Saliva RSID Saliva RSID Saliva
MQAX88	Negative KM	Positive/ Indicative RSID PSA	Not Tested
MY39JM	Negative KM	Positive/ Indicative AP	Not Tested *saliva testing not performed in this lab
MZHQHM	Negative Hemastix <u>Other Tests Performed</u> Contact DNA - Sent for DNA analysis	Positive/ Indicative Fast Blue (AP) Micro	Not Tested Current protocols do not allow for Saliva testing on swabs.
N3ZCLU	Negative KM	Positive/ Indicative AP p30 KPIC	Inconclusive Phadebas

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
N9PM3J	Not Tested No red-brown discoloration	Positive/ Indicative AP Micro	Negative Amylase Diffusion
N9PRB2	Not Tested stains without any visible red/brown color are not routinely tested for blood	Positive/ Indicative AP Micro (sperm search)	Not Tested semen positive stains are not routinely tested for saliva
NQ6YNH	Not Tested	Positive/ Indicative AP Micro	Negative RSID-Saliva
P8P4F9	Not Tested No reddish/brown stain observed during visual inspection.	Positive/ Indicative AP Microscopy	Negative RSID Saliva
PCGF7Y	Not Tested no red-brown staining	Positive/ Indicative AP Micro	Negative Amylase with Radial Diffusion
PDCWKG	Negative KM O-tol ABA Card	Positive/ Indicative AP	Negative SALlgAE
PJ2EYR	Not Tested Not tested based on staining.	Positive/ Indicative ALS AP Micro	Not Tested Testing not offered.
PJ7RQU	Negative KM	Positive/ Indicative AP Seratec Micro	Negative RSID
PKWW3E	Negative KM Visual	Positive/ Indicative AP Slide made (Micro)	Inconclusive Phadebas
PMLJF3	Not Tested Item 4 was not tested for blood because of the location and color of the stain.	Positive/ Indicative AP Micro ALS	Positive/ Indicative RSID

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
PNF7UR	Negative TMB HemaTrace	Positive/ Indicative AP PSA Micro	Negative Phadebas (amylase)
QADXPB	Positive/ Indicative Hemastix (TMB)	Positive/ Indicative Microscopy	Negative RSID a-amylase
QP6VRD	Not Tested	Positive/ Indicative BCIP Microscopic Sperm Search	Positive/ Indicative Amylase Diffusion
R9VBPE	Negative LMG	Positive/ Indicative Microscopy AP	Inconclusive Starch Iodine test
RBKG8Q	Negative KM	Positive/ Indicative AP P-30 KPIC	Negative Phadebas
RDTPGB	Negative O-tol	Positive/ Indicative Micro	Negative RSID
RNACNQ	Positive/ Indicative TMB HemaTrace	Positive/ Indicative AP Micro	Negative Amylase
RNQTHJ	Negative Visually	Positive/ Indicative AP Micro	Not Tested
RVQ877	Not Tested No apparent blood staining	Positive/ Indicative AP Micro	Negative Amylase diffusion
T2V4RR	Negative KM	Positive/ Indicative AP Seratec Micro	Negative RSID

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
T7LPCY	Not Tested white swab, test for blood not applicable to the item in scenario	Positive/ Indicative Microscopic	Not Tested No currently validated procedures for saliva test.
TJBFHW	Not Tested NT (Per Policy) - no reddish-brown color - not required	Positive/ Indicative AP Micro	Not Tested NT (Per Policy) - Item not eligible for saliva testing
TKLJQE	Negative LMG	Positive/ Indicative AP Micro PSA	Negative Amylase Diffusion
TX4EYX	Negative KM	Positive/ Indicative Micro	Not Tested
U3J3HU	Not Tested No reddish brown stains seen on swabs; not required to test	Positive/ Indicative AP Micro	Not Tested Per Policy, not required to test
UAYML8	Negative KM Test	Positive/ Indicative ALS AP Test RSID-Semen Microscopy	Negative ALS RSID-Saliva
UAZKRV	Not Tested Swabs were not stained	Positive/ Indicative Micro	Not Tested
UF9YH8	Not Tested No red-brown staining	Positive/ Indicative ALS AP Micro	Positive/ Indicative ALS RSID
UHA7HB	Negative KM	Positive/ Indicative AP	Negative Phadebas
UPT6GW	Negative Visual examination	Positive/ Indicative AP PSA	Positive/ Indicative Phadebas amylase

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
UPU6B6	Not Tested	Positive/ Indicative AP MICRO	Not Tested
V4BH9W	Not Tested No visual indication of blood being present	Positive/ Indicative ALS AP Micro	Not Tested Case scenario
V7F8B6	Negative O-Tol	Positive/ Indicative AP Micro	Negative RSID
V9Z287	Negative KM	Positive/ Indicative AP	Negative PHADEBAS
VKV7HR	Negative Tetrabase	Positive/ Indicative Micro	Negative RSID
VM2TTB	Negative LMG	Positive/ Indicative AP Micro	Negative Amylase Diffusion
W4LNZR	Negative KM	Positive/ Indicative AP Micro exam	Negative Amylase with radial diffusion
WBNQ8D	Negative KM	Positive/ Indicative AP Micro P30	Negative Phadebas
WF2Q7P	Negative KM Hematrace <u>Other Tests Performed</u> Rsid urine	Positive/ Indicative AP Micro	Negative RSID Saliva
WM2X2R	Negative PTMB Slight yellowish color on swabs. <u>Other Tests Performed</u> Urine - Creatinine	Positive/ Indicative Acid Phosphatase (AP) Micro	Negative Phadebas (amylase)

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
WXGZ8Q	Not Tested No visibled[sic] reddish-brown staining - not tested as per policy - therefore not required to test.	Positive/ Indicative AP Micro	Not Tested Per policy, item does not meet criteria to be tested for this body fluid - therefore not required to test.
XD83L9	Negative KM	Positive/ Indicative AP	Not Tested Not competent in saliva testing
XHZBMU	Not Tested no red stain observed	Positive/ Indicative Micro	Not Tested
XJC9UK	Negative KM	Positive/ Indicative AP Micro	
XXUUMQ	Not Tested No visible reddish-brown staining - not required to test	Positive/ Indicative AP Micro	Not Tested Saliva test not performed on this item per laboratory policy - not required to test
Y28ZFH	Not Tested Testing not performed due to the probative value of the sample based on the scenario given.	Positive/ Indicative AP Micro PSA ALS	Not Tested This laboratory does not offer saliva testing.
Y8GAQK	Negative KM	Positive/ Indicative AP Micro - Xmas tree stain	Inconclusive Phadebas Amylase
Y8WHUU	Negative LMG	Positive/ Indicative AP PSA Micro	Not Tested
Y9VLM6	Negative Phenolphthalein	Negative RSID	Negative Saligae
YBDZQ6	Negative Visual KM	Positive/ Indicative AP Microscopic exam - spermatazoa[sic] identified	Not Tested Phadebas tube test not accredited

Table 1 - Item 4

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
YLWGGQ	Not Tested Not brown stains	Positive/ Indicative Microscopic	Not Tested
YMCGYU	Negative KM	Positive/ Indicative RSID-Semen PSA	Not Tested
YN3XWA	Negative Hemastix -ve @ 60 secs	Positive/ Indicative Micro	Negative RSID -ve @ 8½ mins
YQT4CZ	Not Tested No red brown staining	Positive/ Indicative AP Micro	Not Tested No indication in Scenario of saliva on victims stomach just sexual assault
YTUAPK	Negative LMG	Positive/ Indicative AP PSA Micro	Not Tested
YZDAA4	Negative Phenolphthalein	Negative RSID	Negative SaligaE
ZCTE3K	Negative Visual	Positive/ Indicative AP P30 Microscopic	Negative RSID
ZKJZ9J	Negative KM	Positive/ Indicative AP Seratec Micro	Negative RSID
ZRYNXE	Negative KM <u>Other Tests Performed</u> URINE - CREATININE - FAINT	Positive/ Indicative AP MICRO	Positive/ Indicative PHADEBAS (TUBE METHOD)
ZYFE47	Not Tested No red discoloration	Positive/ Indicative AP Micro	Negative Amylase Diffusion

Response Summary for Table 1 - Item 4

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	5	146	11
Negative	93	3	79
Inconclusive	1	0	9
Not Tested (NT)	50	1	50

Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
27CG2C	KM	AP	PHADEBAS
28BLGK	Positive/ Indicative LMG HemaTrace Ouchterlony <u>Other Tests Performed</u> Biological Fluids - ALS	Negative AP	Not Tested
2EMTZU	Positive/ Indicative KM	Negative AP	Not Tested Not required per policy
2FHERZ	Positive/ Indicative O-tol ABACARD Hematrace	Negative ALS	Negative ALS
2PQVGK	Positive/ Indicative KM	Negative PSA	Negative RSID
2VEVGY	Positive/ Indicative O-tol HemaTrace	Negative ALS AP Micro RSID	Negative RSID
3FF4BN	Positive/ Indicative KM Hematrace	Negative AP Micro	Negative Phadebas
3MXG82	Positive/ Indicative KM	Not Tested Would not expect semen to be on headwear	Negative Phadebas paper test
3P4YCN	Positive/ Indicative KM	Negative ALS AP Micro	Negative ALS RSID-Saliva
4APZYB	Positive/ Indicative KM HemaTrace	Negative AP Micro	Negative Amylase (radial diffusion)

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
4BJLQG	Positive/ Indicative LMG HemaTrace Ouchterlony <u>Other Tests Performed</u> Biological Fluids - ALS	Negative AP	Not Tested
4HXUT4	Positive/ Indicative LMG	Negative AP	Not Tested
4JVXRK	Positive/ Indicative KM	Negative AP Micro RSID	Not Tested Not tested per policy *Not Required
68YCXV	Positive/ Indicative OT HemaTrace	Not Tested ALS	Not Tested ALS
6HVCDK	Positive/ Indicative Visual KM HemaTrace	Negative ALS AP	Negative ALS RSID-Saliva
6PZ4XB	Positive/ Indicative KM	Not Tested Case scenario did not require testing for semen	Not Tested Case scenario did not require testing for saliva
6UWLGW	Positive/ Indicative O-tol Hema Trace	Negative ALS AP Micro RSID	Negative ALS RSID
6ZGMFL	Positive/ Indicative Kastle-Meyer	Negative Micro P30 ALS	Not Tested
7M6VAR	Positive/ Indicative LMG <u>Other Tests Performed</u> ALS	Negative ACP	Negative Amylase

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
7MNLCK	Positive/ Indicative Kastle-Meyer Test	Negative Microscopic P30 ALS	Not Tested Per Agency's Policy
7XNELB	Positive/ Indicative KM	Negative AP P30 KPIC	Negative Phadebas
8C6JGB	Positive/ Indicative KM Modified Takayama	Negative AP	Negative Phadebas
8NUMDA	Positive/ Indicative LMG HemaTrace <u>Other Tests Performed</u> ALS	Not Tested This testing was not performed due to case scenario.	Not Tested This laboratory does not perform testing for saliva.
8PP7HW	Positive/ Indicative LMG spot test	Not Tested	Negative Phadebas press test
8XLW8D	Positive/ Indicative LMG HemaTrace Visual	Negative ALS	Negative RSID-Saliva
8YD3JW	Positive/ Indicative LMG: Visible blood-like staining observed on inside front of item, LMG presumptive test performed, LMG positive	Not Tested No requirement to test item (headpiece) for semen	Not Tested Item did not appear to have been worn over the face, therefore item not presumptively tested for saliva
9B677N	Positive/ Indicative KM	Negative AP Micro	Not Tested Not tested per policy.
9GQWTH	Positive/ Indicative KM	Not Tested Not applicable to case scenario	Not Tested Not applicable to case scenario

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
9HLHJN	Positive/ Indicative KM	Negative AP	Not Tested Per laboratory policy, only intimate swabs from the vaginal area are saliva tested.
9JH363	Positive/ Indicative KM ABACard HemaTrace	Negative AP	Negative Phadebas Press
9Q3XCX	Positive/ Indicative KM	Negative ALS	Negative ALS
9UE4NY	Positive/ Indicative LMG	Negative AP	Negative Amylase Diffusion
9W3CXZ	Positive/ Indicative KM ABAcard Hematrace	Negative AP	Negative Phadebas press
AWHZRY	Positive/ Indicative KM <u>Other Tests Performed</u> Semen - ALS	Negative AP PSA Micro	Negative Phadebas
B2DAF7	Positive/ Indicative Hemastix - pos 5 secs	Not Tested Semen not tested as headpiece from suspect.	Negative Polylight[sic] - neg Phadebas Press Test - neg
B33FWW	Positive/ Indicative KM RSID	Negative AP	Not Tested
BA8DNA	Positive/ Indicative Visual LMG HemaTrace	Negative AP P30 Spermatozoa (Micro) ALS	Negative RSID
BHXXRM	Positive/ Indicative Item 5A: O-tol	Negative Item 5: ALS Item 5A: AP, Micro	Not Tested Not tested due to item type.

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
BMUKKB	Positive/ Indicative KM Hematrace	Not Tested	Not Tested
BNKQXU	Positive/ Indicative KM HEMATRACE	Negative PSA AP MICRO	Negative RSID
BU8H89	Positive/ Indicative LMG HemaTrace Ouchterlony Species Determination <u>Other Tests Performed</u> Bodily Fluid - ALS	Negative AP	Not Tested
BYM4XV	Positive/ Indicative LMG HemaTrace	Negative AP Micro ALS	Negative Amylase diffusion
C4UYB4	Positive/ Indicative LMG HemaTrace	Negative ALS Micro P30	Not Tested Saliva testing is not offered as a service in this laboratory.
CGLJZ8	Positive/ Indicative KM Modified Takayama	Not Tested Not tested based on case scenario	Not Tested Not tested based on case scenario
DL73AZ	Positive/ Indicative KM ABACard HemaTrace	Negative AP	Negative Phadebas Press
EB38QV	Positive/ Indicative KM	Negative PSA MICRO	Negative PHADEBAS
EB9GR9	Positive/ Indicative TMB	Inconclusive ALS	Inconclusive ALS
EFF3F6	Positive/ Indicative KM Takayama	Negative AP (Seri AP)	Negative Amylase (Phadebas)

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
ERDV3B	Positive/ Indicative KM	Negative AP ALS	Not Tested Not routinely performed on KM positive stains/ ALS negative
EVPZDB	Positive/ Indicative KM	Negative ALS Micro P30	Not Tested Saliva testing not being done in this lab at this time
EXGQ74	Positive/ Indicative KM Mod Taka	Negative AP (SAP)	Negative Phadebas Amylase
F42C9M	Positive/ Indicative Phenolphthalein	Negative RSID	Negative SaligaE
FVMNEL	O-Tol Hemtrace[sic]	Negative ALS	Negative ALS
FXB6NT	Positive/ Indicative KM	Negative AP PSA Micro	Negative Phadebas
G74RZB	Positive/ Indicative O-tol HemaTrace	Negative ALS AP	Not Tested No information to indicate oral contact and ALS was negative
G8XFEZ	Positive/ Indicative TMB HemaTrace Luminol	Negative AP	Not Tested
G9NLVQ	Positive/ Indicative LMG	Positive/ Indicative AP	Not Tested
GFC6FU	Positive/ Indicative KM ABA	Negative AP	Negative Phadebas
GPVJEP	Positive/ Indicative Visual Exam KM	Negative AP	Negative Phadebas Paper

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
GV43GN	Positive/ Indicative LMG	Negative ALS AP Micro	Negative Amylase diffusion
GZXZM8	Positive/ Indicative KM	Negative Micro P30 ALS	Not Tested
H7J77Z	Positive/ Indicative KM Takayama	Not Tested Analyst's discretion (based on scenario provided).	Not Tested Analyst's discretion (based on scenario provided).
H9ADNB	Positive/ Indicative O-tol HemaTrace	Negative ALS AP	Not Tested Bloodstain was very heavy and concentrated and did not appear diluted
HDYUYF	Positive/ Indicative KM Abacard Hematrace	Negative AP Micro Abacard p30 ALS	Not Tested
HEDJ2T	Positive/ Indicative KM ABAcad HemaTrace	Negative AP	Negative Phadebas test
HHBZEM	Positive/ Indicative LMG HemaTrace	Negative AP Micro PSA ALS	Negative Amylase diffusion
HQ63HJ	Positive/ Indicative O-tol Hematrace	Negative AP Micro RSID ALS	Negative RSID ALS
HWRXDW	Positive/ Indicative KM HemaTrace	Negative AP Micro	Negative Radial Diffusion

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
J23ATA	Positive/ Indicative KM Hexagon OBTI	Not Tested Based on scenario, semen test not conducted.	Not Tested
J94HFP	Positive/ Indicative TMB	Not Tested	Negative Phadebas
JHVYAZ	Positive/ Indicative KM Modified Takayama	Negative AP	Negative Phadebas
JJM2UU	Positive/ Indicative KM	Negative AP blot	Negative Phadebas blot
JNFWZ9	Positive/ Indicative K-M	Not Tested No scenario information that would make it likely to find semen	Not Tested Same reason as above [Item 5 - Semen Tests Performed]
JPQ3UN	Positive/ Indicative LMG	Positive/ Indicative AP	Not Tested
JRWLWQ	Positive/ Indicative Hemochromogen (+ve) <u>Other Tests Performed</u> Hair exam - Microscopic exam Contact - DNA sent for DNA analysis	Not Tested	Not Tested Current protocol requires a specific request for saliva as it is not a routine test.
JTMRDG	Positive/ Indicative KM RSID-Blood	Negative ALS AP	Negative ALS RSID-Saliva
JYWP6A	Positive/ Indicative KM Hexagon OBTI	Not Tested Based on the scenario information provided, semen test was not conducted.	Not Tested
K6XJEV	Positive/ Indicative KM	Not Tested Case scenario did not require testing for semen	Not Tested Case scenario did not require testing for saliva

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
K7UYQR	Positive/ Indicative Visual ALS KM	Negative AP	Negative Phadebas
KQ3XGF	Positive/ Indicative KM Takayama ALS Infrared	Negative AP Micro ABAcad p30 ALS	Inconclusive radial diffusion
KUE3RF	Positive/ Indicative O-TOL HEMATRACE	Negative AP MICRO RSID ALS	Negative RSID ALS
KXYWNG	Positive/ Indicative HemaTrace	Negative ABAcad P30	Negative Saligae
KZ3R7Y	Positive/ Indicative LMG HemaTrace Ouchterlony <u>Other Tests Performed</u> Biological fluid - ALS	Negative AP	Not Tested
LXURT6	Positive/ Indicative Visual exam KM HemaTrace	Negative AP ALS	Negative RSID ALS
M9CTN6	Positive/ Indicative Visual KM	Negative ALS	Not Tested Assumed headpiece was worn as designed.
MED8A7	Positive/ Indicative O-Tol Hematrace	Negative ALS AP Micro RSID	Negative ALS RSID
MEWRED	Positive/ Indicative O-tol <u>Other Tests Performed</u> Semen, Saliva, Sweat - ALS (See additional comments page [Table 3])	Negative AP, Micro	Not Tested Did not test for saliva based on case scenario.

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
MQAX88	Positive/ Indicative KM Hexagon OBTI	Not Tested Based on scenario information, semen tests not conducted.	Not Tested
MY39JM	Positive/ Indicative KM	Negative ALS AP	Not Tested *saliva testing not performed in this lab
MZHQHM	Positive/ Indicative Hemochromogen <u>Other Tests Performed</u> Contact DNA - sent for DNA analysis Hair Exam - Microscopic Exam	Not Tested	Not Tested Current protocol requires a specific request for saliva as it is not a routine test.
N3ZCLU	Positive/ Indicative KM	Negative AP p30 KPIC ALS	Negative Phadebas ALS
N9PM3J	Positive/ Indicative LMG HemaTrace	Negative AP Micro ALS	Negative Amylase Diffusion
N9PRB2	Positive/ Indicative KM	Negative AP ALS	Not Tested KM positive stains are not routinely tested for saliva
NQ6YNH	Positive/ Indicative KM RSID-Blood	Negative AP RSID-Semen	Not Tested
P8P4F9	Positive/ Indicative O-tol	Negative AP Microscopy	Negative RSID Saliva
PCGF7Y	Positive/ Indicative KM	Not Tested Scenario didn't indicate the need for this testing.	Not Tested Scenario didn't indicate the need for this testing.

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
PDCWKG	Positive/ Indicative KM O-tol ABA Card	Negative AP	Negative SALlgAE
PJ2EYR	Positive/ Indicative LMG Hematrace	Not Tested ALS No further testing performed based on case circumstances and ALS results	Not Tested Test not offered.
PJ7RQU	Positive/ Indicative KM Hematrace	Negative AP Seratec Micro	Negative RSID
PKWW3E	Positive/ Indicative Visual KM	Not Tested Examining for 'wearer' therefore item did not necessitate	Not Tested Blood found therefore no further examination as headwear
PMLJF3	Positive/ Indicative O-tol HemaTrace	Negative AP ALS	Negative RSID ALS
PNF7UR	Positive/ Indicative TMB Luminol HemaTrace Haemochromogen	Negative AP	Negative Phadebas (amylase)
QADXPB	Positive/ Indicative Hemastix (TMB)	Not Tested Testing for semen not done on this item - not applicable according to case scenario (suspect's clothing)	Negative Phadebas Press Test
QP6VRD	Positive/ Indicative KM <u>Other Tests Performed</u> ALS	Not Tested	Not Tested
R9VBPE	Positive/ Indicative HemaTrace LMG	Negative AP	Not Tested

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
RBKG8Q	Positive/ Indicative KM <u>Other Tests Performed</u> ALS	Negative AP P-30 KPIC	Negative Phadebas
RDTPGB	Positive/ Indicative O-tol HemaTrace	Negative AP ALS	Negative RSID ALS
RNACNQ	Positive/ Indicative TMB HemaTrace Luminol	Negative AP	Negative Amylase
RNQTHJ	Positive/ Indicative LMG OBTI <u>Other Tests Performed</u> Hair: Micro Forensic light source	Negative AP Micro	Not Tested
RVQ877	Positive/ Indicative LMG HemaTrace	Negative ALS AP Micro	Negative Amylase diffusion
T2V4RR	Positive/ Indicative KM Hematrace	Not Tested No biological fluorescence observed	Not Tested No biological fluorescence observed
T7LPCY	Positive/ Indicative KM	Not Tested The item was allegedly worn by suspect, test for semen would not connect suspect to victim in this scenario.	Not Tested No currently validated procedures for saliva test.
TJBFHW	Positive/ Indicative KM	Negative AP	Not Tested NT (Per Policy) - Item not eligible for saliva testing
TKLJQE	Positive/ Indicative LMG HemaTrace	Negative ALS AP Micro	Negative Amylase Diffusion

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
TX4EYX	Positive/ Indicative KM	Negative ALS Micro P30	Not Tested
U3J3HU	Positive/ Indicative KM	Negative AP	Not Tested Per Policy, not required to test
UAYML8	Positive/ Indicative KM Test RSID-Blood	Negative ALS AP Test	Negative ALS RSID-Saliva
UAZKRV	Positive/ Indicative KM	Negative Micro P30 ALS	Not Tested
UF9YH8	Positive/ Indicative O-tol Hematrace	Not Tested ALS	Not Tested ALS
UHA7HB	Positive/ Indicative KM		
UPT6GW	Positive/ Indicative KM	Negative AP	Negative Phadebas amylase
UPU6B6	Positive/ Indicative KM	Not Tested	Not Tested
V4BH9W	Positive/ Indicative KM	Negative ALS	Not Tested Case scenario
V7F8B6	Positive/ Indicative O-Tol RSID	Negative ALS RSID AP Micro	Negative RSID
V9Z287	Positive/ Indicative KM	Negative AP	Not Tested

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
VKV7HR	Positive/ Indicative Tetrabase RSID	Negative Micro PSA RSID AP	Negative RSID
VM2TTB	Positive/ Indicative LMG Hematrace	Negative AP Micro ALS	Negative Amylase Diffusion
W4LNZR	Positive/ Indicative KM HemaTrace	Negative AP Micro exam	Negative Amylase with radial diffusion
WBNQ8D	Positive/ Indicative KM HemaTrace	Negative AP	Negative Phadebas
WF2Q7P	Positive/ Indicative KM Hematrace <u>Other Tests Performed</u> Microscopic[sic] examination for hairs Minitaped	Negative AP	Negative Phadebas
WM2X2R	Positive/ Indicative PTMB HemaTrace	Negative Acid Phosphatase (AP)	Negative Phadebas (amylase)
WXGZ8Q	Positive/ Indicative K-M	Negative AP	Not Tested Per policy, item does not meet criteria to be tested for this body fluid - therefore not required to test.
XD83L9	Positive/ Indicative KM	Negative ALS	Not Tested Not competent in saliva testing
XHZBMU	Positive/ Indicative KM	Negative ALS Micro P30	Not Tested
XJC9UK	Positive/ Indicative KM	Negative AP	

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
XXUUMQ	Positive/ Indicative KM	Negative AP	Not Tested Saliva test not performed on this item per laboratory policy - not required to test
Y28ZFH	Positive/ Indicative LMG HemaTrace	Not Tested Testing not performed due to the probative value of the sample based on the scenario given.	Not Tested This laboratory does not offer saliva testing.
Y8GAQK	Positive/ Indicative KM Modified Takayama	Negative AP	Negative Phadebas Amylase
Y8WHUU	Positive/ Indicative LMG HemaTrace Ouchterlony - Human <u>Other Tests Performed</u> Biological Fluids - ALS	Negative AP	Not Tested
Y9VLM6	Positive/ Indicative Phenolphthalein	Negative RSID	Negative SaligaE
YBDZQ6	Positive/ Indicative Visual KM	Negative AP	Negative Phadebas
YLWGGQ	Positive/ Indicative Kastle-Meyer Test	Negative ALS Microscopic P30	Not Tested
YMCGYU	Positive/ Indicative KM Hexagon OBTI	Not Tested Based on scenario information provided, semen test was not conducted.	Not Tested
YN3XWA	Positive/ Indicative Hemastix (+)ve @ 5 secs	Not Tested Polilight AP Test - Not Tested. Based on information supplied there is no suggestion there would be semen on an item worn on the assailants head.	Negative Phadebas Press Test

Table 1 - Item 5

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
YQT4CZ	Positive/ Indicative KM	Not Tested No other stains observed	Not Tested No other stains observed
YTUAPK	Positive/ Indicative LMG HemaTrace Ouchterlony (human) <u>Other Tests Performed</u> Bodily Fluids - ALS	Negative AP	Not Tested
YZDAA4	Positive/ Indicative Phenolphthalein	Negative RSID	Negative SaligaE
ZCTE3K	Positive/ Indicative LMG HemaTrace Visual	Negative AP ALS	Negative RSID
ZKJZ9J	Positive/ Indicative KM Hematrace	Negative AP Seratec Micro	Negative RSID
ZRYNXE	Positive/ Indicative KM	Negative AP	Negative PHADEBAS (BLOTTING METHOD)
ZYFE47	Positive/ Indicative LMG Hematrace	Negative ALS AP Micro	Negative Amylase Diffusion

Response Summary for Table 1 - Item 5

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	149	2	0
Negative	0	114	77
Inconclusive	0	1	2
Not Tested (NT)	0	32	69

Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
27CG2C	KM	AP MICRO	PHADEBAS
28BLGK	Positive/ Indicative LMG HemaTrace Ouchterlony <u>Other Tests Performed</u> Biological Fluids - ALS	Positive/ Indicative AP P30 Micro	Not Tested
2EMTZU	Positive/ Indicative KM	Positive/ Indicative AP Micro	Not Tested Not required per policy
2FHERZ	Positive/ Indicative O-tol	Positive/ Indicative micro ALS AP	Negative ALS RSID Saliva
2PQVGK	Positive/ Indicative KM <u>Other Tests Performed</u> Sperm Search	Positive/ Indicative PSA	Negative RSID
2VEVGY	Positive/ Indicative O-tol HemaTrace	Positive/ Indicative ALS AP Micro	Negative RSID
3FF4BN	Positive/ Indicative KM Hematrace	Positive/ Indicative AP Micro	Negative Phadebas
3MXG82	Positive/ Indicative KM	Positive/ Indicative AP Micro	Negative Phadebas paper test
3P4YCN	Positive/ Indicative KM	Positive/ Indicative ALS AP Micro	Not Tested Semen was detected and there is no reason to believe saliva would also be present in the underwear.

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
4APZYB	Positive/ Indicative KM HemaTrace	Positive/ Indicative AP Micro	Negative Amylase (radial diffusion)
4BJLQG	Positive/ Indicative LMG HemaTrace Ouchterlony <u>Other Tests Performed</u> Biological Fluids - ALS	Positive/ Indicative AP PSA Micro	Not Tested
4HXUT4	Positive/ Indicative LMG	Positive/ Indicative AP	Not Tested
4JVXRK	Positive/ Indicative KM	Positive/ Indicative AP Micro	Not Tested Not tested per policy *Not Required
68YCXV	Positive/ Indicative OT	Positive/ Indicative Micro ALS AP	Negative RSID
6HVC DK	Positive/ Indicative Visual KM HemaTrace	Positive/ Indicative ALS AP micro	Negative ALS RSID-saliva
6PZ4XB	Positive/ Indicative KM	Positive/ Indicative AP Micro	Not Tested Case scenario did not require testing for saliva
6UWLGW	Positive/ Indicative O-tol Heme Trace[sic]	Positive/ Indicative ALS AP Micro	Negative ALS RISD[sic]
6ZGMFL	Positive/ Indicative Kastle-Meyer	Positive/ Indicative Micro ALS	Not Tested
7M6VAR	Positive/ Indicative LMG <u>Other Tests Performed</u> ALS	Positive/ Indicative ACP, Micro	Negative Amylase

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
7MNLCK	Positive/ Indicative Kastle-Meyer Test	Positive/ Indicative Microscopic ALS	Not Tested Per Agency's Policy
7XNELB	Positive/ Indicative KM	Positive/ Indicative AP P30 KPIC	Negative Phadebas
8C6JGB	Positive/ Indicative KM Modified Takayama	Positive/ Indicative AP Micro	Negative Phadebas
8NUMDA	Positive/ Indicative LMG HemaTrace	Positive/ Indicative Micro PSA ALS	This laboratory does not perform testing for saliva.
8PP7HW	Positive/ Indicative LMG spot test	Not Tested	Negative Phadebas press test
8XLW8D	Positive/ Indicative LMG HemaTrace Visual	Positive/ Indicative ALS AP P30 Micro	Negative RSID-Saliva
8YD3JW	Positive/ Indicative LMG: Reddish staining observed on inside crotch area, LMG presumptive test performed, LMG positive <u>Other Tests Performed</u> ALS: crimelite screen for semen & saliva	Positive/ Indicative AP Back: AP negative AP Front: AP positive Micro: 1+ heads, negative for NECs	Negative Phadebas: Phadebas negative
9B677N	Positive/ Indicative KM	Positive/ Indicative AP Micro	Not Tested Not tested per policy.
9GQWTH	Positive/ Indicative KM	Positive/ Indicative ALS AP Micro	Not Tested Not applicable to case scenario

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
9HLHJN	Positive/ Indicative KM	Positive/ Indicative AP Micro	Not Tested Per laboratory policy, only intimate swabs from the vaginal area are saliva tested.
9JH363	Positive/ Indicative KM ABACard HemaTrace	Positive/ Indicative AP PSA	Negative Phadebas Press
9Q3XCX	Positive/ Indicative KM	Negative AP ALS	Negative AMYLASE DIFFUSION ALS
9UE4NY	Positive/ Indicative LMG	Positive/ Indicative AP Micro	Amylase Diffusion
9W3CXZ	Positive/ Indicative KM ABAcard Hematrace	Positive/ Indicative AP PSA Semiquant	Negative Phadebas press
AWHZRY	Positive/ Indicative KM <u>Other Tests Performed</u> Semen - ALS	Positive/ Indicative AP PSA Micro	Negative Phadebas
B2DAF7	Positive/ Indicative Hemastix - pos 5 secs	Polylight[sic] - v. faint fluorescence AP - pos 10 secs Micro - pos	Not Tested No indication in scenario about oral rape and was semen positive.
B33FWW	Positive/ Indicative KM RSID	Positive/ Indicative AP RSID Micro	Not Tested
BA8DNA	Positive/ Indicative Visual LMG HemaTrace	Positive/ Indicative ALS AP P30 Spermatozoa (Micro)	Negative RSID

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
BHXXRM	Positive/ Indicative Item 6A: O-tol	Positive/ Indicative Item 6: ALS Item 6A: AP, Micro	Negative Item 6A: RSID
BMUKKB	Positive/ Indicative KM Hematrace	Positive/ Indicative AP Seratec Micro	Negative RSID
BNKQXU	Positive/ Indicative KM HEMATRACE	Positive/ Indicative PSA AP MICRO	Negative RSID
BU8H89	Positive/ Indicative LMG HemaTrace Ouchterlony Species Determination <u>Other Tests Performed</u> Bodily Fluid - ALS	Positive/ Indicative AP PSA Micro	Not Tested
BYM4XV	Positive/ Indicative LMG HemaTrace	Positive/ Indicative AP Micro ALS AP mapping	Negative Amylase diffusion Amylase mapping
C4UYB4	Positive/ Indicative LMG HemaTrace	Positive/ Indicative ALS Micro	Not Tested Saliva testing is not offered as a service in this laboratory.
CGLJZ8	Positive/ Indicative KM Modified Takayama	Positive/ Indicative AP Micro	Not Tested Not tested based on case scenario
DL73AZ	Positive/ Indicative KM ABACard HemaTrace	Positive/ Indicative AP PSA Micro	Negative Phadebas Press
EB38QV	Positive/ Indicative KM	Positive/ Indicative PSA MICRO	Negative PHADEBAS

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
EB9GR9	Positive/ Indicative TMB	Positive/ Indicative ALS AP on bloodstain AP near bloodstain	Inconclusive ALS RSID Saliva Field Kit (on bloodstain)
EFF3F6	Positive/ Indicative KM Takayama	Positive/ Indicative AP (Seri AP) Micro (Xmas Stain)	Negative Amylase (Phadebas)
ERDV3B	Positive/ Indicative KM	Positive/ Indicative ALS AP PSA Micro	Negative RSID
EVPZDB	Positive/ Indicative KM	Positive/ Indicative ALS Micro	Not Tested Saliva testing not being done in this lab at this time
EXGQ74	Positive/ Indicative KM Mod Taka	Positive/ Indicative AP (SAP) Micro (Xmas)	Not Tested Stain in crotch, where vaginal secretions may be, (known to contain amylase and react + with phadebas test) would not test crotch for saliva with phadebas test
F42C9M	Positive/ Indicative Phenolphthalein	Positive/ Indicative RSID Christmas Tree Stain	Negative SaligaE
FVMNEL	Not Tested No visible red-brown stains	Positive/ Indicative ALS AP Micro	Negative ALS RSID
FXB6NT	Positive/ Indicative KM	Positive/ Indicative AP PSA Micro	Negative Phadebas
G74RZB	Positive/ Indicative O-tol HemaTrace	Positive/ Indicative AP Micro	Not Tested No information to indicate oral contact by suspect and semen is present in this stain.

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
G8XFEZ	Positive/ Indicative TMB HemaTrace Luminol	Positive/ Indicative AP PSA Microscopic	Not Tested
G9NLVQ	Positive/ Indicative LMG	Positive/ Indicative AP	Not Tested
GFC6FU	Positive/ Indicative KM ABA	Positive/ Indicative AP PSA	Negative Phadebas
GPVJEP	Positive/ Indicative Visual Exam KM	Positive/ Indicative AP H&E Slide	Not Tested
GV43GN	Positive/ Indicative LMG	Positive/ Indicative ALS AP Micro PSA	Negative Amylase diffusion
GZXZM8	Positive/ Indicative KM	Positive/ Indicative Micro ALS	Not Tested
H7J77Z	Positive/ Indicative KM Takayama	Positive/ Indicative AP (Map) Micro (Xmas)	Not Tested Analyst's discretion (based on scenario provided).
H9ADNB	Positive/ Indicative O-tol HemaTrace	Positive/ Indicative ALS AP Micro	Positive/ Indicative RSID
HDYUYF	Positive/ Indicative KM Abacard Hematrace	Positive/ Indicative AP Micro Abacard p30 ALS	Not Tested
HEDJ2T	Positive/ Indicative KM ABAcad HemaTrace	Positive/ Indicative AP ABAcad P30 Micro	Negative Phadebas test

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
HHBZEM	Positive/ Indicative LMG HemaTrace	Positive/ Indicative AP Micro ALS Evidence mapping	Negative Amylase diffusion Evidence mapping
HQ63HJ	Positive/ Indicative O-tol Hematrace	Positive/ Indicative AP Micro ALS	Negative RSID ALS
HWRXDW	Positive/ Indicative KM HemaTrace	Positive/ Indicative AP Micro	Negative Radial Diffusion
J23ATA	Positive/ Indicative KM Hexagon OBTI	Positive/ Indicative AP RSID-Semen PSA	Not Tested
J94HFP	Positive/ Indicative TMB	Positive/ Indicative AP Micro	Not Tested
JHVYAZ	Positive/ Indicative KM Modified Takayama	Positive/ Indicative AP Micro	Negative Phadebas
JJM2UU	Positive/ Indicative KM	Positive/ Indicative AP blot Micro	Negative Phadebas blot
JNFWZ9	Positive/ Indicative K-M	Positive/ Indicative ALS AP Microscopic sperm search	Not Tested No scenario info provided that would indicate saliva as a possibility
JPQ3UN	Positive/ Indicative LMG	Positive/ Indicative AP	Not Tested
JRWLWQ	Positive/ Indicative Hemochromogen (+ve)	Inconclusive Fast Blue (+ve) (AP) Micro (-ve) P30 (+ve)	Not Tested Current protocol requires a specific request for saliva as it is not a routine test.

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
JTMRDG	Positive/ Indicative KM RSID-Blood	Positive/ Indicative ALS AP RSID-Semen Microscopy	Negative ALS RSID-Saliva
JYWP6A	Positive/ Indicative KM Hexagon OBTI	Positive/ Indicative AP RSID PSA	Not Tested
K6XJEV	Positive/ Indicative KM	Positive/ Indicative AP Micro	Not Tested Case scenario did not require testing for saliva
K7UYQR	Positive/ Indicative Visual KM	Positive/ Indicative AP Micro	Negative Phadebas
KQ3XGF	Positive/ Indicative KM Takayama	Positive/ Indicative AP Micro ABAcad p30 ALS	Inconclusive radial diffusion
KUE3RF	Positive/ Indicative O-TOL HEMATRACE	Positive/ Indicative AP MICRO ALS	Negative RSID ALS
KXYWNG	Positive/ Indicative HemaTrace	Positive/ Indicative ABAcad P30	Negative Saligae
KZ3R7Y	Positive/ Indicative LMG HemaTrace Ouchterlony <u>Other Tests Performed</u> Biological Fluids - ALS	Positive/ Indicative AP PSA Microscopic	Not Tested
LXURT6	Positive/ Indicative Visual exam KM HemaTrace	Positive/ Indicative AP Micro ALS	Negative RSID ALS

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
M9CTN6	Positive/ Indicative Visual KM	Positive/ Indicative ALS AP Micro	Not Tested Would normally not test panties for saliva unless oral sodomy was reported.
MED8A7	Positive/ Indicative O-Tol	Positive/ Indicative ALS AP Micro	Negative ALS RSID
MEWRED	Positive/ Indicative O-tol <u>Other Tests Performed</u> Semen, Saliva, Sweat - ALS	Positive/ Indicative AP, Micro	Negative RSID Saliva
MQAX88	Positive/ Indicative KM Hexagon OBTI	Positive/ Indicative AP RSID PSA	Not Tested
MY39JM	Positive/ Indicative KM	Positive/ Indicative ALS AP	Not Tested *saliva testing not performed in this lab
MZHQHM	Positive/ Indicative Hemochromogen Hemastix	Negative Fast Blue (AP)	Not Tested Current protocol requires a specific request for saliva as it is not a routine test.
N3ZCLU	Positive/ Indicative KM	Positive/ Indicative AP p30 KPIC ALS	Negative Phadebas ALS
N9PM3J	Positive/ Indicative LMG HemaTrace	Positive/ Indicative AP Micro ALS	Negative Amylase Diffusion
N9PRB2	Positive/ Indicative KM	Positive/ Indicative AP Micro (sperm search) PSA ALS	Not Tested underwear crotch stains are not routinely tested for saliva unless oral sex is indicated

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
NQ6YNH	Positive/ Indicative KM RSID-Blood	Positive/ Indicative AP Micro	Not Tested
P8P4F9	Positive/ Indicative O-tol	Positive/ Indicative AP Microscopy	Negative RSID Saliva
PCGF7Y	Positive/ Indicative KM	Positive/ Indicative AP Micro	Negative Amylase with Radial Diffusion
PDCWKG	Positive/ Indicative KM O-tol ABA Card	Positive/ Indicative AP	Negative SALlgAE
PJ2EYR	Positive/ Indicative LMG Hematrace	Positive/ Indicative ALS Micro P30	Not Tested Testing not offered.
PJ7RQU	Positive/ Indicative KM Hematrace	Positive/ Indicative AP Seratec Micro	Negative RSID
PKWW3E	Positive/ Indicative Visual KM <u>Other Tests Performed</u> Item also fibre taped	Positive/ Indicative AP Slide made	Negative Phadebas (only inside front and crotch tested)
PMLJF3	Positive/ Indicative O-tol HemaTrace	Positive/ Indicative AP Micro PSA ALS	Not Tested Item 6 was not tested for saliva because the case scenario did not indicate an oral assault.
PNF7UR	Positive/ Indicative TMB Luminol HemaTrace	Positive/ Indicative AP PSA Micro	Not Tested Deemed not relevant to scenario

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
QADXPH	Positive/ Indicative Hemastix (TMB)	Positive/ Indicative AP Test Microscopy	Negative Phadebas Press Test
QP6VRD	Positive/ Indicative KM <u>Other Tests Performed</u> ALS	Positive/ Indicative BCIP P30 Microscopic Sperm Search	Positive/ Indicative Amylase Diffusion
R9VBPE	Positive/ Indicative HemaTrace LMG	Negative Microscopy AP	Not Tested
RBKG8Q	Positive/ Indicative KM <u>Other Tests Performed</u> ALS	Positive/ Indicative AP P-30 KPIC	Negative Phadebas
RDTPGB	Positive/ Indicative O-tol HemaTrace	Positive/ Indicative AP Micro ALS	Negative RSID ALS
RNACNQ	Positive/ Indicative TMB HemaTrace Luminol	Positive/ Indicative AP Micro	Negative Amylase
RNQTHJ	Positive/ Indicative LMG OBTI <u>Other Tests Performed</u> Forensic light source	Positive/ Indicative AP Micro	Not Tested
RVQ877	Positive/ Indicative LMG HemaTrace	Positive/ Indicative ALS AP Micro AP mapping	Negative Amylase diffusion Amylase mapping (Phadebas)
T2V4RR	Positive/ Indicative KM Hematrace	Positive/ Indicative AP Seratec Micro	Negative RSID

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
T7LPCY	Positive/ Indicative KM	Positive/ Indicative Microscopic ALS	Not Tested No currently validated procedures for saliva test.
TJBFHW	Positive/ Indicative KM	Positive/ Indicative AP Micro	Not Tested NT (Per Policy) - Item not eligible for saliva testing
TKLJQE	Positive/ Indicative LMG Hematrace	Positive/ Indicative ALS AP Micro PSA	Negative Amylase Diffusion
TX4EYX	Positive/ Indicative KM	Positive/ Indicative ALS Micro	Not Tested
U3J3HU	Positive/ Indicative KM	Positive/ Indicative AP Micro	Not Tested Per Policy, not required to test
UAYML8	Positive/ Indicative KM Test RSID-Blood	Positive/ Indicative ALS AP RSID-Semen Microscopy	Negative ALS RSID-Saliva
UAZKRV	Positive/ Indicative KM	Positive/ Indicative Micro ALS	Not Tested
UF9YH8	Positive/ Indicative O-tol	Positive/ Indicative ALS AP Micro	Negative ALS RSID
UHA7HB	Positive/ Indicative KM	AP	
UPT6GW	Positive/ Indicative KM	Positive/ Indicative PSA Micro	Negative Phadebas amylase

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
UPU6B6	Not Tested	Positive/ Indicative AP MICRO PSA	Not Tested
V4BH9W	Positive/ Indicative KM	Positive/ Indicative ALS AP Micro	Not Tested Case scenario
V7F8B6	Positive/ Indicative O-Tol RSID	Positive/ Indicative AP ALS Micro	Negative RSID
V9Z287	Positive/ Indicative KM	Positive/ Indicative AP	Not Tested PHADEBAS
VKV7HR	Positive/ Indicative Tetrabase RSID	Positive/ Indicative AP Micro PSA RSID	Negative RSID
VM2TTB	Positive/ Indicative LMG Hematrace	Positive/ Indicative AP Micro ALS	Negative Amylase Diffusion
W4LNZR	Positive/ Indicative KM HemaTrace	Positive/ Indicative AP Micro exam P30	Negative Amylase with radial diffusion
WBNQ8D	Positive/ Indicative KM HemaTrace	Positive/ Indicative AP Micro P30	Negative Phadebas
WF2Q7P	Positive/ Indicative KM Hematrace	Positive/ Indicative AP Micro Rsid semen	Negative Phadebas

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
WM2X2R	Positive/ Indicative PTMB HemaTrace <u>Other Tests Performed</u> Urine - Creatinine Fecal - Micro Fecal - Urobilinogen	Positive/ Indicative Acid Phosphatase (AP) Micro	Negative Phadebas (amylase)
WXGZ8Q	Positive/ Indicative K-M	Positive/ Indicative AP Micro	Not Tested Per policy, item does not meet criteria to be tested for this body fluid - therefore not required to test.
XD83L9	Positive/ Indicative KM	Positive/ Indicative ALS AP	Not Tested Not competent in saliva testing
XHZBMU	Positive/ Indicative KM	Positive/ Indicative ALS Micro	Not Tested
XJC9UK	Positive/ Indicative KM	Positive/ Indicative AP Micro ABA p30	Negative Amylase
XXUUMQ	Positive/ Indicative KM	Positive/ Indicative AP Micro	Not Tested Saliva test not performed on this item per laboratory policy - not required to test
Y28ZFH	Positive/ Indicative LMG HemaTrace	Positive/ Indicative AP Micro PSA ALS	Not Tested This laboratory does not offer saliva testing.
Y8GAQK	Positive/ Indicative KM Modified Takayama	Positive/ Indicative AP Micro - Xmas tree stain	Not Tested As per case scenario, this item is from the victim. Samples from the vaginal area can give a false positive result for amylase, so this test was not performed.

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
Y8WHUU	Positive/ Indicative LMG HemaTrace Ouchterlony-Human <u>Other Tests Performed</u> Biological Fluids - ALS	Positive/ Indicative AP PSA Micro	Not Tested
Y9VLM6	Positive/ Indicative Phenolphthalein	Positive/ Indicative RSID Christmas Tree Stain	Negative SaligaE
YBDZQ6	Positive/ Indicative Visual KM	Positive/ Indicative AP microscopic exam - spermatazoa[sic] identified	Negative Phadebas
YLWGGQ	Positive/ Indicative Kastle-Meyer Test	Positive/ Indicative ALS Microscopic	Not Tested
YMCGYU	Positive/ Indicative KM Hexagon OBTI	Positive/ Indicative AP RSID-Semen PSA	Not Tested
YN3XWA	Positive/ Indicative Hemastix (+)ve @ 5 secs	Positive/ Indicative Polilight AP Test Micro	Not Tested Phadebas - Not Tested. Hemastix and AP positive tests place the assailant and victim together. No suggestion of oral rape.
YQT4CZ	Positive/ Indicative KM	Positive/ Indicative AP Micro PSA	Not Tested blood and semen found no indication in scenario of possible saliva
YTUAPK	Positive/ Indicative LMG HemaTrace Ouchterlony (human) <u>Other Tests Performed</u> Bodily Fluids - ALS	Positive/ Indicative AP PSA Micro	Not Tested
YZDAA4	Positive/ Indicative Phenolphthalein	Positive/ Indicative RSID Christmas Tree Stain	Negative SaligaE

Table 1 - Item 6

Webcode	Blood Conclusion Tests Performed	Semen Conclusion Tests Performed	Saliva Conclusion Tests Performed
ZCTE3K	Positive/ Indicative LMG HemaTrace Visual	Positive/ Indicative AP P30 Microscopic ALS	Negative RSID
ZKJZ9J	Positive/ Indicative KM Hematrace	Positive/ Indicative AP Seratec Micro	Negative RSID
ZRYNXE	Positive/ Indicative KM	Positive/ Indicative AP MICRO	Negative PHADEBAS (BLOTTING METHOD)
ZYFE47	Positive/ Indicative LMG Hematrace	Positive/ Indicative ALS AP Micro	Negative Amylase Diffusion

Response Summary for Table 1 - Item 6

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	148	143	2
Negative	0	3	78
Inconclusive	0	1	2
Not Tested (NT)	2	1	65

Serology Screening Results

Indicate the results of any screening tests performed on each Item.

TABLE 2 - Item 1

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	KM	Negative		Visual	Negative
28BLGK	LMG	Negative	8YD3JW	Yellowish staining observed, but no visible blood-like staining observed, therefore no presumptive test for blood performed	
2EMTZU	KM	Negative			
2FHERZ	O-tol	Negative			
2PQVGK	KM	Negative	9B677N	KM	Negative
2VEVGY	O-tol	Negative	9GQWTH	KM	Negative
3FF4BN	KM	Negative	9HLHJN	KM	Negative
3MXG82	KM	Negative	9JH363	KM	Negative
3P4YCN	KM	Negative	9Q3XCX	KM	Negative
4APZYB	KM	Negative	9UE4NY	Visual exam - no staining observed	
4BJLQG	LMG	Negative	9W3CXZ	KM	Negative
4HXUT4	LMG	Negative	AWHZRY	KM	Negative
4JVXRK	KM	Negative	B2DAF7	Hemastix - wk pos 25 secs	Positive
68YCXV	OT	Negative		Hematrace - pos 3 mins	Positive
6HVCDK	Visual	Negative	B33FWW	Not tested: No red-brown staining observed on the Items.	
6PZ4XB	KM	Negative	BA8DNA	Visual	Positive
6UWLGW	O-tol	Negative		LMG	Negative
6ZGMFL	Swab was not red stained		BHXXRM	Not tested due to lack of reddish-brown color on swabs.	
7M6VAR	No visible stain observed.		BMUKKB	KM	Negative
7MNLCK	Swabs have no red stains		BNKQXU	Hematrace	Positive
7XNELB	KM	Negative	BU8H89	LMG	Negative
8C6JGB	KM	Negative	BYM4XV	no apparent bloodstain observed	
8NUMDA	LMG	Negative	C4UYB4	LMG	Negative
	HemaTrace	Positive		HemaTrace	Positive
8PP7HW	LMG test	Negative			
8XLW8D	LMG	Negative			

TABLE 2 - Item 1

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
CGLJZ8	KM	Negative		ABACard HemaTrace	Negative
DL73AZ	KM	Negative	HHBZEM	A red stain was not observed.	
EB38QV	KM	Negative	HQ63HJ	O-tol	Negative
EB9GR9	TMB	Negative	HWRXDW	KM	Negative
EFF3F6	KM	Negative	J23ATA	KM	Negative
ERDV3B	KM	Negative	JHVYAZ	KM	Negative
EVPZDB	White scene or kit swabs not tested for blood in this lab		JJM2UU	KM	Negative
EXGQ74	KM	Negative	JNFWZ9	K-M	Negative
F42C9M	Phenolphthalein	Negative	JPQ3UN	LMG	Negative
FVMNEL	O-Tol	Negative	JRWLWQ	Visual Exam	Negative
FXB6NT	KM	Negative	JTMRDG	KM	Negative
G74RZB	O-tol	Negative	JYWP6A	KM	Negative
G8XFEZ	TMB	Positive	K6XJEV	KM	Negative
	HemaTrace	Positive	K7UYQR	Visual	Negative
	Haemochromogen	Negative		KM	Negative
G9NLVQ	LMG	Negative	KQ3XGF	KM	Negative
GFC6FU	KM	Negative	KUE3RF	O-TOL	Negative
GPVJEP	KM	Negative	KXYWNG	HemaTrace	Negative
	Visual exam	Negative	KZ3R7Y	LMG	Negative
GV43GN	no staining on swab		LXURT6	Visual exam	Negative
GZXZM8	*swabs were not red stained		M9CTN6	Visual	Negative
H7J77Z	KM	Negative		KM	Negative
H9ADNB	O-tol	Negative	MED8A7	O-Tol	Negative
			MEWRED	Swabs were not tested for blood - does not display the characteristic reddish-brown stain associated with blood.	
HDYUYF	KM	Negative			
	Abacard Hematrace	Positive	MQAX88	KM	Negative
HEDJ2T	KM	Negative	MY39JM	KM	Negative

TABLE 2 - Item 1

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
MZHQHM	Hemastix	Negative	T2V4RR	KM	Negative
N3ZCLU	KM	Negative	T7LPCY	KM	Negative
N9PM3J	No red-brown discoloration		TJBFHW	KM	Negative
N9PRB2	KM	Negative	TKLJQE	LMG	Negative
P8P4F9	No reddish/brown stain observed during visual inspection.		TX4EYX	KM	Negative
PCGF7Y	KM	Negative	U3J3HU	Kastle-Meyer Test	Negative
PDCWKG	KM	Negative	UAYML8	KM test	Negative
	O-tol	Negative	UAZKRV	KM	Negative
	ABA Card	Negative	UF9YH8	o-tol	Negative
PJ2EYR	LMG	Negative	UHA7HB	KM	Negative
	Hematrace	Positive	UPT6GW	Visual examination	Negative
PJ7RQU	KM	Negative	V4BH9W	KM	Negative
PKWW3E	KM	Negative	V7F8B6	O-TOL	Negative
	Visual	Negative	V9Z287	KM	Negative
PMLJF3	*NT because no red brown stain observed.		VKV7HR	Tetabase	Negative
PNF7UR	TMB	Inconclusive	VM2TTB	LMG	Negative
	HemaTrace	Positive	W4LNZR	KM	Negative
QADXPB	Hemastix (TMB)	Negative	WBNQ8D	KM	Negative
QP6VRD	KM	Negative	WF2Q7P	KM TEST	Negative
R9VBPE	LMG	Negative		HemaTrace	Negative
RBKG8Q	KM	Negative	WM2X2R	PTMB	
RDTPGB	O-tol	Negative		Slight yellowish stain present.	
RNACNQ	TMB	Positive	WXGZ8Q	K-M	Negative
	HemaTrace	Positive	XD83L9	KM	Negative
	Haemachromogen[sic]	Negative	XJC9UK	KM	Negative
RNQTHJ	Visually	Negative	XXUUMQ	KM	Negative
RVQ877	No visual apparent blood staining		Y28ZFH	LMG	Negative
				HemaTrace	Positive

TABLE 2 - Item 1

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
Y8GAQK	KM	Negative			
Y8WHUU	LMG	Negative			
Y9VLM6	Phenolphthalein	Negative			
YBDZQ6	Visual	Negative			
	KM	Negative			
YLWGGQ	Not a brown stains[sic]				
YMCGYU	KM	Negative			
YN3XWA	Hemastix -ve @ 60 secs	Negative			
YQT4CZ	No red brown staining				
YTUAPK	LMG	Negative			
YZDAA4	Phenolphthalein	Negative			
ZCTE3K	Visual	Negative			
ZKJZ9J	KM	Negative			
ZRYNXE	KM	Negative			
ZYFE47	No red discoloration				

TABLE 2 - Item 1

Semen						
Webcode	Test Performed	Results	Webcode	Test Performed	Results	
27CG2C	AP	Positive	6PZ4XB	AP	Positive	
	MICRO	Positive		Micro	Positive	
28BLGK	AP	Positive	6UWLGW	Micro	Positive	
	Micro	Positive		6ZGMFL	Micro	Positive
	P30	Positive		7M6VAR	AP, Micro	Positive
2EMTZU	AP	Positive	7MNLCK	Microscopic	Positive	
	Micro	Positive		7XNELB	AP	Positive
2FHERZ	micro	Positive		P30	Positive	
2PQVGK	PSA	Positive		KPIC	Positive	
2VEVGY	AP	Positive	8C6JGB	AP	Positive	
	Micro	Positive		Micro	Positive	
3FF4BN	AP	Positive	8NUMDA	AP	Positive	
	Micro	Positive		Micro	Positive	
3MXG82	AP	Positive		PSA	Positive	
	Micro	Positive		ALS	Positive	
3P4YCN	AP	Positive	8XLW8D	AP	Positive	
	Micro	Positive		P30	Positive	
4APZYB	AP	Positive		Micro	Positive	
	Micro	Positive	8YD3JW	AP: x4 AP positive	Positive	
4BJLQG	AP	Positive		Micro: 2+ heads & tails, negative for NECs	Positive	
	PSA	Positive	9B677N	AP	Positive	
	Micro	Positive		Micro	Positive	
4HXUT4	AP	Positive	9GQWTH	AP	Positive	
4JVXRK	AP	Positive		Micro	Positive	
	Micro	Positive	9HLHJN	AP	Positive	
68YCXV	Micro	Positive		Micro	Positive	
6HVCDK	AP	Positive	9JH363	AP	Positive	
	micro	Positive		PSA	Positive	

TABLE 2 - Item 1

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
9Q3XCX	AP	Positive		Micro	Positive
	MICRO	Positive	BYM4XV	AP	Positive
9UE4NY	AP	Positive		Micro	Positive
	Micro	Positive	C4UYB4	ALS	Positive
	P30	Positive		AP	Positive
9W3CXZ	AP	Positive		Micro	Positive
	PSA Semiquant	Positive		P30	Positive
AWHZRY	AP	Positive	CGLJZ8	AP	Positive
	PSA	Positive		Micro	Positive
	Micro	Positive	DL73AZ	AP	Positive
B2DAF7	AP - pos 8 secs	Positive		PSA	Positive
	Micro - pos	Positive		Micro	Positive
	p30 - pos 1 min	Positive	EB38QV	PSA	Positive
B33FWW	AP	Positive		MICRO	Positive
	RSID	Positive	EB9GR9	AP	Positive
	Micro	Positive	EFF3F6	AP (Seri AP)	Positive
BA8DNA	AP	Positive		Micro (Xmas Stain)	Positive
	P30	Positive	ERDV3B	AP	Positive
	Spermatozoa (Micro)	Positive		PSA	Positive
BHXXRM	Item 1 AnB: AP, Micro	Positive		Micro	Positive
	Item 1 CnD: AP, Micro	Positive	EVPZDB	Micro	Positive
BMUKKB	AP	Positive	EXGQ74	AP (SAP)	Positive
	Seratec	Positive		PSA (ABA)	Positive
	Micro	Positive		Micro (Xmas)	Positive
BNKQXU	PSA	Positive	F42C9M	RSID	Negative
	Micro	Positive	FVMNEL	Micro	Positive
BU8H89	AP	Positive			
	PSA	Positive			

TABLE 2 - Item 1

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
FXB6NT	AP	Positive	HHBZEM	AP	Positive
	P30	Positive		Micro	Positive
	Micro	Positive		PSA	Positive
G74RZB	AP	Positive	HQR63HJ	ALS	Positive
	Micro	Positive		Micro	Positive
G8XFEZ	AP	Positive	HWRXDW	AP	Positive
	PSA	Positive		Micro	Positive
	Microscopic	Positive		J23ATA	RSID-Semen
G9NLVQ	AP	Positive		PSA	Positive
GFC6FU	AP	Positive	J94HFP	AP	Positive
	PSA	Positive		Micro	Positive
GPVJEP	AP	Positive	JHVYAZ	AP	Positive
	H&E Slide	Positive		Micro	Positive
GV43GN	AP	Positive	JIM2UU	AP	Positive
	Micro	Positive		Micro	Positive
GZXZM8	Micro	Positive	JNFWZ9	AP	Positive
H7J77Z	AP (Map)	Positive	JPQ3UN	Sperm search w/ microscope	Positive
	Micro (Xmas)	Positive		AP	Positive
H9ADNB	ALS	Positive	JRWLWQ	Fast blue (+ve) (AP)	Positive
	AP	Positive		Micro (+ve)	Positive
	Micro	Positive	JTMRDG	ALS	Positive
	PSA	Positive		AP	Positive
HDYUYF	AP	Positive	JYWP6A	RSID-Semen	Positive
	Micro	Positive		Microscopy	Positive
	Abacard p30	Positive		RSID	Positive
HEDJ2T	AP	Positive	K6XJEV	PSA	Positive
	ABAcad P30	Positive		AP	Positive
	Micro	Positive		Micro	Positive

TABLE 2 - Item 1

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
K7UYQR	AP	Positive		p30	Positive
	Micro	Positive		KPIC	Positive
KQ3XGF	AP	Positive	N9PM3J	AP	Positive
	Micro	Positive		Micro	Positive
	ABAcad p30	Positive		PSA	Positive
KUE3RF	AP	Positive	N9PRB2	AP	Positive
	MICRO	Positive		Micro (sperm search)	Positive
KXYWNG	ABAcad P30	Positive		PSA	Positive
KZ3R7Y	AP	Positive	NQ6YNH	AP	Positive
	PSA	Positive		Micro	Positive
	Microscopic	Positive		P8P4F9	AP
LXURT6	AP	Positive		Microscopy	Positive
	Micro	Positive	PCGF7Y	AP	Positive
M9CTN6	AP	Positive		Micro	Positive
	Micro	Positive	PDCWKG	AP	Positive
MED8A7	ALS	Positive	PJ2EYR	AP	Positive
	AP	Positive		Micro	Positive
	Micro	Positive		P30	Positive
MEWRED	AP, Micro	Positive		ALS	Positive
	AP, Micro	Positive	PJ7RQU	AP	Positive
	AP, Micro	Positive		Seratec	Positive
	AP, Micro	Positive		Micro	Positive
MQAX88	RSID	Positive	PKWW3E	AP	Positive
	PSA	Positive		Slide made (Micro)	Positive
MY39JM	AP	Positive	PMLJF3	AP	Positive
MZHQHM	Fast Blue (AP)	Positive		Micro	Positive
	Micro	Positive			
N3ZCLU	AP	Positive			

TABLE 2 - Item 1

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
PNF7UR	AP	Positive	TX4EYX	Micro	Positive
	PSA	Positive	U3J3HU	Acid Phosphatase	Positive
	Micro	Positive		Microscopic	Positive
QADXPB	Microscopy	Positive	UAYML8	ALS	Positive
QP6VRD	BCIP	Positive		AP Test	Positive
	Microscopic Sperm Search	Positive		RSID-Semen	Positive
R9VBPE	Microscopy	Positive		Microscopy	Positive
	AP	Positive	UAZKRV	Micro	Positive
RBKG8Q	AP	Positive	UF9YH8	ALS	Positive
	P-30	Positive		AP	Positive
	KPIC	Positive		Micro	Positive
RDTPGB	Micro	Positive	UHA7HB	AP	Positive
RNACNQ	AP	Positive	UPT6GW	AP	Positive
	PSA	Positive		PSA	Positive
	Micro	Positive	UPU6B6	AP	Positive
RNQTHJ	AP	Positive		MICRO	Positive
	Micro	Positive	V4BH9W	ALS	Positive
RVQ877	AP	Positive		AP	Positive
	Micro	Positive		Micro	Positive
	PSA	Positive	V7F8B6	AP	Positive
T2V4RR	AP	Positive		Micro	Positive
	Seratec	Positive	V9Z287	AP	Positive
	Micro	Positive	VKV7HR	Micro	Positive
T7LPCY	Microscopic	Positive	VM2TTB	AP	Positive
TJBFHW	AP	Positive		Micro	Positive
	Micro	Positive		PSA	Positive
TKLJQE	AP	Positive	W4LNZR	AP	Positive
	Micro	Positive		Micro exam	Positive

TABLE 2 - Item 1

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
WBNQ8D	AP	Positive	YMCGYU	RSID-Semen	Positive
	Micro	Positive		PSA	Positive
	P30	Positive	YN3XWA	Micro	Positive
WF2Q7P	AP	Positive	YQT4CZ	AP	Positive
	Micro H+E stain	Positive		Micro	Positive
WM2X2R	Acid Phosphatase (AP)	Positive	YTUAPK	AP	Positive
	microscopic (micro)	Positive		PSA	Positive
WXGZ8Q	AP	Positive		Micro	Positive
	Micro	Positive	YZDAA4	RSID	Negative
XD83L9	AP	Positive	ZCTE3K	AP	Positive
XHZBMU	Micro	Positive		P30	Positive
XJC9UK	AP	Positive		Microscopic	Positive
	Micro	Positive	ZKJZ9J	AP	Positive
XXUUMQ	AP	Positive		Seratec	Positive
	Micro	Positive		Micro	Positive
Y28ZFH	AP	Positive	ZRYNXE	AP	Positive
	Micro	Positive		MICRO	Positive
	PSA	Positive	ZYFE47	AP	Positive
	ALS	Positive		Micro	Positive
Y8GAQK	AP	Positive		PSA	Positive
	Micro - Xmas tree stain	Positive			
Y8WHUU	AP	Positive			
	PSA	Positive			
	Micro	Positive			
Y9VLM6	RSID	Negative			
YBDZQ6	AP	Positive			
	Microscopic exam - spermatzoa[sic] identified	Positive			
YLWGGQ	Microscopic	Positive			

TABLE 2 - Item 1

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	PHADEBAS	Positive	9UE4NY	Amylase Diffusion	Negative
2EMTZU	RSID	Positive	9W3CXZ	Phadebas Press	Negative
2FHERZ	RSID Saliva	Negative	AWHZRY	Phadebas	Negative
2PQVGK	RSID	Negative	B2DAF7	RSID - neg 8.5 min	Negative
2VEVGY	RSID	Positive	BA8DNA	RSID	Negative
3FF4BN	Phadebas	Negative	BHXXRM	Item 1 AnB: RSID	Negative
3MXG82	Phadebas tube test	Positive		Item 1 CnD: RSID	Negative
3P4YCN	RSID-Saliva	Negative	BMUKKB	RSID	Negative
4APZYB	Amylase (radial diffusion)	Negative	BNKQXU	RSID	Negative
4JVXRK	RSID	Positive	BYM4XV	Amylase diffusion	Negative
68YCXV	RSID	Negative	C4UYB4	Saliva testing is not offered as a service in this laboratory.	
6HVCDK	RSID-Saliva	Negative	CGLJZ8	Phadebas Amylase	Negative
6PZ4XB	Phadebas	Positive	DL73AZ	Phadebas Press	Negative
6UWLGW	RSID	Inconclusive	EB38QV	PHADEBAS	Negative
7M6VAR	Amylase	Negative	EB9GR9	RSID Saliva Field Kit	Negative
7MNLCK	Per Agency's Policy		EFF3F6	Amylase (Phadebas)	Negative
7XNELB	Phadebas	Negative	ERDV3B	RSID	Negative
8C6JGB	Phadebas	Negative	EVPZDB	Saliva testing not currently being done in this lab	
8NUMDA	This laboratory does not perform testing for saliva.		EXGQ74	Phadebas Amylase	Negative
8PP7HW	Phadebas tube test	Inconclusive	F42C9M	SaligaE	Negative
8XLW8D	RSID - Saliva	Negative	FVMNEL	RSID	Negative
8YD3JW	Phadebas: x4 Phadebas negative	Negative	FXB6NT	Phadebas	Negative
9B677N	RSID	Positive	G74RZB	RSID	Positive
9GQWTH	RSID	Negative	G9NLVQ	PHADEBAS PAPER TEST	Negative
9HLHJN	RSID	Positive	GFC6FU	Phadebas	Negative
9JH363	Phadebas Press	Negative	GPVJEP	Phadebas paper	Negative
9Q3XCX	AMYLASE DIFFUSION	Positive	GV43GN	Amylase diffusion	Negative

TABLE 2 - Item 1

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
H7J77Z	a-Amylase	Negative	LXURT6	RSID	Negative
H9ADNB	RSID	Positive	M9CTN6	RSID	Negative
HEDJ2T	Phadebas test	Negative	MED8A7	ALS	Positive
HHBZEM	Amylase diffusion	Negative		RSID	Negative
HQ63HJ	RSID	Negative	MEWRED	RSID Saliva	Negative
HWRXDW	Radial Diffusion	Negative		RSID Saliva	Negative
J23ATA	Saliva test is not done if item is tested positive for blood or semen. No saliva test was done for Item 1 as it was tested positive for semen			RSID Saliva	Negative
				RSID Saliva	Negative
J94HFP	Phadebas	Negative	MY39JM	*saliva testing not performed in this lab	
JHVYAZ	Phadebas	Negative	MZHQHM	Current protocol does not allow for Saliva testing on swabs.	
JJM2UU	Phadebas tube test	Positive	N3ZCLU	Phadebas	Inconclusive
	Micro	Negative	N9PM3J	Amylase Diffusion	Negative
JNFWZ9	RSID-Saliva	Negative	N9PRB2	RSID	Negative
JPQ3UN	PHADEBAS PAPER TEST	Negative	NQ6YNH	RSID-Saliva	Negative
JRWLWQ	Current protocol does not allow for Saliva testing on swabs.		P8P4F9	RSID Saliva	Negative
JTMRDG	ALS	Negative	PCGF7Y	Amylase testing by Radial Diffusion	Negative
	RSID-Saliva	Negative	PDCWKG	SALigAE	Negative
JYWP6A	Saliva test is not done if item is tested positive for either blood or semen. Item 1 is tested positive for semen, thus no saliva test was performed.		PJ2EYR	Testing not offered.	
			PJ7RQU	RSID	Negative
K6XJEV	Phadebas	Positive	PKWW3E	Phadebas	Inconclusive
K7UYQR	Phadebas	Positive	PMLJF3	RSID	Positive
KQ3XGF	radial diffusion	Inconclusive	PNF7UR	Phadebas (amylase)	Negative
KUE3RF	RSID	Negative	QADXPB	RSID a-amylase	Negative
KXYWNG	Saligae	Negative	QP6VRD	Amylase Diffusion	Positive
			RBKG8Q	Phadebas	Negative
			RDTPGB	RSID	Negative
			RNACNQ	Amylase	Negative

TABLE 2 - Item 1

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
RVQ877	Amylase diffusion	Negative	Y28ZFH	This laboratory does not currently offer saliva testing	
T2V4RR	RSID	Negative	Y8GAQK	Phadebas Amylase	Positive
T7LPCY	No currently validated procedures for saliva test.		Y9VLM6	SaligaE	Negative
TJBFHW	RSID	Positive	YBDZQ6	Phadebas tube test not accredited	
TKLJQE	Amylase Diffusion	Negative	YMCGYU	Not tested for saliva since item tested positive for semen.	
U3J3HU	RSID Saliva	Positive	YN3XWA	RSID -ve @ 8½ mins	Negative
UAYML8	ALS	Positive	YQT4CZ	Radial Diffusion	Inconclusive
	RSID-Saliva	Negative	YZDAA4	SaligaE	Negative
UF9YH8	ALS	Positive	ZCTE3K	RSID	Negative
	RSID	Negative	ZKJZ9J	RSID	Negative
UHA7HB	Phadebas	Negative	ZRYNXE	PHADEBAS (TUBE METHOD)	Positive
UPT6GW	Phadebas amylase	Positive	ZYFE47	Amylase Diffusion	Negative
UPU6B6	AMYLASE PLATE	Inconclusive			
V4BH9W	RSID	Negative			
	ALS	Positive			
V7F8B6	RSID	Positive			
V9Z287	Phadebas	Negative			
VKV7HR	RSID	Negative			
VM2TTB	Amylase Diffusion	Negative			
W4LNZR	Amylase with radial diffusion	Negative			
WBNQ8D	Phadebas	Negative			
WF2Q7P	Rsid saliva	Negative			
WM2X2R	Phadebas	Negative			
WXGZ8Q	RSID Saliva	Positive			
XD83L9	Not competent in saliva testing				
XXUUMQ	RSID	Positive			

TABLE 2 - Item 1

Other Body Fluid Types					
Webcode	Body Fluid Type - Test Performed	Results	Webcode	Body Fluid Type - Test Performed	Results
2PQVGK	Sperm Search	Positive		Urine - Urease	Negative
3MXG82	Urine - Creatinine	Negative		Fecal Material - Urobilinogen	Negative
	Urine - DMAC	Positive	ZRYNXE	URINE, TEST PERFORMED CREATININE - FAINT	Positive
	Faeces - Edelmanns	Negative			
8C6JGB	Urine - Creatinine	Negative			
	Feces - Urobilinogen	Negative			
CGLJZ8	Urine - Creatinine	Negative			
	Urine - DMAC	Negative			
	Fecal material - Urobilinogen	Negative			
EFF3F6	Urine - Creatinine, DMAC	Negative			
	Fecal Material - Urobilinogen	Negative			
EXGQ74	Urine - Creatinine	Negative			
	Fecal Material - Urobilinogen	Negative			
G8XFEZ	Ouchterlony (Human origin)	Positive			
H7J7Z	Urine - Creatinine	Negative			
	Fecal Material - Urobilinogen	Negative			
JHVYAZ	Urine - Creatinine	Negative			
	Feces - Urobilinogen	Negative			
JJM2UU	Urine - DMAC	Negative			
JNFWZ9	Urine - Visual & DMAC	Negative			
	Feces - Visual & Microscopy	Negative			
PNF7UR	Human - Ouchterlony	Positive			
	Faeces - Price (urobilin)	Negative			
RNACNQ	Faeces - Prices	Negative			
WF2Q7P	Urine RSID	Negative			
WM2X2R	Urine - Creatinine	Negative			
Y8GAQK	Urine - Creatinine	Negative			

TABLE 2 - Item 2

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	KM	Negative	8YD3JW	No visible blood-like staining observed, therefore no presumptive test for blood performed	
28BLGK	LMG	Negative			
2EMTZU	KM	Negative			
2FHERZ	O-tol	Negative	9B677N	KM	Negative
2PQVGK	KM	Negative	9GQWTH	KM	Negative
2VEVGY	O-tol	Negative	9HLHJN	KM	Negative
3FF4BN	KM	Negative	9JH363	KM	Negative
3MXG82	KM	Negative	9Q3XCX	KM	Negative
3P4YCN	KM	Negative	9UE4NY	Visual exam - no staining observed	
4APZYB	KM	Negative	9W3CXZ	KM	Negative
4BJLQG	LMG	Negative	AWHZRY	KM	Negative
4HXUT4	LMG	Negative	B2DAF7	Hemastix - neg 60 secs	Negative
4JVXRK	KM	Negative		Hematrace - neg 10 mins	Negative
68YCXV	OT	Negative	B33FWW	Not tested: No red-brown staining observed on the Items.	
6HVCDK	Visual	Negative	BA8DNA	Visual	Negative
6PZ4XB	KM	Negative	BHXXRM	Not tested due to lack of reddish-brown color on swabs.	
6UWLGW	O-tol	Negative	BMUKKB	KM	Negative
6ZGMFL	Swab was not red stained		BNKQXU	Hematrace	Negative
7M6VAR	No visible stain observed		BU8H89	LMG	Negative
7MNLCK	Swabs have no red stains		BYM4XV	no apparent bloodstain observed	
7XNELB	KM	Negative	C4UYB4	LMG	Negative
8C6JGB	KM	Negative		HemaTrace	Negative
8NUMDA	LMG	Negative	CGLJZ8	KM	Negative
	HemaTrace	Negative	DL73AZ	KM	Negative
8PP7HW	LMG test	Negative	EB38QV	KM	Negative
8XLW8D	Visual	Negative	EB9GR9	TMB	Negative

TABLE 2 - Item 2

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
EFF3F6	KM	Negative	JHVYAZ	KM	Negative
ERDV3B	KM	Negative	JJM2UU	KM	Negative
EVPZDB	White scene or kit swabs are not tested for blood in this lab		JNFWZ9	K-M	Negative
EXGQ74	KM	Negative	JPQ3UN	LMG	Negative
F42C9M	Phenolphthalein	Negative	JRWLWQ	Visual Exam	Negative
FVMNEL	O-Tol	Negative	JTMRDG	KM	Negative
FXB6NT	KM	Negative	JYWP6A	KM	Negative
G74RZB	O-tol	Negative	K6XJEV	KM	Negative
G8XFEZ	TMB	Inconclusive	K7UYQR	Visual	Negative
	HemaTrace	Inconclusive		KM	Negative
	Haemochromogen	Negative	KQ3XGF	KM	Negative
G9NLVQ	LMG	Negative	KUE3RF	O-TOL	Negative
GFC6FU	KM	Negative	KXYWNG	HemaTrace	Negative
GPVJEP	Visual exam	Negative	KZ3R7Y	LMG	Negative
	KM	Negative	LXURT6	Visual exam	Negative
GV43GN	no staining on swab		M9CTN6	Visual	Negative
GZXZM8	*swabs were not red stained		MED8A7	O-Tol	Negative
H7J77Z	KM	Negative	MEWRED	Swabs were not tested for blood - does not display the characteristic reddish-brown stain associated with blood.	
H9ADNB	O-tol	Negative	MQAX88	KM	Negative
HDYUYF	KM	Negative	MZHQHM	Visual Exam	Negative
	Abacard Hematrace	Negative	N3ZCLU	KM	Negative
HEDJ2T	KM	Negative	N9PM3J	No red-brown discoloration	
	ABAcad HemaTrace	Negative	N9PRB2	KM	Negative
HHBZEM	A red stain was not observed.		P8P4F9	No reddish/brown stain observed during visual inspection.	
HQ63HJ	O-tol	Negative	PCGF7Y	KM	Negative
HWRXDW	KM	Negative			
J23ATA	KM	Negative			

TABLE 2 - Item 2

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
PDCWKG	KM	Negative	UAZKRV	Swabs were not stained	
	O-tol	Negative	UF9YH8	o-tol	Negative
	ABA Card	Negative	UHA7HB	KM	Negative
PJ2EYR	LMG	Negative	UPT6GW	Visual examination	Negative
	Hematrace	Negative	V4BH9W	KM	Negative
PJ7RQU	KM	Negative	V7F8B6	O-Tol	Negative
PKWW3E	KM	Negative	V9Z287	KM	Negative
	Visual	Negative	VKV7HR	tetrabase	Negative
PMLJF3	*NT because no red brown stain observed.		VM2TTB	No apparent discoloration to swabs	
PNF7UR	TMB	Negative	W4LNZR	KM	Negative
	HemaTrace	Negative	WBNQ8D	KM	Negative
QADXPB	Hemastix (TMB)	Negative	WF2Q7P	KM	Negative
QP6VRD	KM	Negative		HemaTrace	Negative
R9VBPE	LMG	Negative	WM2X2R	no visible dried red stain present	
RBKG8Q	KM	Negative	WXGZ8Q	K-M	Negative
RDTPGB	O-tol	Negative	XD83L9	KM	Negative
RNACNQ	TMB	Negative	XHZBMU	no red stain observed	
	HaemaTrace[sic]	Negative	XJC9UK	KM	Negative
	Haemachromogen[sic]	Negative	XXUUMQ	KM	Negative
RNQTHJ	Visually	Negative	Y28ZFH	LMG	Negative
RVQ877	No visual apparent blood staining			HemaTrace	Negative
T2V4RR	KM	Negative	Y8GAQK	KM	Negative
T7LPCY	KM	Negative	Y8WHUU	LMG	Negative
TJBFHW	KM	Negative	Y9VLM6	Phenolphthalein	Negative
TKLJQE	No stain consistent with blood		YBDZQ6	Visual	Negative
TX4EYX	KM	Negative		KM	Negative
U3J3HU	KM	Negative	YLWGGQ	Not brown stains	
UAYML8	KM Test	Negative	YMCGYU	KM	Negative

TABLE 2 - Item 2

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
YN3XWA	Hemastix -ve @ 60 secs	Negative			
YQT4CZ	No red brown staining				
YTUAPK	LMG	Negative			
YZDAA4	Phenolphthalein	Negative			
ZCTE3K	Visual	Negative			
ZKJZ9J	KM	Negative			
ZRYNXE	KM	Negative			
ZYFE47	No red discoloration				

TABLE 2 - Item 2

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	AP	Negative	6UWLGW	Micro	Negative
	MICRO	Negative		RSID	Negative
28BLGK	AP	Negative	6ZGMFL	Micro	Negative
2EMTZU	AP	Negative		P-30	Negative
	Micro	Negative	7M6VAR	AP	Negative
2FHERZ	RSID Semen	Negative	7MNLCK	Microscopic	Negative
	Micro	Negative		P30	Negative
2PQVGK	PSA	Negative	7XNELB	AP	Negative
2VEVGY	AP	Negative		P30	Negative
	Micro	Negative		KPIC	Negative
	RSID	Negative	8C6JGB	AP	Negative
3FF4BN	AP	Negative	8NUMDA	AP	Negative
	Micro	Negative		Micro	Negative
3MXG82	AP	Negative		PSA	Negative
	Micro	Negative		ALS	Positive
3P4YCN	AP	Negative	8XLW8D	AP	Negative
	Micro	Negative	8YD3JW	AP: x4 AP negative	Negative
4APZYB	AP	Negative	9B677N	AP	Negative
	Micro	Negative	9GQWTH	AP	Negative
4BJLQG	AP	Negative		Micro	Negative
4HXUT4	AP	Negative	9HLHJN	AP	Negative
4JVXRK	AP	Negative	9JH363	AP	Negative
	Micro	Negative		Micro	Negative
68YCXV	Micro	Negative	9Q3XCX	AP	Negative
	RSID		9UE4NY	AP	Negative
6HVCDK	AP	Negative		Micro	Negative
6PZ4XB	AP	Negative	9W3CXZ	AP	Negative
	Micro	Negative		Micro	Negative

TABLE 2 - Item 2

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
AWHZRY	AP	Negative	EB9GR9	AP	Negative
	PSA	Negative	EFF3F6	AP (Seri AP)	Negative
	Micro	Negative	ERDV3B	AP	Negative
B2DAF7	AP - neg 60 secs	Negative		PSA	Negative
	Micro - neg	Negative	EVPZDB	Micro (one (1) head observed)	Positive
	p30 - neg 10 mins	Negative		P30	Positive
B33FWW	AP	Negative	EXGQ74	AP (SAP)	Negative
BA8DNA	AP	Negative		PSA (ABA)	Negative
BHXXRM	Item 2 AnB: AP, Micro	Negative		Micro (Xmas)	Negative
	Item 2 CnD: AP, Micro	Negative	F42C9M	RSID	Negative
BMUKKB	AP	Negative	FVMNEL	Micro	Negative
	Seratec	Negative		RSID	Negative
	Micro	Negative	FXB6NT	AP	Inconclusive
BNKQXU	PSA	Negative		PSA	Negative
	Micro	Negative		Micro	Negative
BU8H89	AP	Negative	G74RZB	AP	Negative
BYM4XV	AP	Negative	G8XFEZ	AP	Negative
	Micro	Negative		PSA	Negative
C4UYB4	ALS	Positive		Microscopic	Negative
	AP	Negative	G9NLVQ	AP	Negative
	Micro	Negative	GFC6FU	AP	Negative
	P30	Negative		Micro	Negative
CGLJZ8	AP	Negative	GPVJEP	AP	Negative
	PSA	Negative		H&E Slide	Negative
DL73AZ	AP	Negative	GV43GN	AP	Negative
	Micro	Negative		Micro	Negative
EB38QV	PSA	Negative	GZXZM8	Micro	Negative
	MICRO	Negative		P30	Negative

TABLE 2 - Item 2

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
H7J77Z	AP (Map)	Negative	JRWLWQ	Fast Blue (Negative) (AP)	Negative
H9ADNB	ALS	Positive		Micro (Negative)	Negative
	AP	Negative		P30 (Negative) (PSA)	Negative
	Micro	Negative	JTMRDG	ALS	Positive
	PSA	Negative		AP	Negative
HDYUYF	AP	Negative	JYWP6A	RSID	Negative
	Micro	Negative		PSA	Negative
	Abacard p30	Negative	K6XJEV	AP	Negative
HEDJ2T	AP	Negative		Micro	Negative
	ABAcad P30	Negative	K7UYQR	AP	Negative
	Micro	Negative		Micro	Negative
HHBZEM	AP	Negative	KQ3XGF	AP	Negative
	PSA	Negative		Micro	Negative
	Micro	Negative		ABAcad p30	Negative
	ALS	Negative	KUE3RF	AP	Negative
HQ63HJ	Micro	Negative		MICRO	Negative
	RSID	Negative		RSID	Negative
HWRXDW	AP	Negative	KXYWNG	ABAcad P30	Negative
	Micro	Negative	KZ3R7Y	AP	Negative
J23ATA	RSID-Semen	Negative	LXURT6	AP	Negative
	PSA	Negative	M9CTN6	ALS	Negative
J94HFP	AP	Negative		AP	Negative
JHVYAZ	AP	Negative		Micro	Negative
JJM2UU	AP	Negative	MED8A7	ALS	Positive
	Micro	Negative		AP	Negative
JNFWZ9	AP	Negative		Micro	Negative
	Microscopic sperm search	Negative		RSID	Negative
JPQ3UN	AP	Negative			

TABLE 2 - Item 2

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
MEWRED	AP, Micro	Negative	PJ7RQU	AP	Negative
	AP, Micro	Negative		Seratec	Negative
	AP, Micro	Negative		Micro	Negative
	AP, Micro	Negative			
MQAX88	RSID	Negative	PKWW3E	AP	Negative
	PSA	Negative		Slide made (micro)	Negative
MY39JM	AP	Negative	PMLJF3	AP	Negative
MZHQHM	Fast Blue (AP)	Negative	PNF7UR	AP	Negative
	Micro	Negative		PSA	Positive
	p30 (PSA)	Negative		Micro	Negative
N3ZCLU	AP	Negative	QADXPB	AP negative @ 60 secs	Negative
	p30	Negative		Microscopy	Negative
	KPIC	Negative		P30	Negative
N9PM3J	AP	Negative	QP6VRD	BCIP	Negative
	Micro	Negative	R9VBPE	AP	Negative
N9PRB2	AP	Negative	RBKG8Q	AP	Negative
	PSA	Negative		P-30	Negative
NQ6YNH	AP	Negative	RDTPGB	Micro	Negative
	RSID-Semen	Negative		RSID	Negative
P8P4F9	AP	Negative	RNACNQ	AP	Negative
	Microscopy	Negative		PSA	Negative
PCGF7Y	AP	Negative	RNQTHJ	Micro	Negative
	Micro	Negative		AP	Negative
PDCWKG	AP	Negative	RVQ877	AP	Negative
PJ2EYR	ALS	Positive	T2V4RR	Micro	Negative
	AP	Negative		AP	Negative
	micro	Negative		Seratec	Negative
	P30	Negative			

TABLE 2 - Item 2

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
T7LPCY	Microscopic	Negative	VKV7HR	Micro	Negative
	P30	Negative		PSA	Negative
TJBFHW	AP	Negative		RSID	Negative
TKLJQE	AP	Inconclusive	VM2TTB	AP	Negative
	Micro	Negative		Micro	Negative
	PSA	Negative	W4LNZR	AP	Negative
TX4EYX	Micro	Negative		Micro exam	Negative
	P30	Negative	WBNQ8D	AP	Negative
U3J3HU	AP	Negative		Micro	Negative
UAYML8	ALS	Positive	WF2Q7P	AP	Negative
	AP Test	Negative		Micro	Negative
UAZKRV	Micro	Negative		Rsid semen	Negative
	P30	Negative	WM2X2R	Acid Phosphatase (AP)	Negative
UF9YH8	ALS	Negative		P30/PSA	Negative
	AP	Negative		Micro	Negative
	Micro	Negative	WXGZ8Q	AP	Negative
	RSID-semenogelin	Negative	XD83L9	AP	Negative
UHA7HB	AP	Negative	XHZBMU	Micro	Negative
UPT6GW	AP	Negative		P30	Negative
UPU6B6	AP	Negative	XJC9UK	AP	Negative
	MICRO	Negative	XXUUMQ	AP	Negative
V4BH9W	ALS	Positive	Y28ZFH	AP	Negative
	AP	Negative		Micro	Negative
	Micro	Negative		PSA	Negative
V7F8B6	AP	Negative		ALS	Positive
	RSID	Negative	Y8GAQK	AP	Negative
	Micro	Negative	Y8WHUU	AP	Negative
V9Z287	AP	Negative	Y9VLM6	RSID	Negative

TABLE 2 - Item 2

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
YBDZQ6	AP	Negative			
	Microscopic exam - no spermatazoa[sic]	Negative			
YLWGGQ	Microscopic	Negative			
	P30	Negative			
YMCGYU	RSID-Semen	Negative			
	PSA	Negative			
YN3XWA	Micro	Negative			
	AP -ve @ 60 secs	Negative			
YQT4CZ	AP	Negative			
	Micro	Negative			
YTUAPK	AP	Negative			
YZDAA4	RSID	Negative			
ZCTE3K	AP	Negative			
ZKJZ9J	AP	Negative			
	Seratec	Negative			
	Micro	Negative			
ZRYNXE	AP	Negative			
	MICRO	Negative			
ZYFE47	AP	Negative			
	Micro	Negative			

TABLE 2 - Item 2

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	PHADEBAS	Positive	9UE4NY	Amylase Diffusion	Positive
2EMTZU	RSID	Positive	9W3CXZ	Phadebas press	Positive
2FHERZ	RSID Saliva	Positive	AWHZRY	Phadebas	Positive
2PQVGK	RSID	Positive	B2DAF7	RSID - pos 1 min	Positive
2VEVGY	RSID	Positive	B33FWW	RSID	Positive
3FF4BN	Phadebas	Positive	BA8DNA	RSID	Positive
3MXG82	Phadebas	Positive	BHXXRM	Item 2 AnB: RSID	Positive
3P4YCN	RSID-Saliva	Positive		Item 2 CnD: RSID	Positive
4APZYB	Amylase (radial diffusion)	Positive	BMUKKB	RSID	Positive
4JVXRK	RSID	Positive	BNKQXU	RSID	Positive
68YCXV	RSID	Positive	BYM4XV	Amylase diffusion	Positive
6HVCDK	RSID-Saliva	Positive	C4UYB4	Saliva testing is not offered as a service in this laboratory.	
6PZ4XB	Phadebas	Positive	CGLJZ8	Phadebas Amylase	Positive
6UWLGW	RSID	Positive	DL73AZ	Phadebas Press	Positive
7M6VAR	Amylase	Positive	EB38QV	PHADEBAS	Positive
7MNLCK	Per Agency's Policy		EB9GR9	RSID Saliva Field Kit	Positive
7XNELB	Phadebas	Positive	EFF3F6	Amylase (Phadebas)	Positive
8C6JGB	Phadebas	Positive	ERDV3B	RSID	Positive
8NUMDA	This laboratory does not perform testing for saliva.		EVPZDB	Saliva testing not being done in this lab at this time	
8PP7HW	Phadebas tube test	Positive	EXGQ74	Phadebas Amylase	Positive
8XLW8D	RSID-Saliva	Positive	F42C9M	SaligaE	Positive
8YD3JW	Phadebas: x4 Phadebas positive	Positive	FVMNEL	RSID	Positive
9B677N	RSID	Positive	FXB6NT	Phadebas	Positive
9GQWTH	RSID	Positive	G74RZB	RSID	Positive
9HLHJN	RSID	Positive	G9NLVQ	PHADEBAS PAPER TEST	Positive
9JH363	Phadebas Press	Positive	GFC6FU	Phadebas	Positive
9Q3XCX	AMYLASE DIFFUSION	Positive			

TABLE 2 - Item 2

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
GPVJEP	Phadebas Paper	Positive	M9CTN6	RSID	Positive
	H&E slide	Positive	MED8A7	ALS	Positive
GV43GN	Amylase diffusion	Positive		RSID	Positive
H7J77Z	α-Amylase	Positive	MEWRED	RSID Saliva	Positive
H9ADNB	RSID	Positive		RSID Saliva	Positive
HEDJ2T	Phadebas test	Positive		RSID Saliva	Positive
HHBZEM	Amylase diffusion	Positive		RSID Saliva	Positive
HQ63HJ	RSID	Positive	MQAX88	Starch-iodine	Positive
HWRXDW	Radial Diffusion	Positive	MY39JM	*saliva testing not performed in this lab	
J23ATA	Starch-Iodine Assay	Positive	MZHQHM	Current protocol does not allow for Saliva testing on swabs.	
J94HFP	Phadebas	Positive	N3ZCLU	Phadebas	Positive
JHVYAZ	Phadebas	Positive	N9PM3J	Amylase Diffusion	Positive
JJM2UU	Phadebas tube test	Positive	N9PRB2	RSID	Positive
	Micro	Positive	NQ6YNH	RSID-Saliva	Positive
JNFWZ9	RSID-Saliva	Positive	P8P4F9	RSID Saliva	Positive
JPQ3UN	PHADEBAS PAPER TEST	Positive	PCGF7Y	Amylase Testing by Radial Diffusion	Positive
JRWLWQ	Current protocol does not allow for Saliva testing on swabs.		PDCWKG	SALlgAE	Positive
JTMRDG	ALS	Positive	PJ2EYR	Testing not offered.	
	RSID-Saliva	Positive	PJ7RQU	RSID	Positive
JYWP6A	Starch-iodine assay	Positive	PKWW3E	Phadebas	Positive
K6XJEV	Phadebas	Positive	PMLJF3	RSID	Positive
K7UYQR	Phadebas	Positive	PNF7UR	Phadebas (amylase)	Positive
KQ3XGF	radial diffusion	Positive	QADXPB	RSID (+)ve @ 36 secs	Positive
	Micro (for nucleated epithelial cells)	Positive	QP6VRD	Amylase Diffusion	Positive
KUE3RF	RSID	Positive	RBKG8Q	Phadebas	Positive
KXYWNG	Saligae	Positive	RDTPGB	RSID	Positive
LXURT6	RSID	Positive			

TABLE 2 - Item 2

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
RNACNQ	Amylase	Positive	XJC9UK	Amylase	Positive
	Micro	Positive	XXUUMQ	RSID	Positive
RVQ877	Amylase diffusion	Positive	Y28ZFH	This laboratory does not currently offer saliva testing	
T2V4RR	RSID	Positive	Y8GAQK	Phadebas Amylase	Positive
T7LPCY	No currently validated procedures for saliva test.		Y9VLM6	SaligaE	Positive
TJBFHW	RSID	Positive	YBDZQ6	Phadebas tube test not accredited	
TKLJQE	Amylase Diffusion	Positive	YMCGYU	Starch-iodine test	Positive
U3J3HU	RSID	Positive	YN3XWA	RSID (+)ve @ 1 min	Positive
UAYML8	ALS	Positive	YQT4CZ	Radial Diffusion	Positive
	RSID-Saliva	Positive	YZDAA4	SaligaE	Positive
UF9YH8	ALS	Negative	ZCTE3K	RSID	Positive
	RSID	Positive	ZKJZ9J	RSID	Positive
UHA7HB	Phadebas	Positive	ZRYNXX	PHADEBAS (TUBE METHOD)	Positive
UPT6GW	Phadebas amylase	Positive	ZYFE47	Amylase Diffusion	Positive
UPU6B6	AMYLASE PLATE	Positive			
V4BH9W	ALS	Positive			
	RSID	Positive			
V7F8B6	RSID	Positive			
V9Z287	PHADEBAS	Positive			
VKV7HR	RSID	Positive			
VM2TTB	Amylase Diffusion	Positive			
W4LNZR	Amylase with radial diffusion	Positive			
WBNQ8D	Phadebas	Positive			
WF2Q7P	Rsid saliva	Positive			
WM2X2R	Phadebas (amylase)	Positive			
WXGZ8Q	RSID Saliva	Positive			
XD83L9	Not competent in saliva testing				

TABLE 2 - Item 2

Other Body Fluid Types					
Webcode	Body Fluid Type - Test Performed	Results	Webcode	Body Fluid Type - Test Performed	Results
3MXG82	Urine - Creatinine	Negative		Fecal Material - Urobilinogen	Negative
	Urine - DMAC	Negative			
	Faeces - Edelmans	Negative			
8C6JGB	Urine - Creatinine	Negative			
	Feces - Urobilinogen	Negative			
CGLJZ8	Urine - Creatinine	Negative			
	Urine - DMAC	Negative			
	Fecal Material - Urobilinogen	Negative			
EFF3F6	Urine - Creatinine, DMAC	Negative			
	Fecal Material - Urobilinogen	Negative			
EXGQ74	Urine - Creatinine	Negative			
	Fecal Material - Urobilinogen	Negative			
G8XFEZ	Ouchterlony (Human origin)	Negative			
H7J77Z	Urine - Creatinine	Negative			
	Fecal Material - Urobilinogen	Negative			
JHVYAZ	Urine - Creatinine	Negative			
	Feces - Urobilinogen	Negative			
JJM2UU	Urine - DMAC	Negative			
JNFWZ9	Urine - Visual & DMAC	Negative			
	Feces - Visual & Microscopy	Negative			
PNF7UR	Human - Ouchterlony	Negative			
	Faeces - Price (urobilin)	Negative			
RNACNQ	Faeces - Prices	Negative			
RNQTHJ	DNA-analysis	Positive			
WF2Q7P	Rsid Urine	Negative			
WM2X2R	Urine - Creatinine	Negative			
Y8GAQK	Urine - Creatinine	Negative			

TABLE 2 - Item 3

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	KM	Positive		Hema Trace	Positive
28BLGK	LMG	Positive	6ZGMFL	Kastle-Meyer	Positive
	HemaTrace	Positive	7M6VAR	LMG	Positive
	Ouchterlony	Positive	7MNLCK	Kastle-Meyer Test	Positive
2EMTZU	KM	Positive	7XNELB	KM	Positive
2FHERZ	O-tol	Positive	8C6JGB	KM	Positive
	ABACARD Hematrace	Positive		Modified Takayama	Positive
2PQVGK	KM	Positive	8NUMDA	LMG	Positive
2VEVGY	O-tol	Positive		HemaTrace	Positive
	HemaTrace	Positive	8PP7HW	LMG test	Positive
3FF4BN	KM	Positive	8XLW8D	LMG	Positive
	Hematrace	Positive		HemaTrace	Positive
3MXG82	KM	Positive		Visual	Positive
3P4YCN	KM	Positive	8YD3JW	LMG: Heavy blood-like staining observed, LMG presumptive test performed, x4 LMG positive	Positive
4APZYB	KM	Positive	9B677N	KM	Positive
	HemaTrace	Positive	9GQWTH	KM	Positive
4BJLQG	LMG	Positive	9HLHJN	KM	Positive
	HemaTrace	Positive	9JH363	KM	Positive
	Ouchterlony	Positive		ABACard HemaTrace	Positive
4HXUT4	LMG	Positive	9Q3XCX	KM	Positive
4JVXRK	KM	Positive		ABACARD HEMATrace	Positive
68YCXV	OT		9UE4NY	LMG	Positive
	HemaTrace	Positive	9W3CXZ	KM	Positive
6HVCDK	Visual	Positive		ABAcad Hematrace	Positive
	KM	Positive	AWHZRY	KM	Positive
	HemaTrace	Positive	B2DAF7	Hemastix - pos 5 secs	Positive
6PZ4XB	KM	Positive		Hematrace - pos 2 mins	Positive
6UWLGW	O-tol	Positive			

TABLE 2 - Item 3

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
B33FWW	KM	Positive		Mod Taka	Positive
	RSID	Positive	F42C9M	Phenolphthalein	Negative
BA8DNA	LMG	Positive	FVMNEL	O-Tol	Positive
	HemaTrace	Positive		Hematrace	Positive
	Visual	Positive	FXB6NT	KM	Positive
BHXXRM	Item 3 AnB: O-tol	Positive	G74RZB	O-tol	Positive
	Item 3 CnD: O-tol	Positive		HemaTrace	Positive
BMUKKB	KM	Positive	G8XFEZ	TMB	Positive
	Hematrace	Positive		HemaTrace	Positive
BNKQXU	Hematrace	Positive		Haemochromogen	Positive
BU8H89	LMG	Positive	G9NLVQ	LMG	Positive
	HemaTrace	Positive	GFC6FU	KM	Positive
	Ouchterlony Species Determination	Positive		ABA	Positive
BYM4XV	LMG	Positive	GPVJEP	Visual exam	Positive
	HemaTrace	Positive		KM	Positive
C4UYB4	LMG	Positive	GV43GN	LMG	Positive
	HemaTrace	Positive	GZXZM8	KM	Positive
CGLJZ8	KM	Positive	H7J77Z	KM	Positive
	Modified Takayama	Positive		Takayama	Positive
DL73AZ	KM	Positive	H9ADNB	O-tol	Positive
	ABACard HemaTrace	Positive		HemaTrace	Positive
EB38QV	KM	Positive	HDYUYF	KM	Positive
EB9GR9	TMB	Positive		Abacard Hematrace	Positive
EFF3F6	KM	Positive	HEDJ2T	KM	Positive
	Takayama	Positive		ABACard HemaTrace	Positive
ERDV3B	KM	Positive	HHBZEM	LMG	Positive
EVPZDB	KM	Positive		Hematrace	Positive
EXGQ74	KM	Positive			

TABLE 2 - Item 3

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
HQ63HJ	O-tol	Positive	LXURT6	Visual exam	Positive
	Hematrace	Positive		KM	Positive
HWRXDW	KM	Positive		HemaTrace	Positive
	HemaTrace	Positive	M9CTN6	Visual	Positive
J23ATA	KM	Positive		KM	Positive
	HO	Positive		HemaTrace	Positive
J94HFP	TMB	Positive	MED8A7	O-Tol	Positive
JHVYAZ	KM	Positive		Hematrace	Positive
	Modified Takayama	Positive	MEWRED	O-tol	Positive
JJM2UU	KM	Positive		O-tol	Positive
JNFWZ9	K-M	Positive		O-tol	Positive
JPQ3UN	LMG	Positive		O-tol	Positive
JRWLWQ	Hemochromogen (Positive)	Positive	MQAX88	KM	Positive
JTMRDG	KM	Positive		Hexagon OBTI	Positive
	RSID-Blood	Positive	MY39JM	KM	Positive
JYWP6A	KM	Positive	MZHQHM	Hemochromogen	Positive
	Hexagon OBTI	Positive	N3ZCLU	KM	Positive
K6XJEV	KM	Positive	N9PM3J	LMG	Positive
K7UYQR	Visual	Positive		HemaTrace	Positive
	KM	Positive	N9PRB2	KM	Positive
KQ3XGF	KM	Positive	NQ6YNH	KM	Positive
	Takayama	Positive		RSID-Blood	Positive
KUE3RF	O-TOL	Positive	P8P4F9	O-tol	Positive
	HEMATRACE	Positive	PCGF7Y	KM	Positive
KXYWNG	HemaTrace	Positive		HemaTrace	Positive
KZ3R7Y	LMG	Positive	PDCWKG	KM	Positive
	HemaTrace	Positive		O-tol	Positive
	Ouchterlony	Positive		ABA Card	Positive

TABLE 2 - Item 3

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
PJ2EYR	LMG	Positive	T7LPCY	KM	Positive
	HemaTrace	Positive	TJBFHW	KM	Positive
PJ7RQU	KM	Positive	TKLJQE	LMG	Positive
	Hematrace	Positive		HemaTrace	Positive
PKWW3E	KM	Positive	TX4EYX	KM	Positive
	Visual	Positive	U3J3HU	KM	Positive
PMLJF3	O-tol	Positive	UAYML8	KM Test	Positive
	HemaTrace	Positive		RSID-Blood	Positive
PNF7UR	TMB	Positive	UAZKRV	KM	Positive
	HemaTrace	Positive	UF9YH8	ALS	Negative
	Haemochromogen	Positive		O-tol	Positive
QADXPB	Hemastix (TMB)	Positive		Hematrace	Positive
	HemaTrace	Positive	UHA7HB	KM	Positive
QP6VRD	KM	Positive	UPT6GW	KM	Positive
	Hematrace	Positive	UPU6B6	KM	
R9VBPE	HemaTrace	Positive	V4BH9W	KM	Positive
	LMG	Positive	V7F8B6	O-Tol	Positive
RBKG8Q	KM	Positive		RSID	Positive
RDTPGB	O-tol	Positive	V9Z287	KM	Positive
	HemaTrace	Positive	VKV7HR	Tetrabase	Positive
RNACNQ	TMB	Positive		RSID	Positive
	HemaTrace	Positive	VM2TTB	LMG	Positive
RNQTHJ	LMG	Positive		Hematrace	Positive
	OBTI	Positive	W4LNZR	KM	Positive
RVQ877	LMG	Positive		HemaTrace	Positive
	HemaTrace	Positive	WBNQ8D	KM	Positive
T2V4RR	KM	Positive		HemaTrace	Positive
	Hematrace	Positive			

TABLE 2 - Item 3

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
WF2Q7P	KM	Positive	YZDAA4	Phenolphthalein	Negative
	Hematrace	Positive	ZCTE3K	LMG	Positive
WM2X2R	PTMB	Positive		HemaTrace	Positive
	HemaTrace	Positive		Visual	Positive
WXGZ8Q	K-M	Positive	ZKJZ9J	KM	Positive
XD83L9	KM	Positive		Hematrace	Positive
XHZBMU	KM	Positive	ZRYNXE	KM	Positive
XJC9UK	KM	Positive	ZYFE47	LMG	Positive
XXUUMQ	KM	Positive		Hematrace	Positive
Y28ZFH	LMG	Positive			
	HemaTrace	Positive			
Y8GAQK	KM	Positive			
	Modified Takayama	Positive			
Y8WHUU	LMG	Positive			
	HemaTrace	Positive			
	Ouchterlony - Human	Positive			
Y9VLM6	Phenolphthalein	Negative			
YBDZQ6	Visual	Positive			
	KM	Positive			
YLWGGQ	Kastle-Meyer Test	Positive			
YMCGYU	KM	Positive			
	Hexagon OBTI	Positive			
YN3XWA	Hemastix (+)ve @ 5 secs	Positive			
	Hematrace (+)ve @ 20 secs	Positive			
YQT4CZ	KM	Positive			
YTUAPK	LMG	Positive			
	HemaTrace	Positive			
	Ouchterlony (human)	Positive			

TABLE 2 - Item 3

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	AP	Negative	6ZGMFL	Micro	Negative
	MICRO	Negative		P30	Negative
28BLGK	AP	Negative	7M6VAR	AP, Micro, P30	Negative
2EMTZU	AP	Negative	7MNLCK	Microscopic	Negative
2FHERZ	micro	Negative		P30	Negative
	RSID Semen	Negative	7XNELB	AP	Negative
2PQVGK	PSA	Negative		P30	Negative
2VEVGY	AP	Negative		KPIC	Negative
	Micro	Negative	8C6JGB	AP	Negative
	RSID	Negative	8NUMDA	Testing was not performed due to case scenario.	
3FF4BN	AP	Negative	8XLW8D	AP	Negative
	Micro	Negative			
3MXG82	AP	Negative	8YD3JW	As visible blood-like staining that gave an LMG positive result was observed, no presumptive test for semen performed	
	Micro	Negative			
3P4YCN	The swabs were heavily stained reddish-brown and there is no reason to believe semen would also be present on these swabs.		9B677N	AP	Negative
			9GQWTH	Not applicable to case scenario	
4APZYB	AP	Negative	9HLHJN	AP	Negative
	Micro	Negative	9JH363	AP	Negative
4BJLQG	AP	Negative		Micro	Negative
4HXUT4	AP	Negative	9Q3XCX	NOT TESTED DUE TO SCENARIO, LOCATION OF STAIN, AND POSITIVE KM TEST.	
4JVXRK	AP	Negative			
68YCXV	Micro	Negative	9UE4NY	AP	Negative
	RSID			Micro	Negative
6HVCDK	AP	Negative		P30	Negative
6PZ4XB	AP	Negative			
6UWLGW	Micro	Negative	9W3CXZ	AP	Negative
	RSID	Negative		Micro	Negative

TABLE 2 - Item 3

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
AWHZRY	AP	Inconclusive	EFF3F6	AP (Seri AP)	Negative
	P30	Negative	ERDV3B	AP	Negative
	Micro	Negative	EVPZDB	Brown or red stained scene swabs are not tested for semen in this lab	
B2DAF7	AP - neg 60 secs	Negative	EXGQ74	AP (SAP)	Negative
	Micro - neg	Negative		PSA (ABA)	Negative
B33FWW	AP	Negative		Micro (Xmas)	Negative
BA8DNA	AP	Negative	F42C9M	RSID	Negative
	P30	Negative	FVMNEL	Micro	Negative
	Spermatozoa (Micro)	Negative		RSID	Negative
BHXXRM	Item 3 AnB: AP, Micro	Negative	FXB6NT	AP	Inconclusive
	Item 3 CnD: AP, Micro	Negative		PSA	Negative
BMUKKB	AP	Negative		Micro	Negative
	Seratec	Negative	G74RZB	No information to indicate that the human bloodstain on the wall would contain semen and no specific request for this testing	
	Micro	Negative	G8XFEZ	AP	Negative
BNKQXU	PSA	Negative	G9NLVQ	AP	Negative
	MICRO	Negative	GFC6FU	AP	Negative
BU8H89	AP	Negative		Micro	Negative
BYM4XV	AP	Negative	GPVJEP	AP	Negative
	Micro	Negative		H&E Slide	Negative
C4UYB4	ALS	Negative	GV43GN	AP	Negative
	Micro	Negative		Micro	Negative
	P30	Negative	GZXZM8	Micro	Negative
CGLJZ8	Not tested based on case scenario			P30	Negative
DL73AZ	AP	Negative	H7J77Z	Analyst's discretion (based on scenario provided).	
	Micro	Negative			
EB38QV	PSA	Negative			
	MICRO	Negative			
EB9GR9	AP	Negative			

TABLE 2 - Item 3

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
H9ADNB	ALS	Negative		PSA	Negative
	AP	Negative	K6XJEV	AP	Negative
HDYUYF	AP	Negative		Micro	Negative
	Micro	Negative	K7UYQR	AP	Negative
	Abacard p30	Negative		Micro	Negative
HEDJ2T	AP	Negative	KQ3XGF	AP	Negative
	Micro	Negative		Micro	Negative
HHBZEM	AP	Negative		ABAcad p30	Negative
	PSA	Negative	KUE3RF	AP	Negative
	Micro	Negative		MICRO	Negative
	ALS	Negative		RSID	Negative
HQ63HJ	Micro	Negative	KXYWNG	ABAcad P30	Negative
	RSID	Negative	KZ3R7Y	AP	Negative
HWRXDW	AP	Negative		PSA	Negative
	Micro	Negative		Microscopic	Negative
J23ATA	RSID-Semen	Negative	LXURT6	AP	Negative
	PSA	Negative	M9CTN6	AP	Negative
J94HFP	AP	Negative		Micro	Negative
JHVYAZ	AP	Negative	MED8A7	ALS	Negative
JJM2UU	AP	Negative		Micro	Negative
	Micro	Negative		RSID	Negative
JNFWZ9	No scene information was provided that made it seem likely this stain would contain more than one body fluid		MEWRED	AP, Micro	Negative
				AP, Micro	Negative
				AP, Micro	Negative
JPQ3UN	AP	Negative		AP, Micro	Negative
JRWLWQ	Fast Blue (-ve) (AP)	Negative	MQAX88	RSID	Negative
JTMRDG	ALS	Negative		PSA	Negative
	AP	Negative	MY39JM	AP	Negative
JYWP6A	RSID	Negative			

TABLE 2 - Item 3

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
MZHQHM	Fast Blue (AP)	Negative	R9VBPE	AP	Negative
N3ZCLU	AP	Inconclusive	RBKG8Q	AP	Negative
	p30	Negative		P-30	Negative
	KPIC	Negative		KPIC	Negative
N9PM3J	AP	Negative	RDTPGB	Micro	Negative
	Micro	Negative		RSID	Negative
N9PRB2	AP	Negative	RNACNQ	AP	Negative
NQ6YNH	AP	Negative	RNQTHJ	AP	Negative
	RSID-Semen	Negative		Micro	Negative
P8P4F9	AP	Negative	RVQ877	AP	Negative
	Microscopy	Negative		Micro	Negative
PCGF7Y	Scenario didn't indicate the need for this testing.		T2V4RR	AP	Negative
PDCWKG	AP	Negative		Seratec	Negative
PJ2EYR	AP	Negative	T7LPCY	Microscopic	Negative
	ALS	Positive		P30	Negative
PJ7RQU	AP	Negative	TJBFHW	AP	Negative
	Seratec	Negative	TKLJQE	AP	Negative
	Micro	Negative		Micro	Negative
PKWW3E	Area recovered - wasn't described as sex swab		TX4EYX	Micro	Negative
				P30	Negative
PMLJF3	Item 3 was not tested for semen because of the location of the stain.		U3J3HU	AP	Negative
PNF7UR	AP	Negative	UAYML8	ALS	Negative
	PSA	Negative		AP Test	Negative
	Micro	Negative	UAZKRV	Micro	Negative
QADXPB	AP Test	Negative		P-30	Negative
	Microscopy	Negative	UF9YH8	ALS	Negative
	P30	Negative		Micro	Negative
QP6VRD	BCIP	Negative		RSID-semenogelin	Negative

TABLE 2 - Item 3

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
UHA7HB	AP	Negative	Y28ZFH	Testing not performed due to the probative value of the sample based on the scenario given.	
UPT6GW	AP	Negative			
V4BH9W	ALS	Inconclusive			
	AP	Negative	Y8GAQK	AP	Negative
V7F8B6	AP	Negative	Y8WHUU	AP	Negative
	RSID	Negative	Y9VLM6	RSID	Negative
	Micro	Negative	YBDZQ6	AP	Negative
V9Z287	AP	Negative		Microscopic exam - no spermatazoa[sic]	Negative
VKV7HR	Micro	Negative	YLGWGGQ	Microscopic	Negative
	PSA	Negative		P30	Negative
	RSID	Negative	YMCGYU	RSID-Semen	Negative
VM2TTB	AP	Negative		PSA	Negative
	Micro	Negative	YN3XWA	Micro	Negative
W4LNZR	AP	Negative		AP -ve @ 60 secs	Negative
	Micro exam	Negative	YQT4CZ	No senario[sic] to indicate other bodily fluids present	
WBNQ8D	AP	Negative	YTUAPK	AP	Negative
	Micro	Negative	YZDAA4	RSID	Negative
WF2Q7P	AP	Negative	ZCTE3K	AP	Negative
	Micro	Negative		P30	Negative
	Rsid Semen	Negative		Microscopic	Negative
WM2X2R	Acid Phosphatase (AP)	Negative	ZKJZ9J	AP	Negative
WXGZ8Q	AP	Negative		Seratec	Negative
XD83L9	AP	Negative		Micro	Negative
XHZBMU	Micro	Negative	ZRYNXXE	AP	Negative
	P30	Negative		MICRO	Negative
XJC9UK	AP	Negative		CHOLINE	Negative
XXUUMQ	AP	Negative	ZYFE47	AP	Negative
				Micro	Negative

TABLE 2 - Item 3

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	PHADEBAS	Inconclusive	9B677N	Not tested per policy.	
2EMTZU	Not required per policy		9GQWTH	Not applicable to case scenario	
2FHERZ	RSID Saliva	Negative	9HLHJN	Per laboratory policy, only intimate samples are saliva tested.	
2PQVGK	RSID	Negative	9JH363	Phadebas Press	Negative
2VEVGY	RSID	Negative	9Q3XCX	NOT TESTED DUE TO SCENARIO, LOCATION OF STAIN, AND POSITIVE KM TEST.	
3FF4BN	Phadebas	Negative	9UE4NY	Amylase Diffusion	Negative
3MXG82	Phadebas tube test	Inconclusive	9W3CXZ	Phadebas press	Negative
3P4YCN	The swabs were heavily stained reddish-brown and there is no reason to believe saliva would also be present on these swabs.		AWHZRY	Phadebas	Negative
4APZYB	Amylase (radial diffusion)	Negative	B2DAF7	No indication in scenario of possibility of saliva	
4JVXRK	Not tested per policy *Not Required		BA8DNA	RSID	Negative
68YCXV	RSID	Negative	BHXXRM	Not tested due to swab type.	
6HVCDK	RSID-Saliva	Negative	BMUKKB	RSID	Negative
6PZ4XB	Case scenario did not require testing for saliva		BNKQXU	RSID	Negative
6UWLGW	RSID	Negative	BYM4XV	Amylase diffusion	Negative
7M6VAR	Amylase	Negative	C4UYB4	Saliva testing is not offered as a service in this laboratory.	
7MNLCK	Per Agency's Policy		CGLJZ8	Not tested based on case scenario	
7XNELB	Phadebas	Negative	DL73AZ	Phadebas Press	Negative
8C6JGB	Phadebas	Negative	EB38QV	PHADEBAS	Negative
8NUMDA	This laboratory does not perform testing for saliva.		EB9GR9	RSID Saliva Field Kit	Negative
8PP7HW	Phadebas tube test	Inconclusive	EFF3F6	Amylase (Phadebas)	Negative
8XLW8D	RSID-Saliva	Negative	ERDV3B	Not routinely performed on KM positive stains	
8YD3JW	As visible blood-like staining that gave an LMG positive result was observed, no presumptive test for saliva performed		EVPZDB	Saliva testing not being done in this lab at this time	

TABLE 2 - Item 3

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
EXGQ74	Phadebas Amylase	Negative	K6XJEV	Phadebas	Inconclusive
F42C9M	SaligaE	Negative	K7UYQR	Phadebas	Inconclusive
FVMNEL	RSID	Negative	KQ3XGF	radial diffusion	Inconclusive
FXB6NT	Phadebas	Negative	KUE3RF	RSID	Negative
G74RZB	No information to indicate oral contact and no specific request for this testing		KXYWNG	Saligae	Negative
GFC6FU	Phadebas	Negative	LXURT6	RSID	Negative
GPVJEP	Phadebas Paper	Negative	M9CTN6	RSID	Negative
GV43GN	Amylase diffusion	Negative	MED8A7	ALS	Negative
H7J77Z	Analyst's discretion (based on scenario provided).			RSID	Negative
H9ADNB	The swabs were very dark and there was no evidence of dilution		MEWRED	Did not test for saliva based on case scenario.	
HEDJ2T	Phadebas test	Negative	MY39JM	*saliva testing not performed in this lab	
HHBZEM	Amylase diffusion	Negative	MZHQHM	Current protocols do not allow for Saliva testing on swabs.	
HQ63HJ	RSID	Negative	N3ZCLU	Phadebas	Negative
HWRXDW	Radial Diffusion	Negative	N9PM3J	Amylase Diffusion	Negative
J94HFP	Phadebas	Negative	N9PRB2	KM positive stains are not routinely tested for saliva	
JHVYAZ	Phadebas	Negative	P8P4F9	RSID Saliva	Negative
JJM2UU	Phadebas tube test	Negative	PCGF7Y	Scenario didn't indicate the need for this testing.	
	Micro	Negative	PDCWKG	SALigAE	Negative
JNFWZ9	No scene information was provided that made it seem likely this stain would contain more than one body fluid		PJ2EYR	Testing not offered	
JRWLWQ	Current protocols do not allow for Saliva testing on swabs.		PJ7RQU	RSID	Negative
JTMRDG	ALS	Negative	PKWW3E	Area recovered from - was deemed unnecessary (& not dilute looking)	
	RSID-Saliva	Negative	PMLJF3	Item 3 was not tested for saliva because of the location of the stain.	
JYWP6A			PNF7UR	Phadebas (amylase)	Negative
			QADXPB	RSID a-amylase	Negative

TABLE 2 - Item 3

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
QP6VRD	Amylase Diffusion	Positive	WXGZ8Q	Per policy, item does not meet criteria to be tested for this body fluid - therefore not required to test.	
R9VBPE	Scenario specific - case details does not declare biting nor oral contact				
RBKG8Q	Phadebas	Negative	XD83L9	Not competent in saliva testing	
RDTPGB	RSID	Negative	XXUUMQ	Saliva test not performed on this item per laboratory policy - not required to test	
RNACNQ	Amylase	Negative			
RVQ877	Amylase diffusion	Negative	Y28ZFH	This laboratory does not offer saliva testing.	
T2V4RR	RSID	Negative	Y8GAQK	Phadebas Amylase	Negative
T7LPCY	No currently validated procedures for saliva test.		Y9VLM6	SaligaE	Negative
TJBFHW	NT (Per Policy) - Item not eligible for saliva testing		YBDZQ6	Phadebas tube test not accredited	
TKLJQE	Amylase Diffusion	Negative	YN3XWA	RSID -ve @ 8½ mins	Negative
U3J3HU	Per Policy; not required to test		YQT4CZ	No scenario to indicate other bodily fluids present	
UAYML8	ALS	Negative	YZDAA4	SaligaE	Negative
	RSID-Saliva	Negative	ZCTE3K	RSID	Negative
UF9YH8	ALS	Negative	ZKJZ9J	RSID	Negative
	RSID	Negative	ZRYNXE	PHADEBAS (TUBE METHOD)	Inconclusive
UPT6GW	Phadebas amylase	Negative	ZYFE47	Amylase Diffusion	Negative
V4BH9W	Case scenario & lab policy/procedures				
V7F8B6	RSID	Negative			
V9Z287	PHADEBAS	Negative			
VKV7HR	RSID	Negative			
VM2TTB	Amylase Diffusion	Negative			
W4LNZR	Amylase with radial diffusion	Negative			
WBNQ8D	Phadebas	Negative			
WF2Q7P	Rsid Saliva	Negative			
WM2X2R	Phadebas (Amylase)	Negative			

TABLE 2 - Item 3

Other Body Fluid Types					
Webcode	Body Fluid Type - Test Performed	Results	Webcode	Body Fluid Type - Test Performed	Results
2PQVGK	Sperm Search	Negative			
JJM2UU	Urine - DMAC	Negative			
WF2Q7P	Rsid Urine	Positive			
WM2X2R	Urine - Creatinine	Negative			
	Fecal - Micro	Negative			
	Fecal - Urobilinogen	Negative			

TABLE 2 - Item 4

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	KM	Negative	8PP7HW	LMG test	Negative
28BLGK	LMG	Negative	8XLW8D	LMG	Negative
2EMTZU	Not required per policy			Visual	Negative
2FHERZ	Swabs were not stained red-brown		8YD3JW	No visible blood-like staining observed, therefore no presumptive test for blood performed	
2PQVGK	KM	Negative	9B677N	KM	Negative
2VEVGY	O-tol	Negative	9GQWTH	No staining consistent with blood observed	
3FF4BN	KM	Negative	9HLHJN	KM	Negative
3MXG82	KM	Negative	9JH363	KM	Negative
3P4YCN	Light yellow staining was noted on each swab; no reason to test for blood.		9Q3XCX	NOT TESTED DUE TO LOCATION OF SWABBING, NO VISIBLE STAINING, AND POSITIVE SEMEN TEST.	
4APZYB	No characteristic staining observed, would not test based on staining & scenario		9UE4NY	Visual exam - no staining observed	
4BJLQG	LMG	Negative	9W3CXZ	KM	Negative
4HXUT4	LMG	Negative	AWHZRY	KM	Negative
4JVXRK	No reddish-brown stain. Not tested per policy *Not Required		B2DAF7	Hemastix - wk pos 25 secs. Did not confirm presence of blood as semen positive on victim.	Positive
68YCXV	OT	Negative	B33FWW	Not tested: No red-brown staining observed on the Items.	
6HVCDK	Visual	Negative	BA8DNA	Visual	Positive
6PZ4XB	Case scenario did not require testing for blood, absence of red-brown staining			LMG	Negative
6UWLGW	O-tol	Negative	BHXXRM	Not tested due to lack of reddish-brown color on swabs.	
6ZGMFL	Swab not red stained		BMUKKB	KM	Negative
7M6VAR	no visible stain observed		BNKQXU	HEMATRACE	Positive
7MNLCK	Swabs have no red stains		BU8H89	LMG	Negative
7XNELB	KM	Negative			
8C6JGB	KM	Negative			
8NUMDA	Testing was not performed due to case scenario.				

TABLE 2 - Item 4

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
BYM4XV	no apparent bloodstain observed		HDYUYF	KM	Negative
C4UYB4	No red/brown staining was visualized on the samples; therefore, blood testing was not performed.			Abacard Hematrace	Positive
			HEDJ2T	KM	Negative
			HHBZEM	A red stain was not observed.	
CGLJZ8	KM	Negative	HQ63HJ	O-tol	Negative
DL73AZ	KM	Negative	HWRXDW	KM	Negative
EB38QV	KM	Negative	J23ATA	KM	Negative
EB9GR9	TMB	Negative	JHVYAZ	KM	Negative
EFF3F6	KM	Negative	JJM2UU	KM	Negative
ERDV3B	Not routinely performed on DRS-negative stains		JNFWZ9	No visual indication that blood was present	
EVPZDB	White scene or kit swabs are not tested for blood in this lab		JPQ3UN	LMG	Negative
			JRWLWQ	Visual Exam	Negative
EXGQ74	KM	Negative	JTMRDG	KM	
F42C9M	Phenolphthalein	Negative	JYWP6A	KM	Negative
FVMNEL	Not applicable to case history		K6XJEV	KM	Negative
FXB6NT	KM	Negative	K7UYQR	Visual	Negative
G74RZB	O-tol	Negative		KM	Negative
G8XFEZ	TMB	Positive	KQ3XGF	KM	Negative
	HemaTrace	Positive	KUE3RF	O-TOL	Negative
G9NLVQ	LMG	Negative	KXYWNG	HemaTrace	Negative
GFC6FU	KM	Negative	KZ3R7Y	LMG	Negative
GPVJEP	Visual Exam	Negative	LXURT6	Visual exam	Negative
	KM	Negative	M9CTN6	Visual	Negative
GV43GN	no staining on swab		MED8A7	O-Tol	Negative
GZXZM8	*swabs were not red stained		MEWRED	Swabs were not tested for blood - does not display the characteristic reddish-brown stain associated with blood.	
H7J77Z	Analyst's discretion (based on scenario provided).				
H9ADNB	O-tol	Negative	MQAX88	KM	Negative

TABLE 2 - Item 4

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
MY39JM	KM	Negative	RVQ877	No apparent blood staining	
MZHQHM	Hemastix	Negative	T2V4RR	KM	Negative
N3ZCLU	KM	Negative	T7LPCY	white swab, test for blood not applicable to the item in scenario	
N9PM3J	No red-brown discoloration		TJBFHW	NT (Per Policy) - no reddish-brown color - not required	
N9PRB2	stains without any visible red/brown color are not routinely tested for blood		TKLJQE	LMG	Negative
P8P4F9	No reddish/brown stain observed during visual inspection.		TX4EYX	KM	Negative
PCGF7Y	no red-brown staining		U3J3HU	No reddish brown stains seen on swabs; not required to test	
PDCWKG	KM	Negative	UAYML8	KM Test	Negative
	O-tol	Negative	UAZKRV	Swabs were not stained	
	ABA Card	Negative	UF9YH8	No red-brown staining	
PJ2EYR	Not tested based on staining.		UHA7HB	KM	Negative
PJ7RQU	KM	Negative	UPT6GW	Visual examination	Negative
PKWW3E	KM	Negative	V4BH9W	No visual indication of blood being present	
	Visual	Negative	V7F8B6	O-Tol	Negative
PMLJF3	Item 4 was not tested for blood because of the location and color of the stain.		V9Z287	KM	Negative
PNF7UR	TMB	Inconclusive	VKV7HR	Tetrabase	Negative
	HemaTrace	Negative	VM2TTB	LMG	Negative
QADXPB	Hemastix (TMB)	Negative	W4LNZR	KM	Negative
R9VBPE	LMG	Negative	WBNQ8D	KM	Negative
RBKG8Q	KM	Negative	WF2Q7P	KM	Negative
RDTPGB	O-tol	Negative		Hematrace	Negative
RNACNQ	TMB	Positive	WM2X2R	PTMB	Negative
	HemaTrace	Positive		Slight yellowish color on swabs.	
RNQTHJ	Visually	Negative	WXGZ8Q	No visibled[sic] reddish-brown staining - not tested as per policy - therefore not required to test.	

TABLE 2 - Item 4

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
XD83L9	KM	Negative			
XHZBMU	no red stain observed				
XJC9UK	KM	Negative			
XXUUMQ	No visible reddish-brown staining - not required to test				
Y28ZFH	Testing not performed due to the probative value of the sample based on the scenario given.				
Y8GAQK	KM	Negative			
Y8WHUU	LMG	Negative			
Y9VLM6	Phenolphthalein	Negative			
YBDZQ6	Visual	Negative			
	KM	Negative			
YLWGGQ	Not brown stains				
YMCGYU	KM	Negative			
YN3XWA	Hemastix -ve @ 60 secs	Negative			
YQT4CZ	No red brown staining				
YTUAPK	LMG	Negative			
YZDAA4	Phenolphthalein	Negative			
ZCTE3K	Visual	Negative			
ZKJZ9J	KM	Negative			
ZRYNXE	KM	Negative			
ZYFE47	No red discoloration				

TABLE 2 - Item 4

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	AP	Positive	6PZ4XB	AP	Positive
	MICRO	Positive		Micro	Positive
28BLGK	AP	Positive	6UWLGW	Micro	Positive
	P30	Positive	6ZGMFL	Micro	Positive
	Micro	Positive	7M6VAR	ACP, Micro	Positive
2EMTZU	AP	Positive	7MNLCK	Microscopic	Positive
	Micro	Positive	7XNELB	AP	Positive
2FHERZ	micro	Positive		P30	Positive
2PQVGK	PSA	Positive		KPIC	Positive
2VEVGY	AP	Positive	8C6JGB	AP	Positive
	Micro	Positive		Micro	Positive
3FF4BN	AP	Positive	8NUMDA	AP	Positive
	Micro	Positive		Micro	Positive
3MXG82	AP	Positive		ALS	Positive
	Micro	Positive	8XLW8D	AP	Positive
3P4YCN	AP	Positive		P30	Positive
	Micro	Positive		Micro	Positive
4APZYB	AP	Positive	8YD3JW	AP: x4 AP positive	Positive
	Micro	Positive		Micro: 2+ heads & tails, negative for NECs	Positive
4BJLQG	AP	Positive	9B677N	AP	Positive
	PSA	Positive		Micro	Positive
	Micro	Positive	9GQWTH	AP	Positive
4HXUT4	AP	Positive		Micro	Positive
4JVXRK	AP	Positive	9HLHJN	AP	Positive
	Micro	Positive		Micro	Positive
68YCXV	Micro	Positive	9JH363	AP	Positive
6HVCDK	AP	Positive		PSA	Positive
	micro	Positive			

TABLE 2 - Item 4

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
9Q3XCX	AP	Positive		Micro	Positive
	MICRO	Positive		PSA	Positive
9UE4NY	AP	Positive	C4UYB4	ALS	Positive
	Micro	Positive		AP	Positive
9W3CXZ	AP	Positive		Micro	Positive
	PSA Semiquant	Positive	CGLJZ8	AP	Positive
AWHZRY	AP	Positive		Micro	Positive
	PSA	Positive	DL73AZ	AP	Positive
	Micro	Positive		PSA	Positive
B2DAF7	AP - pos 8 secs	Positive		Micro	Positive
	Micro - pos	Positive	EB38QV	PSA	Positive
B33FWW	AP	Positive		MICRO	Positive
	RSID	Positive	EB9GR9	AP	Positive
	Micro	Positive	EFF3F6	AP (Seri AP)	Positive
BA8DNA	AP	Positive		Micro (Xmas Stain)	Positive
	P30	Positive	ERDV3B	PSA	Positive
	Spermatozoa (Micro)	Positive		AP	Positive
BHXRM	Item 4 AnB: AP, Micro	Positive		Micro	Positive
	Item 4 CnD: AP, Micro	Positive	EVPZDB	Micro	Positive
BMUKKB	AP	Positive	EXGQ74	AP (SAP)	Positive
	Seratec	Positive		PSA (ABA)	Positive
	Micro	Positive		Micro (Xmas)	Positive
BNKQXU	PSA	Positive	F42C9M	RSID	Negative
	MICRO	Positive	FVMNEL	Micro	Positive
BU8H89	AP	Positive	FXB6NT	AP	Positive
	PSA	Positive		PSA	Positive
	Micro	Positive		Micro	Positive
BYM4XV	AP	Positive			

TABLE 2 - Item 4

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
G74RZB	AP	Positive		Micro	Positive
	Micro	Positive	J23ATA	RSID-Semen	Positive
G8XFEZ	AP	Positive		PSA	Positive
	PSA	Positive	J94HFP	AP	Positive
	Microscopic	Positive		Micro	Positive
G9NLVQ	AP	Positive	JHVYAZ	AP	Positive
GFC6FU	AP	Positive		Micro	Positive
	PSA	Positive	JJM2UU	AP	Positive
GPVJEP	AP	Positive		Micro	Positive
	H&E Slide	Positive	JNFWZ9	AP	Positive
GV43GN	AP	Positive		Microscopic sperm search	Positive
	Micro	Positive	JPQ3UN	AP	Positive
GZXZM8	Micro	Positive	JRWLWQ	Fast Blue (+ve) (AP)	Positive
H7J77Z	AP (Map)	Positive		Micro (+ve)	Positive
	Micro (Xmas)	Positive	JTMRDG	ALS	Positive
H9ADNB	ALS	Positive		AP	Positive
	AP	Positive		RSID-Semen	Positive
	Micro	Positive		Microscopy	Positive
HDYUYF	AP	Positive	JYWP6A	RSID	Positive
	Micro	Positive		PSA	Positive
	Abacard p30	Positive	K6XJEV	AP	Positive
HEDJ2T	AP	Positive		Micro	Positive
	Micro	Positive	K7UYQR	AP	Positive
HHBZEM	AP	Positive		Micro	Positive
	Micro	Positive	KQ3XGF	AP	Positive
	ALS	Positive		Micro	Positive
HQ63HJ	Micro	Positive		ABAcad p30	Positive
HWRXDW	AP	Positive			

TABLE 2 - Item 4

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
KUE3RF	AP	Positive		Micro (sperm search)	Positive
	MICRO	Positive	NQ6YNH	AP	Positive
KXYWNG	ABAcad P30	Positive		Micro	Positive
KZ3R7Y	AP	Positive	P8P4F9	AP	Positive
	PSA	Positive		Microscopy	Positive
	Microscopic	Positive	PCGF7Y	AP	Positive
LXURT6	AP	Positive		Micro	Positive
	Micro	Positive	PDCWKG	AP	Positive
M9CTN6	ALS	Positive	PJ2EYR	ALS	Positive
	AP	Positive		AP	Positive
	Micro	Positive		Micro	Positive
MED8A7	ALS	Positive	PJ7RQU	AP	Positive
	Micro	Positive		Seratec	Positive
MEWRED	AP, Micro	Positive		Micro	Positive
	AP, Micro	Positive	PKWW3E	AP	Positive
	AP, Micro	Positive		Slide made (Micro)	Positive
	AP, Micro	Positive	PMLJF3	AP	Positive
MQAX88	RSID	Positive		Micro	Positive
	PSA	Positive		ALS	Positive
MY39JM	AP	Positive	PNF7UR	AP	Positive
MZHQHM	Fast Blue (AP)	Positive		PSA	Positive
	Micro	Positive		Micro	Positive
N3ZCLU	AP	Positive	QADXPB	Microscopy	Positive
	p30	Positive	QP6VRD	BCIP	Positive
	KPIC	Positive		Microscopic Sperm Search	Positive
N9PM3J	AP	Positive	R9VBPE	Microscopy	Positive
	Micro	Positive		AP	Positive
N9PRB2	AP	Positive			

TABLE 2 - Item 4

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
RBKG8Q	AP	Positive		AP	Positive
	P-30	Positive		Micro	Positive
	KPIC	Positive	UHA7HB	AP	Positive
RDTPGB	Micro	Positive	UPT6GW	AP	Positive
RNACNQ	AP	Positive		PSA	Positive
	Micro	Positive	UPU6B6	AP	Positive
RNQTHJ	AP	Positive		MICRO	Positive
	Micro	Positive	V4BH9W	ALS	Positive
RVQ877	AP	Positive		AP	Positive
	Micro	Positive		Micro	Positive
T2V4RR	AP	Positive	V7F8B6	AP	Positive
	Seratec	Positive		Micro	Positive
	Micro	Positive	V9Z287	AP	Positive
T7LPCY	Microscopic	Positive	VKV7HR	Micro	Positive
TJBFHW	AP	Positive	VM2TTB	AP	Positive
	Micro	Positive		Micro	Positive
TKLJQE	AP	Positive	W4LNZR	AP	Positive
	Micro	Positive		Micro exam	Positive
	PSA	Positive	WBNQ8D	AP	Positive
TX4EYX	Micro	Positive		Micro	Positive
U3J3HU	AP	Positive		P30	Positive
	Micro	Positive	WF2Q7P	AP	Positive
UAYML8	ALS	Positive		Micro	Positive
	AP Test	Positive	WM2X2R	Acid Phosphatase (AP)	Positive
	RSID-Semen	Positive		Micro	Positive
	Microscopy	Positive	WXGZ8Q	AP	Positive
UAZKRV	Micro	Positive		Micro	Positive
UF9YH8	ALS	Positive	XD83L9	AP	Positive

TABLE 2 - Item 4

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
XHZBMU	Micro	Positive		P30	Positive
XJC9UK	AP	Positive		Microscopic	Positive
	Micro	Positive	ZKJZ9J	AP	Positive
XXUUMQ	AP	Positive		Seratec	Positive
	Micro	Positive		Micro	Positive
Y28ZFH	AP	Positive	ZRYNXE	AP	Positive
	Micro	Positive		MICRO	Positive
	PSA	Positive	ZYFE47	AP	Positive
	ALS	Positive		Micro	Positive
Y8GAQK	AP	Positive			
	Micro - Xmas tree stain	Positive			
Y8WHUU	AP	Positive			
	PSA	Positive			
	Micro	Positive			
Y9VLM6	RSID	Negative			
YBDZQ6	AP	Positive			
	Microscopic exam - spermatzoa[sic] identified	Positive			
YLWGGQ	Microscopic	Positive			
YMCGYU	RSID-Semen	Positive			
	PSA	Positive			
YN3XWA	Micro	Positive			
YQT4CZ	AP	Positive			
	Micro	Positive			
YTUAPK	AP	Positive			
	PSA	Positive			
	Micro	Positive			
YZDAA4	RSID	Negative			
ZCTE3K	AP	Positive			

TABLE 2 - Item 4

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	PHADEBAS	Positive	9HLHJN	Per laboratory policy, only intimate swabs from the vaginal area are saliva tested.	
2EMTZU	Not required per policy				
2FHERZ	RSID Saliva	Negative	9JH363	Phadebas Press	Negative
2PQVGK	RSID	Negative	9Q3XCX	NO TESTING DUE TO POS INDICATION OF SEMEN.	
2VEVGY	RSID	Positive	9UE4NY	Amylase Diffusion	Negative
3FF4BN	Phadebas	Negative	9W3CXZ	Phadebas press	Negative
3MXG82	Phadebas	Positive	AWHZRY	Phadebas	Negative
3P4YCN	Semen was detected and there is no reason to believe saliva would also be present on these swabs.		B2DAF7	RSID - neg 8.5 mins	Negative
4APZYB	Amylase (radial diffusion)	Negative	B33FWW	RSID	Negative
4JVXRK	Not tested per policy *Not Required		BA8DNA	RSID	Negative
68YCXV	RSID	Negative	BHXXRM	Item 4 AnB: RSID	Negative
6HVCDK	RSID-Saliva	Negative		Item 4 CnD: RSID	Negative
6PZ4XB	Case scenario did not require testing for saliva		BMUKKB	RSID	Negative
6UWLGW	RSID	Positive	BNKQXU	RSID	Negative
7M6VAR	Amylase	Negative	BYM4XV	Amylase diffusion	Negative
7MNLCK	Per Agency's Policy		C4UYB4	Saliva testing is not offered as a service in this laboratory.	
7XNELB	Phadebas	Negative	CGLJZ8	Phadebas Amylase	Negative
8C6JGB	Phadebas	Negative	DL73AZ	Phadebas Press	Negative
8NUMDA	This laboratory does not perform testing for saliva.		EB38QV	PHADEBAS	Negative
8PP7HW	Phadebas tube test	Inconclusive	EB9GR9	RSID Saliva Field Kit	Negative
8XLW8D	RSID-Saliva	Negative	EFF3F6	Amylase (Phadebas)	Negative
8YD3JW	Phadebas: x4 Phadebas negative	Negative	ERDV3B	RSID	Negative
9B677N	Not tested per policy.		EVPZDB	Saliva testing not being done in this lab at this time	
9GQWTH	Not applicable to case scenario		EXGQ74	Phadebas Amylase	Positive
			F42C9M	SaligaE	Negative
			FVMNEL	RSID	Negative

TABLE 2 - Item 4

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
FXB6NT	Phadebas	Negative	KUE3RF	RSID	Negative
G74RZB	No information to indicate oral contact by suspect and semen is present in this stain.		KXYWNG	Saligae	Negative
G9NLVQ	PHADEBAS PAPER TEST	Negative	LXURT6	RSID	Negative
GFC6FU	Phadebas	Negative	M9CTN6	RSID	Negative
GPVJEP	Phadebas Paper	Negative	MED8A7	ALS	Positive
GV43GN	Amylase diffusion	Negative		RSID	Negative
H7J7ZZ	α -Amylase	Negative	MEWRED	RSID Saliva	Negative
H9ADNB	RSID	Positive		RSID Saliva	Negative
HEDJ2T	Phadebas test	Negative		RSID Saliva	Negative
HHBZEM	Amylase diffusion	Negative	MY39JM	*saliva testing not performed in this lab	
HQ63HJ	RSID	Positive	MZHQHM	Current protocols do not allow for Saliva testing on swabs.	
HWRXDW	Radial Diffusion	Negative	N3ZCLU	Phadebas	Inconclusive
J94HFP	Phadebas	Negative	N9PM3J	Amylase Diffusion	Negative
JHVYAZ	Phadebas	Negative	N9PRB2	semen positive stains are not routinely tested for saliva	
JJM2UU	Phadebas tube test	Positive	NQ6YNH	RSID-Saliva	Negative
	Micro	Negative	P8P4F9	RSID Saliva	Negative
JNFWZ9	No senario[sic] information provided that supports likelihood of a mixed semen/saliva stain		PCGF7Y	Amylase with Radial Diffusion	Negative
JPQ3UN	PHADEBAS PAPER TEST	Negative	PDCWKG	SALlgAE	Negative
JRWLWQ	Current protocols do not allow for Saliva testing on swabs.		PJ2EYR	Testing not offered.	
JTMRDG	ALS	Positive	PJ7RQU	RSID	Negative
	RSID-Saliva	Negative	PKWW3E	Phadebas	Inconclusive
JYWP6A			PMLJF3	RSID	Positive
K6XJEV	Phadebas	Positive	PNF7UR	Phadebas (amylase)	Negative
K7UYQR	Phadebas	Positive	QADXPB	RSID α -amylase	Negative
KQ3XGF	radial diffusion	Inconclusive	QP6VRD	Amylase Diffusion	Positive

TABLE 2 - Item 4

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
R9VBPE	Starch Iodine test	Inconclusive	WXGZ8Q	Per policy, item does not meet criteria to be tested for this body fluid - therefore not required to test.	
RBKG8Q	Phadebas	Negative			
RDTPGB	RSID	Negative			
RNACNQ	Amylase	Negative	XD83L9	Not competent in saliva testing	
RVQ877	Amylase diffusion	Negative	XXUUMQ	Saliva test not performed on this item per laboratory policy - not required to test	
T2V4RR	RSID	Negative			
T7LPCY	No currently validated procedures for saliva test.		Y28ZFH	This laboratory does not offer saliva testing.	
TJBFHW	NT (Per Policy) - Item not eligible for saliva testing		Y8GAQK	Phadebas Amylase	Positive
TKLJQE	Amylase Diffusion	Negative	Y9VLM6	SaligaE	Negative
U3J3HU	Per Policy, not required to test		YBDZQ6	Phadebas tube test not accredited	
UAYML8	ALS	Positive	YN3XWA	RSID -ve @ 8½ mins	Negative
	RSID-Saliva	Negative	YQT4CZ	No indication in Scenario of saliva on victims stomach just sexual assault	
UF9YH8	ALS	Positive	YZDAA4	SaligaE	Negative
	RSID	Positive	ZCTE3K	RSID	Negative
UHA7HB	Phadebas	Negative	ZKJZ9J	RSID	Negative
UPT6GW	Phadebas amylase	Positive	ZRYNXE	PHADEBAS (TUBE METHOD)	Positive
V4BH9W	Case scenario		ZYFE47	Amylase Diffusion	Negative
V7F8B6	RSID	Negative			
V9Z287	PHADEBAS	Negative			
VKV7HR	RSID	Negative			
VM2TTB	Amylase Diffusion	Negative			
W4LNZR	Amylase with radial diffusion	Negative			
WBNQ8D	Phadebas	Negative			
WF2Q7P	RSID Saliva	Negative			
WM2X2R	Phadebas (amylase)	Negative			

TABLE 2 - Item 4

Other Body Fluid Types					
Webcode	Body Fluid Type - Test Performed	Results	Webcode	Body Fluid Type - Test Performed	Results
2PQVGK	Sperm Search	Positive			
JJM2UU	Urine - DMAC	Negative			
MZHQHM	Contact DNA - Sent for DNA analysis				
WF2Q7P	Rsid urine	Negative			
WM2X2R	Urine - Creatinine	Negative			
ZRYNXE	URINE - CREATININE - FAINT	Positive			

TABLE 2 - Item 5

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	KM	Positive		Hema Trace	Positive
28BLGK	LMG	Positive	6ZGMFL	Kastle-Meyer	Positive
	HemaTrace	Positive	7M6VAR	LMG	Positive
	Ouchterlony	Positive	7MNLCK	Kastle-Meyer Test	Positive
2EMTZU	KM	Positive	7XNELB	KM	Positive
2FHERZ	O-tol	Positive	8C6JGB	KM	Positive
	ABACARD Hematrace	Positive		Modified Takayama	Positive
2PQVGK	KM	Positive	8NUMDA	LMG	Positive
2VEVGY	O-tol	Positive		HemaTrace	Positive
	HemaTrace	Positive	8PP7HW	LMG spot test	Positive
3FF4BN	KM	Positive	8XLW8D	LMG	Positive
	Hematrace	Positive		HemaTrace	Positive
3MXG82	KM	Positive		Visual	Positive
3P4YCN	KM	Positive	8YD3JW	LMG: Visible blood-like staining observed on inside front of item, LMG presumptive test performed, LMG positive	Positive
4APZYB	KM	Positive	9B677N	KM	Positive
	HemaTrace	Positive	9GQWTH	KM	Positive
4BJLQG	LMG	Positive	9HLHJN	KM	Positive
	HemaTrace	Positive	9JH363	KM	Positive
	Ouchterlony	Positive		ABACard HemaTrace	Positive
4HXUT4	LMG	Positive	9Q3XCX	KM	Positive
4JVXRK	KM	Positive	9UE4NY	LMG	Positive
68YCXV	OT		9W3CXZ	KM	Positive
	HemaTrace	Positive		ABAcad Hematrace	Positive
6HVCDK	Visual	Positive	AWHZRY	KM	Positive
	KM	Positive	B2DAF7	Hemastix - pos 5 secs	Positive
	HemaTrace	Positive			
6PZ4XB	KM	Positive			
6UWLGW	O-tol	Positive			

TABLE 2 - Item 5

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
B33FWW	KM	Positive		Mod Taka	Positive
	RSID	Positive	F42C9M	Phenolphthalein	Positive
BA8DNA	Visual	Positive	FVMNEL	O-Tol	Positive
	LMG	Positive		Hemtrace[sic]	Positive
	HemaTrace	Positive	FXB6NT	KM	Positive
BHXXRM	Item 5A: O-tol	Positive	G74RZB	O-tol	Positive
BMUKKB	KM	Positive		HemaTrace	Positive
	Hematrace	Positive	G8XFEZ	TMB	Positive
BNKQXU	KM	Positive		HemaTrace	Positive
	HEMATRACE	Positive		Luminol	Positive
BU8H89	LMG	Positive	G9NLVQ	LMG	Positive
	HemaTrace	Positive	GFC6FU	KM	Positive
	Ouchterlony Species Determination	Positive		ABA	Positive
BYM4XV	LMG	Positive	GPVJEP	Visual Exam	Positive
	HemaTrace	Positive		KM	Positive
C4UYB4	LMG	Positive	GV43GN	LMG	Positive
	HemaTrace	Positive	GZXZM8	KM	Positive
CGLJZ8	KM	Positive	H7J77Z	KM	Positive
	Modified Takayama	Positive		Takayama	Positive
DL73AZ	KM	Positive	H9ADNB	O-tol	Positive
	ABACard HemaTrace	Positive		HemaTrace	Positive
EB38QV	KM	Positive	HDYUYF	KM	Positive
EB9GR9	TMB	Positive		Abacard Hematrace	Positive
EFF3F6	KM	Positive	HEDJ2T	KM	Positive
	Takayama	Positive		ABACard HemaTrace	Positive
ERDV3B	KM	Positive	HHBZEM	LMG	Positive
EVPZDB	KM	Positive		HemaTrace	Positive
EXGQ74	KM	Positive			

TABLE 2 - Item 5

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
HQ63HJ	O-tol	Positive		HemaTrace	Positive
	Hematrace	Positive		Ouchterlony	Positive
HWRXDW	KM	Positive	LXURT6	Visual exam	Positive
	HemaTrace	Positive		KM	Positive
J23ATA	KM	Positive		HemaTrace	Positive
	Hexagon OBTI	Positive	M9CTN6	Visual	Positive
J94HFP	TMB	Positive		KM	Positive
JHVYAZ	KM	Positive	MED8A7	O-Tol	Positive
	Modified Takayama	Positive		Hematrace	Positive
JJM2UU	KM	Positive	MEWRED	O-tol	Positive
JNFWZ9	K-M	Positive	MQAX88	KM	Positive
JPQ3UN	LMG	Positive		Hexagon OBTI	Positive
JRWLWQ	Hemochromogen (+ve)	Positive	MY39JM	KM	Positive
JTMRDG	KM	Positive	MZHQHM	Hemochromogen	Positive
	RSID-Blood	Positive	N3ZCLU	KM	Positive
JYWP6A	KM	Positive	N9PM3J	LMG	Positive
	Hexagon OBTI	Positive		HemaTrace	Positive
K6XJEV	KM	Positive	N9PRB2	KM	Positive
K7UYQR	Visual	Positive	NQ6YNH	KM	Positive
	ALS	Positive		RSID-Blood	Positive
	KM	Positive	P8P4F9	O-tol	Positive
KQ3XGF	KM	Positive	PCGF7Y	KM	Positive
	Takayama	Positive	PDCWKG	KM	Positive
	ALS Infrared	Positive		O-tol	Positive
KUE3RF	O-TOL	Positive		ABA Card	Positive
	HEMATRACE	Positive	PJ2EYR	LMG	Positive
KXYWNG	HemaTrace	Positive		Hematrace	Positive
KZ3R7Y	LMG	Positive			

TABLE 2 - Item 5

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
PJ7RQU	KM	Positive	TKLJQE	LMG	Positive
	Hematrace	Positive		HemaTrace	Positive
PKWW3E	Visual	Positive	TX4EYX	KM	Positive
	KM	Positive	U3J3HU	KM	Positive
PMLJF3	O-tol	Positive	UAYML8	KM Test	Positive
	HemaTrace	Positive		RSID-Blood	Positive
PNF7UR	TMB	Positive	UAZKRV	KM	Positive
	Luminol	Positive	UF9YH8	O-tol	Positive
	HemaTrace	Positive		Hematrace	Positive
	Haemochromogen	Positive	UHA7HB	KM	Positive
QADXPB	Hemastix (TMB)	Positive	UPT6GW	KM	Positive
QP6VRD	KM	Positive	UPU6B6	KM	Positive
R9VBPE	HemaTrace	Positive	V4BH9W	KM	Positive
	LMG	Positive	V7F8B6	O-Tol	Positive
RBKG8Q	KM	Positive		RSID	Positive
RDTPGB	O-tol	Positive	V9Z287	KM	Positive
	HemaTrace	Positive	VKV7HR	Tetrabase	Positive
RNACNQ	TMB	Positive		RSID	Positive
	HemaTrace	Positive	VM2TTB	LMG	Positive
	Luminol	Positive		Hematrace	Positive
RNQTHJ	LMG	Positive	W4LNZR	KM	Positive
	OBTI	Positive		HemaTrace	Positive
RVQ877	LMG	Positive	WBNQ8D	KM	Positive
	HemaTrace	Positive		HemaTrace	Positive
T2V4RR	KM	Positive	WF2Q7P	KM	Positive
	Hematrace	Positive		Hematrace	Positive
T7LPCY	KM	Positive	WM2X2R	PTMB	Positive
TJBFHW	KM	Positive		HemaTrace	Positive

TABLE 2 - Item 5

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
WXGZ8Q	K-M	Positive		Hematrace	Positive
XD83L9	KM	Positive	ZRYNXE	KM	Positive
XHZBMU	KM	Positive	ZYFE47	LMG	Positive
XJC9UK	KM	Positive		Hematrace	Positive
XXUUMQ	KM	Positive			
Y28ZFH	LMG	Positive			
	HemaTrace	Positive			
Y8GAQK	KM	Positive			
	Modified Takayama	Positive			
Y8WHUU	LMG	Positive			
	HemaTrace	Positive			
	Ouchterlony - Human	Positive			
Y9VLM6	Phenolphthalein	Positive			
YBDZQ6	Visual	Positive			
	KM	Positive			
YLWGGQ	Kastle-Meyer Test	Positive			
YMCGYU	KM	Positive			
	Hexagon OBTI	Positive			
YN3XWA	Hemastix (+)ve @ 5 secs	Positive			
YQT4CZ	KM	Positive			
YTUAPK	LMG	Positive			
	HemaTrace	Positive			
	Ouchterlony (human)	Positive			
YZDAA4	Phenolphthalein	Positive			
ZCTE3K	LMG	Positive			
	HemaTrace	Positive			
	Visual	Positive			
ZKJZ9J	KM	Positive			

TABLE 2 - Item 5

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	AP	Negative		AP	Negative
28BLGK	AP	Negative		Micro	Negative
2EMTZU	AP	Negative		RSID	Negative
2FHERZ	ALS	Negative	6ZGMFL	Micro	Negative
2PQVGK	PSA	Negative		P30	Negative
2VEVGY	ALS	Negative		ALS	Negative
	AP	Negative	7M6VAR	ACP	Negative
	Micro	Negative	7MNLCK	Microscopic	Negative
	RSID	Negative		P30	Negative
3FF4BN	AP	Negative		ALS	Negative
	Micro	Negative	7XNELB	AP	Negative
3MXG82	Would not expect semen to be on headwear			P30	Negative
				KPIC	Negative
3P4YCN	ALS	Positive	8C6JGB	AP	Negative
	AP	Negative	8NUMDA	This testing was not performed due to case scenario.	
	Micro	Negative			
4APZYB	AP	Positive	8XLW8D	ALS	Negative
	Micro	Positive	8YD3JW	No requirement to test item (headpiece) for semen	
4BJLQG	AP	Negative			
4HXUT4	AP	Negative	9B677N	AP	Negative
4JVXRK	AP	Inconclusive		Micro	Negative
	Micro	Negative	9GQWTH	Not applicable to case scenario	
	RSID	Negative	9HLHJN	AP	Negative
68YCXV	ALS	Negative	9JH363	AP	Negative
6HVCDK	ALS	Negative	9Q3XCX	ALS	Negative
	AP	Negative	9UE4NY	AP	Negative
6PZ4XB	Case scenario did not require testing for semen		9W3CXZ	AP	Negative
6UWLGW	ALS	Negative			

TABLE 2 - Item 5

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
AWHZRY	AP	Negative		ALS	Negative
	PSA	Negative	EVPZDB	ALS	Negative
	Micro	Negative		Micro	Negative
B2DAF7	Semen not tested as headpiece from suspect.			P30	Negative
B33FWW	AP	Negative	EXGQ74	AP (SAP)	Negative
BA8DNA	AP	Negative	F42C9M	RSID	Negative
	P30	Negative	FVMNEL	ALS	Negative
	Spermatozoa (Micro)	Negative	FXB6NT	AP	Negative
	ALS	Negative		PSA	Negative
BHXXRM	Item 5: ALS	Negative		Micro	Negative
	Item 5A: AP, Micro	Negative	G74RZB	ALS	Negative
BNKQXU	PSA	Negative		AP	Negative
	AP	Negative	G8XFEZ	AP	Negative
	MICRO	Negative	G9NLVQ	AP	Positive
BU8H89	AP	Negative	GFC6FU	AP	Negative
BYM4XV	AP	Negative	GPVJEP	AP	Negative
	Micro	Negative	GV43GN	ALS	Inconclusive
	ALS	Negative		AP	Negative
C4UYB4	ALS	Negative		Micro	Negative
	Micro	Negative	GZXZM8	Micro	Negative
	P30	Negative		P30	Negative
CGLJZ8	Not tested based on case scenario			ALS	Negative
DL73AZ	AP	Negative	H7J77Z	Analyst's discretion (based on scenario provided).	
EB38QV	PSA	Negative	H9ADNB	ALS	Negative
	MICRO	Negative		AP	Negative
EB9GR9	ALS	Negative			
EFF3F6	AP (Seri AP)	Negative			
ERDV3B	AP	Negative			

TABLE 2 - Item 5

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
HDYUYF	AP	Negative	KQ3XGF	AP	Negative
	Micro	Negative		Micro	Negative
	Abacard p30	Negative		ABAcad p30	Negative
	ALS	Negative		ALS	Negative
HEDJ2T	AP	Negative	KUE3RF	AP	Negative
HHBZEM	AP	Negative		MICRO	Negative
	Micro	Negative		RSID	Negative
	PSA	Negative		ALS	Negative
	ALS	Negative	KXYWNG	ABAcad P30	Negative
HQ63HJ	AP	Negative	KZ3R7Y	AP	Negative
	Micro	Negative	LXURT6	AP	Negative
	RSID	Negative		ALS	Negative
	ALS	Negative	M9CTN6	ALS	Negative
HWRXDW	AP	Negative	MED8A7	ALS	Negative
	Micro	Negative		AP	Negative
J23ATA	Based on scenario, semen test not conducted.			Micro	Negative
				RSID	Negative
JHVYAZ	AP	Negative	MEWRED	AP, Micro	Negative
JJM2UU	AP blot	Negative	MQAX88	Based on scenario information, semen tests not conducted.	
JNFWZ9	No scenario information that would make it likely to find semen		MY39JM	ALS	Negative
JPQ3UN	AP	Positive		AP	Negative
JTMRDG	ALS	Negative	N3ZCLU	AP	Inconclusive
	AP	Negative		p30	Negative
JYWP6A	Based on the scenario information provided, semen test was not conducted.			KPIC	Negative
				ALS	Negative
K6XJEV	Case scenario did not require testing for semen		N9PM3J	AP	Negative
				Micro	Negative
K7UYQR	AP	Negative		ALS	Negative

TABLE 2 - Item 5

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
N9PRB2	AP	Negative	RNACNQ	AP	Negative
	ALS	Negative	RNQTHJ	AP	Negative
NQ6YNH	AP	Negative		Micro	Negative
	RSID-Semen	Negative	RVQ877	ALS	Negative
P8P4F9	AP	Negative		AP	Negative
	Microscopy	Negative		Micro	Negative
PCGF7Y	Scenario didn't indicate the need for this testing.		T2V4RR	No biological fluorescence observed	
PDCWKG	AP	Negative	T7LPCY	The item was allegedly worn by suspect, test for semen would not connect suspect to victim in this scenario.	
PJ2EYR	ALS	Negative			
	No further testing performed based on case circumstances and ALS results		TJBFHW	AP	Negative
PJ7RQU	AP	Negative	TKLJQE	ALS	Negative
	Seratec	Negative		AP	Negative
	Micro	Negative		Micro	Negative
PKWW3E	Examining for 'wearer' therefore item did not necessitate		TX4EYX	ALS	Negative
				Micro	Negative
PMLJF3	AP	Negative		P30	Negative
	ALS	Positive	U3J3HU	AP	Negative
PNF7UR	AP	Negative	UAYML8	ALS	Negative
				AP Test	Negative
QADXPB	Testing for semen not done on this item - not applicable according to case scenario (suspect's clothing)		UAZKRV	Micro	Negative
				P30	Negative
R9VBPE	AP	Negative		ALS	Negative
RBKG8Q	AP	Negative	UF9YH8	ALS	Negative
	P-30	Negative	UPT6GW	AP	Negative
	KPIC	Negative	V4BH9W	ALS	Negative
RDTPGB	AP	Negative			
	ALS	Negative			

TABLE 2 - Item 5

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
V7F8B6	ALS	Negative	Y9VLM6	RSID	Negative
	RSID	Negative	YBDZQ6	AP	Negative
	AP	Negative	YLWGGQ	ALS	Negative
	Micro	Negative		Microscopic	Negative
V9Z287	AP	Negative		P30	Negative
VKV7HR	Micro	Negative	YMCGYU	Based on scenario information provided, semen test was not conducted.	
	PSA	Negative	YN3XWA	Polilight	Negative
	RSID	Negative		AP Test - Not Tested. Based on information supplied there is no suggestion there would be semen on an item worn on the assailants head.	
	AP	Negative			
VM2TTB	AP	Negative			
	Micro	Negative			
	ALS	Negative			
W4LNZR	AP	Negative	YQT4CZ	No other stains observed	
	Micro exam	Negative	YTUAPK	AP	Negative
WBNQ8D	AP	Negative	YZDAA4	RSID	Negative
WF2Q7P	AP	Negative	ZCTE3K	AP	Negative
WM2X2R	Acid Phosphatase (AP)	Negative		ALS	Negative
WXGZ8Q	AP	Negative	ZKJZ9J	AP	Negative
XD83L9	ALS	Negative		Seratec	Negative
XHZBMU	ALS	Positive		Micro	Negative
	Micro	Negative	ZRYNXE	AP	Negative
	P30	Negative	ZYFE47	ALS	Negative
XJC9UK	AP	Negative		AP	Negative
XXUUMQ	AP	Negative		Micro	Negative
Y28ZFH	Testing not performed due to the probative value of the sample based on the scenario given.				
Y8GAQK	AP	Negative			
Y8WHUU	AP	Negative			

TABLE 2 - Item 5

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	PHADEBAS	Negative	9B677N	Not tested per policy.	
2EMTZU	Not required per policy		9GQWTH	Not applicable to case scenario	
2FHERZ	ALS	Negative	9HLHJN	Per laboratory policy, only intimate swabs from the vaginal area are saliva tested.	
2PQVGK	RSID	Negative			
2VEVGY	RSID	Negative	9JH363	Phadebas Press	Negative
3FF4BN	Phadebas	Negative	9Q3XCX	ALS	Negative
3MXG82	Phadebas paper test	Negative	9UE4NY	Amylase Diffusion	Negative
3P4YCN	ALS	Positive	9W3CXZ	Phadebas press	Negative
	RSID-Saliva	Negative	AWHZRY	Phadebas	Negative
4APZYB	Amylase (radial diffusion)	Negative	B2DAF7	Polylight[sic] - neg	Negative
4JVXRK	Not tested per policy *Not Required			Phadebas Press Test - neg	Negative
68YCXV	ALS	Negative	BA8DNA	RSID	Negative
6HVCDK	ALS	Negative	BHXXRM	Not tested due to item type.	
	RSID-Saliva	Negative	BNKQXU	RSID	Negative
6PZ4XB	Case scenario did not require testing for saliva		BYM4XV	Amylase diffusion	Negative
6UWLGW	ALS	Negative	C4UYB4	Saliva testing is not offered as a service in this laboratory.	
	RSID	Negative	CGLJZ8	Not tested based on case scenario	
7M6VAR	Amylase	Negative	DL73AZ	Phadebas Press	Negative
7MNLCK	Per Agency's Policy		EB38QV	PHADEBAS	Negative
7XNELB	Phadebas	Negative	EB9GR9	ALS	Negative
8C6JGB	Phadebas	Negative	EFF3F6	Amylase (Phadebas)	Negative
8NUMDA	This laboratory does not perform testing for saliva.		ERDV3B	Not routinely performed on KM positive stains/ ALS negative	
8PP7HW	Phadebas press test	Negative	EVPZDB	Saliva testing not being done in this lab at this time	
8XLW8D	RSID-Saliva	Negative	EXGQ74	Phadebas Amylase	Negative
8YD3JW	Item did not appear to have been worn over the face, therefore item not presumptively tested for saliva		F42C9M	SaligaE	Negative

TABLE 2 - Item 5

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
FVMNEL	ALS	Negative	KQ3XGF	radial diffusion	Inconclusive
FXB6NT	Phadebas	Negative	KUE3RF	RSID	Negative
G74RZB	No information to indicate oral contact and ALS was negative			ALS	Negative
GFC6FU	Phadebas	Negative	KXYWNG	Saligae	Negative
GPVJEP	Phadebas Paper	Negative	LXURT6	RSID	Negative
GV43GN	Amylase diffusion	Negative		ALS	Negative
H7J77Z	Analyst's discretion (based on scenario provided).		M9CTN6	Assumed headpiece was worn as designed.	
H9ADNB	Bloodstain was very heavy and concentrated and did not appear diluted		MED8A7	ALS	Negative
HEDJ2T	Phadebas test	Negative		RSID	Negative
HHBZEM	Amylase diffusion	Negative	MEWRED	Did not test for saliva based on case scenario.	
HQ63HJ	RSID	Negative	MY39JM	*saliva testing not performed in this lab	
	ALS	Negative	MZHQHM	Current protocol requires a specific request for saliva as it is not a routine test.	
HWRXDW	Radial Diffusion	Negative	N3ZCLU	Phadebas	Negative
J94HFP	Phadebas	Negative		ALS	Negative
JHVYAZ	Phadebas	Negative	N9PM3J	Amylase Diffusion	Negative
JJM2UU	Phadebas blot	Negative	N9PRB2	KM positive stains are not routinely tested for saliva	
JNFWZ9	Same reason as above [Item 5 - Semen Tests Performed]		P8P4F9	RSID Saliva	Negative
JRWLWQ	Current protocol requires a specific request for saliva as it is not a routine test.		PCGF7Y	Scenario didn't indicate the need for this testing.	
JTMRDG	ALS	Negative	PDCWKG	SALlgAE	Negative
	RSID-Saliva	Negative	PJ2EYR	Test not offered.	
JYWP6A			PJ7RQU	RSID	Negative
K6XJEV	Case scenario did not require testing for saliva		PKWW3E	Blood found therefore no further examination as headwear	
K7UYQR	Phadebas	Negative	PMLJF3	RSID	Negative
				ALS	Positive
			PNF7UR	Phadebas (amylase)	Negative

TABLE 2 - Item 5

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
QADXPB	Phadebas Press Test	Negative	XD83L9	Not competent in saliva testing	
RBKG8Q	Phadebas	Negative	XXUUMQ	Saliva test not performed on this item per laboratory policy - not required to test	
RDTPGB	RSID	Negative			
	ALS	Negative	Y28ZFH	This laboratory does not offer saliva testing.	
RNACNQ	Amylase	Negative	Y8GAQK	Phadebas Amylase	Negative
RVQ877	Amylase diffusion	Negative	Y9VLM6	SaligaE	Negative
T2V4RR	No biological fluorescence observed		YBDZQ6	Phadebas	Negative
T7LPCY	No currently validated procedures for saliva test.		YN3XWA	Phadebas Press Test	Negative
TJBFHW	NT (Per Policy) - Item not eligible for saliva testing		YQT4CZ	No other stains observed	
TKLJQE	Amylase Diffusion	Negative	YZDAA4	SaligaE	Negative
U3J3HU	Per Policy, not required to test		ZCTE3K	RSID	Negative
UAYML8	ALS	Negative	ZKJZ9J	RSID	Negative
	RSID-Saliva	Negative	ZRYNXE	PHADEBAS (BLOTTING METHOD)	Negative
UF9YH8	ALS	Negative	ZYFE47	Amylase Diffusion	Negative
UPT6GW	Phadebas amylase	Negative			
V4BH9W	Case scenario				
V7F8B6	RSID	Negative			
VKV7HR	RSID	Negative			
VM2TTB	Amylase Diffusion	Negative			
W4LNZR	Amylase with radial diffusion	Negative			
WBNQ8D	Phadebas	Negative			
WF2Q7P	Phadebas	Negative			
WM2X2R	Phadebas (amylase)	Negative			
WXGZ8Q	Per policy, item does not meet criteria to be tested for this body fluid - therefore not required to test.				

TABLE 2 - Item 5

Other Body Fluid Types					
Webcode	Body Fluid Type - Test Performed	Results	Webcode	Body Fluid Type - Test Performed	Results
28BLGK	Biological Fluids - ALS	Negative			
4BJLQG	Biological Fluids - ALS	Negative			
7M6VAR	ALS	Negative			
8NUMDA	ALS	Negative			
AWHZRY	Semen - ALS	Negative			
BU8H89	Bodily Fluid - ALS	Negative			
JRWLWQ	Hair exam - Microscopic exam	Negative - (unsuitable) Animal hair			
	Contact - DNA sent for DNA analysis				
KZ3R7Y	Biological fluid - ALS	Negative			
MEWRED	Semen, Saliva, Sweat - ALS (See additional comments page [Table 3])	Negative			
MZHQHM	Contact DNA - sent for DNA analysis				
	Hair Exam - Microscopic Exam	Negative - Unsuitable Animal Hair			
QP6VRD	ALS	Negative			
RBKG8Q	ALS	Negative			
RNQTHJ	Hair: Micro	Negative			
	Forensic light source	Inconclusive			
WF2Q7P	Micrsopic[sic] examination for hairs Minitaped	Negative			
Y8WHUU	Biological Fluids - ALS	Negative			
YTUAPK	Bodily Fluids - ALS	Negative			

TABLE 2 - Item 6

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	KM	Positive	7M6VAR	LMG	Positive
28BLGK	LMG	Positive	7MNLCK	Kastle-Meyer Test	Positive
	HemaTrace	Positive	7XNELB	KM	Positive
	Ouchterlony	Positive	8C6JGB	KM	Positive
2EMTZU	KM	Positive		Modified Takayama	Positive
2FHERZ	O-tol	Positive	8NUMDA	LMG	Positive
2PQVGK	KM	Positive		HemaTrace	Positive
2VEVGY	O-tol	Positive	8PP7HW	LMG spot test	Positive
	HemaTrace	Positive	8XLW8D	LMG	Positive
3FF4BN	KM	Positive		HemaTrace	Positive
	Hematrace	Positive		Visual	Positive
3MXG82	KM	Positive	8YD3JW	LMG: Reddish staining observed on inside crotch area, LMG presumptive test performed, LMG positive	Positive
3P4YCN	KM	Positive			
4APZYP	KM	Positive	9B677N	KM	Positive
	HemaTrace	Positive	9GQWTH	KM	Positive
4BJLQG	LMG	Positive	9HLHJN	KM	Positive
	HemaTrace	Positive	9JH363	KM	Positive
	Ouchterlony	Positive		ABACard HemaTrace	Positive
4HXUT4	LMG	Positive	9Q3XCX	KM	Positive
4JVXRK	KM	Positive	9UE4NY	LMG	Positive
68YCXV	OT	Positive	9W3CXZ	KM	Positive
6HVCDK	Visual	Positive		ABACard Hematrace	Positive
	KM	Positive	AWHZRY	KM	Positive
	HemaTrace	Positive	B2DAF7	Hemastix - pos 5 secs	Positive
6PZ4XB	KM	Positive	B33FWW	KM	Positive
6UWLGW	O-tol	Positive		RSID	Positive
	Heme Trace[sic]	Positive			
6ZGMFL	Kastle-Meyer	Positive			

TABLE 2 - Item 6

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
BA8DNA	Visual	Positive	FVMNEL	No visible red-brown stains	
	LMG	Positive	FXB6NT	KM	Positive
	HemaTrace	Positive	G74RZB	O-tol	Positive
BHXXRM	Item 6A: O-tol	Positive		HemaTrace	Positive
BMUKKB	KM	Positive	G8XFEZ	TMB	Positive
	Hematrace	Positive		HemaTrace	Positive
BNKQXU	KM	Positive		Luminol	Positive
	HEMATRACE	Positive	G9NLVQ	LMG	Positive
BU8H89	LMG	Positive	GFC6FU	KM	Positive
	HemaTrace	Positive		ABA	Positive
	Ouchterlony Species Determination	Positive	GPVJEP	Visual Exam	Positive
BYM4XV	LMG	Positive		KM	Positive
	HemaTrace	Positive	GV43GN	LMG	Positive
C4UYB4	LMG	Positive	GZXZM8	KM	Positive
	HemaTrace	Positive	H7J77Z	KM	Positive
CGLJZ8	KM	Positive		Takayama	Positive
	Modified Takayama	Positive	H9ADNB	O-tol	Positive
DL73AZ	KM	Positive		HemaTrace	Positive
	ABACard HemaTrace	Positive	HDYUYF	KM	Positive
EB38QV	KM	Positive		Abacard Hematrace	Positive
EB9GR9	TMB	Positive	HEDJ2T	KM	Positive
EFF3F6	KM	Positive		ABACard HemaTrace	Positive
	Takayama	Positive	HHBZEM	LMG	Positive
ERDV3B	KM	Positive		HemaTrace	Positive
EVPZDB	KM	Positive	HQ63HJ	O-tol	Positive
EXGQ74	KM	Positive		Hematrace	Positive
	Mod Taka	Positive	HWRXDW	KM	Positive
F42C9M	Phenolphthalein	Positive		HemaTrace	Positive

TABLE 2 - Item 6

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
J23ATA	KM	Positive		KM	Positive
	Hexagon OBTI	Positive	MED8A7	O-Tol	Positive
J94HFP	TMB	Positive	MEWRED	O-tol	Positive
JHVYAZ	KM	Positive	MQAX88	KM	Positive
	Modified Takayama	Positive		Hexagon OBTI	Positive
JJM2UU	KM	Positive	MY39JM	KM	Positive
JNFWZ9	K-M	Positive	MZHQHM	Hemochromogen	Positive
JPQ3UN	LMG	Positive		Hemastix	Positive
JRWLWQ	Hemochromogen (+ve)	Positive	N3ZCLU	KM	Positive
JTMRDG	KM	Positive	N9PM3J	LMG	Positive
	RSID-Blood	Positive		HemaTrace	Positive
JYWP6A	KM	Positive	N9PRB2	KM	Positive
	Hexagon OBTI	Positive	NQ6YNH	KM	Positive
K6XJEV	KM	Positive		RSID-Blood	Positive
K7UYQR	Visual	Positive	P8P4F9	O-tol	Positive
	KM	Positive	PCGF7Y	KM	Positive
KQ3XGF	KM	Positive	PDCWKG	KM	Positive
	Takayama	Positive		O-tol	Positive
KUE3RF	O-TOL	Positive		ABA Card	Positive
	HEMATRACE	Positive	PJ2EYR	LMG	Positive
KXYWNG	HemaTrace	Positive		Hematrace	Positive
KZ3R7Y	LMG	Positive	PJ7RQU	KM	Positive
	HemaTrace	Positive		Hematrace	Positive
	Ouchterlony	Positive	PKWW3E	Visual	Positive
LXURT6	Visual exam	Positive		KM	Positive
	KM	Positive	PMLJF3	O-tol	Positive
	HemaTrace	Positive		HemaTrace	Positive
M9CTN6	Visual	Positive			

TABLE 2 - Item 6

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
PNF7UR	TMB	Positive	UF9YH8	O-tol	Positive
	Luminol	Positive	UHA7HB	KM	Positive
	HemaTrace	Positive	UPT6GW	KM	Positive
QADXPB	Hemastix (TMB)	Positive	V4BH9W	KM	Positive
QP6VRD	KM	Positive	V7F8B6	O-Tol	Positive
R9VBPE	HemaTrace	Positive		RSID	Positive
	LMG	Positive	V9Z287	KM	Positive
RBKG8Q	KM	Positive	VKV7HR	Tetrabase	Positive
RDTPGB	O-tol	Positive		RSID	Positive
	HemaTrace	Positive	VM2TTB	LMG	Positive
RNACNQ	TMB	Positive		Hematrace	Positive
	HemaTrace	Positive	W4LNZR	KM	Positive
	Luminol	Positive		HemaTrace	Positive
RNQTHJ	LMG	Positive	WBNQ8D	KM	Positive
	OBTI	Positive		HemaTrace	Positive
RVQ877	LMG	Positive	WF2Q7P	KM	Positive
	HemaTrace	Positive		Hematrace	Positive
T2V4RR	KM	Positive	WM2X2R	PTMB	Positive
	Hematrace	Positive		HemaTrace	Positive
T7LPCY	KM	Positive	WXGZ8Q	K-M	Positive
TJBFHW	KM	Positive	XD83L9	KM	Positive
TKLJQE	LMG	Positive	XHZBMU	KM	Positive
	Hematrace	Positive	XJC9UK	KM	Positive
TX4EYX	KM	Positive	XXUUMQ	KM	Positive
U3J3HU	KM	Positive	Y28ZFH	LMG	Positive
UAYML8	KM Test	Positive		HemaTrace	Positive
	RSID-Blood	Positive	Y8GAQK	KM	Positive
UAZKRV	KM	Positive		Modified Takayama	Positive

TABLE 2 - Item 6

Blood					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
Y8WHUU	LMG	Positive			
	HemaTrace	Positive			
	Ouchterlony-Human	Positive			
Y9VLM6	Phenolphthalein	Positive			
YBDZQ6	Visual	Positive			
	KM	Positive			
YLWGGQ	Kastle-Meyer Test	Positive			
YMCGYU	KM	Positive			
	Hexagon OBTI	Positive			
YN3XWA	Hemastix (+)ve @ 5 secs	Positive			
YQT4CZ	KM	Positive			
YTUAPK	LMG	Positive			
	HemaTrace	Positive			
	Ouchterlony (human)	Positive			
YZDAA4	Phenolphthalein	Positive			
ZCTE3K	LMG	Positive			
	HemaTrace	Positive			
	Visual	Positive			
ZKJZ9J	KM	Positive			
	Hematrace	Positive			
ZRYNXE	KM	Positive			
ZYFE47	LMG	Positive			
	Hematrace	Positive			

TABLE 2 - Item 6

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	AP	Positive		Micro	Positive
	MICRO	Positive	68YCXV	Micro	Positive
28BLGK	AP	Positive		ALS	
	P30	Positive		AP	
	Micro	Positive	6HVC DK	ALS	Negative
2EMT ZU	AP	Positive		AP	Positive
	Micro	Positive		micro	Positive
2FHERZ	micro	Positive	6PZ4XB	AP	Positive
	ALS	Negative		Micro	Positive
	AP	Positive	6UWLGW	ALS	Negative
2PQVGK	PSA	Positive		AP	Positive
2VEVGY	ALS	Negative		Micro	Positive
	AP	Positive	6ZGMFL	Micro	Positive
	Micro	Positive		ALS	Positive
3FF4BN	AP	Positive	7M6VAR	ACP, Micro	Positive
	Micro	Positive	7MNLCK	Microscopic	Positive
3MXG82	AP	Positive		ALS	Negative
	Micro	Positive	7XNELB	AP	Positive
3P4YCN	ALS	Negative		P30	Positive
	AP	Positive		KPIC	Positive
	Micro	Positive	8C6JGB	AP	Positive
4APZYB	AP	Positive		Micro	Positive
	Micro	Positive	8NUMDA	Micro	Positive
4BJLQG	AP	Positive		PSA	Positive
	PSA	Positive		ALS	Negative
	Micro	Negative			
4HXUT4	AP	Positive			
4JVXRK	AP	Positive			

TABLE 2 - Item 6

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
8XLW8D	ALS	Positive	B33FWW	AP	Positive
	AP	Positive		RSID	Positive
	P30	Positive		Micro	Positive
	Micro	Positive		BA8DNA	ALS
8YD3JW	AP Back: AP negative	Negative		AP	Positive
	AP Front: AP positive	Positive		P30	Positive
	Micro: 1+ heads, negative for NECs	Positive		Spermatozoa (Micro)	Positive
9B677N	AP	Positive	BHXXRM	Item 6: ALS	Negative
	Micro	Positive		Item 6A: AP, Micro	Positive
9GQWTH	ALS	Negative	BMUKKB	AP	Positive
	AP	Positive		Seratec	Positive
	Micro	Positive		Micro	Positive
9HLHJN	AP	Positive	BNKQXU	PSA	Positive
	Micro	Positive		AP	Positive
9JH363	AP	Positive		MICRO	Positive
	PSA	Positive		BU8H89	AP
9Q3XCX	AP	Negative		PSA	Positive
	ALS	Negative		Micro	Positive
9UE4NY	AP	Positive	BYM4XV	AP	Positive
	Micro	Positive		Micro	Positive
9W3CXZ	AP	Positive		ALS	Negative
	PSA Semiquant	Positive		AP mapping	Positive
AWHZRY	AP	Positive	C4UYB4	ALS	Negative
	PSA	Positive		Micro	Positive
	Micro	Positive	CGLJZ8	AP	Inconclusive
B2DAF7	Polylight[sic] - v. faint fluorescence	Positive		Micro	Positive
	AP - pos 10 secs	Positive			
	Micro - pos	Positive			

TABLE 2 - Item 6

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
DL73AZ	AP	Positive	G8XFEZ	AP	Positive
	PSA	Positive		PSA	Positive
	Micro	Positive		Microscopic	Positive
EB38QV	PSA	Positive	G9NLVQ	AP	Positive
	MICRO	Positive	GFC6FU	AP	Positive
EB9GR9	ALS	Negative	GPVJEP	PSA	Positive
	AP on bloodstain	Positive		AP	Positive
	AP near bloodstain	Negative		H&E Slide	Positive
EFF3F6	AP (Seri AP)	Positive	GV43GN	ALS	Inconclusive
	Micro (Xmas Stain)	Positive		AP	Positive
ERDV3B	ALS	Positive	GZXZM8	Micro	Positive
	AP	Positive		PSA	Positive
	PSA	Positive		Micro	Positive
	Micro	Positive		ALS	Negative
EVPZDB	ALS	Negative	H7J77Z	AP (Map)	Positive
	Micro	Positive		Micro (Xmas)	Positive
EXGQ74	AP (SAP)	Positive	H9ADNB	ALS	Positive
	Micro (Xmas)	Positive		AP	Positive
F42C9M	RSID	Positive	HDYUYF	Micro	Positive
	Christmas Tree Stain	Positive		AP	Positive
FVMNEL	ALS	Negative	HEDJ2T	Micro	Positive
	AP	Positive		Abacard p30	Positive
	Micro	Positive		ALS	Negative
FXB6NT	AP	Positive	HEDJ2T	AP	Positive
	PSA	Positive		ABAcad P30	Positive
	Micro	Positive		Micro	Positive
G74RZB	AP	Positive			
	Micro	Positive			

TABLE 2 - Item 6

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
HHBZEM	AP	Positive		Microscopy	Positive
	Micro	Positive	JYWP6A	AP	Positive
	ALS	Positive		RSID	Positive
	Evidence mapping	Positive		PSA	Positive
HQ63HJ	AP	Positive	K6XJEV	AP	Positive
	Micro	Positive		Micro	Positive
	ALS	Negative	K7UYQR	AP	Positive
HWRXDW	AP	Positive		Micro	Positive
	Micro	Positive	KQ3XGF	AP	Positive
J23ATA	AP	Positive		Micro	Positive
	RSID-Semen	Positive		ABAcad p30	Positive
	PSA	Positive		ALS	Positive
J94HFP	AP	Positive	KUE3RF	AP	Positive
	Micro	Positive		MICRO	Positive
JHVYAZ	AP	Positive		ALS	Negative
	Micro	Positive	KXYWNG	ABAcad P30	Positive
JJM2UU	AP blot	Positive	KZ3R7Y	AP	Positive
	Micro	Positive		PSA	Positive
JNFWZ9	ALS	Negative		Microscopic	Positive
	AP	Positive	LXURT6	AP	Positive
	Microscopic sperm search	Positive		Micro	Positive
JPQ3UN	AP	Positive		ALS	Negative
JRWLWQ	Fast Blue (+ve) (AP)	Positive	M9CTN6	ALS	Negative
	Micro (-ve)	Negative		AP	Positive
	P30 (+ve)	Positive		Micro	Positive
JTMRDG	ALS	Positive	MED8A7	ALS	Negative
	AP	Positive		AP	Positive
	RSID-Semen	Positive		Micro	Positive

TABLE 2 - Item 6

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
MEWRED	AP, Micro	Positive	PJ7RQU	AP	Positive
MQAX88	AP	Positive		Seratec	Positive
	RSID	Positive		Micro	Positive
	PSA	Positive	PKWW3E	AP	Positive
MY39JM	ALS	Negative		Slide made	Positive
	AP	Positive	PMLJF3	AP	Positive
MZHQHM	Fast Blue (AP)	Negative		Micro	Negative
N3ZCLU	AP	Positive		PSA	Positive
	p30	Positive		ALS	Negative
	KPIC	Positive	PNF7UR	AP	Positive
	ALS	Negative		PSA	Positive
N9PM3J	AP	Positive		Micro	Negative
	Micro	Positive	QADXPB	AP Test	Positive
	ALS	Negative		Microscopy	Positive
N9PRB2	AP	Positive	QP6VRD	BCIP	Positive
	Micro (sperm search)	Negative		P30	Positive
	PSA	Positive		Microscopic Sperm Search	Negative
	ALS	Positive	R9VBPE	Microscopy	Negative
NQ6YNH	AP	Positive		AP	Positive
	Micro	Positive	RBKG8Q	AP	Positive
P8P4F9	AP	Positive		P-30	Positive
	Microscopy	Positive		KPIC	Positive
PCGF7Y	AP	Positive	RDTPGB	AP	Positive
	Micro	Positive		Micro	Positive
PDCWKG	AP	Positive		ALS	Negative
PJ2EYR	ALS	Negative	RNACNQ	AP	Positive
	Micro	Positive		Micro	Positive
	P30	Positive			

TABLE 2 - Item 6

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
RNQTHJ	AP	Positive		AP	Positive
	Micro	Positive		Micro	Positive
RVQ877	ALS	Negative	UHA7HB	AP	Positive
	AP	Positive	UPT6GW	PSA	Positive
	Micro	Positive		Micro	Positive
	AP mapping	Positive	UPU6B6	AP	Positive
T2V4RR	AP	Positive		MICRO	Negative
	Seratec	Positive		PSA	Positive
	Micro	Positive	V4BH9W	ALS	Negative
T7LPCY	Microscopic	Positive		AP	Positive
	ALS	Negative		Micro	Positive
TJBFHW	AP	Positive	V7F8B6	AP	Positive
	Micro	Positive		ALS	Positive
TKLJQE	ALS	Positive		Micro	Positive
	AP	Positive	V9Z287	AP	Positive
	Micro	Positive	VKV7HR	AP	Positive
	PSA	Positive		Micro	Positive
TX4EYX	ALS	Negative		PSA	Positive
	Micro	Positive		RSID	Positive
U3J3HU	AP	Positive	VM2TTB	AP	Positive
	Micro	Positive		Micro	Positive
UAYML8	ALS	Positive		ALS	Negative
	AP	Positive	W4LNZR	AP	Positive
	RSID-Semen	Positive		Micro exam	Negative
	Microscopy	Positive		P30	Positive
UAZKRV	Micro	Positive	WBNQ8D	AP	Positive
	ALS	Negative		Micro	Positive
UF9YH8	ALS	Negative		P30	Positive

TABLE 2 - Item 6

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
WF2Q7P	AP	Positive	YLVGGQ	microscopic exam - spermatazoa[sic] identified	Positive
	Micro	Negative		ALS	Positive
	Rsid semen	Positive		Microscopic	Positive
WM2X2R	Acid Phosphatase (AP)	Positive	YMCGYU	AP	Positive
	Micro	Positive		RSID-Semen	Positive
WXGZ8Q	AP	Positive	YQ44CZ	PSA	Positive
	Micro	Positive		Polilight	Negative
XD83L9	ALS	Negative	YQ44CZ	AP Test	Positive
	AP	Positive		Micro	Positive
XHZBMU	ALS	Positive	YQ44CZ	AP	Positive
	Micro	Positive		Micro	Negative
XJC9UK	AP	Positive	YTUAPK	PSA	Positive
	Micro	Negative		AP	Positive
	ABA p30	Positive		PSA	Positive
XXUUMQ	AP	Positive	YZDAA4	Micro	Positive
	Micro	Positive		RSID	Positive
Y28ZFH	AP	Positive	ZCTE3K	Christmas Tree Stain	Positive
	Micro	Positive		AP	Positive
	PSA	Positive		P30	Positive
	ALS	Negative		Microscopic	Positive
Y8GAQK	AP	Positive	ZKJZ9J	ALS	Positive
	Micro - Xmas tree stain	Positive		AP	Positive
Y8WHUU	AP	Positive	ZRYNXX	Seratec	Positive
	PSA	Positive		Micro	Positive
	Micro	Positive		AP	Positive
Y9VLM6	RSID	Positive	ZRYNXX	MICRO	Positive
	Christmas Tree Stain	Positive			
YBDZQ6	AP	Positive			

TABLE 2 - Item 6

Semen					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
ZYFE47	ALS	Negative			
	AP	Positive			
	Micro	Positive			

TABLE 2 - Item 6

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
27CG2C	PHADEBAS	Negative	9B677N	Not tested per policy.	
2EMTZU	Not required per policy		9GQWTH	Not applicable to case scenario	
2FHERZ	ALS	Negative	9HLHJN	Per laboratory policy, only intimate swabs from the vaginal area are saliva tested.	
	RSID Saliva	Negative			
2PQVGK	RSID	Negative	9JH363	Phadebas Press	Negative
2VEVGY	RSID	Negative	9Q3XCX	AMYLASE DIFFUSION	Negative
3FF4BN	Phadebas	Negative		ALS	Negative
3MXG82	Phadebas paper test	Negative	9UE4NY	Amylase Diffusion	Negative
3P4YCN	Semen was detected and there is no reason to believe saliva would also be present in the underwear.		9W3CXZ	Phadebas press	Negative
4APZYB	Amylase (radial diffusion)	Negative	AWHZRY	Phadebas	Negative
4JVXK	Not tested per policy *Not Required		B2DAF7	No indication in scenario about oral rape and was semen positive.	
68YCXV	RSID	Negative	BA8DNA	RSID	Negative
6HVCDK	ALS	Negative	BHXXRM	Item 6A: RSID	Negative
	RSID-saliva	Negative	BMUKKB	RSID	Negative
6PZ4XB	Case scenario did not require testing for saliva		BNKQXU	RSID	Negative
6UWLGW	ALS	Negative	BYM4XV	Amylase diffusion	Negative
	RISD[sic]	Negative		Amylase mapping	Negative
7M6VAR	Amylase	Negative	C4UYB4	Saliva testing is not offered as a service in this laboratory.	
7MNLCK	Per Agency's Policy		CGLJZ8	Not tested based on case scenario	
7XNELB	Phadebas	Negative	DL73AZ	Phadebas Press	Negative
8C6JGB	Phadebas	Negative	EB38QV	PHADEBAS	Negative
8NUMDA	This laboratory does not perform testing for saliva.		EB9GR9	ALS	Negative
8PP7HW	Phadebas press test	Negative		RSID Saliva Field Kit (on bloodstain)	Negative
8XLW8D	RSID-Saliva	Negative	EFF3F6	Amylase (Phadebas)	Negative
8YD3JW	Phadebas: Phadebas negative	Negative	ERDV3B	RSID	Negative

TABLE 2 - Item 6

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
EVPZDB	Saliva testing not being done in this lab at this time		JRWLWQ	Current protocol requires a specific request for saliva as it is not a routine test.	
EXGQ74	Stain in crotch, where vaginal secretions may be, (known to contain amylase and react + with phadebas test) would not test crotch for saliva with phadebas test		JTMRDG	ALS	Negative
				RSID-Saliva	Negative
			JYWP6A		
			K6XJEV	Case scenario did not require testing for saliva	
F42C9M	SaligaE	Negative	K7UYQR	Phadebas	Negative
FVMNEL	ALS	Negative	KQ3XGF	radial diffusion	Inconclusive
	RSID	Negative	KUE3RF	RSID	Negative
FXB6NT	Phadebas	Negative		ALS	Negative
G74RZB	No information to indicate oral contact by suspect and semen is present in this stain.		KXYWNG	Saligae	Negative
			LXURT6	RSID	Negative
GFC6FU	Phadebas	Negative		ALS	Negative
GV43GN	Amylase diffusion	Negative	M9CTN6	Would normally not test panties for saliva unless oral sodomy was reported.	
H7J77Z	Analyst's discretion (based on scenario provided).		MED8A7	ALS	Negative
H9ADNB	RSID	Positive		RSID	Negative
HEDJ2T	Phadebas test	Negative	MEWRED	RSID Saliva	Negative
HHBZEM	Amylase diffusion	Negative	MY39JM	*saliva testing not performed in this lab	
	Evidence mapping	Negative	MZHQHM	Current protocol requires a specific request for saliva as it is not a routine test.	
HQ63HJ	RSID	Negative			
	ALS	Negative	N3ZCLU	Phadebas	Negative
HWRXDW	Radial Diffusion	Negative		ALS	Negative
JHVYAZ	Phadebas	Negative	N9PM3J	Amylase Diffusion	Negative
JJM2UU	Phadebas blot	Negative	N9PRB2	underwear crotch stains are not routinely tested for saliva unless oral sex is indicated	
JNFWZ9	No scenario info provided that would indicate saliva as a possibility		P8P4F9	RSID Saliva	Negative

TABLE 2 - Item 6

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
PCGF7Y	Amylase with Radial Diffusion	Negative	UPT6GW	Phadebas amylase	Negative
PDCWKG	SALlgAE	Negative	V4BH9W	Case scenario	
PJ2EYR	Testing not offered.		V7F8B6	RSID	Negative
PJ7RQU	RSID	Negative	V9Z287	PHADEBAS	
PKWW3E	Phadebas (only inside front and crotch tested)	Negative	VKV7HR	RSID	Negative
PMLJF3	Item 6 was not tested for saliva because the case scenario did not indicate an oral assault.		VM2TTB	Amylase Diffusion	Negative
PNF7UR	Deemed not relevant to scenario		W4LNZR	Amylase with radial diffusion	Negative
QADXPB	Phadebas Press Test	Negative	WBNQ8D	Phadebas	Negative
QP6VRD	Amylase Diffusion	Positive	WF2Q7P	Phadebas	Negative
RBKG8Q	Phadebas	Negative	WM2X2R	Phadebas (amylase)	Negative
RDTPGB	RSID	Negative	WXGZ8Q	Per policy, item does not meet criteria to be tested for this body fluid - therefore not required to test.	
	ALS	Negative	XD83L9	Not competent in saliva testing	
RNACNQ	Amylase	Negative	XJC9UK	Amylase	Negative
RVQ877	Amylase diffusion	Negative	XXUUMQ	Saliva test not performed on this item per laboratory policy - not required to test	
	Amylase mapping (Phadebas)	Negative	Y28ZFH	This laboratory does not offer saliva testing.	
T2V4RR	RSID	Negative	Y8GAQK	As per case scenario, this item is from the victim. Samples from the vaginal area can give a false positive result for amylase, so this test was not performed.	
T7LPCY	No currently validated procedures for saliva test.		Y9VLM6	SaligaE	Negative
TJBFHW	NT (Per Policy) - Item not eligible for saliva testing		YBDZQ6	Phadebas	Negative
TKLJQE	Amylase Diffusion	Negative	YN3XWA	Phadebas - Not Tested. Hemastix and AP positive tests place the assailant and victim together. No suggestion of oral rape.	
U3J3HU	Per Policy, not required to test				
UAYML8	ALS	Positive			
	RSID-Saliva	Negative			
UF9YH8	ALS	Negative			
	RSID	Negative			

TABLE 2 - Item 6

Saliva					
Webcode	Test Performed	Results	Webcode	Test Performed	Results
YQT4CZ	blood and semen found no indication in scenario of possible saliva				
YZDAA4	SaligaE	Negative			
ZCTE3K	RSID	Negative			
ZKJZ9J	RSID	Negative			
ZRYNXE	PHADEBAS (BLOTTING METHOD)	Negative			
ZYFE47	Amylase Diffusion	Negative			

TABLE 2 - Item 6

Other Body Fluid Types					
Webcode	Body Fluid Type - Test Performed	Results	Webcode	Body Fluid Type - Test Performed	Results
28BLGK	Biological Fluids - ALS	Negative			
2PQVGK	Sperm Search	Positive			
4BJLQG	Biological Fluids - ALS	Positive			
7M6VAR	ALS	Positive			
8YD3JW	ALS: crimelite screen for semen & saliva	Inconclusive			
AWHZRY	Semen - ALS	Negative			
BU8H89	Bodily Fluid - ALS	Negative			
KZ3R7Y	Biological Fluids - ALS	Negative			
MEWRED	Semen, Saliva, Sweat - ALS	Positive			
PKWW3E	Item also fibre taped				
QP6VRD	ALS	Negative			
RBKG8Q	ALS	Positive			
RNQTHJ	Forensic light source	Negative			
WM2X2R	Urine - Creatinine	Negative			
	Fecal - Micro	Negative			
	Fecal - Urobilinogen	Negative			
Y8WHUU	Biological Fluids - ALS	Negative			
YTUAPK	Bodily Fluids - ALS	Negative			

Additional Comments

TABLE 3

WebCode	Additional Comments
27CG2C	Items 1 and 4 - both gave positive presumptive tests for the presence of amylase and also positive for the presence of semen. Therefore the amylase could be salivary or pancreatic from the semen. Item 3 - phadebas tube test inconclusive due to heavy blood staining. Item 6 - two areas of outside surface of the gusset gave AP positive results but were negative for spermatozoa under microscopic exam. Associated area of AP positive staining on the inside surface of the gusset was positive for spermatozoa under microscopic exam and was bloodstained. None of these items smelled of urine or were visibly stained with urine/faeces therefore we did not carry out presumptive tests for urine/faeces as per our laboratory SOPs.
28BLGK	Saliva testing was not performed on any of the items because that test is not one performed by the laboratory.
2FHERZ	Human test was not performed on brown-stain from item 1-6 because according to our protocols it does not have to be run on a stain located in an intimate area
3P4YCN	This proficiency test consisted of too many items and was extremely time consuming to work. Proficiency tests should not consist of this many items when multiple tests need to be performed on most items. I feel strongly that the number of samples in this proficiency test should be cut in half. In addition to making the number of items less, I also feel that it would be better to simply ask analysts to test all items for all body fluids as the scenarios for these proficiency tests never provide enough information for analysts to make the best testing decisions possible. Because this particular scenario did not give any specifics other than that it was a violent sexual assault, an analyst could easily interpret this to mean you need to examine all items for all body fluids since you have no way of knowing what actually happened (e.g., Was the victim and/or suspect bleeding? Did the suspect lick the victim's abdominal area so saliva on swabs from the victim's stomach would be important?). Regarding Item 5 (the headpiece allegedly worn by the suspect) in particular, it would have been helpful to have more information to determine what testing should be performed. Are we trying to determine that this headpiece did indeed belong to the suspect or are we trying to find the victim on this item or a combination of both? A positive result with a presumptive test for blood was obtained from a reddish-brown stain on this item. I also tested this stain for semen and saliva based on the appearance of the stain and I did not swab for wearer DNA due to the stain being present (e.g., if the suspect's DNA profile is obtained from the stain, no additional examination may be necessary). I also did not collect swabs of the reddish-brown stain from Item 5 for DNA analysis due to this being a proficiency test. If this was an actual case, I would have contacted the agency to determine whether DNA analysis of the reddish-brown stain was necessary and whether we also need to determine the wearer of the headpiece via DNA analysis. I also would have asked for further information to determine if the examinations for semen and saliva were necessary. Due to the limited scenario provided, many analysts may be performing more testing on some items that is not necessary.
4BJLQG	We do not test for saliva @ [Laboratory]
4JVXRK	Item 5: The AP test was "inconclusive" due to the substrate color and reddish-brown stain affecting the ability to see a potential color change when the reagent was applied.
68YCXV	Following the case scenario and our lab's procedures for Item 1-5, after using an alternate light source and not detecting additional stains, I determined that further semen and saliva testing was not necessary. Following the case scenario and our lab's procedures, a confirmatory test for blood was not done on Item 1-6. The portal would not allow me to select the answers for multiple tests performed. So for certain body fluids requiring multiple tests, since I was only able give results for one of the tests, I selected the most confirmatory test.
6UWLGW	Item 1: RSID Amylase results are being considered inconclusive because the line was too light. I was unable to determine if the line was pink or grey. I ran the sample twice with the same result.
6ZGMFL	Saliva tests were not performed as the lab does not currently have a procedure or validated method.
7MNLCK	Agency has no approved testing for saliva.
8NUMDA	For items 1 and 2 this laboratory would not accept swabs for which the laboratory did not know the

TABLE 3

WebCode	Additional Comments
	source or type of body fluid needing to be tested nor would this laboratory perform all tests, presumptive and confirmatory, for every body fluid for which we have a screening test. All tests, presumptive and confirmatory, were performed on all the samples provided even though this is not done in casework typically to conserve samples or because the laboratory has knowledge of source of the body fluid being tested. For semen testing in this laboratory, we report seminal fluid and sperm results separately. No overall conclusion about semen is provided. This laboratory does not perform testing for saliva.
8PP7HW	Items have not been screened for semen - examination for semen is not carried out at the laboratory at which these items were examined ([Laboratory]). Items 5 and 6 were examined visually for the presence of blood-like staining - any blood-like stains were spot tested with LMG. No further examination for blood on these items was conducted.
9HLHJN	It may be more consistent for results between labs to remove the portion of the instructions that allows labs to process items by their policies. Analyzing all samples with all tests may be more accessible for both analyzing and grading the items.
AWHZRY	Items 5 and 6 were taped for removal of apparent hairs, fibers and debris. Clear tape into original packaging
B33FWW	Item 1: Saliva is not routinely tested in lab casework but tested only when the case circumstances specify. There is no case scenario for this item to suggest the presence of saliva. Items 3, 5 & 6: Not tested for Saliva: Scenario (location of recovery) does not suggest/indicate that saliva could be present on these items. These items are also non-intimate samples. Saliva is not routinely tested in sexual assault cases, only when specified in the circumstances.
BU8H89	This laboratory does not perform tests for saliva.
C4UYB4	For items 1 and 2: this laboratory would not accept swabs for which the source or type of body fluid needing to be tested was unknown nor would this laboratory perform all tests, presumptive and confirmatory, for every body fluid for which we have a screening test. All tests, both presumptive and confirmatory, were performed on the samples provided for items 1 and 2 even though this is typically not done in casework to conserve the sample or because the laboratory has knowledge of the source of the body fluid being tested. For semen testing in this laboratory, seminal fluid and sperm results are reported separately. No overall conclusion for semen is provided. For items 2, 3, and 5: results for seminal fluid and sperm were both negative. For items 1 and 4: results for seminal fluid and sperm were both positive. For item 6: results for sperm were positive. Results for seminal fluid would be reported as follows: Presumptive testing for seminal fluid was not performed due to the observed staining. Confirmatory testing for seminal fluid was not performed due to microscopic results. In this laboratory all samples are visually examined prior to any testing and all observations are documented. All samples included in this test were visually examined in addition to the testing performed on them. Saliva testing was not performed on any of the samples included in this test because it is not a service provided by this laboratory.
CGLJZ8	Item 6: Inconclusive result for AP due to red-brown color change from stain sample.
EB9GR9	As forensic experts at crime scene investigation and bloodstain pattern analysts of the Mobile Forensic Team of the [Laboratory], we are on-call for 24/7 and we act upon request by the [Country] police forces in [Country] on scenes of crime (homicide and murder cases). As a result, we don't carry out tests typically designed for the laboratory circumstances like the Phadebas (amylase) screening test or the AP screening test. The tests we do carry out are the TMB test, AP spot test (phosphatesmo KM) and the RSID Field Kits Semen & Saliva. If needed, we also have the RSID Blood field kit available, but because of both cost effectiveness and speed we still prefer the TMB test for initial testing for blood. Of course, we also use an ALS for screening purposes in the first place. This is the reason why we have not screened (except for ALS screening) both clothing items (5 and 6) for saliva and semen, because we don't use screening tests like the a-amylase and/or the PSA press test. These types of investigation in casework is normally carried out in a laboratory environment at the [Laboratory]'s department of human biological traces.
EFF3F6	Because Item 5 is black in color, the stain on Item 5 was photographed with an infrared camera to better visualize and observe the location of the stain. Weak fluorescence was also observed around the

TABLE 3

WebCode	Additional Comments
	stain on Item 5 when viewed with an alternate light source. No fluorescence was observed in or around the stain in Item 6 when viewed with an alternate light source.
ERDV3B	For Item 6 - A few sperm heads were seen microscopically. This was verified by a second examiner as per laboratory policy.
EXGQ74	Item 4 Saliva Phadebas Amylase Inconclusive Result Explanation: Amylase is an enzyme that is present in many body fluids in addition to saliva; however, it is present in saliva at higher concentrations. One body fluid that also contains amylase is semen, which was identified on the swabs in Item 4. Phadebas Amylase testing was weak positive for Item 4; therefore, saliva could be present in Item 4 based on the scenario. However, based on the strong semen stains identified in Item 4 and because amylase is present in semen, the Phadebas Amylase weak positive results could be due to the semen. Therefore, the overall conclusion for saliva in Item 4 is inconclusive.
FXB6NT	Item 1: No blood indicated, Semen identified, No saliva indicated. Item 2: No blood indicated, No semen identified, Saliva indicated. Item 3: Blood indicated, No semen identified, No saliva indicated. Item 4: No blood indicated, Semen identified, No saliva indicated. Item 5: Blood indicated, No semen identified, No saliva indicated. Item 6: Blood indicated, Semen identified, No saliva indicated.
G74RZB	For items 3-6, where there is no information to indicate there may be saliva present and no specific request for saliva testing by the submitting agency, I would not test for it. Also, if an oral assault had occurred, the most likely areas to find saliva would be on Items 4 and 6, and both of those are positive for semen, so we would presumably have male DNA present to compare to the suspects DNA.
G8XFEZ	Saliva testing not undertaken by examiner as not authorised in this type of testing. Item 2 - initial testing of swab material gave very weak positive result for both TMB and HT. Not able to be reproduced with additional testing. Assumed that possibly trace amount of blood present in whatever material may be present. Unable to undertake saliva testing.
GPVJEP	On item 2 there were NECs on the H&E slide - possible indication but not confirmatory. On item 6 the bloodstain was coincident with the semen stain in the gusset.
GV43GN	Item 1: Semen with spermatozoa. Item 2: Possible saliva (elevated amylase). Item 3: Apparent blood. Item 4: Semen with spermatozoa. Item 5: Apparent blood. Item 6: Semen with spermatozoa mixed with apparent blood.
GZXZM8	No saliva tests were performed, as the lab does not currently have a procedure or validation method.
H9ADNB	Items 1, 4, and 6 were all positive for seminal fluid/semen. These samples were also tested for saliva and the RSID saliva test kit revealed a weak positive signal for each of these semen positive samples. Item 2 was very clearly a strong rapid positive[sic] saliva result on the test kit. Based on the weak positive results for 1, 4 and 6, it suggests trace amounts of saliva in these samples. Also, it would have been helpful to know the nature of the assault (vaginal or oral or both) in determining what types of testing to do.
HDYUYF	Saliva tests are not performed in our laboratory. The product insert for the ABACard Hematrace test states that, "Due to the extreme sensitivity of the test, trace levels of hemoglobin might be detected occasionally in body fluid samples other than blood (e.e. urine, semen[,] stool, saliva, vaginal fluid, perspiration). Our lab protocol does not require confirmatory testing with Hematrace if the presumptive phenolphthalein test gives a negative result. However, for the purposes of the proficiency test, every available presumptive and confirmatory test was performed on each item provided. Our current lab report wording for a phenolphthalein negative, Hematrace positive sample would be as follows: "The presumptive test for blood provided negative results. The confirmatory test for blood provided positive results." Item 1 and Item 4 were both phenolphthalein negative and Hematrace positive. Therefore, according to our current method of reporting, both Item 1 and Item 4 have been reported on this proficiency test as being positive for blood.
J23ATA	Item 1: General Screening- Saliva test not done. Saliva will only be tested if item tested negative for both blood and semen. Item 3 to 6: Saliva tests not done. Saliva not routinely tested for in sexual assault cases. It is only performed when specifically requested for by the client.
JJM2UU	Items 1 and 4: Although a moderate positive Phadebas reaction was obtained from these swabs, the

TABLE 3

WebCode	Additional Comments
	presence of semen could provide some amylase activity. We cannot say that the amylase activity is necessarily due to the presence of saliva. Also no cells were found. Therefore, our overall conclusion for the presence of saliva on Items 1 & 4 is inconclusive.
JNFWZ9	Scenario needs to be more detailed: 1. What sex acts are alleged? 2. Was Item 5 collected from the suspect? Or was it left behind at scene? 3. How was Item 5 worn - as a hair cover? Or as a face mask? These details would assist in determining what types of body fluids to look for on the different items. For Item 5, I made the following assumptions/guesses - it was worn as a hair cover, not a mask - therefore, not examined for saliva. It was left behind at scene - therefore, swabs were collected for possible wearer DNA. For Items 4 & 6 - I assumed/guessed that oral sodomy was not involved - therefore these items were not examined for saliva.
JRWLWQ	Saliva tests are not routinely performed unless specifically requested by the investigating agency.
JYWP6A	For Item 1 (General screening): Saliva Test is not done if the item is tested positive for either blood or semen. Saliva test will be performed for the item if both blood and semen tests are negative. For Item 3 to Item 6 (Case-specific screening): Saliva not routinely tested for in sexual assault cases. Saliva test is only performed when specifically requested for by the client.
K6XJEV	Item 3: Saliva test was inconclusive because water was red/brown in color and could not tell if water turned blue or stayed clear.
K7UYQR	Note* all samples stored at room temperature when received for examination. ITEM 3 SWAB - Inconclusive for saliva due to sample colour as heavily bloodstained.
KZ3R7Y	Saliva tests not performed at our laboratory.
LXURT6	A positive overall conclusion for semen was obtained from the stain on Item 6. It is likely that this stain did not exhibit fluorescence with exposure to an ALS because of quenching from the blood that was also present.
M9CTN6	Too many samples. Reported scenario was too vague and poorly worded - Item 3 - was "outside wall of a washing machine" meant to be "exterior side of a washing machine"? Item 4 - was "victim's stomach" meant to be "victim's abdomen" or "victim's torso"? Item 5 - was it recovered at the scene? Scenario of "violent sexual assault of a woman" can mean many different acts - it is a very vague description.
MED8A7	Scenario didn't provide enough information in order to make testing decisions.
MEWRED	Item 5 (black headpiece): The inside of this item was swabbed using 2 sterile wet/dry swabs for possible wearer's DNA. Approximately half of each swab was sampled and combined into a sterile tube, and stored for further testing if needed.
MQAX88	For item 1 : General screening - Saliva test will only be conducted if the item tested negative for both blood and semen. For items 3 to 6: Case-specific screening. Saliva test not routinely conducted in sexual assault cases. It is only conducted when specifically requested for by the client.
MY39JM	*saliva testing not performed in this lab.
MZHQHM	Saliva testing not routinely performed unless specifically requested by investigating agency.
N3ZCLU	Note: Phadebas test results for Items 1 and 4 were inconclusive. A picture is available of these test results which shows blue color in the supernatant of each sample. In our laboratory these items would not have been tested for saliva as they were already found to be positive for semen.
PJ2EYR	Item 5 was also swabbed for the possible presence of epithelial cells. Swabs were collected from the interior surface of the black head piece. Microscopic examination was not performed in order to conserve the sample for future analysis. Item 6 was also swabbed for the possible presence of epithelial cells. Swabs were collected from the exterior waistband. Microscopic examination was not performed in order to conserve the sample for future analysis. Item 7 labeled "Control Swabs" was not opened or analyzed in accordance with this laboratory's protocols. For items 1 and 2 this laboratory would not accept swabs from which the laboratory did not know the source or type of body fluid needing to be tested nor would this laboratory perform all tests, presumptive and confirmatory, for every body fluid for which we have a screening test. All tests, presumptive and confirmatory, were performed on all the

TABLE 3

WebCode	Additional Comments
	samples provided even though this is not done in casework typically to conserve sample or because the laboratory has knowledge of source of the body fluid being tested. For semen testing in this laboratory, we report seminal fluid and sperm separately. No overall conclusion about semen is provided.
PMLJF3	Item 5: Examination with an alternate light source revealed three small stains with very pale fluorescence on the left front side of the headpiece near the tie. Neither semen nor saliva was detected in these stains.
PNF7UR	Item 1 was weakly positive with a HemaTrace test. The result was negative when repeated. Item 2 was weakly positive with a PSA test. The result was negative when repeated. The Phadebas (amylase) test is used as a screening test and a positive result is indicative of the presence of saliva only. Two samples of the stain on item 6 were examined microscopically for the presence of sperm and none were seen.
QADXPB	Item 3 was hematrax tested - however in actual casework this may not be done unless there is a requirement to confirm that the source of the blood is of human origin eg. DNA profiling results are not as expected or animals are known to have bled at the scene. Therefore hematrax testing was not carried out on all hemastix positive items. For swabs microscopy is carried out first to search for semen then AP testing and P30 testing if the result is negative. For items of clothing etc. AP testing is done first when searching for semen to "map" areas of interest for further examination/testing.
R9VBPE	The starch-iodine test is not routinely performed at the [Laboratory] unless the case background specifically declares that oral contact, or biting etc. occurred. Thus all submitted samples simply underwent AP and LMG testing for screening with application of the starch iodine solely to item 4 where the act is described as violent but not sufficiently specific. Microscopy at the lab is used as an alternative to p30 testing, which in this case p30 kit was not presently available so AP results were confirmed via microscope at 1000x (LEICA DM). The lab's present protocol requires documentation of 6 heads or a single intact spermatozoon in order to call a test positive; no results were obtained via microscopy for item 6.
RNACNQ	Amylase testing is considered a screening test for the presence of saliva in this laboratory. A positive result is reported as 'an indication of the presence of saliva'. The presence of nucleated squamous epithelial cells together with the detection of amylase activity is reported as 'highly indicative of the presence of saliva'.
RNQTHJ	Hexagon OBT1 is a rapid immunographic test for confirming the presence of human bloodtraces. Items 1, 2 and 4: No blood-like stains were observed when visually examined, therefore no further testing. Item 2: A DNA-profile was obtained which indicates human secretion. Item 5: Routine examination of headpieces at the laboratory include search for hairs. Tests for the presence of saliva are not routinely used at the laboratory.
T7LPCY	No currently validated procedures for saliva test.
TX4EYX	Saliva Tests - my lab has no current validated tests or methods for saliva identification, so none of the items were tested for saliva.
UAZKRV	Saliva Tests are unavailable at this time for Items 1-6.
UF9YH8	Control swabs not tested. Could not remove 'Inc.' mark from Item 2-Other Tests.
V4BH9W	Item 5: Our procedures allows us to report "No stains consistent with semen were observed" when we get a negative result w/ ALS. That is what I'm reporting on my internal report. However, I'm limited by your answer key and will be answering "Neg" for semen even though that is a stronger conclusion than I would like to report. Based on the limited scenario, I wouldn't test the red-brown stain w/ additional semen tests.
V9Z287	CLOTHING ITEMS NOT TESTED WITH PHADEBAS PAPER AS INSTRUCTED BY UNIT MANAGER
WF2Q7P	The blood staining observed on the panties (item 6) was very light. The KM reaction was faint.
WM2X2R	#1 Swab - Seminal fluid (presumptive) and spermatozoa detected; no blood, saliva or urine identified. #2 Swab - Saliva detected (presumptive) detected[sic]; no semen or urine identified. #3 Swab (washing machine) - Blood (human) detected (presumptive); no semen, saliva, urine or fecal material identified. #4 Swab (victim's stomach) - seminal fluid (presumptive) and spermatozoa detected; no

TABLE 3

WebCode	Additional Comments
	blood, saliva or urine identified. #5 Headpiece (do-rag)-suspect - Human blood detected (presumptive); no seminal fluid or saliva detected. #6 Underwear (victim) - Human blood detected (presumptive) and spermatozoa detected; no saliva, urine or feces detected. Contact the laboratory if further analysis is needed (i.e. submit victim/suspect buccal swabs for DNA comparison).
XHZBMU	No saliva tests were performed as the laboratory does not currently have a procedure or validated method. Items 1, 2, and 4 were not tested for blood because a red stain was not observed on these items.
XJC9UK	Control swabs submitted within the packaging are added as Item #7 in LIMS & on evidence receipt. See communication log. [Communication log was not provided by participant.]
Y28ZFH	1. For items 1 and 2 this laboratory would not accept swabs for which the laboratory did not know the source or type of body fluid needing to be tested nor would this laboratory perform all tests, presumptive and confirmatory, for every body fluid for which we have a screening test. 2. For items 1 and 2 all tests, presumptive and confirmatory, were performed on all the samples provided even though this is not done in casework typically to conserve the samples or because the laboratory has knowledge of source of the body fluid being tested. 3. For semen testing in this laboratory, we report seminal fluid and sperm results separately. No overall conclusion about semen is provided.
Y8GAQK	Testing for urine and fecal material was not performed on Items 3, 4, 5 and 6 due to the case scenario given for these samples. Items 1 and 4: Chemical testing for the presence of saliva was inconclusive due to very weak light blue colored results found during Phadebas amylase testing, giving a very weak positive result for this test. Spermatozoa were identified in these two samples (Items 1 and 4) and literature has shown that seminal material contains low levels of amylase, possibly explaining the very weak results obtained for the Phadebas amylase test.
Y8WHUU	Saliva testing is not performed at this laboratory.
YBDZQ6	Semen and blood on item 6 - mixed / coincidental staining. Swabs not tested for saliva - but supernatant from extraction could have been - Phadebas tube test not accredited at this laboratory.
YLWGGQ	Saliva test was not performed because the laboratory does not have a validated procedure for saliva.
YMCGYU	Item 1: General Screening- Saliva test not done. Saliva will only be tested if item tested negative for both blood and semen. Item 3 to 6: Case-specific Screening- Saliva not routinely tested for in sexual assault cases. It is only performed when specifically requested for by the client.
YTUAPK	Saliva tests not performed by lab.
ZRYNXE	All 6 items were visually examined for the presence of urine and faeces. Nothing of significance was found in relation to faeces for all 6 items. There was a faint yellow discolouration present on the swabs from items 1 and 4. These swabs were tested for creatinine, a constituent present in large amounts in urine but also present in other body fluids such as sweat, saliva and semen. Both swabs gave a faint positive reaction. These two swabs (items 1 and 4) were also positive for semen. There could have been traces of urine present in the semen or the test could be detecting creatinine in the semen. Item 3 gave a positive reaction to a chemical test for blood (KM). The swab was heavily blood stained discolouring the contents of the tube. The result for the saliva tube method for this swab was inconclusive.

Appendix: Data Sheet

Collaborative Testing Services ~ Forensic Testing Program

Test No. 14-579: Body Fluid Identification Test

DATA MUST BE RECEIVED BY November 10, 2014 TO BE INCLUDED IN THE REPORT

Participant Code:

WebCode:

Accreditation Release Statement

CTS submits external proficiency test data directly to ASCLD/LAB and ANSI-ASQ NAB/FQS. Please select one of the following statements to ensure your data is handled appropriately.

This participant's data is intended for submission to ASCLD/LAB and/or ANSI-ASQ NAB/FQS. (Accreditation Release section on the last page must be completed and submitted.)

This participant's data is **NOT** intended for submission to ASCLD/LAB or ANSI-ASQ NAB/FQS.

Scenario (this scenario applies only to Items 3 - 6):

Police are investigating the violent sexual assault of a woman who was attacked in her apartment complex's laundry room. The following items are being submitted to your department for examination. One envelope containing four duplicate swabs (Item 3) which were taken from a stain found on the outside wall of one of the washing machines. One envelope containing four duplicate swabs (Item 4) which were taken from the victim's stomach. A black headpiece (Item 5) allegedly worn by the suspect during the incident. A pair of red women's underwear (Item 6) which had been worn by the victim during the incident.

Items Submitted (Sample Pack BF12):

- Item 1 Swab (provided as 4 duplicates) - intended for general screening for body fluids tested in normal laboratory casework; scenario not applicable
- Item 2 Swab (provided as 4 duplicates) - intended for general screening for body fluids tested in normal laboratory casework; scenario not applicable
- Item 3 Swab (provided as 4 duplicates) - taken from a stain found on the outside wall of washing machine
- Item 4 Swab (provided as 4 duplicates) - taken from the victim's stomach
- Item 5 One black headpiece
- Item 6 One pair of red women's underwear
- 1 package of control swabs (unused)

Online Data Entry

Visit www.cts-portal.com to enter your proficiency test results online. If you have any questions please do not hesitate to contact CTS.

Please return all pages of this data sheet.

Page 1 of 10

Instructions:

Please indicate the overall conclusion for each body fluid using the check boxes provided. List the tests performed and the results of each test. Item descriptions (as applicable) are provided at the top of each Item page.

For the purposes of this test, all swabs within an Item have been prepared as duplicates in order to provide sufficient sample quantity. A package of unused swabs is also provided for control purposes.

Please use the abbreviations listed in this response key to fill in the Screening Test tables. This is not an all inclusive list of tests, and should not be used to determine what tests should be performed. TESTS NOT ON THIS LIST MAY BE USED FOR SCREENING.

Test	Abbreviation	Test	Abbreviation
Acid Phosphatase	AP	Ortho-tolidine	O-tol
Alternate Light Source	ALS	Phenolphthalein-Tetramethyl benzidine	PTMB
Kastle Meyer	KM	Prostate Specific Antigen	PSA
Leucomalachite Green	LMG	Rapid Stain Identification	RSID
Microscopic	Micro	Tetramethyl benzidine	TMB

Example		Pos./Ind.	Neg.	Inc.	NT*
Blood Tests					
What is the overall conclusion for blood?		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s)</u> Performed	HemaTrace	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	KM	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Pos./Ind. - Positive/Indicative; Neg. - Negative; Inc. - Inconclusive; NT - Not Tested

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Item 1: Cotton tipped swabs with white plastic handle

Scenario: Scenario is not applicable for this Item. Utilizing your laboratory policies, please examine this item for body fluids tested in normal laboratory casework.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for blood?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s)</u> <u>Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for semen?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s)</u> <u>Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for saliva?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s)</u> <u>Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Test Results		
Body Fluid	Test(s) Performed	Pos./Ind.	Neg.	Inc.
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Please return all pages of this data sheet.

Item 2: Cotton tipped swabs with pink plastic handle

Scenario: Scenario is not applicable for this Item. Utilizing your laboratory policies, please examine this item for body fluids tested in normal laboratory casework.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for blood?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s) Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for semen?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s) Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for saliva?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s) Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Test Results		
Body Fluid	Test(s) Performed	Pos./Ind.	Neg.	Inc.
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Please return all pages of this data sheet.

Item 3: Cotton tipped swabs with blue plastic handle

Scenario: Four duplicate blue handled swabs taken from a stain found on the outside wall of a washing machine.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for blood?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s)</u> <u>Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for semen?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s)</u> <u>Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for saliva?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s)</u> <u>Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Test Results		
Body Fluid	Test(s) Performed	Pos./Ind.	Neg.	Inc.
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Please return all pages of this data sheet.

Item 4: Cotton tipped swabs with yellow plastic handle

Scenario: Four duplicate yellow handled swabs taken from the victim's stomach.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for blood?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s) Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for semen?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s) Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for saliva?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s) Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Test Results		
Body Fluid	Test(s) Performed	Pos./Ind.	Neg.	Inc.
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Please return all pages of this data sheet.

Item 5: One black headpiece

Scenario: The item of evidence is a black headpiece allegedly worn by the suspect during the incident.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for blood?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s) Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for semen?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s) Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for saliva?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s) Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Test Results		
Body Fluid	Test(s) Performed	Pos./Ind.	Neg.	Inc.
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Please return all pages of this data sheet.

Item 6: One pair of red women's underwear

Scenario: The item of evidence is a pair of red women's underwear worn by the victim during the attack.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for blood?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s) Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for semen?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s) Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for saliva?		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<u>Test(s) Performed</u>	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Test Results		
Body Fluid	Test(s) Performed	Pos./Ind.	Neg.	Inc.
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Please return all pages of this data sheet.

RELEASE OF DATA TO ACCREDITATION BODIES

The following Accreditation Releases will apply only to:

Participant Code:

WebCode:

for Test No. **14-579: Body Fluid Identification**

This release page must be completed and received by **November 10, 2014** to have this participant's submitted data included in the reports forwarded to the respective Accreditation Bodies.

ASCLD/LAB RELEASE

If your lab has been accredited by ASCLD/LAB and you are submitting this data as part of their external proficiency test requirements, have the laboratory's designated individual complete the following.

The information below must be completed in its entirety for the results to be submitted to ASCLD/LAB.

ASCLD/LAB Legacy Certificate No. _____ ASCLD/LAB International Certificate No. _____

Signature _____ Date _____

Laboratory Name _____

Location (City/State) _____

ANSI-ASQ NAB/FQS RELEASE

If your laboratory maintains its accreditation through ANSI-ASQ NAB/FQS, please complete the following form in its entirety to have your results forwarded.

ANSI-ASQ NAB/FQS Certificate No. _____

Signature and Title _____ Date _____

Laboratory Name _____

Location (City/State) _____

Accreditation Release

Return Instructions

Please submit the completed Accreditation Release at the same time as your full data sheet. See Data Sheet Return Instructions on the previous page.

*Questions? Contact us 8 am-4:30 pm EST
Telephone: +1-571-434-1925
email: forensics@cts-interlab.com*

Please return all pages of this data sheet.

Page 10 of 10