



Body Fluid Identification

Test No. 23-5781 Summary Report

Each participant received a sample set consisting of six items for examination. Participants were asked to examine each item and report the body fluid(s) identified. Data were returned from 224 participants and are compiled into the following tables:

Report Contents:	Page:
<u>Manufacturer's Information</u>	2
<u>Summary Comments</u>	3
<u>Key for Screening Tests Used</u>	4
<u>Table 1: Screening Conclusions & Test(s) Performed</u>	5
<u>Table 2: Additional Comments</u>	144
<u>Appendix: Data Sheet</u>	

This report contains the data received from the participants in this test. Since these participants are located in many countries around the world, and it is their option how the samples are to be used (e.g., training exercise, known or blind proficiency testing, research and development of new techniques, etc.), the results compiled in the Summary Report are not intended to be an overview of the quality of work performed in the profession and cannot be interpreted as such. The Summary Comments are included for the benefit of participants to assist with maintaining or enhancing the quality of their results. These comments are not intended to reflect the general state of the art within the profession.

Participant results are reported using a randomly assigned "WebCode". This code maintains participant's anonymity, provides linking of the various report sections, and will change with every report.

Manufacturer's Information

Each sample set consisted of six separate items for testing. Two items contained non-scenario related swabs (Items 1 and 2), two items contained swabs that related to the provided scenario (Items 3 and 4), and two items (Items 5 and 6) contained complete items of evidence that also related to the provided scenario. Participants were requested to examine Items 1 and 2 for body fluids that they typically encounter in routine casework examinations. For Items 3 through 6, the participants were requested to examine the items using the information provided by the scenario for guidance.

SAMPLE PREPARATION: The stain on Item 1 (Swabs: white plastic handle) was prepared using human saliva. The stain on Item 2 (Swabs: pink plastic handle) was prepared using human blood. The stain on Item 3 (Swabs: blue plastic handle) was prepared using human blood. The stain on Item 4 (Swabs: yellow plastic handle) was prepared using a suspension of human semen and saliva. The stain on Item 5 (tube top) was prepared using a suspension of human blood and semen. The stain on Item 6 (tampon) was prepared using human semen. For the swab items where semen was used, the semen was mixed in a 1:1 ratio with TAE buffer. For the evidence items where semen was used, the semen was not mixed with TAE buffer. In all items where a mixture is present, the two body fluids were mixed in a 1:1 ratio. Stains from different sources were prepared at separate times and were packaged once thoroughly dried.

SAMPLE SET ASSEMBLY: For each sample set, all Items (1-6) were packaged into separate envelopes and then placed together in a pre-labeled sample set envelope and sealed. The sealed sample set envelopes were then packaged in pre-labeled heat seal envelopes and sealed. This process was repeated until all of the sample sets were prepared.

VERIFICATION: All predistribution laboratories confirmed the presence of the body fluid(s) that comprised each item.

Sample Preparation Schedule

<u>Item</u>	<u>Substrate</u>	<u>Body Fluid</u>	<u>Volume</u>
1	4 Swabs: white plastic handle	Saliva	100µL per swab
2	4 Swabs: pink plastic handle	Blood	100µL per swab
3	4 Swabs: blue plastic handle	Blood	100µL per swab
4	4 Swabs: yellow plastic handle	Semen & Saliva Mixture	100µL per swab
5	Tube Top	Blood & Semen Mixture	200µL
6	Tampon	Semen	200µL

Completed sample sets were stored at -20°C until shipment on March 20, 2023, following completion of the verification stage.

Summary Comments

The Body Fluid Identification test was designed to allow participants to assess their proficiency in the identification of stains that may contain dried body fluids. Each participant received six items for testing; four items were swabs and two were whole evidence items. A scenario was not provided for Items 1 and 2; participants were requested to test for body fluids that they typically encounter in casework. A scenario was provided for Items 3 - 6, and participants were asked to treat these items as they would normal casework, testing for the fluids they deemed necessary after reviewing the scenario.

Item 1 consisted of four swabs that were prepared using saliva. The majority of responding participants reported the expected results: "Negative" for blood (100%), "Negative" for semen (99%) and "Positive/Indicative" for saliva (100%).

Item 2 consisted of four swabs that were prepared using blood. The majority of responding participants reported the expected results: "Positive/Indicative" for blood (100%), "Negative" for semen (99%), and "Negative" for saliva (97%).

Item 3 consisted of four swabs that were prepared using blood. The majority of responding participants reported the expected results: "Positive/Indicative" for blood (100%), "Negative" for semen (99%), and "Negative" for saliva (99%).

Item 4 consisted of four swabs that were prepared using a mixture of semen and saliva. The majority of responding participants reported the expected results: "Negative" for blood (96%), "Positive/Indicative" for semen (100%), and "Positive/Indicative" for saliva (99%).

Item 5 consisted of a tube top that was prepared using a mixture of blood and semen. The majority of responding participants reported the expected results: "Positive/Indicative" for blood (100%), "Positive/Indicative" for semen (97%), and "Negative" for saliva (98%).

Item 6 consisted of a tampon that was prepared using semen. The majority of responding participants reported the expected results: "Negative" for blood (96%), "Positive/Indicative" for semen (100%), and "Negative" for saliva (98%).

Key for Screening Tests Used

Participants were asked to use, where possible, the following chart of abbreviated screening test names. This was not an all inclusive list and was not designed to determine what tests should be performed. Participants were advised that tests not on this list may be used for screening.

<u>Test</u>	<u>Abbreviation</u>
Acid Phosphatase	AP
Alternate Light Source	ALS
Kastle Meyer	KM
Leucomalachite Green	LMG
Microscopic	Micro
Ortho-tolidine	O-tol
Phenolphthalein-Tetramethyl benzidine	PTMB
Prostate Specific Antigen	PSA
Rapid Stain Identification	RSID
Tetramethyl benzidine	TMB

Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
2289QN	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
22BP8V	KM <i>Negative</i> HemaTrace <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Microscopic <i>Negative</i> p30-ABACard <i>Negative</i> Conclusion: Negative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
2G2TAX	Not tested due to absence of red-brown staining. Conclusion: Not Tested	Acid Phosphatase <i>Negative</i> p30 <i>Negative</i> Microscopy <i>Negative</i> Conclusion: Negative	SALigAE <i>Positive/Indicative</i> Microscopy <i>Negative</i> Conclusion: Positive/Indicative
<u>Other Tests Performed</u> Urine - Urea Test <i>Negative</i>			
2Q8KTZ	LMG <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
2WBMVF	Visual Examination <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	RSID Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
2Y2QL2	LMG <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
34RKYT	physical exam = no red-brown = no testing Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
38HFLL	BLUESTAR (LUMINOL) <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
<u>Other Tests Performed</u> CRIME LITE <i>Inconclusive</i>			
3C2TKZ	No visible staining - not tested Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
3GUVNZ	No red/brown staining on swabs Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
3VAJDJ	KM - Hemophan (paper test for Blood erythrocytes, haemoglobin) by Erba Lachema - non-specific test for Blood <i>Negative</i> FOB - cassette test for determination of faecal occult Blood by NADAL - specific test for Human Blood <i>Negative</i> Conclusion: Negative	AP - Phosphatesmo (paper test for the determination of acid phosphatase) by Macherey-Nagel <i>Negative</i> PSA - cassette test for prostate specific antigen by SERATEC <i>Negative</i> Conclusion: Negative	RSID - Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
47AY3T	O-tol <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
4DB9ED	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
4DP4V9	Molecular Serology <i>Negative</i> Conclusion: Negative	Molecular Serology <i>Negative</i> Conclusion: Negative	Molecular Serology <i>Positive/Indicative</i> Conclusion: Positive/Indicative
4GQCZH	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> RSID <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
4R82UG	Conclusion: Not Tested	ALS <i>Positive/Indicative</i> PSA <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
4W3QGF	NO VISIBLE STAINING OBSERVED Conclusion: Not Tested	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: No Response	Phadebas <i>Positive/Indicative</i> Conclusion: No Response
	<u>Other Tests Performed</u> Visual / ODOUR TEST FOR POSSIBLE URINE / FAECES <i>Negative</i>		
4XR299	Swabs were not red in color Conclusion: Not Tested	AP <i>Negative</i> Micro Sperm <i>Negative</i> p30 <i>Negative</i> Conclusion: Negative	Test no longer performed Conclusion: Not Tested
62AAYQ	O-tol <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
62P4BF	Conclusion: Not Tested	ALS <i>Positive/Indicative</i> PSA <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
6AVWHW	TMB	Negative	AP	Negative	Phadebas tube/paper tests (Amylase)	Positive/Indicative
			Micro	Negative	Micro	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
6MWPPY	KM	Negative	AP	Negative	RSID	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
6QXKKV	KM	Negative	AP	Negative	RSID	Positive/Indicative
			Micro	Negative		
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
6R22DE	KM	Negative	AP	Negative	Phadebas	Positive/Indicative
			Micro	Negative		
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
6VNAP7	Visual	Negative	ALS	Negative	ALS	Negative
	KM	Negative	AP	Negative	Amylase Diffusion	Positive/Indicative
			Micro	Negative		
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
6ZW78F	Molecular Serology	Negative	Molecular Serology	Negative	Molecular Serology	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
6ZW9RD	Colorless		Acid Phosphatase (AP)	Negative	Lab does not test for Saliva	
			Sperm Search	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
7ARU7X	Micro	Negative	AP	Negative	RSID	Positive/Indicative
	KM	Negative				
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
7LHP24	Visual	Negative	ALS (positive only at the bottom of the swabs)	Positive/Indicative	Amylase Diffusion	Positive/Indicative
	KM	Negative	AP (top/bottom)	Negative		
			Micro	Negative		
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
7N8RBX	KM	Negative	Micro	Negative	[Laboratory] does not conduct Saliva tests.	
			PSA	Positive/Indicative		
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
7ULBZK	OT	Negative	Serated PSA	Negative	Amylase Diffusion	Positive/Indicative
	Seratec HemDirect	Negative	AP	Negative		
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
7YKAAU	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> PSA <i>Negative</i> RSID-Semen <i>Negative</i> Sperm Hy-Liter PI <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
82LQ3H	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
83UJJM	HemaTrace <i>Negative</i> Conclusion: Negative	PSA <i>Negative</i> Conclusion: Negative	SALigAE <i>Positive/Indicative</i> Conclusion: Positive/Indicative
8ERLLK	KM <i>Negative</i> HemaTrace <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	The [Laboratory] does not conduct forensic serological tests for Saliva. Conclusion: Not Tested
8RNLK4	Visual Examination <i>Negative</i> KM <i>Negative</i> HemaTrace <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Amylase Diffusion <i>Positive/Indicative</i> Conclusion: Positive/Indicative
8RT27A	KM <i>Negative</i> Conclusion: Negative	Brentamine- AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
8TTFRM	No reddish-brown staining observed Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	RSID Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
	<u>Other Tests Performed</u> Body Fluid - Test Performed: ALS <i>Negative</i>		
8VWJK6	No Blood like stains could be observed Conclusion: Not Tested	Kind's reagent <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
8Y4DL2	Blood testing not conducted based on observations (no Visual indication of possible Blood) Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
8YL7UQ	KM <i>Negative</i> HemaTrace <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Micro <i>Negative</i> p30-ABAcad <i>Negative</i> Conclusion: Negative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
8YVTV7	Rapid test for the detection of Human hemoglobin for use in legal medicine (Hexagon OBTI) <i>Negative</i>	<i>Negative</i>	Sperm identification by Contrast Staining Christmas tree <i>Negative</i> Rapid Stain Identification of Human Semen (RSID) <i>Negative</i>	<i>Negative</i>	Seratec – Amylase Test <i>Positive/Indicative</i>	<i>Positive/Indicative</i>
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
946QJ6	KM <i>Negative</i>	<i>Negative</i>	AP <i>Negative</i>	<i>Negative</i>	Amylase <i>Positive/Indicative</i>	<i>Positive/Indicative</i>
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
963BVP	KM <i>Negative</i>	<i>Negative</i>	AP <i>Negative</i>	<i>Negative</i>	Phadebas <i>Positive/Indicative</i>	<i>Positive/Indicative</i>
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
	<u>Other Tests Performed</u>					
	Urine - Creatinine <i>Negative</i>					
98U8JJ	Micro <i>Negative</i> KM <i>Negative</i>	<i>Negative</i>	AP <i>Negative</i> Micro <i>Negative</i>	<i>Negative</i>	RSID <i>Positive/Indicative</i>	<i>Positive/Indicative</i>
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
9GXGNN	Micro <i>Negative</i> KM <i>Negative</i>	<i>Negative</i>	FBB <i>Negative</i> Micro <i>Negative</i>	<i>Negative</i>	RSID <i>Positive/Indicative</i>	<i>Positive/Indicative</i>
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
9H2PM2	KM <i>Negative</i>	<i>Negative</i>	AP <i>Negative</i>	<i>Negative</i>		
	Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested	
9H9FT6	Visual <i>Negative</i> O-tol <i>Negative</i>	<i>Negative</i>	Visual <i>Negative</i> ALS <i>Inconclusive</i> AP <i>Negative</i>	<i>Negative</i>		
	Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested	
9HBZPK	No RBS observed.		PSA <i>Negative</i> Micro <i>Negative</i>	<i>Negative</i>	[Laboratory] does not routinely test for Saliva. Item was preserved.	
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
9KDWA9			AP <i>Negative</i>	<i>Negative</i>		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
9KH9X6	Bluestar Obti <i>Negative</i> LMG <i>Negative</i>	<i>Negative</i>	Seratec PSA Semiquant <i>Negative</i>	<i>Negative</i>	Phadebas® Amylase Test Tablets <i>Positive/Indicative</i>	<i>Positive/Indicative</i>
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
9T4NET	KM <i>Negative</i>	<i>Negative</i>				
	Conclusion: Negative		Conclusion: Not Tested		Conclusion: Not Tested	

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
9UXQ9H	NO RBS		ALS Negative AP Negative		Amylase PLATE Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Positive/Indicative	
A4FUMV	KM Negative		not in our scope		not in our scope	
	Conclusion: Negative		Conclusion: Not Tested		Conclusion: Not Tested	
A4QTAY	Visual Negative		AP Negative		RSID-Saliva Positive/Indicative	
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
AAEF3J	KM Negative		AP Negative PSA Negative Micro Negative		not testing for Saliva	
	Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested	
AANRDT	No red/brown discoloration		Acid Phosphatase Negative		RSID Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Positive/Indicative	
AD6FYL	Not tested due to no red-brown staining present.					
	Conclusion: Not Tested		Conclusion: Not Tested		Conclusion: Not Tested	
ADC2PF	Visual Negative KM Negative		ALS Negative p30 Negative			
	Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested	
ADPV82			AP Negative			
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
AHKH3B	KM Negative		AP Negative PSA Negative Micro Negative			
	Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested	
AJHZYK	not tested due to no visible RBS on swabs		AP Negative		RSID Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Positive/Indicative	
APAZLT	KM Negative		AP Negative Micro Negative		RSID Positive/Indicative	
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
AWH3H9	HemaTrace (hHb) Negative		PSA Negative		is not within the scope of this laboratory, we only perform identification of human blood and semen (P30 enzyme)	
	Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested	

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
BAYGU7	No red-brown staining Conclusion: Not Tested	AP Negative p30 Negative Microscopy Negative Conclusion: Negative	Amylase radial Diffusion Positive/Indicative Microscopy Positive/Indicative Conclusion: Positive/Indicative
BKPKDT	LMG Negative Conclusion: Negative	AP Negative Conclusion: Negative	Conclusion: Not Tested
BPVRLH	LMG Negative Conclusion: Negative	AP Negative Conclusion: Negative	Conclusion: Not Tested
BTGR6Q	No rbs noted to swabs Conclusion: Not Tested	AP Negative Conclusion: Negative	RSID Positive/Indicative Conclusion: Positive/Indicative
BXGUFK	Swabs were white in color Conclusion: Not Tested	AP Negative Micro Negative Conclusion: Negative	Amylase Diffusion Positive/Indicative Conclusion: Positive/Indicative
BXPWFJ	Phenolphthalein Negative Conclusion: Negative	p30 Negative Micro Negative Conclusion: Negative	Phadebas Positive/Indicative Conclusion: Positive/Indicative
CFMNR6	Conclusion: Not Tested	ALS Positive/Indicative PSA Negative Conclusion: Negative	Conclusion: Not Tested
CHEUWG	KM Negative Conclusion: Negative	AP Negative Conclusion: Negative	Conclusion: Not Tested
CLAAZR	FOB Negative Hemophan Negative Conclusion: Negative	RSID Semen Negative PSA Negative Conclusion: Negative	RSID Saliva Positive/Indicative Conclusion: Positive/Indicative
CLWA6M	LMG Negative Conclusion: Negative	AP Negative Conclusion: Negative	Conclusion: Not Tested
CQ7C8H	Conclusion: Not Tested	Micro Negative p30 Negative Conclusion: Negative	Conclusion: Not Tested
CTEPNK	no dry red stain observed Conclusion: Not Tested	AP Negative PSA Negative Micro Negative Conclusion: Negative	Phadebas Positive/Indicative Conclusion: Positive/Indicative
	<u>Other Tests Performed</u> Urine - RSID Negative		

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
CU4ETT	Kastle Meyer <i>Negative</i> Hexagon OBTI <i>Negative</i> Conclusion: Negative	not in our scope Conclusion: Not Tested	not in our scope Conclusion: Not Tested
CXB6NT	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
	<u>Other Tests Performed</u> Urine - Creatinine <i>Negative</i>		
CZGMMP	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
D7LP7V	no Blood testing conducted due to no observation of red/brown staining on any of the swabs Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
D8PPVU	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> RSID-Semen <i>Negative</i> PSA <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
D9NLU6	Visual Examination <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
D9VJDG	KM <i>Negative</i> Conclusion: Negative	Acid Phosphatase <i>Negative</i> Conclusion: Negative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
DGJNZM	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
	<u>Other Tests Performed</u> Urine - Creatinine <i>Negative</i>		
DHUHPV	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> RSID-Semen <i>Negative</i> PSA <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
DKDWTV	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> RSID-Semen <i>Negative</i> PSA <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
DRK8DB	Micro KM	Negative Negative	AP Micro	Negative Negative	RSID	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
END392	Visual O-tol	Negative Negative	ALS AP PSA Micro	Negative Negative Negative Negative		
	Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested	
EYWP69	KM TMB HemaTrace	Negative Negative Negative	ALS AP PSA	Negative Negative Negative		This analyst does not perform Saliva testing.
	Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested	
F36YQ6			AP	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
	<u>Other Tests Performed</u>					
	ALS Negative					
FN7W9X	HemDirect Hemoglobin Seratec	Negative	Christmas Tree Stain (Micro) PSA Seratec	Negative Negative	Amylase Test Seratec	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
FUZFGT	KM	Negative	AP	Negative		Not in analyst's protocol to test for Saliva.
	Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested	
FW4ZGL	not tested as no obvious rbs discoloration		AP Micro	Negative Negative	RSID	Positive/Indicative
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Positive/Indicative	
FWLA6G	Blood testing not conducted due to no visible RBS		AP	Negative	RSID	Positive/Indicative
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Positive/Indicative	
FYRKAW	KM	Negative	AP Micro	Negative Negative	RSID	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
FZGMQZ	HemDirect Hemoglobin (SERATEC)	Negative	Christmas Tree Stain Slide (Micro) PSA Semiquant (SERATEC)	Negative Negative	Amylase Test (SERATEC)	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	

Table 1 - Item 1
 Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
G2DR8Y	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
G6XGVH	There did not appear to be reddish staining on swab Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	Our Crime Scene unit does not test for Saliva Conclusion: Not Tested
GC99WV	Molecular Serology <i>Negative</i> Conclusion: Negative	Molecular Serology <i>Negative</i> Conclusion: Negative	Molecular Serology <i>Positive/Indicative</i> Conclusion: Positive/Indicative
GCA8Q3	Phenol <i>Negative</i> Conclusion: Negative	PSA <i>Negative</i> Conclusion: Negative	Amylase <i>Positive/Indicative</i> Conclusion: Positive/Indicative
GD34A6	Conclusion: Not Tested	ALS <i>Positive/Indicative</i> PSA <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
GXYAE	O-tol <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
H24RD4	Conclusion: Not Tested	ALS <i>Positive/Indicative</i> PSA <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
H8YVMF	No DRS observed Conclusion: Not Tested	PSA <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
	<u>Other Tests Performed</u> Urine-RSID <i>Negative</i>		
H9DDT8	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
HJFLMB	Hemastix <i>Negative</i> RSID-Blood <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> RSID-Semen <i>Negative</i> Microscope slide <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
HK3CJA	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
HQ68XW	SERATEC HemDirect <i>Negative</i> Conclusion: Negative	Christmas Tree <i>Negative</i> Immunochromatographic Test Rapid Stain Identification of Human Semen (RSID-Semen) <i>Negative</i> Conclusion: Negative	SERATEC Amylase Test <i>Positive/Indicative</i> Conclusion: Positive/Indicative
J74M6N	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	Phadebas Amylase <i>Positive/Indicative</i> Conclusion: Positive/Indicative
	<u>Other Tests Performed</u> Urine - Creatinine <i>Negative</i>		
J9AY8J	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
JBD4EC	Visual, white swab/no red color Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	Amylase Diffusion <i>Positive/Indicative</i> Conclusion: Positive/Indicative
JF76CK	Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
JGETMY	Conclusion: Not Tested	ALS <i>Positive/Indicative</i> PSA <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
JKBG7R	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
JRX62W	No stains consistent with Blood were observed Conclusion: Not Tested	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
JU8UHJ	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
JUFNLR	KM <i>Negative</i> Hemident <i>Negative</i> ALS <i>Negative</i> Conclusion: Negative	PSA <i>Negative</i> ALS <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
JWCGWJ	KM	Negative	AP	Negative	RSID-Saliva	Positive/Indicative
		Conclusion: Negative	RSID-Semen	Negative		
			PSA	Negative		
			Conclusion: Negative		Conclusion: Positive/Indicative	
JXBCVT	KM	Negative	AP	Negative	Phadebas Paper	Positive/Indicative
		Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative
JXJDAU	Molecular Serology	Negative	Molecular Serology	Negative	Molecular Serology	Positive/Indicative
		Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative
JYAL9W	LMG	Negative	AP	Negative	Amylase Radial Diffusion	Positive/Indicative
		Conclusion: Negative	Micro	Negative		
			PSA	Negative		
			Conclusion: Negative		Conclusion: Positive/Indicative	
JZDH48	KM	Negative	AP	Negative	Amylase	Positive/Indicative
		Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative
	<u>Other Tests Performed</u>					
	Urine - Creatinine Negative					
	Feces - Urobilinogen Negative					
K62L7A	The item is not brown stained.		Micro and p30		The laboratory has no validated procedure for Saliva testing.	
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
K68HAF			Micro	Negative	Phadebas	Positive/Indicative
		Conclusion: Not Tested	PSA (p30)	Negative		Conclusion: Positive/Indicative
			Conclusion: Negative			
K93MAD	No red-brown staining		Acid Phosphatase		RSID Saliva	
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Positive/Indicative	
KDVHHD	KM	Negative	AP	Negative	Our lab doesn't test for Saliva, testing for Amylase is conducted.	
		Conclusion: Negative	Micro	Negative	Conclusion: Not Tested	
	<u>Other Tests Performed</u>					
	Amylase Diffusion Positive/Indicative					
KG2YYG	No visible staining		PSA (p30)		RSID	
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Positive/Indicative	
			Micro	Negative		
KL4PPV	KM	Negative	AP	Negative	Not qualified	
		Conclusion: Negative		Conclusion: Negative	Conclusion: Not Tested	

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
KPFR9R	PTMB	Negative Conclusion: Negative	CS Proficiency for Blood screening only	Conclusion: Not Tested	CS Proficiency for Blood screening only	Conclusion: Not Tested
KPU8K8		Conclusion: Not Tested	Acid Phosphatase	Negative Conclusion: Negative	Phadebas paper	Positive/Indicative Conclusion: Positive/Indicative
KUY28N		no red/brown staining observed, no Blood testing performed Conclusion: Not Tested	AP	Negative Conclusion: Negative	RSID	Positive/Indicative Conclusion: Positive/Indicative
KX88VB	LMG	Negative Conclusion: Negative	AP	Negative Conclusion: Negative	Amylase Diffusion	Positive/Indicative Conclusion: Positive/Indicative
KXNE3D	KM	Negative Conclusion: Negative	AP Micro PSA	Negative Negative Negative Conclusion: Negative	Phadebas	Positive/Indicative Conclusion: Positive/Indicative
KYUEVT		Molecular Serology Conclusion: Negative	Molecular Serology	Negative Conclusion: Negative	Molecular Serology	Positive/Indicative Conclusion: Positive/Indicative
KZMBDV	Hexagon	Negative Conclusion: Negative	Micro Christmas tree PSA	Negative Negative Conclusion: Negative	Seratec	Positive/Indicative Conclusion: Positive/Indicative
LNCNP2	KM	Negative Conclusion: Negative	Acid Phosphatase PSA Microscopy	Negative Negative Negative Conclusion: Negative		Conclusion: Not Tested
LUH87E		No visible staining observed, swabs white in color Conclusion: Not Tested	AP	Negative Conclusion: Negative	RSID	Positive/Indicative Conclusion: Positive/Indicative
LVVYWA		No staining observed Conclusion: Not Tested	AP	Negative Conclusion: Negative	RSID	Positive/Indicative Conclusion: Positive/Indicative
LXYP8V	Visual	Negative Conclusion: Negative	AP	Negative Conclusion: Negative	RSID	Positive/Indicative Conclusion: Positive/Indicative
M3BX8F	LMG	Negative Conclusion: Negative	AP	Negative Conclusion: Negative		Conclusion: Not Tested
MEKR4Q		No Visual indications of Blood Conclusion: Not Tested	AP Micro	Negative Negative Conclusion: Negative	RSID Saliva	Positive/Indicative Conclusion: Positive/Indicative

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
MKKFX7	ABACard HemaTrace	Negative	ABACard p30	Negative	RSID Saliva	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
MKTUPU	Visual	Negative	AP	Negative	RSID	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
MRPREH	KM	Negative	AP	Negative	Phad (Micro tech)	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
	<u>Other Tests Performed</u>					
	Urine- Creatinine Negative					
MTZTZG	KM	Negative	AP	Negative	Phadebas Amylase	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
	<u>Other Tests Performed</u>					
	Urine-creatinine Negative					
MU7NDM	LMG	Negative	AP	Negative	Amylase Radial Diffusion	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
N4E6QJ	KM	Negative	AP	Negative	Phadebas	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
	<u>Other Tests Performed</u>					
	Urine - Creatinine Negative					
N7HDRH	LM (Leuco-Malachite)	Negative	AP	Negative	CSI does not test for Saliva	
	Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested	
NBP26J	LMG	Negative	AP	Negative	Amylase Radial Diffusion	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: No Response	
NE3WGH	PTMB	Negative	AP	Negative	Phadebas Tablet	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
NJPXAF	No red-brown staining observed on item.		AP	Negative	RSID	Positive/Indicative
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Positive/Indicative	
NKEQ4T	HemaTrace (hHb)	Negative	PSA	Negative	is not within the scope of this laboratory, we only perform identification of human blood and semen (P30 enzyme)	
	Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested	

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
P3Y9WY	KM	Negative	AP	Negative		
			PSA	Negative		
			Micro	Negative		
	Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested	
P7VMXN	Visual Exam- No staining Visually consistent with Blood.		AP		RSID	
			Micro		Visual Exam	
			Visual Exam			
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Positive/Indicative	
PF6RCZ	ALS	Negative	ALS	Inconclusive	Our lab does not test for Saliva	
			AP	Negative		
	Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested	
PJ6CF6	Not Tested Due to Visual Appearance		PSA	Negative	The [Laboratory] does not perform tests for Saliva	
			Micro	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
PJZK2N	No stains Visually consistent with Blood were observed.		AP	Negative		
			PSA	Negative		
			Micro	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
PKKRVR	KM	Negative	AP	Negative	RSID	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
PLPWZU			ALS	Positive/Indicative		
			PSA	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
PP7CBL	HemDirect Hemoglobin (SERATEC) Negative		Christmas Tree Coloring Negative Immunochromatographic Test (Semen-RSID) Negative		Amylase Test (SERATEC) Positive/Indicative	
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
PQ6UUG	no red brown stains observed		AP	Negative	RSID Saliva	Positive/Indicative
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Positive/Indicative	
PVZNRH	no reddish brown staining observed		AP	Negative	RSID-Saliva	Positive/Indicative
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Positive/Indicative	
	<u>Other Tests Performed</u>					
	ALS Positive/Indicative					

Table 1 - Item 1
 Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
PWAG6Q	KM	Negative	AP	Negative	Phadebas	Positive/Indicative
		Conclusion: Negative	Micro	Negative		Conclusion: Positive/Indicative
PWE34X	KM	Negative	AP	Negative	RSID	Positive/Indicative
		Conclusion: Negative	Micro	Negative		Conclusion: Positive/Indicative
PXZMHV	KM	Negative	AP	Negative		Conclusion: Not Tested
		Conclusion: Negative				
QDBYT8	no staining consistent with Blood		Acid Phosphatase	Negative	RSID Saliva	Positive/Indicative
		Conclusion: Not Tested		Conclusion: Negative		Conclusion: Positive/Indicative
QGEQ4T	KM	Negative	AP	Negative		Conclusion: Not Tested
		Conclusion: No Response	PSA	Negative		
			Micro	Negative		
				Conclusion: No Response		
QH6HUY	KM Test	Negative	Brentamine Test	Negative	Phadebas Forensic Press Test	Positive/Indicative
		Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative
QKUQYU			AP	Negative		Conclusion: Not Tested
		Conclusion: Not Tested		Conclusion: Negative		
QLAK2U	KM	Negative	AP	Negative	No testing for Saliva on commercial PTs.	
		Conclusion: Negative	PSA	Negative		Conclusion: Not Tested
			Micro	Negative		
				Conclusion: Negative		
QWM64W	KM	Negative	AP SPOT	Negative		Conclusion: Not Tested
		Conclusion: Negative		Conclusion: Negative		
R4TNA8	O-tol	Negative	ABACard p30	Negative	Test not performed as method verification for RSID field kit is currently being conducted.	
		Conclusion: Negative		Conclusion: Negative		Conclusion: Not Tested
R6CW3F	Not tested due to absence of red brown staining		AP	Negative	RSID	Positive/Indicative
		Conclusion: Not Tested	PSA	Negative		Conclusion: Positive/Indicative
				Conclusion: Negative		
RBLTBV			AP	Negative		Conclusion: Not Tested
		Conclusion: Not Tested		Conclusion: Negative		

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
RKF2NZ	Visual examination <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
RLN9YE	no reddish-brown staining/possible Blood observed Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
RPRUWD	PTMB <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> PSA <i>Negative</i> Conclusion: Negative	Phadebas Tablet <i>Positive/Indicative</i> Conclusion: Positive/Indicative
RTMQHX	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
RWNZ6K	Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
T8MP4A	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
TFHUWB	Visual <i>Negative</i> KM <i>Negative</i> Conclusion: Negative	Visual <i>Negative</i> ALS <i>Inconclusive</i> AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
THH6LJ	Conclusion: Not Tested	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
TKQCQW	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
TLYEQA	KM <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested	Conclusion: Not Tested
TUKQGY	LMG <i>Negative</i> RSID-Blood <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> RSID-Semen <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
TY7HZB		Conclusion: Not Tested	AP <i>Negative</i>	Conclusion: Negative		Conclusion: Not Tested
U77Y9V	lab policy, case scenario	Conclusion: Not Tested	AP	Conclusion: Negative	Phadebas	Conclusion: Positive/Indicative
U8QUAN	KM <i>Negative</i>	Conclusion: No Response	AP <i>Negative</i> Micro <i>Negative</i>	Conclusion: No Response	Phadebas tube test <i>Positive/Indicative</i>	Conclusion: No Response
UBCY2C	KM <i>Negative</i>	Conclusion: Negative	AP <i>Negative</i> RSID-Semen <i>Negative</i> PSA <i>Negative</i> Sperm Hy-Liter PI <i>Negative</i>	Conclusion: Negative	RSID-Saliva <i>Positive/Indicative</i>	Conclusion: Positive/Indicative
UETNJJ	LMG <i>Negative</i>	Conclusion: Negative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i>	Conclusion: Negative	This laboratory does not test for Saliva.	Conclusion: Not Tested
UF93XK	Phenolphthalein KM <i>Negative</i>	Conclusion: No Response	PSA and Acid Phosphatase <i>Negative</i>	Conclusion: No Response		Conclusion: Not Tested
	<u>Other Tests Performed</u> ALS using CSS range an orange filter- visible fluorescence <i>Negative</i>					
UGXDND	LMG <i>Negative</i> HemaTrace <i>Negative</i>	Conclusion: Negative	AP <i>Negative</i> p30 <i>Negative</i> Microscopic examination for Sperm cells <i>Negative</i>	Conclusion: Negative	Amylase <i>Positive/Indicative</i>	Conclusion: Positive/Indicative
UH277Z	LMG <i>Negative</i>	Conclusion: Negative	AP <i>Negative</i>	Conclusion: Negative		Conclusion: Not Tested
UZEWY4	Not tested due to no dried red staining observed on swabs	Conclusion: Not Tested	PSA <i>Negative</i>	Conclusion: Negative	RSID-Saliva <i>Positive/Indicative</i>	Conclusion: Positive/Indicative
	<u>Other Tests Performed</u> Urine-RSID <i>Negative</i>					
V3FG72	no RBS on swabs observed	Conclusion: Not Tested	AP <i>Negative</i>	Conclusion: Negative	RSID <i>Positive/Indicative</i>	Conclusion: Positive/Indicative

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
V9QJFB	PTMB <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> PSA <i>Negative</i> Conclusion: Negative	Phadebus Tablet <i>Positive/Indicative</i> Conclusion: Positive/Indicative
VHPP6X	LMG Conclusion: Negative	AP Conclusion: Negative	RSID Conclusion: Positive/Indicative
VKJRCT	no notable staining Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
VHHLQB	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
	<u>Other Tests Performed</u> Urine - Creatinine <i>Negative</i>		
W238Q9	no apparent discoloration, no rbs observed Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
WUV8A7	KM <i>Negative</i> Conclusion: Negative	ALS <i>Negative</i> AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
WV79W6	ALS <i>Negative</i> KM <i>Negative</i> RSID <i>Negative</i> Conclusion: Negative	ALS <i>Negative</i> AP <i>Negative</i> RSID <i>Negative</i> Conclusion: Negative	ALS <i>Negative</i> RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
WVXB33	LMG <i>Negative</i> Conclusion: Negative	ALS <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
WWQCWW	LMG <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Micro <i>Negative</i> PSA <i>Negative</i> Conclusion: Negative	Amylase Diffusion <i>Positive/Indicative</i> Conclusion: Positive/Indicative
WZB9GG	Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
X824XC	Phenolphthalein <i>Negative</i> Conclusion: Negative	Semen ABA Card p30 <i>Negative</i> Conclusion: Negative	Saliva RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
XAEG39	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
	<u>Other Tests Performed</u> Urine-Creatinine <i>Negative</i>		

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
XE3WZ9	LMG <i>Negative</i> Conclusion: No Response	Conclusion: Not Tested	Conclusion: Not Tested
XLJL39	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> PSA <i>Negative</i> RSID-Semen <i>Negative</i> Sperm HY-Liter PI <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
XMEKXM	Conclusion: Not Tested	Conclusion: Not Tested	Conclusion: Not Tested
XXDLD7	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> H&E slide <i>Negative</i> Conclusion: Negative	Phadebas tube <i>Positive/Indicative</i> Conclusion: Positive/Indicative
Y74ZXA	None done due to no discoloration on swabs Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
Y9D7F4	LMG <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
YA8RBG	Visual <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
YBF4NC	Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
YD44NM	Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
YFVAQD	KM <i>Negative</i> Conclusion: Negative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
YQRBGG	LC HR/MS <i>Negative</i> Conclusion: Negative	LC HR/MS <i>Negative</i> Conclusion: Negative	LC HR/MS <i>Positive/Indicative</i> Conclusion: Positive/Indicative
Z6BHHU	No red/brown color observed Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	Amylase Diffusion <i>Positive/Indicative</i> Conclusion: Positive/Indicative

Table 1 - Item 1
Evidence - Four white plastic swabs containing Saliva

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
ZC9797	KM	Negative	AP	Negative	RSID-Saliva	Positive/Indicative
			RSID-Semen	Negative		
			PSA	Negative		
			Sperm Hy-Liter PI	Negative		
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
ZD4X6C	KM	Negative	AP	Negative	RSID-Saliva	Positive/Indicative
			Micro	Negative		
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
ZFMHRM	No visible R/B		A/P	Negative	No Saliva test on screening PT worksheet	
			Control	Positive/Indicative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
ZFR8GY	KM	Negative	AP	Negative	Phadebas	Positive/Indicative
	Conclusion: Negative		Conclusion: Negative		Conclusion: Positive/Indicative	
	<u>Other Tests Performed</u>					
	Urine-Creatinine Negative					
ZUEEDZ	not red brown stained		AP, UV	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	

Response Summary for Table 1 - Item 1: Four white plastic swabs containing Saliva

	Blood	Semen	Saliva
Positive/Indicative	0	1	139
Negative	146	211	0
Inconclusive	0	0	0
Not Tested (NT)	74	8	82
No Response	4	4	3
Total	224	224	224

Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 2

Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
2289QN	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
22BP8V	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Microscopic <i>Negative</i> p30-ABACard <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
2G2TAX	Kastle Meyer <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Acid Phosphatase <i>Negative</i> p30 <i>Negative</i> Microscopy <i>Negative</i> Conclusion: Negative	SALigAE <i>Negative</i> Microscopy <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Urine - Urea Test <i>Negative</i> Feces - Urobilinogen Test <i>Negative</i> Feces - Microscopy <i>Negative</i>		
2Q8KTZ	LMG <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
2WBMVF	Visual Examination <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID Saliva <i>Negative</i> Conclusion: Negative
2Y2QL2	LMG <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
34RKYT	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
38HFFL	BLUESTAR (LUMINOL) <i>Positive/Indicative</i> HEXAGON OBTI TEST <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> CRIME LITE <i>Negative</i>		

Table 1 - Item 2
Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
3C2TKZ	Hemastix <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Luminol <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Negative</i> Conclusion: Negative
3GUVNZ	Phenolphthalein <i>Positive/Indicative</i> HEMASTIX <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Luminol <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
3VAJDJ	KM - Hemophan (paper test for Blood erythrocytes, haemoglobin) by Erba Lachema - non-specific test for Blood <i>Positive/Indicative</i> FOB - cassette test for determination of faecal occult Blood by NADAL - specific test for Human Blood <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP - Phosphatesmo (paper test for the determination of acid phosphatase) by Macherey-Nagel <i>Negative</i> PSA - cassette test for prostate specific antigen by SERATEC <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Negative</i> Conclusion: Negative
47AY3T	O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
4DB9ED	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
4DP4V9	Molecular Serology <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Molecular Serology <i>Negative</i> Conclusion: Negative	Molecular Serology <i>Negative</i> Conclusion: Negative
4GQCZH	KM <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> RSID <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
4R82UG	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
4W3QGF	KM <i>Positive/Indicative</i> Conclusion: No Response	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: No Response	Phadebas <i>Negative</i> Conclusion: No Response
	<u>Other Tests Performed</u> Visual / ODOUR TEST FOR POSSIBLE URINE / FAECES <i>Negative</i>		

Table 1 - Item 2
Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
4XR299	KM	Positive/Indicative	AP	Negative	Test no longer performed	
	HemaTrace	Positive/Indicative	Micro Sperm	Negative		
			p30	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
62AAYQ	O-tol	Positive/Indicative	AP	Negative	RSID Negative	
	HemaTrace	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
62P4BF	KM	Positive/Indicative	ALS	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
6AVWHW	TMB	Positive/Indicative	AP	Negative	Phadebas paper/tube test (Amylase)	
	HemaTrace	Positive/Indicative	Micro	Negative	Negative	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
6MWPPY	KM	Positive/Indicative	AP	Negative	RSID Negative	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
6QXKKV	KM	Positive/Indicative	AP	Negative	RSID Negative	
	HemaTrace	Positive/Indicative	Micro	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
6R22DE	KM	Positive/Indicative	AP	Negative	Phadebas Negative	
	Conclusion: Positive/Indicative		Micro	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
6VNAP7	Visual	Positive/Indicative	ALS	Negative	ALS Negative	
	KM	Positive/Indicative	AP	Negative	Amylase Diffusion Negative	
	Hemastix	Positive/Indicative	Micro	Negative		
	HemaTrace (HT)	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
6ZW78F	Molecular Serology Positive/Indicative		Molecular Serology Negative		Molecular Serology Negative	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
6ZW9RD	Kastle Meyer (KM) Positive/Indicative		AP	Negative	Lab does not test for Saliva	
	Conclusion: Positive/Indicative		Sperm search	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
7ARU7X	Micro	Positive/Indicative	AP	Negative	RSID Negative	
	KM	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	

Table 1 - Item 2
 Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
7LHP24	Visual KM HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	ALS AP Micro	Negative Negative Negative	Amylase Diffusion	Inconclusive
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
7N8RBX	KM Takayama	Positive/Indicative Positive/Indicative	Micro PSA	Negative Positive/Indicative	[Laboratory] does not conduct Saliva tests.	
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
7ULBZK	Seratec HemDirect OT	Positive/Indicative Positive/Indicative	Seratec PSA AP	Negative Negative	Amylase Diffusion	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
7YKAAU	KM Bluestar OBTI	Positive/Indicative Positive/Indicative	AP PSA RSID-Semen Sperm Hy-Liter PI	Negative Negative Negative Negative	RSID-Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
82LQ3H	KM	Positive/Indicative	AP	Negative	[Laboratory] does not conduct forensic serological tests for Saliva.	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
83UJJM	HemaTrace	Positive/Indicative	PSA	Negative	SALigAE	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
8ERLLK	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative	The [Laboratory] does not conduct forensic serological tests for Saliva.	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
8RNLK4	Visual Examination KM HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	AP Micro	Negative Negative	Amylase Diffusion	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
8RT27A	KM	Positive/Indicative	Brentamine - AP	Negative	Phadebas	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
8TTFRM	KM Hemastix HemaTrace Luminol	Positive/Indicative Positive/Indicative Positive/Indicative Positive/Indicative	AP	Negative	RSID Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	Body Fluid - Test Performed: ALS Negative					

Table 1 - Item 2
Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
8VWJK6	Leucomalachite Green <i>Positive/Indicative</i>		Kind's reagent <i>Negative</i>			
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
8Y4DL2	KM <i>Positive/Indicative</i>		AP <i>Negative</i>		RSID-Saliva <i>Negative</i>	
	Hemastix <i>Positive/Indicative</i>		PSA <i>Negative</i>			
	HemaTrace <i>Positive/Indicative</i>		Micro <i>Negative</i>			
	Luminol <i>Positive/Indicative</i>					
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
8YL7UQ	KM <i>Positive/Indicative</i>		AP <i>Negative</i>		Phadebas <i>Negative</i>	
	HemaTrace <i>Positive/Indicative</i>		Micro <i>Negative</i>			
			p30-ABAcad <i>Negative</i>			
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
8YVTV7	Rapid test for the detection of Human hemoglobin for use in legal medicine (Hexagon OBTI) <i>Positive/Indicative</i>		Sperm identification by Contrast Staining Christmas tree <i>Negative</i>		Seratec – Amylase Test <i>Negative</i>	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
946QJ6	KM <i>Positive/Indicative</i>		AP <i>Negative</i>		Amylase <i>Negative</i>	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
963BVP	KM <i>Positive/Indicative</i>		AP <i>Negative</i>		Phadebas <i>Negative</i>	
	Takayama <i>Positive/Indicative</i>					
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	Species Origin - HemaTrace <i>Positive/Indicative</i>					
	Urine - Creatinine <i>Negative</i>					
98U8JJ	Micro <i>Positive/Indicative</i>		AP <i>Negative</i>		RSID <i>Negative</i>	
	KM <i>Positive/Indicative</i>		Micro <i>Negative</i>			
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
9GXGNN	Micro <i>Positive/Indicative</i>		FBB <i>Negative</i>		RSID <i>Negative</i>	
	KM <i>Positive/Indicative</i>		Micro <i>Negative</i>			
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
9H2PM2	KM <i>Positive/Indicative</i>		AP <i>Negative</i>			
	HemaTrace <i>Positive/Indicative</i>					
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	

Table 1 - Item 2
 Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
9H9FT6	Visual O-tol HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	Visual ALS AP	Negative Negative Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
9HBZPK	KM HemaTrace	Positive/Indicative Positive/Indicative	PSA Micro	Negative Negative	[Laboratory] does not routinely test for Saliva.	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
9KDWA9			AP	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
9KH9X6	Bluestar Obti LMG	Positive/Indicative Positive/Indicative	Seratec PSA Semiquant	Negative	Phadebas® Amylase Test Tablets	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
9T4NET	KM HemaTrace	Positive/Indicative Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	
9UXQ9H	KM	Positive/Indicative	ALS AP	Negative Negative	Amylase PLATE	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
A4FUMV	KM Hexagon OBTI	Positive/Indicative Positive/Indicative	not in our scope		not in our scope	
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	
A4QTAY	Visual LMG HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	AP	Negative	RSID-Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
AAEF3J	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative	not testing for Saliva	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
AANRDT	Phenolphthalein HemaTrace HemaStix	Positive/Indicative Positive/Indicative Positive/Indicative	Acid Phosphatase	Negative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	

Table 1 - Item 2
 Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
AD6FYL	O-tol HemaTrace	Positive/Indicative Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	
ADC2PF	Visual KM HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	ALS p30	Negative Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
ADPV82			AP	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
AHKH3B	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
AJHZYK	Phenolphthalein Hemastix HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	AP	Negative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	CrimeLite 82S IR Positive/Indicative					
APAZLT	KM HemaTrace	Positive/Indicative Positive/Indicative	AP Micro	Negative Negative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
AWH3H9	HemaTrace (hHb)	Positive/Indicative	PSA	Negative	is not within the scope of this laboratory, we only perform identification of human blood and semen (P30 enzyme)	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
BAYGU7	KM LMG HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	AP Microscopy	Negative Negative	Amylase radial Diffusion Microscopy	Inconclusive Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
BKPKDT	LMG	Positive/Indicative	AP	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
BPVRLH	LMG, Item D	Positive/Indicative	AP	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	

Table 1 - Item 2
 Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
BTGR6Q	Phenolphthalein HemaTrace Conclusion: Positive/Indicative	Positive/Indicative Positive/Indicative	AP	Negative	RSID	Negative
BXGUFK	LMG Conclusion: Positive/Indicative	Positive/Indicative	AP Micro PSA	Negative Negative Negative	Amylase Diffusion	Negative
BXPWFJ	Phenolphthalein Conclusion: Positive/Indicative	Positive/Indicative	p30 Micro	Negative Negative	Phadebas	Negative
CFMNR6	KM Conclusion: Positive/Indicative	Positive/Indicative	ALS	Negative		Conclusion: Not Tested
CHEUWG	KM Conclusion: Positive/Indicative	Positive/Indicative	AP	Negative		Conclusion: Not Tested
CLAAZR	FOB Hemophan Conclusion: Positive/Indicative	Positive/Indicative Positive/Indicative	PSA RSID Semen	Negative Negative	RSID Saliva	Negative
CLWA6M	LMG Conclusion: Positive/Indicative	Positive/Indicative	AP	Negative		Conclusion: Not Tested
CQ7C8H	KM Conclusion: Positive/Indicative	Positive/Indicative	Micro p30	Negative Negative		Conclusion: Not Tested
CTEPNK	PTMB Conclusion: Positive/Indicative	Positive/Indicative	AP PSA Micro	Negative Negative Negative	Phadebas	Negative
CU4ETT	Kastle Meyer Hexagon OBTI Conclusion: Positive/Indicative	Positive/Indicative Positive/Indicative	not in our scope		not in our scope	
CXB6NT	KM Takayama Conclusion: Positive/Indicative	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Negative
	<u>Other Tests Performed</u> Species Origin - HemaTrace Positive/Indicative Urine - Creatinine Negative					

Table 1 - Item 2
Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
CZGMMP	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
D7LP7V	KM <i>Positive/Indicative</i> ABACard HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
D8PPVU	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> RSID-Semen <i>Negative</i> PSA <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Negative</i> Conclusion: Negative
D9NLU6	Visual Examination <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
D9VJDG	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Acid Phosphatase <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
DGJNZM	KM <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Species - HemaTrace <i>Positive/Indicative</i> Urine - Creatinine <i>Negative</i>		
DHUHPV	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> RSID-Semen <i>Negative</i> PSA <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Negative</i> Conclusion: Negative
DKDWTV	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> RSID-Semen <i>Negative</i> PSA <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Negative</i> Conclusion: Negative
DRK8DB	Micro <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative

Table 1 - Item 2
Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
END392	Visual <i>Positive/Indicative</i> O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Positive/Indicative</i> AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
EYWP69	KM <i>Positive/Indicative</i> TMB <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> PSA <i>Negative</i> Conclusion: Negative	This analyst does not perform Saliva testing. Conclusion: Not Tested
F36YQ6	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested	Conclusion: Not Tested
FN7W9X	HemDirect Hemoglobin Seratec <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Christmas Tree Stain (Micro) <i>Negative</i> PSA Seratec <i>Negative</i> Conclusion: Negative	Amylase test Seratec <i>Negative</i> Conclusion: Negative
FUZFGT	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Not in analyst's protocol to test for Saliva. Conclusion: Not Tested
FW4ZGL	Phenolphthalein <i>Positive/Indicative</i> Hemastix <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
FWLA6G	Hemastix <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Luminol <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
FYRKAW	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
FZGMQZ	HemDirect Hemoglobin (SERATEC) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Christmas Tree Stain Slide (Micro) <i>Negative</i> PSA Semiquant (SERATEC) <i>Negative</i> Conclusion: Negative	Amylase Test (SERATEC) <i>Negative</i> Conclusion: Negative
G2DR8Y	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative

Table 1 - Item 2
Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
G6XGVH	KM <i>Positive/Indicative</i> HemDirect <i>Positive/Indicative</i> Conclusion: Positive/Indicative	There was no case information indicating presence of Semen Conclusion: Not Tested	Our Crime Scene unit does not test for Saliva Conclusion: Not Tested
GC99WV	Molecular Serology <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Molecular Serology <i>Negative</i> Conclusion: Negative	Molecular Serology <i>Negative</i> Conclusion: Negative
GCA8Q3	Phenol <i>Positive/Indicative</i> Conclusion: Positive/Indicative	PSA <i>Negative</i> Conclusion: Negative	Amylase <i>Negative</i> Conclusion: Negative
GD34A6	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
GXYAE	O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
H24RD4	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
H8YVMF	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	PSA <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Urine-RSID <i>Negative</i>		
H9DDT8	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
HJFLMB	Hemastix <i>Positive/Indicative</i> RSID-Blood <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> RSID-Semen <i>Negative</i> Microscope slide <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Negative</i> Conclusion: Negative
HK3CJA	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
HQ68XW	SERATEC HemDirect <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Christmas Tree <i>Negative</i> Immunochromatographic Test Rapid Stain Identification of Human Semen (RSID-Semen) <i>Negative</i> Conclusion: Negative	SERATEC Amylase Test <i>Negative</i> Conclusion: Negative

Table 1 - Item 2
Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
J74M6N	KM	Positive/Indicative	AP	Negative	Phadebas Amylase	Negative
	Takayama	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	Species Origin - HemaTrace		Positive/Indicative			
	Urine - Creatinine		Negative			
J9AY8J	KM	Positive/Indicative	AP	Negative	RSID	Negative
			Micro	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
JBD4EC	LMG	Positive/Indicative	AP	Negative	Amylase Diffusion	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
JF76CK			AP	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
JGETMY	KM	Positive/Indicative	ALS	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
JKBG7R	KM	Positive/Indicative	AP	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
JRX62W	LMG	Positive/Indicative	AP	Negative		
	HemaTrace	Positive/Indicative	PSA	Negative		
			Micro	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
JU8UJH	KM	Positive/Indicative	AP	Negative	RSID	Negative
			Micro	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
JUFNLR	KM	Positive/Indicative	PSA	Negative		
	Hemident	Positive/Indicative	ALS	Negative		
	ALS	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
JWCGWJ	KM	Positive/Indicative	AP	Negative	RSID-Saliva	Negative
	Bluestar OBTI	Positive/Indicative	RSID-Semen	Negative		
			PSA	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
JXBCVT	KM	Positive/Indicative	AP	Negative	RSID	Negative
	HemaTrace	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	

Table 1 - Item 2
 Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
JXJDAU	Molecular Serology <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Molecular Serology <i>Negative</i> Conclusion: Negative	Molecular Serology <i>Negative</i> Conclusion: Negative
JYAL9W	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Amylase Radial Diffusion <i>Negative</i> Conclusion: Negative
JZDH48	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Amylase <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Urine - Creatinine <i>Negative</i> Feces - Uroblinogen <i>Negative</i>		
K62L7A	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Micro and p30 <i>Negative</i> Conclusion: Negative	The laboratory has no validated procedure for Saliva testing. Conclusion: Not Tested
K68HAF	TMB <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Micro <i>Negative</i> PSA (p30) <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
K93MAD	Phenolphthalein <i>Positive/Indicative</i> Hemastix <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Luminol <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Acid Phosphatase <i>Negative</i> Conclusion: Negative	RSID Saliva <i>Negative</i> Conclusion: Negative
KDVHHD	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Our lab doesn't test for Saliva, testing for Amylase is conducted Conclusion: Not Tested
	<u>Other Tests Performed</u> Amylase Diffusion <i>Inconclusive</i>		
KG2YYG	O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	PSA (p30) <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
KL4PPV	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Not qualified Conclusion: Not Tested
KPFR9R	PTMB <i>Positive/Indicative</i> Conclusion: Positive/Indicative	CS Proficiency for Blood screening only Conclusion: Not Tested	CS Proficiency for Blood screening only Conclusion: Not Tested

Table 1 - Item 2
Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
KPU8K8	Kastle Meyer	Positive/Indicative	Acid Phosphatase	Negative	Phadebas paper	Negative
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
KUY28N	KM	Positive/Indicative	AP	Negative	RSID	Negative
	HemaTrace	Positive/Indicative				
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
KX88VB	LMG	Positive/Indicative	AP	Negative	Amylase Diffusion	Negative
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
KXNE3D	KM	Positive/Indicative	AP	Negative	Phadebas	Negative
	HemaTrace	Positive/Indicative	Micro	Negative		
			PSA	Negative		
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
KYUEVT	Molecular Serology	Positive/Indicative	Molecular Serology	Negative	Molecular Serology	Negative
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
KZMBDV	Hexagon	Positive/Indicative	Micro Christmas tree	Negative	Seratec	Negative
			PSA	Negative		
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
LNCNP2	KM	Positive/Indicative	Acid Phosphatase	Negative		
	HemaTrace	Positive/Indicative	PSA	Negative		
			Microscopy	Negative		
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Not Tested
LUH87E	KM	Positive/Indicative	AP	Negative	RSID	Negative
	HemaTrace	Positive/Indicative	Micro	Negative		
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
LWYWA	KM	Positive/Indicative	AP	Negative	RSID	Negative
	Hemastix	Positive/Indicative				
	HemaTrace	Positive/Indicative				
	Luminol	Positive/Indicative				
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
LXYP8V	Visual	Positive/Indicative	AP	Negative	RSID	Negative
	LMG	Positive/Indicative				
	HemaTrace	Positive/Indicative				
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
M3BX8F	LMG	Positive/Indicative	AP	Negative		
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Not Tested

Table 1 - Item 2
 Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
MEKR4Q	KM HemDirect Conclusion: Positive/Indicative	Positive/Indicative Positive/Indicative	AP Micro Conclusion: Negative	Negative Negative	RSID Saliva Conclusion: Negative	Negative
MKKFX7	ABACard HemaTrace Conclusion: Positive/Indicative	Positive/Indicative	ABACard p30 Conclusion: Negative	Negative	RSID Saliva Conclusion: Negative	Negative
MKTUPU	Visual LMG HemaTrace Conclusion: Positive/Indicative	Positive/Indicative Positive/Indicative Positive/Indicative	AP Conclusion: Negative	Negative	RSID Conclusion: Negative	Negative
MRPREH	KM Mod Tak Conclusion: Positive/Indicative	Positive/Indicative Positive/Indicative	AP Conclusion: Negative	Negative	Phad (Micro tech) Conclusion: Negative	Negative
	<u>Other Tests Performed</u> Species- HemaTrace Positive/Indicative Urine- Creatinine Negative					
MTZTZG	KM Modified Takayama Conclusion: Positive/Indicative	Positive/Indicative Positive/Indicative	AP Conclusion: Negative	Negative	Phadebas Amylase Conclusion: Negative	Negative
	<u>Other Tests Performed</u> Species-HemaTrace Positive/Indicative Urine-creatinine Negative					
MU7NDM	LMG HemaTrace Conclusion: Positive/Indicative	Positive/Indicative Positive/Indicative	AP Conclusion: Negative	Negative	Amylase Radial Diffusion Conclusion: Inconclusive	Inconclusive
N4E6QJ	KM Modified Takayama Conclusion: Positive/Indicative	Positive/Indicative Positive/Indicative	AP Conclusion: Negative	Negative	Phadebas Conclusion: Negative	Negative
	<u>Other Tests Performed</u> Species Origin - HemaTrace Positive/Indicative Urine - Creatinine Negative					
N7HDRH	LM (Leuco-Malachite) Conclusion: Positive/Indicative	Positive/Indicative	AP Conclusion: Negative	Negative	CSI does not test for Saliva Conclusion: Not Tested	
NBP26J	LMG HemaTrace Conclusion: Positive/Indicative	Positive/Indicative Positive/Indicative	AP Conclusion: Negative	Negative	Amylase Radial Diffusion Conclusion: Inconclusive	Inconclusive

Table 1 - Item 2
Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
NE3WGH	PTMB HemaTrace	Positive/Indicative Positive/Indicative	AP PSA	Negative Negative	Phadebas Tablet	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
NJPXAF	Hemastix ABACard HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative	RSID	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
NKEQ4T	HemaTrace (hHb)	Positive/Indicative	PSA	Negative	is not within the scope of this laboratory, we only perform identification of human blood and semen (P30 enzyme)	
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
P3Y9WY	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative		
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
P7VMXN	KM Visual Exam		No case scenario to indicate possible Semen. Visual Exam- staining consistent with possible Blood.		No case scenario to indicate possible Saliva. Visual Exam- staining consistent with possible Blood.	
		Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested
PF6RCZ	ALS O-tolidine OBTI	Positive/Indicative Positive/Indicative Positive/Indicative	ALS AP	Inconclusive Negative	Our lab does not test for Saliva	
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
PJ6CF6	KM HemaTrace	Positive/Indicative Positive/Indicative	PSA Micro	Negative Negative	The[Laboratory] does not perform tests for Saliva	
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
PJZK2N	LMG HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative		
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
PKKRVR	KM	Positive/Indicative	AP Micro	Negative Negative	RSID	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
PLPWZU	KM	Positive/Indicative	ALS	Negative		
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested

Table 1 - Item 2
Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
PP7CBL	HemDirect Hemoglobin (SERATEC) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Christmas Tree Coloring <i>Negative</i> Immunochromatographic Test (Semen-RSID) <i>Negative</i> Conclusion: Negative	Amylase Test (SERATEC) <i>Negative</i> Conclusion: Negative
PQ6UUG	Hemastix <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> ABAcad p30 <i>Negative</i> Conclusion: Negative	RSID Saliva <i>Negative</i> Conclusion: Negative
PVZNRH	KM <i>Positive/Indicative</i> Hemastix <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative <u>Other Tests Performed</u> ALS <i>Positive/Indicative</i>	AP <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Negative</i> Conclusion: Negative
PWAG6Q	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
PWE34X	KM <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
PXZMHV	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
QDBYT8	Phenolphthalein <i>Positive/Indicative</i> Hemastix <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Luminol <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Acid Phosphatase <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID Saliva <i>Negative</i> Conclusion: Negative
QGEG4T	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: No Response	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: No Response	Conclusion: Not Tested
QH6HUY	KM Test <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Brentamine Test <i>Negative</i> Conclusion: Negative	Phadebas Forensic Press Test <i>Negative</i> Conclusion: Negative
QKUQYU	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested

Table 1 - Item 2
 Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
QLAK2U	KM	Positive/Indicative	AP	Negative	No Saliva testing on commercial PTs.	
	HemaTrace	Positive/Indicative	PSA	Negative		
			Micro	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
QWM64W	KM	Positive/Indicative	AP	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
R4TNA8	O-tol	Positive/Indicative	ABACard p30	Negative	Test not performed as method verification for RSID field kit is currently being conducted.	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
R6CW3F	KM	Positive/Indicative	AP	Negative	RSID Negative	
	HemaTrace	Positive/Indicative	PSA	Negative		
	Hemastix	Positive/Indicative				
	Luminol	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
RBLTBV			AP	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
RKF2NZ	Visual examination	Positive/Indicative	AP	Negative	Phadebas Negative	
	KM	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
RLN9YE	KM	Positive/Indicative	AP	Negative	RSID Negative	
	Hemastix	Positive/Indicative	Micro	Negative		
	HemaTrace	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
RPRUWD	PTMB	Positive/Indicative	AP	Negative	Phadebas Tablet Negative	
	HemaTrace	Positive/Indicative	PSA	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
RTMQHX	KM	Positive/Indicative	AP	Negative	RSID Negative	
			Micro	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
RWNZ6K			AP	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
T8MP4A	KM	Positive/Indicative	AP	Negative	RSID Negative	
	HemaTrace	Positive/Indicative	Micro	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	

Table 1 - Item 2
 Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
TFHUWB	Visual KM HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	ALS AP Micro	Negative Negative Negative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
THH6LJ			AP Micro	Negative Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
TKQCQW	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
TLYEQA	KM HemaTrace	Positive/Indicative Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	
TUKQGY	LMG RSID-Blood	Positive/Indicative Positive/Indicative	AP RSID-Semen	Negative Negative	RSID-Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
TY7HZB			AP	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
U77Y9V	KM Hexagon OBTI		AP		Phadebas	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
U8QUAN	KM	Positive/Indicative	AP Micro	Negative Negative	Phadebas tube test	Negative
	Conclusion: No Response		Conclusion: No Response		Conclusion: No Response	
UBCY2C	KM Bluestar OBTI	Positive/Indicative Positive/Indicative	AP RSID-Semen PSA Sperm Hy-Liter PI	Negative Negative Negative Negative	RSID-Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
UETNJJ	LMG HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative	This laboratory does not test for Saliva.	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	

Table 1 - Item 2
 Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
UF93XK	Phenolphthalein KM <i>Positive/Indicative</i> Conclusion: No Response	PSA and Acid Phosphatase <i>Negative</i> Conclusion: No Response	Conclusion: Not Tested
	<u>Other Tests Performed</u> no reaction w/ ALS using 415 wavelength and yellow filter as well as CSS w/ orange filter <i>Negative</i>		
UGXDND	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> p30 <i>Negative</i> Microscopic examination for Sperm cells <i>Negative</i> Conclusion: Negative	Amylase <i>Inconclusive</i> Conclusion: Inconclusive
	UH277Z	LMG <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative
UZEWY4	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	PSA <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Urine-RSID <i>Negative</i>		
V3FG72	HemaTrace <i>Positive/Indicative</i> hemastix <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
	V9QJFB	PTMB <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Conclusion: Negative
VHPP6X	LMG Conclusion: Positive/Indicative	AP Conclusion: Negative	RSID Conclusion: Negative
	<u>Other Tests Performed</u> HEMDIRECT +VE <i>Positive/Indicative</i>		
VKJRCT	Hemastix <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
	VHLQB	KM <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative
<u>Other Tests Performed</u> Species - HemaTrace <i>Positive/Indicative</i> Urine - Creatinine <i>Negative</i>			

Table 1 - Item 2
 Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
W238Q9	Hemastix <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	No case details Conclusion: Not Tested	No case details and per [Lab manual] caution should be exercised when when testing for Saliva on a stain that is mixed with Blood and/or Semen. Conclusion: Not Tested
WUV8A7	KM <i>Positive/Indicative</i> BLUE STAR IDENTI-HEM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
WV79W6	ALS <i>Negative</i> KM <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Negative</i> AP <i>Negative</i> RSID <i>Negative</i> Conclusion: Negative	ALS <i>Negative</i> RSID <i>Negative</i> Conclusion: Negative
WVXB33	LMG <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
WWQCWW	LMG <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> PSA <i>Negative</i> Conclusion: Negative	Amylase Diffusion <i>Negative</i> Conclusion: Negative
WZB9GG	Conclusion: Not Tested	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
X824XC	Phenolphthalein <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Semen ABA Card p30 <i>Negative</i> Micro-Sperm search Christmas Tree Staining <i>Negative</i> Conclusion: Negative	Saliva RSID <i>Negative</i> Conclusion: Negative
XAEG39	KM <i>Positive/Indicative</i> Takayama <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Species Origin-HemaTrace <i>Positive/Indicative</i> Urine-Creatinine <i>Negative</i>		
XE3WZ9	LMG <i>Positive/Indicative</i> Conclusion: No Response	Conclusion: Not Tested	Conclusion: Not Tested

Table 1 - Item 2
Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
XLJL39	KM Bluestar OBTI	Positive/Indicative Positive/Indicative	AP PSA RSID-Semen Sperm HY-Liter PI	Negative Negative Negative Negative	RSID-Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
XMEKXM	LMG	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	
XXDL7	KM	Positive/Indicative	AP H&E slide	Negative Negative	Phadebas tube	Positive/Indicative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Positive/Indicative	
Y74ZXA	KM Hemastix HemaTrace Luminol	Positive/Indicative Positive/Indicative Positive/Indicative Positive/Indicative	AP	Negative	RSID-Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
Y9D7F4	LMG	Positive/Indicative	AP	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
YA8RBG	Visual LMG HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	AP	Negative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
YBF4NC			AP	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
YD44NM	KM	Positive/Indicative	AP	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
YFVAQD	KM HT	Positive/Indicative Positive/Indicative	AP	Negative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
YQRBGG	LC HR/MS	Positive/Indicative	LC HR/MS	Negative	LC HR/MS	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
Z6BHHU	LMG	Positive/Indicative	AP	Negative	Amylase Diffusion	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	

Table 1 - Item 2
Evidence - Four pink plastic handle swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
ZC9797	KM Bluestar OBTI	Positive/Indicative Positive/Indicative	AP RSID-Semen PSA Sperm Hy-Liter PI	Negative Negative Negative Negative	RSID-Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
ZD4X6C	KM HemDirect	Positive/Indicative Positive/Indicative	AP Micro	Negative Negative	RSID-Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
ZFMHRM	O-tolidine Hex-OBTI Control	Positive/Indicative Positive/Indicative Positive/Indicative	A/P Control	Negative Positive/Indicative	No Saliva test on screening PT worksheet	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
ZFR8GY	KM Modified Takayama	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	Species-HemaTrace Positive/Indicative					
	Urine-Creatinine Negative					
ZUEEDZ	hemastix	Positive/Indicative	AP, UV	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	

Response Summary for Table 1 - Item 2: Four pink plastic handle swabs containing Blood

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	210	1	1
Negative	0	207	134
Inconclusive	0	0	3
Not Tested (NT)	9	12	84
No Response	5	4	2
Total	224	224	224

Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 3

Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
2289QN	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
22BP8V	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
2G2TAX	Kastle Meyer <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Acid Phosphatase <i>Negative</i> p30 <i>Negative</i> Microscopy <i>Negative</i> Conclusion: Negative	SALigAE <i>Negative</i> Microscopy <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Urine - Urea Test <i>Negative</i> Feces - Urobilinogen <i>Negative</i> Feces - Microscopy <i>Negative</i>		
2Q8KTZ	LMG <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
2WBMVF	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Visual Examination <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID Saliva <i>Negative</i> Conclusion: Negative
2Y2QL2	LMG <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
34RKYT	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	NT based on scenario Conclusion: Not Tested	NT based on scenario Conclusion: Not Tested
38HFFL	BLUESTAR (LUMINOL) <i>Positive/Indicative</i> HEXAGON OBTI TEST <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> CRIME LITE <i>Negative</i>		
3C2TKZ	pheno <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Negative</i> Conclusion: Negative

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
3GUVNZ	Phenolphthalein	Positive/Indicative	AP	Negative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
3VAJDJ	KM - Hemophan (paper test for Blood erythrocytes, haemoglobin) by Erba Lachema - non-specific test for Blood	Positive/Indicative	AP - Phosphatesmo (paper test for the determination of acid phosphatase) by Macherey-Nagel	Negative	RSID-Saliva	Negative
	FOB - cassette test for determination of faecal occult Blood by NADAL - specific test for Human Blood	Positive/Indicative	PSA - cassette test for prostate specific antigen by SERATEC	Negative	Conclusion: Negative	
47AY3T	O-tol	Positive/Indicative	Scenario details do not suggest need for this testing on this item.		Scenario details do not suggest need for this testing on this item.	
	HemaTrace	Positive/Indicative	Conclusion: Not Tested		Conclusion: Not Tested	
4DB9ED	KM	Positive/Indicative	AP	Negative	RSID	Negative
	Conclusion: Positive/Indicative		Micro	Negative	Conclusion: Negative	
4DP4V9	Molecular Serology	Positive/Indicative	Molecular Serology	Negative	Molecular Serology	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
4GQCZH	KM	Positive/Indicative	AP	Negative	RSID	Negative
	RSID	Positive/Indicative	RSID	Negative	Conclusion: Negative	
4R82UG	KM	Positive/Indicative	ALS	Negative	Conclusion: Not Tested	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
4W3QGF	KM	Positive/Indicative	Conclusion: Not Tested		Conclusion: Not Tested	
	Conclusion: No Response		Conclusion: Not Tested		Conclusion: Not Tested	
4XR299	KM	Positive/Indicative	AP	Negative	Test no longer performed	
	HemaTrace	Positive/Indicative	Micro Sperm	Negative	Conclusion: Not Tested	
	Conclusion: Positive/Indicative		p30	Negative	Conclusion: Not Tested	
62AAYQ	O-tol	Positive/Indicative	AP	Negative	Scenario does not suggest Saliva present on suspect hands.	
	HemaTrace	Positive/Indicative	Conclusion: Negative		Conclusion: Not Tested	
62P4BF	KM	Positive/Indicative	ALS	Negative	Conclusion: Not Tested	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
6AVWHW	TMB HemaTrace	Positive/Indicative Positive/Indicative	AP Micro	Negative Negative	Phadebas paper/tube test (Amylase) Micro	Negative Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
6MWPPY	KM	Positive/Indicative	AP	Negative	RSID	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
6QXKKV	KM	Positive/Indicative	scenario does not indicate Semen for this sample		RSID	Negative
		Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Negative
6R22DE	KM	Positive/Indicative				
		Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested
6VNAP7	Visual KM HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	ALS AP	Inconclusive Negative	Not probative	
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
6ZW78F	Molecular Serology	Positive/Indicative	Molecular Serology	Negative	Molecular Serology	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
6ZW9RD	KM	Positive/Indicative	AP Sperm Search	Negative Negative	Lab does not test for Saliva	
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
7ARU7X	Micro KM	Positive/Indicative Positive/Indicative	AP	Negative	RSID	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
7LHP24	Visual KM Hemastix HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative Positive/Indicative	ALS AP	Negative Negative	Not relevant	
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
7N8RBX	KM HemaTrace	Positive/Indicative Positive/Indicative	Micro PSA	Negative Positive/Indicative	[Laboratory] does not conduct Saliva tests.	
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested
7ULBZK	Seratec HemDirect OT	Positive/Indicative Positive/Indicative	Seratec PSA AP	Negative Negative	Amylase Diffusion	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative

Table 1 - Item 3
 Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
7YKAAU	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Not tested based on case scenario Conclusion: Not Tested	Not tested based on case scenario Conclusion: Not Tested
82LQ3H	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
83UJJM	HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	PSA <i>Negative</i> Conclusion: Negative	SALigAE <i>Negative</i> Conclusion: Negative
8ERLLK	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	The [Laboratory] does not conduct forensic serological tests for Saliva. Conclusion: Not Tested
8RNLK4	Visual Examination <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Amylase Diffusion <i>Negative</i> Conclusion: Negative
8RT27A	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Brentamine- AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
8TTFRM	KM <i>Positive/Indicative</i> Hemastix <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Luminol <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID Saliva <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Body Fluid - Test Performed: ALS <i>Negative</i>		
8VWJK6	Leucomalachite Green <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Kind's reagent <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
8Y4DL2	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	No Semen testing conducted based on case information. Focus would be to ID a body fluid from the victim on the suspect's hand. Conclusion: Not Tested	RSID-Saliva <i>Negative</i> Conclusion: Negative
8YL7UQ	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
8YVTV7	Rapid test for the detection of Human hemoglobin for use in legal medicine (Hexagon OBTI) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Sperm identification by Contrast Staining Christmas tree <i>Negative</i> Rapid Stain Identification of Human Semen (RSID) <i>Negative</i> Conclusion: Negative	Seratec – Amylase Test <i>Negative</i> Conclusion: Negative
946QJ6	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested	Conclusion: Not Tested
963BVP	KM <i>Positive/Indicative</i> Takayama <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
98U8JJ	Micro <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
9GXGNN	Micro <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	FBB <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
9H2PM2	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
9H9FT6	Visual <i>Positive/Indicative</i> O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Visual <i>Negative</i> ALS <i>Positive/Indicative</i> AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
9HBZPK	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Scenario did not warrant testing. Conclusion: Not Tested	[Laboratory] does not routinely test for Saliva. Conclusion: Not Tested
9KDWA9	Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
9KH9X6	Bluestar Obti <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Seratec PSA Semiquant <i>Negative</i> Conclusion: Negative	Phadebas® Amylase Test Tablets <i>Negative</i> Conclusion: Negative
9T4NET	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested	Conclusion: Not Tested

Table 1 - Item 3
 Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
9UXQ9H	KM	Positive/Indicative	ALS	Negative	Amylase PLATE	Negative
		Conclusion: Positive/Indicative	AP	Negative		Conclusion: Negative
A4FUMV	KM	Positive/Indicative	not in our scope		not in our scope	
	Hexagon OBTI	Positive/Indicative	Conclusion: Not Tested		Conclusion: Not Tested	
A4QTAY	Visual	Positive/Indicative	AP	Negative	RSID-Saliva	Negative
	LMG	Positive/Indicative	Conclusion: Negative		Conclusion: Negative	
	HemaTrace	Positive/Indicative				
AAEF3J	KM	Positive/Indicative	scenario		not testing for Saliva	
	HemaTrace	Positive/Indicative	Conclusion: Not Tested		Conclusion: Not Tested	
AANRDT	Phenolphthalein	Positive/Indicative	Acid Phosphatase	Negative	RSID	Negative
		Conclusion: Positive/Indicative	Conclusion: Negative		Conclusion: Negative	
AD6FYL	O-tol	Positive/Indicative				
	HemaTrace	Positive/Indicative	Conclusion: Not Tested		Conclusion: Not Tested	
ADC2PF	Visual	Positive/Indicative	ALS	Inconclusive		
	KM	Positive/Indicative	p30	Negative	Conclusion: Not Tested	
	HemaTrace	Positive/Indicative	Conclusion: Negative			
ADPV82			AP	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
AHKH3B	KM	Positive/Indicative				
	HemaTrace	Positive/Indicative	Conclusion: Not Tested		Conclusion: Not Tested	
AJHZYK	Phenolphthalein	Positive/Indicative	AP	Negative	RSID	Negative
	Luminol	Positive/Indicative	Conclusion: Negative		Conclusion: Negative	
APAZLT	KM	Positive/Indicative	Scenario did not provide info to support testing for Semen.		Scenario did not provide info to support testing for Saliva.	
		Conclusion: Positive/Indicative	Conclusion: Not Tested		Conclusion: Not Tested	

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
AWH3H9	HemaTrace (hHb)	Positive/Indicative	PSA	Negative	is not within the scope of this laboratory, we only perform identification of human blood and semen (P30 enzyme)	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
BAYGU7	KM	Positive/Indicative	Microscopy	Negative	Amylase radial Diffusion Inconclusive	
	LMG	Positive/Indicative	AP	Negative	Microscopy Negative	
	HemaTrace	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
BKPKDT	LMG	Positive/Indicative	AP	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
BPVRLH	LMG, Item D	Positive/Indicative	AP	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
BTGR6Q	Hemastix	Positive/Indicative	AP	Negative	RSID Negative	
	HemaTrace	Positive/Indicative				
	Luminol	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
BXGUFK	LMG	Positive/Indicative	AP	Negative	Amylase Diffusion Negative	
	HemaTrace	Positive/Indicative	Micro	Negative		
			PSA	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
BXPWFJ	Phenolphthalein	Positive/Indicative	p30	Negative	Phadebas Negative	
			Micro	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
CFMNR6	KM	Positive/Indicative	ALS	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
CHEUWG	KM	Positive/Indicative	AP	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
CLAAZR	FOB	Positive/Indicative	PSA	Negative	RSID Saliva Negative	
	Hemophan	Positive/Indicative	RSID Semen	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
CLWA6M	LMG	Positive/Indicative	AP	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
CQ7C8H	KM	Positive/Indicative	Micro	Negative		
			p30	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
CTEPNK	PTMB <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Saliva not suspected on suspect's hand based on scenario Conclusion: Not Tested
CU4ETT	Kastle Meyer <i>Positive/Indicative</i> Hexagon OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	not in our scope Conclusion: Not Tested	not in our scope Conclusion: Not Tested
CXB6NT	KM <i>Positive/Indicative</i> Takayama <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
CZGMMP	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
D7LP7V	HEMASTIX <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
D8PPVU	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Semen not tested based on case scenario Conclusion: Not Tested	Saliva not tested based on case scenario Conclusion: Not Tested
D9NLU6	Visual Examination <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
D9VJDG	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Acid Phosphatase <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
DGJNZM	KM <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative
DHUHPV	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Semen not tested based on case scenario. Conclusion: Not Tested	Saliva not tested based on case scenario. Conclusion: Not Tested
DKDWTV	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Semen not tested based on case scenario. Conclusion: Not Tested	Saliva not tested based on case scenario. Conclusion: Not Tested

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
DRK8DB	Micro KM	Positive/Indicative Positive/Indicative	AP Micro	Negative Negative	RSID	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
END392	Visual O-tol HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	ALS AP	Positive/Indicative Negative		
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
EYWP69	KM HemaTrace	Positive/Indicative Positive/Indicative				This analyst does not perform Saliva testing.
		Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested
F36YQ6	KM HemaTrace	Positive/Indicative Positive/Indicative				
		Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested
FN7W9X	HemDirect Hemoglobin Seratec	Positive/Indicative				
		Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested
FUZFGT	KM HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative		Not in analyst's protocol to test for Saliva.
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
FW4ZGL	Phenolphthalein HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative	RSID	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
FWLA6G	KM	Positive/Indicative	AP	Negative	RSID	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
FYRKAW	KM	Positive/Indicative	AP Micro	Negative Negative	RSID	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
FZGMQZ	HemDirect Hemoglobin (SERATEC)	Positive/Indicative	Christmas Tree Stain (Micro) PSA Semiquant (SERATEC)	Negative Negative	Amylase Test (SERATEC)	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
G2DR8Y	KM	Positive/Indicative	AP Micro	Negative Negative	RSID	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative

Table 1 - Item 3
 Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
G6XGVH	KM	Positive/Indicative	AP	Negative	Our Crime Scene unit does not test for Saliva	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
GC99WV	Molecular Serology	Positive/Indicative	Molecular Serology	Negative	Molecular Serology Negative	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
GCA8Q3	Phenol	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	
GD34A6	KM	Positive/Indicative	ALS	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
GXYAE	O-tol	Positive/Indicative	Synopsis doesn't indicate a reason to test for this		Synopsis doesn't indicate a reason this test for this	
	HemaTrace	Positive/Indicative	Conclusion: Not Tested		Conclusion: Not Tested	
H24RD4	KM	Positive/Indicative	ALS	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
H8YVMF	KM	Positive/Indicative	PSA	Negative	RSID-Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
H9DDT8	KM	Positive/Indicative	AP	Negative		
	HemaTrace	Positive/Indicative	PSA	Negative		
			Micro	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
HJFLMB	Hemastix	Positive/Indicative	AP	Negative	RSID-Saliva	Negative
	RSID-Blood	Positive/Indicative	Conclusion: Negative		Conclusion: Negative	
HK3CJA	KM	Positive/Indicative				
	HemaTrace	Positive/Indicative	Conclusion: Not Tested		Conclusion: Not Tested	
HQ68XW	SERATEC HemDirect	Positive/Indicative	Christmas Tree	Negative	SERATEC Amylase Test Negative	
			Immunochromatographic Test Rapid Stain Identification of Human Semen (RSID-Semen)	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
J74M6N	KM	Positive/Indicative	AP	Negative	Phadebas Amylase	Negative
	Takayama	Positive/Indicative	Conclusion: Negative		Conclusion: Negative	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
J9AY8J	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
JBD4EC	LMG <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Amylase Diffusion <i>Negative</i> Conclusion: Negative
JF76CK	Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
JGETMY	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
JKBG7R	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
JRX62W	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
JU8UJH	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
JUFNLR	KM <i>Positive/Indicative</i> Hemident <i>Positive/Indicative</i> ALS <i>Positive/Indicative</i> Conclusion: Positive/Indicative	PSA <i>Negative</i> ALS <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
JWCGWJ	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Semen not tested based on case scenario Conclusion: Not Tested	Saliva not tested based on case scenario Conclusion: Not Tested
JXBCVT	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
JXJDAU	Molecular Serology <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Molecular Serology <i>Negative</i> Conclusion: Negative	Molecular Serology <i>Negative</i> Conclusion: Negative
JYAL9W	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	Amylase Radial Diffusion <i>Negative</i> Conclusion: Negative

Table 1 - Item 3
 Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
JZDH48	KM HemaTrace	Positive/Indicative Positive/Indicative				
	Conclusion:	Positive/Indicative	Conclusion:	Not Tested	Conclusion:	Not Tested
K62L7A	KM	Positive/Indicative	Micro and p30	Negative	The laboratory has no validated procedure for Saliva testing.	
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Not Tested
K68HAF	TMB	Positive/Indicative				
	Conclusion:	Positive/Indicative	Conclusion:	Not Tested	Conclusion:	Not Tested
K93MAD	Phenolphthalein	Positive/Indicative	Acid Phosphatase	Negative	RSID Saliva	Negative
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
KDVHHD	KM	Positive/Indicative	per scenario		per scenario	
	Conclusion:	Positive/Indicative	Conclusion:	Not Tested	Conclusion:	Not Tested
KG2YYG	O-tol	Positive/Indicative	No evidentiary value in finding his own Semen on hands.		Not relevant based on scenario	
	Conclusion:	Positive/Indicative	Conclusion:	Not Tested	Conclusion:	Not Tested
KL4PPV	KM HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative	Not qualified	
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Not Tested
KPFR9R	PTMB	Positive/Indicative	CS Proficiency for Blood screening only		CS Proficiency for Blood screening only	
	Conclusion:	Positive/Indicative	Conclusion:	Not Tested	Conclusion:	Not Tested
KPU8K8	Kastle Meyer	Positive/Indicative	Acid Phosphatase	Negative	Phadebas paper	Negative
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
KUY28N	Hemastix	Positive/Indicative	AP	Negative	RSID	Negative
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
KX88VB	LMG	Positive/Indicative	AP	Negative	Amylase Diffusion	Negative
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
KXNE3D	KM HemaTrace	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Negative
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative
KYUEVT	Molecular Serology	Positive/Indicative	Molecular Serology	Negative	Molecular Serology Negative	
	Conclusion:	Positive/Indicative	Conclusion:	Negative	Conclusion:	Negative

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
KZMBDV	Hexagon	Positive/Indicative	Micro Christmas tree	Negative	Seratec	Negative
		Conclusion: Positive/Indicative	PSA	Negative		Conclusion: Negative
LNCNP2	KM	Positive/Indicative				
	HemaTrace	Positive/Indicative				
		Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested
LUH87E	Hemastix	Positive/Indicative	AP	Negative	RSID	Negative
	HemaTrace	Positive/Indicative	Micro	Negative		
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
LWVWA	KM	Positive/Indicative	AP	Negative	RSID	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
LXYP8V	Visual	Positive/Indicative	AP	Negative	RSID	Negative
	LMG	Positive/Indicative				
	HemaTrace	Positive/Indicative				
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
M3BX8F	LMG	Positive/Indicative	AP	Negative		
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
MEKR4Q	KM	Positive/Indicative	AP	Negative	RSID Saliva	Negative
	HemDirect	Positive/Indicative	RSID - Semen	Negative		
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
MKKFX7	ABACard HemaTrace	Positive/Indicative	ABACard p30	Negative	RSID Saliva	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
MKTUPU	Visual	Positive/Indicative	AP	Negative	RSID	Negative
	LMG	Positive/Indicative				
	HemaTrace	Positive/Indicative				
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
MRPREH	KM	Positive/Indicative	AP	Negative	Phad (Micro tech)	Negative
	Mod Tak	Positive/Indicative				
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
MTZTZG	KM	Positive/Indicative	AP	Negative	Phadebas Amylase	Negative
	Modified Takayama	Positive/Indicative				
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
MU7NDM	LMG HemaTrace	Positive/Indicative Positive/Indicative		scenario did not indicate need for testing		scenario did not indicate need for testing
		Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested
N4E6QJ	KM Modified Takayama	Positive/Indicative Positive/Indicative	AP	Negative	Phadebas	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
N7HDRH	LM (Leuco-Malachite)	Positive/Indicative	AP	Negative	CSI does not test for Saliva	
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
NBP26J	LMG HemaTrace	Positive/Indicative Positive/Indicative		Scenario does not indicate the need for this test.		Scenario does not indicate the need for this test.
		Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested
NE3WGH	PTMB HemaTrace	Positive/Indicative Positive/Indicative	AP PSA	Negative Negative	NT per case scenario.	
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
NJPXAF	KM	Positive/Indicative	AP	Negative	RSID	Negative
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative
NKEQ4T	HemaTrace (hHb)	Positive/Indicative	PSA	Negative	is not within the scope of this laboratory, we only perform identification of human blood and semen (P30 enzyme)	
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
P3Y9WY	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative		
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
P7VMXN	KM Visual Exam			Based on case scenario not needed.		Based on case scenario, not needed.
		Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: No Response
PF6RCZ	ALS O-tolidine OBTI	Positive/Indicative Positive/Indicative Positive/Indicative	ALS AP	Inconclusive Negative	Our lab does not test for Saliva	
		Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested
PJ6CF6	KM HemaTrace	Positive/Indicative Positive/Indicative		Not tested based on scenario of incident/location of collection		The [Laboratory] does not perform tests for Saliva
		Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
PJZK2N	LMG HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Negative Negative Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
PKKRVR	KM	Positive/Indicative	AP Micro	Negative Negative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
PLPWZU	KM	Positive/Indicative	ALS	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
PP7CBL	HemDirect Hemoglobin (SERATEC)	Positive/Indicative	Christmas Tree Coloring Immunochromatographic Test (Semen-RSID)	Negative Negative	Amylase Test (SERATEC)	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
PQ6UUG	KM	Positive/Indicative	AP Micro ABAcad p30	Negative Negative Negative	RSID Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
PVZNRH	KM	Positive/Indicative	AP	Negative	RSID-Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	ALS Positive/Indicative					
PWAG6Q	KM	Positive/Indicative	AP Micro	Negative Negative	Phadebas	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
PWE34X	KM Micro	Positive/Indicative Positive/Indicative	AP	Negative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
PXZMHV	KM	Positive/Indicative	AP	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
QDBYT8	Phenolphthalein	Positive/Indicative	Acid Phosphatase Micro	Negative Negative	RSID Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
QGEQ4T	KM Positive/Indicative HemaTrace Positive/Indicative		AP Negative PSA Negative Micro Negative			
	Conclusion: No Response		Conclusion: No Response		Conclusion: Not Tested	
QH6HUY	KM Test Positive/Indicative		Brentamine Test Negative		Phadebas Forensic Press Test Negative	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
QKUQYU	KM Positive/Indicative		AP Negative			
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
QLAK2U	KM Positive/Indicative HemaTrace Positive/Indicative		Swabs came from male suspect's hand		No Saliva testing on commercial PTs.	
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	
QWM64W	KM Positive/Indicative		AP Negative			
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
R4TNA8	O-tol Positive/Indicative		Swab not tested for Semen due to to synosis of incident and where the swab was collected from.		Test not performed as method verification for RSID field kit is currently being conducted.	
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	
R6CW3F	KM Positive/Indicative HemaTrace Positive/Indicative Hemastix Positive/Indicative		AP Negative PSA Negative		RSID Negative	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
RBLTBV			AP Negative			
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
RKF2NZ	Visual examination Positive/Indicative KM Positive/Indicative		AP Negative		Phadebas Negative	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
RLN9YE	KM Positive/Indicative				RSID Negative	
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Negative	
RPRUWD	PTMB Positive/Indicative HemaTrace Positive/Indicative		AP Negative PSA Negative		not tested because case scinario does not indicate victim licked assailants hand	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
RTMQHX	KM Positive/Indicative		AP Negative Micro Negative		RSID Negative	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	

Table 1 - Item 3
 Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
RWNZ6K		Conclusion: Not Tested	AP Negative	Conclusion: Negative		Conclusion: Not Tested
T8MP4A	KM Positive/Indicative	Conclusion: Positive/Indicative	not supported by scenario	Conclusion: Not Tested	not supported by scenario	Conclusion: Not Tested
TFHUWB	Visual Positive/Indicative KM Positive/Indicative	Conclusion: Positive/Indicative	no documentation that Semen was on (s) hand	Conclusion: Not Tested	RSID Negative	Conclusion: Negative
THH6LJ		Conclusion: Not Tested	AP Negative	Conclusion: Negative		Conclusion: Not Tested
TKQCQW	KM Positive/Indicative HemaTrace Positive/Indicative	Conclusion: Positive/Indicative	AP Negative PSA Negative Micro Negative	Conclusion: Negative		Conclusion: Not Tested
TLYEQA	KM Positive/Indicative HemaTrace Positive/Indicative	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested
TUKQGY	LMG Positive/Indicative RSID-Blood Positive/Indicative	Conclusion: Positive/Indicative	AP Negative RSID-Semen Negative	Conclusion: Negative	RSID-Saliva Negative	Conclusion: Negative
TY7HZB		Conclusion: Not Tested	AP Negative	Conclusion: Negative		Conclusion: Not Tested
U77Y9V	KM Hexagon OBTI	Conclusion: Positive/Indicative	LAB POLICY, case scenario	Conclusion: Not Tested	lab policy, case scenario	Conclusion: Not Tested
U8QUAN	KM Positive/Indicative	Conclusion: No Response	Finding S own Semen on his hand wouldn't assist and no expectation of detecting Semen from another source ie no digital penetration alleged	Conclusion: Not Tested	Finding S and/or IP Saliva on his hand would not assist addressing allegation.	Conclusion: Not Tested
UBCY2C	KM Positive/Indicative Bluestar OBTI Positive/Indicative	Conclusion: Positive/Indicative	Semen not tested based on case scenario.	Conclusion: Not Tested	Saliva not tested based on case scenario.	Conclusion: Not Tested

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
UETNJJ	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> PSA <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	This laboratory does not test for Saliva. Conclusion: Not Tested
UF93XK	Phenolphthalein KM <i>Positive/Indicative</i> Conclusion: No Response	PSA and Acid Phosphatase Conclusion: Negative	Conclusion: Not Tested
	<u>Other Tests Performed</u> ALS using 415 wavelength and yellow filter as well as CSS w/ orange filter <i>Positive/Indicative</i>		
UGXDND	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> p30 <i>Negative</i> Microscopic examination for Sperm cells <i>Negative</i> Conclusion: Negative	Amylase <i>Inconclusive</i> Conclusion: Inconclusive
UH277Z	LMG <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
UZEWY4	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Not applicable based on case scenario Conclusion: Not Tested	Not applicable based on case scenario Conclusion: Not Tested
V3FG72	KM <i>Positive/Indicative</i> luminol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> hemastix <i>Positive/Indicative</i> Conclusion: Positive/Indicative	per scenario details, item only tested for Blood Conclusion: Not Tested	per scenario details, item only tested for Blood Conclusion: Not Tested
V9QJFB	PTMB <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Not indicated Conclusion: Not Tested	Not indicated Conclusion: Not Tested
VHPP6X	LMG Conclusion: Positive/Indicative	AP Conclusion: Negative	RSID Conclusion: Negative
	<u>Other Tests Performed</u> HEMDIRECT +VE <i>Positive/Indicative</i>		
VKJRCT	Hemastix <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
WHLQB	KM <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Phadebas <i>Negative</i> Conclusion: Negative

Table 1 - Item 3
 Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
W238Q9	KM	Positive/Indicative	Per case details, Semen testing not performed.		Per case details, Saliva testing not performed.	
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	
WUV8A7	KM	Positive/Indicative	ALS	Negative		
	BLUE STAR IDENTI-HEM	Positive/Indicative	AP	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
WV79W6	ALS	Negative	ALS	Negative	ALS	Negative
	KM	Positive/Indicative	AP	Negative	RSID	Negative
	RSID	Positive/Indicative	RSID	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
WVXB33	LMG	Positive/Indicative	ALS	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
WWQCWW	LMG	Positive/Indicative	AP	Negative	Amylase Diffusion	Negative
			Micro	Negative		
			PSA	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
WZB9GG			AP	Negative		
			Micro	Negative		
	Conclusion: Not Tested		Conclusion: Negative		Conclusion: Not Tested	
X824XC	Phenolphthalein	Positive/Indicative	Semen ABA Card p30	Negative	Saliva RSID	Negative
	HemaTrace	Positive/Indicative	Micro-Sperm search Christmas Tree staining	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
XAEG39	KM	Positive/Indicative	AP	Negative	Phadebas	Negative
	Takayama	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
XE3WZ9	LMG	Positive/Indicative				
	Conclusion: No Response		Conclusion: Not Tested		Conclusion: Not Tested	
XLJL39	KM	Positive/Indicative	Semen not tested based on case scenario.		Saliva not tested based on case scenario.	
	Bluestar OBTI	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	
XMEKXM	LMG	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
XXDLD7	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Not tested given case circumstances and location swabs taken from Conclusion: Not Tested	Phadebas tube <i>Negative</i> Conclusion: Negative
Y74ZXA	KM <i>Positive/Indicative</i> Hemastix <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Negative</i> Conclusion: Negative
Y9D7F4	LMG <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
YA8RBG	Visual <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
YBF4NC	Conclusion: Not Tested	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
YD44NM	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
YFVAQD	KM <i>Positive/Indicative</i> HT <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Conclusion: Negative	RSID <i>Negative</i> Conclusion: Negative
YQRBGG	LC HR/MS <i>Positive/Indicative</i> Conclusion: Positive/Indicative	LC HR/MS <i>Negative</i> Conclusion: Negative	LC HR/MS <i>Negative</i> Conclusion: Negative
Z6BHHU	LMG <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Negative</i> Conclusion: Negative	Amylase Diffusion <i>Negative</i> Conclusion: Negative
ZC9797	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Semen not tested based on case scenario Conclusion: Not Tested	Saliva not tested based on case scenario Conclusion: Not Tested
ZD4X6C	KM <i>Positive/Indicative</i> HemDirect <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Negative</i> Micro <i>Negative</i> Conclusion: Negative	RSID-Saliva <i>Negative</i> Conclusion: Negative
ZFMHRM	O-tolidine <i>Positive/Indicative</i> Hex-OBTI <i>Positive/Indicative</i> Control <i>Positive/Indicative</i> Conclusion: Positive/Indicative	A/P <i>Negative</i> Control <i>Positive/Indicative</i> Conclusion: Negative	No Saliva test on screening PT worksheet Conclusion: Not Tested

Table 1 - Item 3
Evidence - Four blue plastic swabs containing Blood

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
ZFR8GY	KM	Positive/Indicative	AP	Negative	Phadebas	Negative
	Modified Takayama	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
ZUEEDZ	hemastix	Positive/Indicative	AP, UV	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	

Response Summary for Table 1 - Item 3: Four blue plastic swabs containing Blood

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	210	1	0
Negative	0	168	105
Inconclusive	0	0	1
Not Tested (NT)	9	54	117
No Response	5	1	1
Total	224	224	224

Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 4

Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
2289QN	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
22BP8V	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Microscopic <i>Positive/Indicative</i> p30-ABACard <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
2G2TAX	Not tested due to absence of red-brown staining Conclusion: Not Tested	Acid Phosphatase <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> Microscopy <i>Positive/Indicative</i> Conclusion: Positive/Indicative	SALigAE <i>Positive/Indicative</i> Microscopy <i>Negative</i> Conclusion: Positive/Indicative
	<u>Other Tests Performed</u> Urine - Urea Test <i>Inconclusive</i>		
2Q8KTZ	LMG <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
2WBMVF	Visual Examination <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> RSID Semen <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
2Y2QL2	LMG <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
34RKYT	physical exam = no red-brown = no testing Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
38HFFL	BLUESTAR (LUMINOL) <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> SERATEC PSA SEMI QUANT <i>Positive/Indicative</i> RSID Semen <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
	<u>Other Tests Performed</u> CRIME LITE <i>Positive/Indicative</i>		

Table 1 - Item 4
 Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
3C2TKZ	No visible staining - not tested Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	RSID-Saliva Positive/Indicative Conclusion: Positive/Indicative
3GUVNZ	No red/brown staining on swabs Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
3VAJDJ	KM - Hemophan (paper test for Blood erythrocytes, haemoglobin) by Erba Lachema - non-specific test for Blood Negative FOB - cassette test for determination of faecal occult Blood by NADAL - specific test for Human Blood Negative Conclusion: Negative	AP - Phosphatesmo (paper test for the determination of acid phosphatase) by Macherey-Nagel Positive/Indicative PSA - cassette test for prostate specific antigen by SERATEC Positive/Indicative Conclusion: Positive/Indicative	RSID-Saliva Positive/Indicative Conclusion: Positive/Indicative
47AY3T	No red brown staining observed on item. Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
4DB9ED	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
4DP4V9	Molecular Serology Negative Conclusion: Negative	Molecular Serology Positive/Indicative Conclusion: Positive/Indicative	Molecular Serology Positive/Indicative Conclusion: Positive/Indicative
4GQCZH	 Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative RSID Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
4R82UG	 Conclusion: Not Tested	ALS Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
4W3QGF	KM Negative Conclusion: No Response	AP Positive/Indicative Micro Positive/Indicative Conclusion: No Response	Phadebas TUBE TEST Positive/Indicative Conclusion: No Response

Table 1 - Item 4

Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
4XR299	Swabs were not red in color Conclusion: Not Tested	AP Positive/Indicative Micro Sperm Positive/Indicative p30 Positive/Indicative Conclusion: Positive/Indicative	Test no longer performed Conclusion: Not Tested
62AAYQ	O-tol Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
62P4BF	 Conclusion: Not Tested	ALS Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
6AVWHW	TMB Negative Conclusion: Negative	AP Positive/Indicative RSID Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas tube/paper test (Amylase) Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative
6MWPPY	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
6QXKKV	scenario does not indicate Blood for this sample Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
6R22DE	 Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas Positive/Indicative Conclusion: Positive/Indicative
6VNAP7	Visual Negative KM Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	ALS Positive/Indicative Amylase Diffusion Positive/Indicative Conclusion: Positive/Indicative
6ZW78F	Molecular Serology Negative Conclusion: Negative	Molecular Serology Positive/Indicative Conclusion: Positive/Indicative	Molecular Serology Positive/Indicative Conclusion: Positive/Indicative
6ZW9RD	Colorless Conclusion: Not Tested	AP Positive/Indicative Sperm search Positive/Indicative Conclusion: Positive/Indicative	Lab does not test for Saliva Conclusion: Not Tested

Table 1 - Item 4
Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
7ARU7X	Micro Negative KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
7LHP24	Visual Negative KM Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	Amylase Diffusion Positive/Indicative Conclusion: Positive/Indicative
7N8RBX	KM Negative Conclusion: Negative	Micro Positive/Indicative Conclusion: Positive/Indicative	[Laboratory] does not conduct Saliva tests. Conclusion: Not Tested
7ULBZK	Seratec HemDirect Negative OT Negative Conclusion: Negative	Seratec PSA Positive/Indicative AP Positive/Indicative Conclusion: Positive/Indicative	Amylase Diffusion Positive/Indicative Conclusion: Positive/Indicative
7YKAAU	Not tested based on case scenario and no stain observed Conclusion: Not Tested	PSA Positive/Indicative RSID-Semen Positive/Indicative Conclusion: Positive/Indicative	RSID-Saliva Positive/Indicative Conclusion: Positive/Indicative
82LQ3H	KM Negative Conclusion: Negative	AP Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
83UJJM	HemaTrace Positive/Indicative Conclusion: Positive/Indicative	PSA Positive/Indicative Conclusion: Positive/Indicative	SALigAE Positive/Indicative Conclusion: Positive/Indicative
8ERLLK	KM Negative Conclusion: Negative	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	The [Laboratory] does not conduct forensic serological tests for Saliva. Conclusion: Not Tested
8RNLK4	Visual Examination Negative KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Amylase Diffusion Positive/Indicative Conclusion: Positive/Indicative
8RT27A	KM Negative Conclusion: Negative	Brentamine - AP Positive/Indicative Conclusion: Positive/Indicative	Phadebas Positive/Indicative Conclusion: Positive/Indicative
8TTFRM	No reddish-brown staining observed Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative ABA card PSA Positive/Indicative Seratec PSA Positive/Indicative Conclusion: Positive/Indicative	RSID Saliva Positive/Indicative Conclusion: Positive/Indicative
	Other Tests Performed Body Fluid - Test Performed: ALS Positive/Indicative		

Table 1 - Item 4

Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
8VWJK6	No Blood like stains could be observed		Kind's reagent	Positive/Indicative		
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
8Y4DL2	Not conducted based on observations, no Visual indication of possible Blood.		AP	Positive/Indicative	RSID-Saliva	Positive/Indicative
	Conclusion: Not Tested		PSA	Positive/Indicative	Conclusion: Positive/Indicative	
			Micro	Positive/Indicative		
8YL7UQ	KM	Negative	AP	Positive/Indicative	Phadebas	Positive/Indicative
	Conclusion: Negative		Micro	Positive/Indicative	Conclusion: Positive/Indicative	
			p30-ABAcad	Positive/Indicative		
8YVTV7	Rapid test for the detection of Human hemoglobin for use in legal medicine (Hexagon OBTI)	Negative	Sperm identification by Contrast Staining	Christmas tree	Seratec – Amylase Test	Positive/Indicative
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
946QJ6			AP	Positive/Indicative	Amylase	Positive/Indicative
	Conclusion: Not Tested		Spermatozoa Microscopic ID	Positive/Indicative	Conclusion: Positive/Indicative	
963BVP	KM	Negative	AP	Positive/Indicative	Phadebas	Positive/Indicative
	Conclusion: Negative		Micro	Positive/Indicative	Conclusion: Positive/Indicative	
98U8JJ	Micro	Negative	AP	Positive/Indicative	RSID	Positive/Indicative
	KM	Negative	Micro	Positive/Indicative	Conclusion: Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative			
9GXGNN	Micro	Negative	FBB	Positive/Indicative	RSID	Positive/Indicative
	KM	Negative	Micro	Positive/Indicative	Conclusion: Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative			
9H2PM2	KM	Negative	AP	Positive/Indicative		
	Conclusion: Negative		Micro	Positive/Indicative	Conclusion: Not Tested	
9H9FT6	Visual	Negative	Visual	Negative		
	O-tol	Negative	ALS	Positive/Indicative		
	Conclusion: Negative		AP	Positive/Indicative		
			PSA	Positive/Indicative		
			Micro	Positive/Indicative	Conclusion: Not Tested	
			Conclusion: Positive/Indicative			

Table 1 - Item 4

Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
9HBZPK	No RBS observed.		PSA Positive/Indicative Micro Positive/Indicative		[Laboratory] does not routinely test for Saliva.	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
9KDWA9			AP Positive/Indicative Micro Positive/Indicative			
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
9KH9X6	Bluestar Obti Positive/Indicative LMG Negative		Seratec PSA Semiquant Positive/Indicative		Phadebas® Amylase Test Tablets Positive/Indicative	
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
9T4NET	KM Negative					
	Conclusion: Negative		Conclusion: Not Tested		Conclusion: Not Tested	
9UXQ9H	NO RBS		ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative		Amylase PLATE Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
A4FUMV	KM Negative		not in our scope		not in our scope	
	Conclusion: Negative		Conclusion: Not Tested		Conclusion: Not Tested	
A4QTAY	Visual Negative		AP Positive/Indicative RSID-Semen Positive/Indicative Micro Positive/Indicative		RSID-Saliva Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
AAEF3J	KM Negative		AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative		not testing for Saliva	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
AANRDT	No red/brown discoloration		Acid phosphatase Positive/Indicative Microscopic Sperm Search Positive/Indicative PSA Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
AD6FYL	Not tested due to no red-brown staining present.					
	Conclusion: Not Tested		Conclusion: Not Tested		Conclusion: Not Tested	
ADC2PF	Visual Negative KM Negative		ALS Positive/Indicative p30 Positive/Indicative			
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Not Tested	

Table 1 - Item 4
Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
ADPV82			AP Positive/Indicative Micro Positive/Indicative			
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
AHKH3B			AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative			
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
AJHZYK	not tested due to no visible RBS on swabs		AP Positive/Indicative Microscopic exam Positive/Indicative Seratec PSA Semiquant Positive/Indicative ABAcad p30 Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
APAZLT	No staining observed on the swabs.		AP Positive/Indicative Micro Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
AWH3H9	HemaTrace (hHb) Positive/Indicative		PSA Positive/Indicative		is not within the scope of this laboratory, we only perform identification of human blood and semen (P30 enzyme)	
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
BAYGU7	No red-brown staining		AP Positive/Indicative p30 Positive/Indicative Microscopy Positive/Indicative		Amylase radial Diffusion Positive/Indicative Microscopy Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
BKPKDT	LMG Negative		AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative			
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
BPVRLH	LMG Negative		AP, PSA Positive/Indicative			
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
BTGR6Q	No rbs noted to swabs		AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	

Table 1 - Item 4
Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
BXGUFK	Swabs were white in color Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	Amylase Diffusion Positive/Indicative Conclusion: Positive/Indicative
BXPWFJ	Phenolphthalein Negative Conclusion: Negative	p30 Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas Positive/Indicative Conclusion: Positive/Indicative
CFMNR6	 Conclusion: Not Tested	ALS Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
CHEUWG	KM Negative Conclusion: Negative	AP Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
CLAAZR	FOB Negative Hemophan Negative Conclusion: Negative	PSA Positive/Indicative RSID Semen Positive/Indicative Conclusion: Positive/Indicative	RSID Saliva Positive/Indicative Conclusion: Positive/Indicative
CLWA6M	LMG Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
CQ7C8H	 Conclusion: Not Tested	Micro Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
CTEPNK	no dry red stain present Conclusion: Not Tested	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas Positive/Indicative Conclusion: Positive/Indicative
CU4ETT	Kastle Meyer Negative Hexagon OBTI Negative Conclusion: Negative	not in our scope Conclusion: Not Tested	not in our scope Conclusion: Not Tested
CXB6NT	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas Positive/Indicative Conclusion: Positive/Indicative
CZGMMP	Visual exam did not support chemical testing for Blood Negative Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative

Table 1 - Item 4
 Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
D7LP7V	no Blood testing conducted due to no observation of red/brown staining on any of the swabs Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative ABAcad p30 Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
D8PPVU	Blood not tested based on case scenario and no stain observed Conclusion: Not Tested	RSID-Semen Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	RSID-Saliva Positive/Indicative Conclusion: Positive/Indicative
D9NLU6	Visual Examination Negative Conclusion: Negative	AP Positive/Indicative RSID Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
D9VJDG	KM Negative Conclusion: Negative	Acid Phosphatase Positive/Indicative Sperm Microscopy Positive/Indicative Conclusion: Positive/Indicative	Phadebas Positive/Indicative Conclusion: Positive/Indicative
DGJNZM	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas Positive/Indicative Conclusion: Positive/Indicative
DHUHPV	Blood not tested based on case scenario and no stain observed. Conclusion: Not Tested	RSID-Semen Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	RSID-Saliva Positive/Indicative Conclusion: Positive/Indicative
DKDWTV	Blood not tested based on case scenario and no stain observed. Conclusion: Not Tested	RSID-Semen Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	RSID-Saliva Positive/Indicative Conclusion: Positive/Indicative
DRK8DB	Micro Negative KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
END392	Visual Negative O-tol Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
EYWP69	 Conclusion: Not Tested	AP Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	This analyst does not perform Saliva testing. Conclusion: No Response

Table 1 - Item 4

Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
F36YQ6	<p>Conclusion: Not Tested</p> <p><u>Other Tests Performed</u> ALS <i>Positive/Indicative</i></p>	<p>AP <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>	<p>Conclusion: Not Tested</p>
FN7W9X	<p>Conclusion: Not Tested</p>	<p>Christmas Tree Stain (Micro) <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>	<p>Amylase Test Seratec <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>
FUZFGT	<p>KM <i>Negative</i></p> <p>Conclusion: Negative</p>	<p>AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>	<p>Not in analyst's protocol to test for Saliva.</p> <p>Conclusion: Not Tested</p>
FW4ZGL	<p>not tested as no obvious rbs discoloration</p> <p>Conclusion: Not Tested</p>	<p>AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>	<p>RSID <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>
FWLA6G	<p>Blood testing not conducted due to no visible RBS</p> <p>Conclusion: Not Tested</p>	<p>AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>	<p>RSID <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>
FYRKAW	<p>KM <i>Negative</i></p> <p>Conclusion: Negative</p>	<p>AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>	<p>RSID <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>
FZGMQZ	<p>HemDirect Hemoglobin (SERATEC) <i>Negative</i></p> <p>Conclusion: Negative</p>	<p>Christmas Tree Stain Slide (Micro) <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>	<p>Amylase Test (SERATEC) <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>
G2DR8Y	<p>KM <i>Negative</i></p> <p>Conclusion: Negative</p>	<p>AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>	<p>RSID <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>
G6XGVH	<p>There did not appear to be any reddish staining on swab</p> <p>Conclusion: Not Tested</p>	<p>AP <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>	<p>Our Crime Scene unit does not test for Saliva</p> <p>Conclusion: Not Tested</p>
GC99WW	<p>Molecular Serology <i>Negative</i></p> <p>Conclusion: Negative</p>	<p>Molecular Serology <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>	<p>Molecular Serology <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>
GCA8Q3	<p>Phenol <i>Negative</i></p> <p>Conclusion: Negative</p>	<p>PSA <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>	<p>Amylase <i>Positive/Indicative</i></p> <p>Conclusion: Positive/Indicative</p>

Table 1 - Item 4
Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
GD34A6			ALS Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative			
		Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested
GXYAE	Swabs are unstained		AP Positive/Indicative Micro Positive/Indicative		RSID Positive/Indicative	
		Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative
H24RD4			ALS Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative			
		Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested
H8YVMF	No DRS observed		PSA Positive/Indicative Micro Positive/Indicative		RSID-Saliva Positive/Indicative	
		Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative
H9DDT8	KM Negative		AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative			
		Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Not Tested
HJFLMB	Hemastix Negative		AP Positive/Indicative RSID-Semen Positive/Indicative Microscope slide Positive/Indicative		RSID-Saliva Positive/Indicative	
		Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative
HK3CJA			AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative			
		Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested
HQ68XW	SERATEC HemDirect Negative		Christmas Tree Positive/Indicative		SERATEC Amylase Test Positive/Indicative	
		Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative
J74M6N	KM Negative		AP Positive/Indicative Micro Positive/Indicative		Phadebas Amylase Positive/Indicative	
		Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative
J9AY8J	KM Negative		AP Positive/Indicative Micro Positive/Indicative		RSID Positive/Indicative	
		Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative

Table 1 - Item 4
 Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
JBD4EC	Visual, white swab/not red in color		AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative		Amylase Diffusion Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
JF76CK			AP Positive/Indicative Micro Positive/Indicative			
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
JGETMY			ALS Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative			
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
JKBG7R	KM Negative		AP Positive/Indicative Micro Positive/Indicative			
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
JRX62W	No stains consistent with Blood were observed		AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative			
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
JU8UJH	KM Negative		AP Positive/Indicative Micro Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
JUFNLR	KM Negative Hemident Negative ALS Negative		PSA Positive/Indicative ALS Positive/Indicative			
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
JWCGWJ	Blood not tested based on case scenario and no stain observed		RSID-Semen Positive/Indicative PSA Positive/Indicative		RSID-Saliva Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
JXBCVT			AP Positive/Indicative p30 Positive/Indicative Micro Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
JXJDAU	Molecular Serology Negative		Molecular Serology Positive/Indicative		Molecular Serology Inconclusive	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Inconclusive	

Table 1 - Item 4
 Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
JYAL9W	LMG <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase Radial Diffusion <i>Positive/Indicative</i> Conclusion: Positive/Indicative
JZDH48	 Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase <i>Positive/Indicative</i> Conclusion: Positive/Indicative
K62L7A	The item is not brown stained. Conclusion: Not Tested	Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	The laboratory has no validated procedure for Saliva testing. Conclusion: Not Tested
K68HAF	 Conclusion: Not Tested	Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
K93MAD	No red-brown staining Conclusion: Not Tested	Acid Phosphatase <i>Positive/Indicative</i> Microscopic Sperm <i>Positive/Indicative</i> Seratec PSA <i>Positive/Indicative</i> ABAcad p30 <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
KDVHHD	no discoloration, and per scenario. Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Our lab doesn't test for Saliva, testing for Amylase is conducted Conclusion: Not Tested
	<u>Other Tests Performed</u> Amylase Diffusion <i>Positive/Indicative</i>		
KG2YYG	No visible staining Conclusion: Not Tested	PSA (p30) <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
KL4PPV	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Not qualified Conclusion: Not Tested
KPFR9R	No red stains visible Conclusion: Not Tested	CS Proficiency for Blood screening only Conclusion: Not Tested	CS Proficiency for Blood screening only Conclusion: Not Tested
KPU8K8	 Conclusion: Not Tested	Semen Microscopy <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas paper <i>Positive/Indicative</i> Conclusion: Positive/Indicative

Table 1 - Item 4
Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
KUY28N	no red/brown staining observed, no Blood testing performed		AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
KX88VB	LMG Negative		AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative		Amylase Diffusion Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
KXNE3D	KM Negative		AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative		Phadebas Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
KYUEVT	Molecular Serology Negative		Molecular Serology Positive/Indicative		Molecular Serology Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
KZMBDV	Hexagon Negative		Micro Christmas tree Positive/Indicative PSA Positive/Indicative		SERATEC Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
LNCNP2			Acid Phosphatase Positive/Indicative PSA Positive/Indicative Microscopy Positive/Indicative			
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
LUH87E	KM Negative		AP Positive/Indicative Micro Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
LWYWA	No staining observed		AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
LXYP8V	Visual Negative		AP Positive/Indicative RSID Positive/Indicative Micro Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
M3BX8F	LMG Negative		AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative			
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Not Tested	

Table 1 - Item 4

Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
MEKR4Q	No Visual indications of Blood Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Saliva Positive/Indicative Conclusion: Positive/Indicative
MKKFX7	ABACard HemaTrace Positive/Indicative Conclusion: Positive/Indicative	ABACard p30 Positive/Indicative Conclusion: Positive/Indicative	RSID Saliva Positive/Indicative Conclusion: Positive/Indicative
MKTUPU	Visual Negative Conclusion: Negative	AP Positive/Indicative RSID Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
MRPREH	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phad (Micro tech) Positive/Indicative Conclusion: Positive/Indicative
MTZTZG	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas Amylase Positive/Indicative Conclusion: Positive/Indicative
MU7NDM	staining w/the Visual appearance of Blood not observed Conclusion: Not Tested	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Amylase Radial Diffusion Positive/Indicative Conclusion: Positive/Indicative
N4E6QJ	KM Negative Conclusion: Negative	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas Positive/Indicative Conclusion: Positive/Indicative
N7HDRH	LM (Leuco-Malachite) Negative Conclusion: Negative	AP Positive/Indicative Conclusion: Positive/Indicative	CSI does not test for Saliva Conclusion: Not Tested
NBP26J	Scenario does not indicate the need for this test. Conclusion: Not Tested	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Amylase Radial Diffusion Positive/Indicative Conclusion: Positive/Indicative
NE3WGH	NT per case scenario. Conclusion: Not Tested	AP Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	Phadebas Tablet Positive/Indicative Conclusion: Positive/Indicative

Table 1 - Item 4
 Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
NJPXAF	No red-brown staining observed on item. Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative ABAcad p30 Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
NKEQ4T	HemaTrace (hHb) Positive/Indicative Conclusion: Positive/Indicative	PSA Positive/Indicative Conclusion: Positive/Indicative	is not within the scope of this laboratory, we only perform identification of human blood and semen (P30 enzyme) Conclusion: Not Tested
P3Y9WY	KM Negative Conclusion: Negative	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
P7VMXN	Visual Exam- no staining consistent with Blood. Conclusion: Not Tested	AP Micro Visual Exam Conclusion: Positive/Indicative	RSID Visual Exam Conclusion: Positive/Indicative
PF6RCZ	ALS Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative Conclusion: Positive/Indicative	Our lab does not test for Saliva Conclusion: Not Tested
PJ6CF6	Not tested based on appearance/ scenario of incident/ location of collection Conclusion: Not Tested	PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	The [Laboratory] does not perform tests for Saliva Conclusion: Not Tested
PJZK2N	No stains Visually consistent with Blood were observed. Conclusion: Not Tested	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
PKKRVR	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
PLPWZU	 Conclusion: Not Tested	ALS Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
PP7CBL	HemDirect Hemoglobin (SERATEC) Negative Conclusion: Negative	Christmas Tree Coloring Positive/Indicative Immunochromatographic Test (Semen-RSID) Positive/Indicative Conclusion: Positive/Indicative	Amylase Test (SERATEC) Positive/Indicative Conclusion: Positive/Indicative

Table 1 - Item 4

Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
PQ6UUG	no red brown stains observed Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative ABAcard p30 Positive/Indicative Conclusion: Positive/Indicative	RSID Saliva Positive/Indicative Conclusion: Positive/Indicative
PVZNRH	no reddish brown staining observed Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID-Saliva Positive/Indicative Conclusion: Positive/Indicative
	<u>Other Tests Performed</u> ALS Positive/Indicative		
PWAG6Q	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas Positive/Indicative Conclusion: Positive/Indicative
PWE34X	KM Negative Micro Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Positive/Indicative Conclusion: Positive/Indicative
PXZMHV	KM Negative Conclusion: Negative	AP Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
	<u>Other Tests Performed</u> Micro Positive/Indicative		
QDBYT8	no staining consistent with Blood Conclusion: Not Tested	Acid Phosphatase Positive/Indicative Micro Positive/Indicative PSA (ABAcard) Positive/Indicative PSA (Seratec) Positive/Indicative Conclusion: Positive/Indicative	RSID Saliva Positive/Indicative Conclusion: Positive/Indicative
QGEQ4T	no visible staining Conclusion: Not Tested	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: No Response	Conclusion: Not Tested
QH6HUY	KM Test Negative Conclusion: Negative	Brentamine Test Positive/Indicative Conclusion: Positive/Indicative	Phadebas Forensic Press Test Positive/Indicative Conclusion: Positive/Indicative
QKUQYU	 Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested

Table 1 - Item 4

Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
QLAK2U	No visible staining on swabs. Scenario description.		AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative		No Saliva testing on commercial PTs.	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
QWM64W	KM Negative		AP Positive/Indicative Micro Positive/Indicative			
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
R4TNA8	No test for Blood due to the absence of brown/red staining on the swab.		ABACard p30 Positive/Indicative		Test not performed as method verification for RSID field kit is currently being conducted.	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
R6CW3F	Not tested due to absence of red brown staining		AP Positive/Indicative PSA Positive/Indicative Microscopic Identification of Spermatozoa Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
RBLTBV			AP Positive/Indicative Micro Positive/Indicative			
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
RKF2NZ	Visual examination Negative		AP Positive/Indicative Micro Positive/Indicative		Phadebas Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
RLN9YE	no reddish-brown staining/possible Blood observed		AP Positive/Indicative Micro Positive/Indicative PSA- ABACard p30 Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
RPRUWD	no red staining on swabs indicating possible Blood		AP Positive/Indicative PSA Positive/Indicative		Phadebas Tablet Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
RTMQHX	KM Negative		AP Positive/Indicative Micro Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
RWNZ6K			AP Positive/Indicative Micro Positive/Indicative			
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	

Table 1 - Item 4
 Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
T8MP4A	Visual exam <i>Negative</i> Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
TFHUWB	Visual <i>Negative</i> no stains consistent with Blood observed Conclusion: Not Tested	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
THH6LJ	 Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
TKQCQW	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
TLYEQA	KM <i>Negative</i> Conclusion: Negative	 Conclusion: Not Tested	 Conclusion: Not Tested
TUKQGY	LMG <i>Negative</i> RSID-Blood <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
TY7HZB	 Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
U77Y9V	LAB POLICY, case scenario Conclusion: Not Tested	AP Conclusion: Positive/Indicative	Phadebas Conclusion: Positive/Indicative
U8QUAN	KM <i>Negative</i> Conclusion: No Response	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: No Response	Phadebas tube test <i>Positive/Indicative</i> Conclusion: No Response
UBCY2C	Blood not tested based on case scenario and no stain observed. Conclusion: Not Tested	RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
UETNJJ	The swabs were unstained. The victim was known to be bleeding, so even if Blood was detected it would not be of probative value. Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	This laboratory does not test for Saliva. Conclusion: Not Tested

Table 1 - Item 4
Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
UF93XK	Phenolphthalein KM <i>Negative</i> Conclusion: No Response	PSA test <i>Positive/Indicative</i> Conclusion: No Response	Conclusion: Not Tested
	<u>Other Tests Performed</u> ALS w/ CSS wavelength and orange filter- fluorescence <i>Positive/Indicative</i>		
UGXDND	LMG <i>Negative</i> HemaTrace <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> Microscopic examination for Sperm cells <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase <i>Positive/Indicative</i> Conclusion: Positive/Indicative
UH277Z	LMG <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
UZEWY4	Not tested due to no dried red staining observed Conclusion: Not Tested	PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
V3FG72	no RBS observed on swabs Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative
V9QJFB	Not indicated Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadabus Tablet <i>Positive/Indicative</i> Conclusion: Positive/Indicative
VHPP6X	LMG Conclusion: Negative	AP Conclusion: Positive/Indicative	RSID Conclusion: Positive/Indicative
	<u>Other Tests Performed</u> H&E SLID : +3H <i>Positive/Indicative</i>		
VKJRCT	no notable staining Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
WHLQB	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Positive/Indicative</i> Conclusion: Positive/Indicative
W238Q9	no apparent discoloration, no rbs observed Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Positive/Indicative</i> Conclusion: Positive/Indicative

Table 1 - Item 4
Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
WUV8A7	KM Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
WV79W6	ALS Positive/Indicative KM Negative RSID Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative RSID Positive/Indicative MIC Positive/Indicative Conclusion: Positive/Indicative	ALS Positive/Indicative RSID Positive/Indicative Conclusion: Positive/Indicative
WVXB33	LMG Negative Conclusion: Negative	ALS Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
WWQCWW	LMG Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	Amylase Diffusion Positive/Indicative Conclusion: Positive/Indicative
WZB9GG	 Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
X824XC	Phenolphthalein Negative Conclusion: Negative	Semen ABA Card p30 Positive/Indicative Microscopic Sperm search Christmas Tree Staining Positive/Indicative Conclusion: Positive/Indicative	Saliva RSID Positive/Indicative Conclusion: Positive/Indicative
XAEG39	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas Positive/Indicative Conclusion: Positive/Indicative
XE3WZ9	LMG Negative Conclusion: No Response	 Conclusion: Not Tested	 Conclusion: Not Tested
XLJL39	Blood not tested based on case scenario and no stain observed. Conclusion: Not Tested	PSA Positive/Indicative RSID-Semen Positive/Indicative Conclusion: Positive/Indicative	RSID-Saliva Positive/Indicative Conclusion: Positive/Indicative
XMEKXM	 Conclusion: Not Tested	 Conclusion: Not Tested	 Conclusion: Not Tested
XXDLD7	KM Negative Conclusion: Negative	AP Positive/Indicative H&E slide Positive/Indicative Conclusion: Positive/Indicative	Phadebas tube Positive/Indicative Conclusion: Positive/Indicative

Table 1 - Item 4
Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
Y74ZXA	None done due to no discoloration on swabs		AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative		RSID-Saliva Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
Y9D7F4	LMG Negative		AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative			
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
YA8RBG	Visual Negative		AP Positive/Indicative RSID Positive/Indicative Microscopic Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
YBF4NC			AP Positive/Indicative Micro Positive/Indicative			
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
YD44NM			AP Positive/Indicative Micro Positive/Indicative			
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
YFVAQD	KM Negative		AP Positive/Indicative p30 Positive/Indicative Micro Positive/Indicative		RSID Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
YQRBGG	LC HR/MS Negative		LC HR/MS Positive/Indicative		LC HR/MS Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
Z6BHHU	No red/brown color observed		AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative		Amylase Diffusion Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
ZC9797	Blood not tested based on case scenario and no stain observed.		RSID-Semen Positive/Indicative PSA Positive/Indicative		RSID-Saliva Positive/Indicative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
ZD4X6C	KM Negative		AP Positive/Indicative Micro Positive/Indicative		RSID-Saliva Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	

Table 1 - Item 4

Evidence - Four yellow plastic swabs containing Semen & Saliva Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
ZFMHRM	No Visible R/B		A/P Positive/Indicative Control Positive/Indicative		No Saliva test on screening PT worksheet	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
ZFR8GY	KM Negative		AP Positive/Indicative Micro Positive/Indicative		Phadebas Positive/Indicative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	
ZUEEDZ	not red brown stained		AP, UV Positive/Indicative			
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	

**Response Summary for Table 1 -
Item 4: Four yellow plastic swabs containing Semen & Saliva Mixture**

	<u>Blood</u>	<u>Semen</u>	<u>Saliva</u>
Positive/Indicative	5	212	139
Negative	107	0	0
Inconclusive	0	0	1
Not Tested (NT)	108	8	81
No Response	4	4	3
Total	224	224	224

Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
2289QN	KM <i>Positive/Indicative</i>	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>	RSID <i>Negative</i>
	Conclusion: Positive/Indicative	Conclusion: Positive/Indicative	Conclusion: Negative
<u>Other Tests Performed</u> Semen/Saliva - ALS <i>Positive/Indicative</i>			
22BP8V	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i>	AP <i>Positive/Indicative</i> Microscopic <i>Positive/Indicative</i> p30-ABAcad <i>Positive/Indicative</i>	Phadebas <i>Negative</i>
	Conclusion: Positive/Indicative	Conclusion: Positive/Indicative	Conclusion: Negative
2G2TAX	Kastle Meyer <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i>	Acid Phosphatase <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> Microscopy <i>Positive/Indicative</i> Alternate Light Source <i>Negative</i>	SALigAE <i>Negative</i> Microscopy <i>Negative</i> Alternate Light Source <i>Negative</i>
	Conclusion: Positive/Indicative	Conclusion: Positive/Indicative	Conclusion: Negative
2Q8KTZ	LMG <i>Positive/Indicative</i>	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i>	
	Conclusion: Positive/Indicative	Conclusion: Positive/Indicative	Conclusion: Not Tested
2WBMVF	Visual Examination <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i>	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> RSID Semen <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>	RSID Saliva <i>Negative</i>
	Conclusion: Positive/Indicative	Conclusion: Positive/Indicative	Conclusion: Negative
2Y2QL2	LMG <i>Positive/Indicative</i>	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i>	
	Conclusion: Positive/Indicative	Conclusion: Positive/Indicative	Conclusion: Not Tested
34RKYT	KM <i>Positive/Indicative</i>	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i>	NT based on scenario
	Conclusion: Positive/Indicative	Conclusion: Positive/Indicative	Conclusion: Not Tested

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
38HFFL	BLUESTAR (LUMINOL) <i>Positive/Indicative</i> HEXAGON OBTI TEST <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> SERATEC PSA SEMI QUANT <i>Positive/Indicative</i> RSID Semen <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> CRIME LITE <i>Positive/Indicative</i>		
3C2TKZ	pheno <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID-Saliva <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Semen/Saliva - ALS <i>Positive/Indicative</i> Blood - IR <i>Positive/Indicative</i>		
3GUVNZ	Phenolphthalein <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Semen/Saliva - ALS <i>Positive/Indicative</i>		
3VAJDJ	KM - Hemophan (paper test for Blood erythrocytes, haemoglobin) by Erba Lachema - non-specific test for Blood <i>Positive/Indicative</i> FOB - cassette test for determination of faecal occult Blood by NADAL - specific test for Human Blood <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP - Phosphatesmo (paper test for the determination of acid phosphatase) by Macherey-Nagel <i>Positive/Indicative</i> PSA - cassette test for prostate specific antigen by SERATEC <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID-Saliva <i>Negative</i> Conclusion: Negative
47AY3T	O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> ALS <i>Negative</i> Conclusion: Positive/Indicative	Scenario details do not suggest need for this testing on this item. Conclusion: Not Tested
4DB9ED	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> ALS (Blood, Saliva, Semen) <i>Positive/Indicative</i>		
4DP4V9	Molecular Serology <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Molecular Serology <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Molecular Serology <i>Negative</i> Conclusion: Negative

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
4GQCZH	KM	Positive/Indicative	AP	Negative	RSID	Negative
	RSID	Positive/Indicative				
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	UV LIGHT Negative					
4R82UG	KM	Positive/Indicative	ALS	Positive/Indicative		
			PSA	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
4W3QGF	KM	Positive/Indicative	AP	Positive/Indicative	Amylase PAPER	Negative
			Micro	Positive/Indicative		
	Conclusion: No Response		Conclusion: No Response		Conclusion: No Response	
4XR299	KM	Positive/Indicative	AP	Positive/Indicative		Test no longer performed
	HemaTrace	Positive/Indicative	Micro Sperm	Positive/Indicative		
			p30	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
62AAYQ	O-tol	Positive/Indicative	ALS	Positive/Indicative	RSID	Negative
	HemaTrace	Positive/Indicative	AP	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
62P4BF	KM	Positive/Indicative	ALS	Positive/Indicative		
			PSA	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
6AVWHW	TMB	Positive/Indicative	AP	Positive/Indicative	Phadebas paper test	Negative
	HemaTrace	Positive/Indicative	RSID	Positive/Indicative	Micro	Negative
	Luminol	Positive/Indicative	Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
6MWPPY	KM	Positive/Indicative	AP	Positive/Indicative	RSID	Negative
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
6QXKKV	KM	Positive/Indicative	AP	Positive/Indicative	RSID	Negative
			Micro	Positive/Indicative		
			ALS	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
6R22DE	KM	Positive/Indicative	AP	Positive/Indicative	Amylase Paper	Negative
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		Conclusion: Negative
6VNAP7	Visual	Positive/Indicative	ALS	Positive/Indicative	Amylase mapping	Negative
	KM	Positive/Indicative	AP mapping	Positive/Indicative	Amylase Diffusion	Negative
	HemaTrace	Positive/Indicative	Micro	Positive/Indicative		Conclusion: Negative
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		
6ZW78F	Molecular Serology	Positive/Indicative	Molecular Serology	Positive/Indicative	Molecular Serology	Negative
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative
6ZW9RD	KM	Positive/Indicative	AP	Positive/Indicative	Lab does not test for Saliva	
		Conclusion: Positive/Indicative	Sperm Search	Positive/Indicative		Conclusion: Not Tested
				Conclusion: Positive/Indicative		
7ARU7X	ALS	Positive/Indicative	ALS	Positive/Indicative	RSID	Negative
	KM	Positive/Indicative	AP	Positive/Indicative		Conclusion: Negative
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		
				Conclusion: Positive/Indicative		
7LHP24	Visual	Positive/Indicative	ALS (more fluorescence inside)	Positive/Indicative	Amylase mapping	Negative
	KM	Positive/Indicative	AP (cutting)	Positive/Indicative	Amylase Diffusion	Negative
	HemaTrace	Positive/Indicative	AP (mapping)	Positive/Indicative		Conclusion: Negative
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		
				Conclusion: Positive/Indicative		
7N8RBX	KM	Positive/Indicative	ALS	Positive/Indicative	[Laboratory] does not conduct Saliva tests.	
	Takayama	Positive/Indicative	AP	Positive/Indicative		Conclusion: Not Tested
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		
				Conclusion: Positive/Indicative		
7ULBZK	Seratec HemDirect	Positive/Indicative	Seratec PSA	Positive/Indicative	Amylase Diffusion	Negative
	OT	Positive/Indicative	AP	Positive/Indicative		Conclusion: Negative
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		
7YKAAU	KM	Positive/Indicative	AP	Positive/Indicative	Generally not tested unless specified	
	Bluestar OBTI	Positive/Indicative	PSA	Positive/Indicative		Conclusion: Not Tested
		Conclusion: Positive/Indicative	RSID-Semen	Positive/Indicative		
				Conclusion: Positive/Indicative		
82LQ3H						Conclusion: Not Tested
		Conclusion: Not Tested		Conclusion: Not Tested		Conclusion: Not Tested

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
83UJJM	HemaTrace	Positive/Indicative	PSA	Positive/Indicative	SALigAE	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
8ERLLK	KM	Positive/Indicative	ALS	Positive/Indicative	The [Laboratory] does not conduct forensic serological tests for Saliva.	
	HemaTrace	Positive/Indicative	AP	Positive/Indicative		
			PSA	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
8RNLK4	Visual Examination	Positive/Indicative	ALS	Negative	Amylase Diffusion	
	KM	Positive/Indicative	AP	Positive/Indicative	Negative	
	HemaTrace	Positive/Indicative	Micro	Positive/Indicative		
			PSA	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
8RT27A	KM	Positive/Indicative	Brentamine- AP	Positive/Indicative	Phadebas	
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Negative	
					Conclusion: Negative	
8TTFRM	KM	Positive/Indicative	AP	Positive/Indicative	RSID Saliva	
	Hemastix	Positive/Indicative	Micro	Positive/Indicative	Negative	
	HemaTrace	Positive/Indicative	ABA card PSA	Positive/Indicative		
	Luminol	Positive/Indicative	Seratec PSA	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	Body Fluid - Test Performed: ALS Positive/Indicative					
	Body Fluid - Test Performed: IR Positive/Indicative					
8VWJK6	Leucomalachite Green	Positive/Indicative	Kind's reagent	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
8Y4DL2	KM	Positive/Indicative	AP	Positive/Indicative	RSID-Saliva	
			PSA	Positive/Indicative	Negative	
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
8YL7UQ	KM	Positive/Indicative	AP	Positive/Indicative	Phadebas	
	HemaTrace	Positive/Indicative	Micro	Positive/Indicative	Negative	
			p30-ABAcad	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
8YVTV7	Rapid test for the detection of Human hemoglobin for use in legal medicine (Hexagon OBTI)		Sperm identification by Contrast Staining Christmas tree		Seratec – Amylase Test	
	Positive/Indicative		Positive/Indicative		Negative	
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
946QJ6	KM	Positive/Indicative	AP	Positive/Indicative		
		Conclusion: Positive/Indicative	Spermatozoa Microscopic ID	Positive/Indicative		Conclusion: Not Tested
963BVP	KM	Positive/Indicative	AP	Positive/Indicative	Phadebas	Negative
	Takayama	Positive/Indicative	Micro	Positive/Indicative		Conclusion: Negative
98U8JJ	Micro	Positive/Indicative	ALS	Positive/Indicative	ALS	Positive/Indicative
	ALS	Positive/Indicative	AP	Positive/Indicative	RSID	Negative
	KM	Positive/Indicative	Micro	Positive/Indicative		Conclusion: Negative
9GXGNN	Micro	Positive/Indicative	ALS	Positive/Indicative	RSID	Negative
	KM	Positive/Indicative	FBB	Positive/Indicative		Conclusion: Negative
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		
9H2PM2	KM	Positive/Indicative	AP	Positive/Indicative		
	HemaTrace	Positive/Indicative	Micro	Positive/Indicative		Conclusion: Not Tested
		Conclusion: Positive/Indicative				
9H9FT6	Visual	Positive/Indicative	Visual	Negative		
	O-tol	Positive/Indicative	ALS	Inconclusive		
	HemaTrace	Positive/Indicative	AP	Positive/Indicative		
		Conclusion: Positive/Indicative	PSA	Positive/Indicative		
			Micro	Positive/Indicative		Conclusion: Not Tested
9HBZPK	KM	Positive/Indicative	ALS	Positive/Indicative	[Laboratory] does not routinely test for Saliva.	
	HemaTrace	Positive/Indicative	AP	Positive/Indicative		Conclusion: Not Tested
		Conclusion: Positive/Indicative	PSA	Positive/Indicative		
			Micro	Positive/Indicative		
9KDWA9			AP	Positive/Indicative		
		Conclusion: Not Tested	Micro	Positive/Indicative		Conclusion: Not Tested
9KH9X6	Bluestar Obti	Positive/Indicative	Seratec PSA Semiquant	Positive/Indicative	Phadebas® Amylase Test Tablets	Negative
	LMG	Positive/Indicative				Conclusion: Negative
		Conclusion: Positive/Indicative				

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
9T4NET	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested	Conclusion: Not Tested
9UXQ9H	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase PLATE <i>Negative</i> Conclusion: Negative
A4FUMV	KM <i>Positive/Indicative</i> Hexagon OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	not in our scope Conclusion: Not Tested	not in our scope Conclusion: Not Tested
A4QTAY	Visual <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Negative</i> AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID-Saliva <i>Negative</i> Conclusion: Negative
AAEF3J	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	not testing for Saliva Conclusion: Not Tested
AANRDT	Phenolphthalein <i>Positive/Indicative</i> Luminol <i>Positive/Indicative</i> IR Camera <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Acid Phosphatase <i>Positive/Indicative</i> Microscopic Sperm Search <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
<u>Other Tests Performed</u> Saliva/Semen - ALS <i>Positive/Indicative</i>			
AD6FYL	O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested	Conclusion: Not Tested
ADC2PF	Visual <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
ADPV82	Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
AHKH3B	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
AJZYZK	Phenolphthalein <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Microscopic exam <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
<u>Other Tests Performed</u> CrimeLite 82S blue (with yellow goggles) <i>Negative</i> CrimeLite 82S IR <i>Positive/Indicative</i>			
APAZLT	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Scenario did not provide info to support testing for Saliva. Conclusion: Not Tested
AWH3H9	HemaTrace (hHb) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	is not within the scope of this laboratory, we only perform identification of human blood and semen (P30 enzyme) Conclusion: Not Tested
BAYGU7	KM <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> Microscopy <i>Positive/Indicative</i> ALS <i>Negative</i> AP overlay <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase radial Diffusion <i>Inconclusive</i> Microscopy <i>Negative</i> Conclusion: Negative
BKPKDT	LMG <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
BPVRLH	LMG, Item D <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP, PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
BTGR6Q	IR <i>Positive/Indicative</i> Phenolphthalein <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Positive/Indicative</i> RSID <i>Negative</i> Conclusion: Negative

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
BXGUFK	LMG	Positive/Indicative	ALS AP Micro Mapping	Positive/Indicative Positive/Indicative Positive/Indicative Positive/Indicative	Amylase Diffusion Mapping	Negative Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
BXPWFJ	Phenolphthalein	Positive/Indicative	ALS p30 Micro	Negative Positive/Indicative Positive/Indicative	ALS Phadebas	Negative Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
CFMNR6	KM	Positive/Indicative	ALS PSA Micro	Positive/Indicative Positive/Indicative Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
CHEUWG						
	Conclusion: Not Tested		Conclusion: Not Tested		Conclusion: Not Tested	
CLAAZR	FOB Hemophan	Positive/Indicative Positive/Indicative	PSA RSID Semen	Positive/Indicative Positive/Indicative	RSID Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
CLWA6M	LMG	Positive/Indicative	AP Micro PSA	Positive/Indicative Positive/Indicative Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
CQ7C8H	KM	Positive/Indicative	ALS Micro	Negative Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
CTEPNK	PTMB	Positive/Indicative	AP PSA Micro ALS	Positive/Indicative Positive/Indicative Positive/Indicative Positive/Indicative	Phadebas	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
CU4ETT	Kastle Meyer Hexagon OBTI	Positive/Indicative Positive/Indicative	not in our scope		not in our scope	
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	
CXB6NT	KM Takayama	Positive/Indicative Positive/Indicative	AP Micro	Positive/Indicative Positive/Indicative	Phadebas	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
CZGMMP	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
D7LP7V	KM <i>Positive/Indicative</i> Luminol <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Infrared Light <i>Positive/Indicative</i> ALS (blue light) <i>Positive/Indicative</i>		
D8PPVU	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Saliva generally not tested unless specified Conclusion: Not Tested
D9NLU6	Visual Examination <i>Positive/Indicative</i> LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Positive/Indicative</i> RSID <i>Negative</i> Conclusion: Negative
D9VJDG	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Acid Phosphatase <i>Positive/Indicative</i> Sperm Microscopy <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Negative</i> Conclusion: Negative
DGJNZM	KM <i>Positive/Indicative</i> Modified Takayama <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Negative</i> Conclusion: Negative
DHUHPV	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Saliva generally not tested unless specified. Conclusion: Not Tested
DKDWTV	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Saliva generally not tested unless specified. Conclusion: Not Tested
DRK8DB	Micro <i>Positive/Indicative</i> KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
END392	Visual <i>Positive/Indicative</i> O-tol <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
EYWP69	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	This analyst does not perform Saliva testing. Conclusion: Not Tested
F36YQ6	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested
FN7W9X	HemDirect Hemoglobin Seratec <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Christmas Tree Stain (Micro) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
FUZFGT	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Not in analyst's protocol to test for Saliva. Conclusion: Not Tested
FW4ZGL	Phenolphthalein <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative <u>Other Tests Performed</u> Semen/Saliva - ALS <i>Positive/Indicative</i> Blood - Luminol <i>Positive/Indicative</i>	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
FWLA6G	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative <u>Other Tests Performed</u> Crime-Lite 82S Blue ALS <i>Positive/Indicative</i> Crime-Lite 82S IR <i>Positive/Indicative</i>	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> p30 (ABAcad) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
FYRKAU	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
FZGMQZ	HemDirect Hemoglobin (SERATEC) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Christmas Tree Stain Slide (Micro) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase Test (SERATEC) <i>Negative</i> Conclusion: Negative

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
G2DR8Y	KM	Positive/Indicative	AP	Positive/Indicative	RSID	Negative
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		Conclusion: Negative
G6XGVH	KM	Positive/Indicative	AP	Positive/Indicative	Our Crime Scene unit does not test for Saliva	
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	Conclusion: Not Tested	
	<u>Other Tests Performed</u>					
	ALS Negative					
GC99WV	Molecular Serology		Molecular Serology		Molecular Serology	
	Positive/Indicative		Positive/Indicative		Negative	
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
GCA8Q3	Phenol	Positive/Indicative	PSA	Positive/Indicative	Amylase	Negative
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative
	<u>Other Tests Performed</u>					
	ALS light - Stains Positive/Indicative					
GD34A6	KM	Positive/Indicative	ALS	Positive/Indicative	Conclusion: Not Tested	
		Conclusion: Positive/Indicative	PSA	Positive/Indicative		
			Micro	Positive/Indicative		
GXYAE	O-tol	Positive/Indicative	ALS	Negative	Synopsis doesn't indicate a reason to test for this	
	HemaTrace	Positive/Indicative	AP	Positive/Indicative	Conclusion: Not Tested	
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		
	<u>Other Tests Performed</u>					
H24RD4	KM	Positive/Indicative	ALS	Positive/Indicative	Conclusion: Not Tested	
		Conclusion: Positive/Indicative	PSA	Positive/Indicative		
			Micro	Positive/Indicative		
H8YVMF	KM	Positive/Indicative	PSA	Positive/Indicative	RSID-Saliva	Negative
		Conclusion: Positive/Indicative	Micro	Positive/Indicative	Conclusion: Negative	
			ALS	Positive/Indicative		
	<u>Other Tests Performed</u>					
H9DDT8	KM	Positive/Indicative	AP	Positive/Indicative	Conclusion: Not Tested	
	HemaTrace	Positive/Indicative	PSA	Positive/Indicative		
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
HJFLMB	Hemastix RSID-Blood	Positive/Indicative Positive/Indicative	ALS AP Microscope slide	Positive/Indicative Positive/Indicative Positive/Indicative	ALS RSID-Saliva	Positive/Indicative Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
HK3CJA	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Positive/Indicative Positive/Indicative Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
HQ68XW	SERATEC HemDirect	Positive/Indicative	Christmas Tree	Positive/Indicative	SERATEC Amylase Test	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
J74M6N	KM Takayama	Positive/Indicative Positive/Indicative	AP Micro	Positive/Indicative Positive/Indicative	Phadebas Amylase	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
J9AY8J	KM	Positive/Indicative	AP Micro	Positive/Indicative Positive/Indicative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
JBD4EC	LMG	Positive/Indicative	ALS AP Micro	Positive/Indicative Positive/Indicative Positive/Indicative	Amylase Diffusion	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
JF76CK			ALS-Blue AP Micro	Positive/Indicative Positive/Indicative Positive/Indicative		
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
JGETMY	KM	Positive/Indicative	ALS PSA Micro	Positive/Indicative Positive/Indicative Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
JKBG7R	KM	Positive/Indicative	AP Micro	Positive/Indicative Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
JRX62W	LMG HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro ALS	Positive/Indicative Positive/Indicative Positive/Indicative Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
JU8UHJ	KM	Positive/Indicative	AP	Positive/Indicative	RSID	Negative
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		Conclusion: Negative
JUFNLR	KM	Positive/Indicative	PSA	Positive/Indicative		
	ALS	Positive/Indicative	ALS	Positive/Indicative		Conclusion: Not Tested
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		
JWCGWJ	KM	Positive/Indicative	AP	Positive/Indicative		Saliva generally not tested unless specified.
	Bluestar OBTI	Positive/Indicative	RSID-Semen	Positive/Indicative		
		Conclusion: Positive/Indicative	PSA	Positive/Indicative		Conclusion: Not Tested
				Conclusion: Positive/Indicative		
JXBCVT	KM	Positive/Indicative	AP	Positive/Indicative	RSID	Negative
	HemaTrace	Positive/Indicative	p30	Positive/Indicative		
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		Conclusion: Negative
				Conclusion: Positive/Indicative		
	<u>Other Tests Performed</u>					
	ALS Negative					
JXJDAU	Molecular Serology Positive/Indicative		Molecular Serology Positive/Indicative		Molecular Serology Negative	
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
JYAL9W	LMG	Positive/Indicative	AP	Positive/Indicative	Amylase Radial Diffusion	Negative
	HemaTrace	Positive/Indicative	PSA	Positive/Indicative		
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		Conclusion: Negative
				Conclusion: Positive/Indicative		
	<u>Other Tests Performed</u>					
	ALS Negative					
JZDH48	KM	Positive/Indicative	AP	Positive/Indicative		
	HemaTrace	Positive/Indicative	PSA	Positive/Indicative		
		Conclusion: Positive/Indicative	ALS	Positive/Indicative		
			Micro	Negative		Conclusion: Not Tested
				Conclusion: Positive/Indicative		
K62L7A	KM	Positive/Indicative	ALS and Micro	Positive/Indicative		The laboratory has no validated procedure for Saliva testing.
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested
K68HAF	TMB	Positive/Indicative	AP	Positive/Indicative	Phadebas	Negative
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		Conclusion: Negative
				Conclusion: Positive/Indicative		

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
K93MAD	Phenolphthalein	Positive/Indicative	Acid Phosphatase	Positive/Indicative	RSID Saliva	Negative
			Microscopic Sperm	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	CrimeLite Blue light		Positive/Indicative			
	CrimeLite IR light		Positive/Indicative			
KDVHHD	KM	Positive/Indicative	AP	Positive/Indicative	Our lab doesn't test for Saliva, testing for Amylase is conducted	
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
	<u>Other Tests Performed</u>					
	Amylase Diffusion		Negative			
	ALS		Negative			
KG2YYG	O-tol	Positive/Indicative	AP	Positive/Indicative	Not relevant based on scenario	
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
KL4PPV	KM	Positive/Indicative	AP	Positive/Indicative	Not qualified	
	HemaTrace	Positive/Indicative	Micro	Positive/Indicative		
			PSA	Positive/Indicative		
			ALS	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
KPF9R	PTMB	Positive/Indicative	CS Proficiency for Blood screening only		CS Proficiency for Blood screening only	
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	
KPU8K8	Kastle Meyer	Positive/Indicative	Acid Phosphatase	Positive/Indicative	Phadebas paper Negative	
			Semen Microscopy	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
KUY28N	Hemastix	Positive/Indicative	AP	Positive/Indicative	RSID Negative	
	Luminol	Positive/Indicative	Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	Saliva/Semen - ALS blue		Positive/Indicative			
	Blood - ALS IR		Positive/Indicative			
KX88VB	LMG	Positive/Indicative	ALS	Positive/Indicative	Amylase Diffusion Negative	
			AP	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
KXNE3D	KM HemaTrace	Positive/Indicative Positive/Indicative	AP Micro PSA	Positive/Indicative Positive/Indicative Positive/Indicative	Phadebas	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
KYUEVT	Molecular Serology	Positive/Indicative	Molecular Serology	Positive/Indicative	Molecular Serology	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
KZMBDV	Hexagon	Positive/Indicative	Micro Christmas tree PSA	Positive/Indicative Positive/Indicative	SERATEC	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
LNCNP2	KM HemaTrace	Positive/Indicative Positive/Indicative	Acid Phosphatase PSA Microscopy	Positive/Indicative Negative Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
LUH87E	KM Luminol HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	AP Micro	Positive/Indicative Positive/Indicative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
LWYWA	KM	Positive/Indicative	AP Micro PSA	Positive/Indicative Positive/Indicative Positive/Indicative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	Blood - Infrared ALS Positive/Indicative					
	Semen/Saliva - ALS Negative					
LXYP8V	Visual LMG HemaTrace	Positive/Indicative Positive/Indicative Positive/Indicative	ALS AP RSID Micro	Negative Positive/Indicative Positive/Indicative Positive/Indicative	ALS RSID	Negative Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
M3BX8F	LMG	Positive/Indicative	AP Micro PSA	Positive/Indicative Positive/Indicative Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
MEKR4Q	KM HemDirect	Positive/Indicative Positive/Indicative	AP Micro ALS	Positive/Indicative Positive/Indicative Positive/Indicative	RSID Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
MKKFX7	ABACard HemaTrace	Positive/Indicative	ABACard p30	Positive/Indicative	RSID Saliva	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
MKTUPU	Visual	Positive/Indicative	ALS	Positive/Indicative	RSID	Negative
	LMG	Positive/Indicative	AP	Positive/Indicative		
	HemaTrace	Positive/Indicative	RSID	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
MRPREH	KM	Positive/Indicative	AP	Positive/Indicative	Phad (Micro tech)	Negative
	Mod Tak	Positive/Indicative	Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
MTZTZG	KM	Positive/Indicative	AP	Positive/Indicative	Phadebas Amylase	Negative
	Modified Takayama	Positive/Indicative	Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
MU7NDM	LMG	Positive/Indicative	AP	Positive/Indicative	scenario did not indicate need for testing	
	HemaTrace	Positive/Indicative	PSA	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
N4E6QJ	KM	Positive/Indicative	AP	Positive/Indicative	Phadebas	Negative
	Modified Takayama	Positive/Indicative	PSA	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
N7HDRH	LM (Leuco-Malachite)	Positive/Indicative	ALS	Negative	CSI does not test for Saliva	
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
NBP26J	LMG	Positive/Indicative	AP	Positive/Indicative	Amylase Radial Diffusion	Negative
	HemaTrace	Positive/Indicative	PSA	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	ALS Negative					
NE3WGH	PTMB	Positive/Indicative	AP	Positive/Indicative	Phadebas Tablet	Negative
	HemaTrace	Positive/Indicative	PSA	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
NJPXAF	KM	Positive/Indicative	AP	Positive/Indicative	RSID	Negative
	Luminol	Positive/Indicative	Micro	Positive/Indicative		
		Conclusion: Positive/Indicative	SERATEC PSA Semiquant	Positive/Indicative	Conclusion: Positive/Indicative	Conclusion: Negative
	<u>Other Tests Performed</u>					
	ALS	Positive/Indicative				
	Infrared ALS	Positive/Indicative				
NKEQ4T	HemaTrace (hHb)	Positive/Indicative	PSA	Positive/Indicative	is not within the scope of this laboratory, we only perform identification of human blood and semen (P30 enzyme)	
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative	Conclusion: Not Tested	
P3Y9WY	KM	Positive/Indicative	AP	Positive/Indicative	Conclusion: Not Tested	
	HemaTrace	Positive/Indicative	PSA	Positive/Indicative		
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		
P7VMXN	KM		AP		Based on case scenario not needed.	
	Visual Exam		Micro		Conclusion: Not Tested	
	ALS -staining darkened under UV		ALS -staining fluoresced			
		Conclusion: Positive/Indicative	Visual			
PF6RCZ	ALS	Positive/Indicative	ALS	Positive/Indicative	Our lab does not test for Saliva	
	O-tolidine	Positive/Indicative	AP	Positive/Indicative	Conclusion: Not Tested	
	OBTI	Positive/Indicative				
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		
PJ6CF6	KM	Positive/Indicative	ALS	Positive/Indicative	The [Laboratory] does not perform tests for Saliva	
	HemaTrace	Positive/Indicative	AP	Positive/Indicative	Conclusion: Not Tested	
		Conclusion: Positive/Indicative	PSA	Positive/Indicative		
			Micro	Positive/Indicative		
				Conclusion: Positive/Indicative		
PJZK2N	LMG	Positive/Indicative	AP	Positive/Indicative	Conclusion: Not Tested	
	HemaTrace	Positive/Indicative	PSA	Positive/Indicative		
		Conclusion: Positive/Indicative	Micro	Positive/Indicative		
			ALS	Positive/Indicative		
				Conclusion: Positive/Indicative		
PKKRVR	KM	Positive/Indicative	AP	Positive/Indicative	RSID	Negative
		Conclusion: Positive/Indicative	Micro	Positive/Indicative	Conclusion: Negative	

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
PLPWZU	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
PP7CBL	HemDirect Hemoglobin (SERATEC) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Christmas Tree Coloring <i>Positive/Indicative</i> Immunochromatographic Test (Semen-RSID) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase Test (SERATEC) <i>Negative</i> Conclusion: Negative
PQ6UUG	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID Saliva <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> ALS Blue <i>Positive/Indicative</i> ALS IR <i>Positive/Indicative</i>		
PVZNRH	KM <i>Positive/Indicative</i> Luminol <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA (Seratec PSA Semiquant) <i>Positive/Indicative</i> PSA (ABAcad p30) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID-Saliva <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> ALS <i>Positive/Indicative</i> IR ALS <i>Positive/Indicative</i>		
PWAG6Q	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase <i>Negative</i> Conclusion: Negative
PWE34X	KM <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> ALS <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
PXZMHV	KM <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
	<u>Other Tests Performed</u> Micro <i>Positive/Indicative</i>		

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
QDBYT8	Phenolphthalein	Positive/Indicative	Acid Phosphatase	Positive/Indicative	RSID Saliva	Negative
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	ALS (Blue) Positive/Indicative					
	ALS (IR) Positive/Indicative					
QGEG4T	KM	Positive/Indicative	AP	Positive/Indicative		
	HemaTrace	Positive/Indicative	PSA	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: No Response		Conclusion: No Response		Conclusion: Not Tested	
QH6HUY	KM Test	Positive/Indicative	Brentamine Test	Positive/Indicative	Phadebas Forensic Press Test	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
QKUQYU	KM	Positive/Indicative	AP	Positive/Indicative		
			Micro			
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
QLAK2U	KM	Positive/Indicative	AP	Positive/Indicative	No Saliva testing on commercial PTs.	
	HemaTrace	Positive/Indicative	PSA	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
QWM64W	KM	Positive/Indicative	ALS	Inconclusive		
			AP	Positive/Indicative		
			MOCRO	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
R4TNA8	O-tol	Positive/Indicative	ABAcad p30	Positive/Indicative	Test not performed as method verification for RSID field kit is currently being conducted.	
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
	<u>Other Tests Performed</u>					
	ALS - 450nm - Orange goggles - for Semen Negative					
R6CW3F	KM	Positive/Indicative	AP	Positive/Indicative	RSID	Negative
	HemaTrace	Positive/Indicative	PSA	Positive/Indicative		
			Microscopic Identification of Spermatozoa	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
RBLTBV			AP	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
RKF2NZ	Visual examination KM	Positive/Indicative Positive/Indicative	AP Micro	Positive/Indicative Positive/Indicative	Phadebas	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
RLN9YE	KM Luminol	Positive/Indicative Positive/Indicative	AP Micro PSA- Seratec PSA Semiquant	Positive/Indicative Positive/Indicative Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
	<u>Other Tests Performed</u>					
	IR Positive/Indicative ALS (blue light) Positive/Indicative					
RPRUWD	PTMB HemaTrace	Positive/Indicative Positive/Indicative	AP PSA	Positive/Indicative Positive/Indicative	does not indicate assailant licked chest area	
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
RTMQHX	KM	Positive/Indicative	AP Micro	Positive/Indicative Positive/Indicative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
RWNZ6K			ALS-blue AP Micro	Positive/Indicative Positive/Indicative Positive/Indicative		
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
T8MP4A	KM	Positive/Indicative	AP Micro ALS	Positive/Indicative Positive/Indicative Positive/Indicative	RSID	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
TFHUWB	Visual KM	Positive/Indicative Positive/Indicative	Visual ALS AP Micro	Negative Negative Positive/Indicative Positive/Indicative	Visual ALS RSID	Negative Negative Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
THH6LJ			AP Micro	Positive/Indicative Positive/Indicative		
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
TKQCQW	KM HemaTrace	Positive/Indicative Positive/Indicative	AP PSA Micro	Positive/Indicative Positive/Indicative Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
TLYEQA	KM <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested	Conclusion: Not Tested
TUKQGY	LMG <i>Positive/Indicative</i> RSID-Blood <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID-Saliva <i>Positive/Indicative</i> Conclusion: Positive/Indicative
TY7HZB	Conclusion: Not Tested	ALS-Blue <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
U77Y9V	KM HEXAGON OBTI Conclusion: Positive/Indicative	AP Conclusion: Positive/Indicative	Phadebas Conclusion: Negative
U8QUAN	KM <i>Positive/Indicative</i> Conclusion: No Response	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: No Response	Amylase paper <i>Negative</i> Conclusion: No Response
UBCY2C	KM <i>Positive/Indicative</i> Bluestar OBTI <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Saliva generally not tested unless specified. Conclusion: Not Tested
UETNJJ	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	This laboratory does not test for Saliva. Conclusion: Not Tested
UF93XK	Phenolphthalein KM Conclusion: Positive/Indicative	Acid Phosphatase Conclusion: Negative	Conclusion: Not Tested
	<u>Other Tests Performed</u> ALS at CSS wavelength w/ orange filter and 415 wavelength w/ yellow filter- interior and exterior searched <i>Negative</i>		
UGXDND	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> Microscopic examination for Sperm cells <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase <i>Inconclusive</i> Conclusion: Inconclusive
	<u>Other Tests Performed</u> ALS <i>Negative</i>		

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
UH277Z	LMG	Positive/Indicative	AP	Positive/Indicative		
			Micro	Positive/Indicative		
			PSA	Positive/Indicative		
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested
UZEWY4	KM	Positive/Indicative	ALS	Positive/Indicative		Not tested based on lab protocol
			PSA	Positive/Indicative		
			Micro	Positive/Indicative		
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested
V3FG72	KM	Positive/Indicative	AP	Positive/Indicative		RSID Negative
	luminol	Positive/Indicative	Micro	Positive/Indicative		
			PSA	Positive/Indicative		
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative
	<u>Other Tests Performed</u>					
	ALS Positive/Indicative					
	IR (infrared alternate light source) Positive/Indicative					
V9QJFB	PTMB	Positive/Indicative	AP	Positive/Indicative		Not indicated
			PSA	Positive/Indicative		
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested
VHPP6X	LMG		AP	Positive/Indicative		Phadebas Negative
		Conclusion: Positive/Indicative		Conclusion: No Response		Conclusion: No Response
	<u>Other Tests Performed</u>					
	H&E SLIDE : +3H Positive/Indicative					
	HEMDIRECT Positive/Indicative					
VKJRCT	Hemastix	Positive/Indicative	AP	Positive/Indicative		
	KM	Positive/Indicative	Micro	Positive/Indicative		
	Luminol	Positive/Indicative				
	HemaTrace	Positive/Indicative				
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested
WHLQB	KM	Positive/Indicative	AP	Positive/Indicative		Phadebas Negative
	Modified Takayama	Positive/Indicative	Micro	Positive/Indicative		
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative
W238Q9	KM	Positive/Indicative	AP	Positive/Indicative		RSID Negative
	Luminol	Positive/Indicative	Micro	Positive/Indicative		
			PSA	Positive/Indicative		
		Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative
	<u>Other Tests Performed</u>					
	IR Positive/Indicative					
	ALS Positive/Indicative					

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
WUV8A7	KM	Positive/Indicative	ALS	Negative		
	BLUE STAR IDENTI-HEM	Positive/Indicative	AP	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
WV79W6	ALS	Negative	ALS	Negative	ALS	Negative
	KM	Positive/Indicative	AP	Positive/Indicative	RSID	Negative
	RSID	Positive/Indicative	RSID	Positive/Indicative		
			MIC	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
WVXB33	LMG	Positive/Indicative	ALS	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
WWQCWW	LMG	Positive/Indicative	ALS	Negative	Amylase Diffusion	Negative
			AP	Positive/Indicative		
			Micro	Positive/Indicative		
			PSA	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
WZB9GG			AP	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
	<u>Other Tests Performed</u>					
	ALS - blue light Positive/Indicative					
X824XC	Phenolphthalein	Positive/Indicative	Acid Phosphatase	Positive/Indicative	Saliva RSID	Negative
	HemaTrace	Positive/Indicative	Microscopic Sperm search christmas tree staining	Positive/Indicative	Phadebas press test	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	ALS Positive/Indicative					
XAEG39	KM	Positive/Indicative	AP	Positive/Indicative	Phadebas	Negative
	Takayama	Positive/Indicative	Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
XE3WZ9	LMG	Positive/Indicative				
	Conclusion: No Response		Conclusion: Not Tested		Conclusion: Not Tested	
XLJL39	KM	Positive/Indicative	AP	Positive/Indicative	Saliva generally not tested unless specified.	
	Bluestar OBTI	Positive/Indicative	PSA	Positive/Indicative		
			RSID-Semen	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
XMEKXM	LMG	Positive/Indicative	ALS	Negative		
	Conclusion: Positive/Indicative		Conclusion: Negative		Conclusion: Not Tested	
XXDLD7	KM	Positive/Indicative	AP	Positive/Indicative	Phadebas paper	Negative
			H&E slide	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
Y74ZXA	KM	Positive/Indicative	AP	Positive/Indicative	RSID-Saliva	Negative
	HemaTrace	Positive/Indicative	Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
	<u>Other Tests Performed</u>					
	Body Fluid- Crime-lite 82S(r) Blue Positive/Indicative					
	Blood - Crime- lite(r) UV-Vis-IR Positive/Indicative					
Y9D7F4	LMG	Positive/Indicative	AP	Positive/Indicative		
			Micro	Positive/Indicative		
			PSA	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
YA8RBG	Visual	Positive/Indicative	ALS	Positive/Indicative	RSID	Negative
	LMG	Positive/Indicative	AP	Positive/Indicative		
	HemaTrace	Positive/Indicative	RSID	Positive/Indicative		
			Microscopic	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
YBF4NC			AP	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
YD44NM	KM	Positive/Indicative	AP	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
	<u>Other Tests Performed</u>					
	FLS Positive/Indicative					
YFVAQD	KM	Positive/Indicative	AP	Positive/Indicative	RSID	Negative
	HT	Positive/Indicative	p30	Positive/Indicative		
			Micro	Positive/Indicative		
			ALS	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
YQRBGG	LC HR/MS	Positive/Indicative	LC HR/MS	Positive/Indicative	LC HR/MS	Negative
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	

Table 1 - Item 5
Evidence - Tube Top containing Blood & Semen Mixture

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
Z6BHHU	LMG	Positive/Indicative	ALS	Positive/Indicative	Amylase Diffusion	Negative
			AP	Positive/Indicative		
			Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
ZC9797	KM	Positive/Indicative	AP	Positive/Indicative	Saliva generally not tested unless specified	
	Bluestar OBTI	Positive/Indicative	RSID-Semen	Positive/Indicative		
			PSA	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
ZD4X6C	KM	Positive/Indicative	AP	Positive/Indicative	RSID-Saliva	Negative
	HemDirect	Positive/Indicative	Micro	Positive/Indicative		
			ALS	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
ZFMHRM	O-tolidine	Positive/Indicative	ALS, 420-470	Positive/Indicative	No Saliva test on screening PT worksheet	
	Hex-OBTI	Positive/Indicative	A/P	Positive/Indicative		
	Control	Positive/Indicative	Control	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Not Tested	
ZFR8GY	KM	Positive/Indicative	AP	Positive/Indicative	Phadebas	Negative
	Modified Takayama	Positive/Indicative	Micro	Positive/Indicative		
	Conclusion: Positive/Indicative		Conclusion: Positive/Indicative		Conclusion: Negative	
ZUEEDZ	Hemastix	Positive/Indicative	UV screening	negative		
	Conclusion: Positive/Indicative		Conclusion: Not Tested		Conclusion: Not Tested	

Response Summary for Table 1 - Item 5: Tube Top containing Blood & Semen Mixture

	Blood	Semen	Saliva
Positive/Indicative	209	204	1
Negative	0	6	115
Inconclusive	0	0	1
Not Tested (NT)	11	10	104
No Response	4	4	3
Total	224	224	224

Screening Conclusions & Test(s) Performed

Indicate the overall conclusion from the results of the screening tests performed.

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
2289QN	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
<u>Other Tests Performed</u> Semen/Saliva - ALS <i>Positive/Indicative</i> Blood - ALS <i>Negative</i>			
22BP8V	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Microscopic <i>Positive/Indicative</i> p30-ABACard <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Negative</i> Conclusion: Negative
2G2TAX	Not tested due to absence of red-brown staining. Conclusion: Not Tested	Acid Phosphatase <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> Microscopy <i>Positive/Indicative</i> Alternate Light Source <i>Positive/Indicative</i> Conclusion: Positive/Indicative	SALigAE <i>Negative</i> Microscopy <i>Negative</i> Alternate Light Source <i>Positive/Indicative</i> Conclusion: Negative
2Q8KTZ	LMG <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
2WBMVF	Visual Examination <i>Negative</i> Conclusion: Negative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> RSID Semen <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID Saliva <i>Negative</i> Conclusion: Negative
2Y2QL2	LMG <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
34RKYT	physical exam = no red-brown = no testing Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA Conclusion: Positive/Indicative	NT based on scenario Conclusion: Not Tested

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
38HFFL	BLUESTAR (LUMINOL) <i>Negative</i>	AP <i>Positive/Indicative</i> SERATEC PSA SEMI QUANT <i>Positive/Indicative</i> RSID Semen <i>Positive/Indicative</i>	Phadebas <i>Negative</i>
	Conclusion: Negative	Conclusion: Positive/Indicative	Conclusion: Negative
<u>Other Tests Performed</u> CRIME LITE <i>Positive/Indicative</i>			
3C2TKZ	No visible staining - Not tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>	RSID Saliva <i>Negative</i>
	Conclusion: Not Tested	Conclusion: Positive/Indicative	Conclusion: Negative
<u>Other Tests Performed</u> Semen/Saliva - ALS <i>Positive/Indicative</i>			
3GUVNZ	ALS - IR Light <i>Negative</i>	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>	RSID <i>Negative</i>
	Conclusion: Negative	Conclusion: Positive/Indicative	Conclusion: Negative
<u>Other Tests Performed</u> Semen/Saliva - ALS <i>Positive/Indicative</i>			
3VAJDJ	KM - Hemophan (paper test for Blood erythrocytes, haemoglobin) by Erba Lachema - non-specific test for Blood <i>Negative</i> FOB - cassette test for determination of faecal occult Blood by NADAL - specific test for Human Blood <i>Negative</i>	AP - Phosphatesmo (paper test for the determination of acid phosphatase) by Macherey-Nagel <i>Positive/Indicative</i> PSA - cassette test for prostate specific antigen by SERATEC <i>Positive/Indicative</i>	
	Conclusion: Negative	Conclusion: Positive/Indicative	Conclusion: Not Tested
47AY3T	No red brown staining observed on item.	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>	Scenario details do not suggest need for this testing on this item.
	Conclusion: Not Tested	Conclusion: Positive/Indicative	Conclusion: Not Tested
4DB9ED	KM <i>Negative</i>	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>	RSID <i>Negative</i>
	Conclusion: Negative	Conclusion: Positive/Indicative	Conclusion: Negative
<u>Other Tests Performed</u> ALS (Blood, Saliva, Semen) <i>Positive/Indicative</i>			
4DP4V9	Molecular Serology <i>Negative</i>	Molecular Serology <i>Positive/Indicative</i>	Molecular Serology <i>Negative</i>
	Conclusion: Negative	Conclusion: Positive/Indicative	Conclusion: Negative
4GQCZH		AP <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i>	RSID
	Conclusion: Not Tested	Conclusion: Positive/Indicative	Conclusion: Negative

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
4R82UG	Conclusion: Not Tested	ALS Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
4W3QGF	KM Negative Conclusion: No Response	AP Positive/Indicative Micro Positive/Indicative Conclusion: No Response	Conclusion: Not Tested
	<u>Other Tests Performed</u> Visual FOR Blood Negative		
4XR299	No red deposits on tampon Conclusion: Not Tested	ALS Positive/Indicative AP Positive/Indicative Micro Sperm Positive/Indicative p30 Positive/Indicative Conclusion: Positive/Indicative	Test no longer performed Conclusion: Not Tested
62AAYQ	O-tol Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Scenario does not indicate oral contact with vaginal area. Conclusion: Not Tested
62P4BF	Conclusion: Not Tested	ALS Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
6AVWHW	TMB Negative Luminol Negative Conclusion: Negative	AP Positive/Indicative ALS Positive/Indicative Micro Positive/Indicative RSID Positive/Indicative Conclusion: Positive/Indicative	Phadebas paper test (Amylase) Negative Conclusion: Negative
6MWPPY	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
6QXKKV	no staining observed on Visual examination Negative Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative ALS Positive/Indicative Conclusion: Positive/Indicative	scenario does not indicate Saliva for this sample Conclusion: Not Tested
6R22DE	Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
6VNAP7	Visual Negative KM Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Not probative Conclusion: Not Tested
6ZW78F	Molecular Serology Negative Conclusion: Negative	Molecular Serology Positive/Indicative Conclusion: Positive/Indicative	Molecular Serology Negative Conclusion: Negative
6ZW9RD	No red-brown stains Conclusion: Not Tested	AP Positive/Indicative Sperm Search Positive/Indicative Conclusion: Positive/Indicative	Lab does not test for Saliva Conclusion: Not Tested
7ARU7X	ALS Negative KM Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
7LHP24	Visual Negative KM (stained area extract) Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	Amylase Diffusion Negative Conclusion: Negative
7N8RBX	KM Negative Conclusion: Negative	AP Positive/Indicative Mirco Positive/Indicative Conclusion: Positive/Indicative	[Laboratory] does not conduct Saliva tests. Conclusion: Not Tested
7ULBZK	Seratec HemDirect Negative OT Negative Conclusion: Negative	Seratec PSA Positive/Indicative AP Positive/Indicative Conclusion: Positive/Indicative	Amylase Diffusion Negative Conclusion: Negative
7YKAAU	Not tested based on case scenario and no stain observed Conclusion: Not Tested	AP Positive/Indicative PSA Positive/Indicative RSID-Semen Positive/Indicative Conclusion: Positive/Indicative	Not tested based on case scenario Conclusion: Not Tested
82LQ3H	Conclusion: Not Tested	Conclusion: Not Tested	Conclusion: Not Tested
83UJJM	HemaTrace Positive/Indicative Conclusion: Positive/Indicative	PSA Positive/Indicative Conclusion: Positive/Indicative	SALlgAE Negative Conclusion: Negative

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
8ERLLK	KM <i>Negative</i> Conclusion: Negative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	The [Laboratory] does not conduct forensic serological tests for Saliva. Conclusion: Not Tested
8RNLK4	Visual Examination <i>Negative</i> Conclusion: Negative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase Diffusion <i>Negative</i> Conclusion: Negative
8RT27A	KM <i>Negative</i> Conclusion: Negative	Brentamine -AP <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Negative</i> Conclusion: Negative
8TTFRM	No reddish-brown staining observed Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> ABA card PSA <i>Positive/Indicative</i> Seratec PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID Saliva <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Body Fluid - Test Performed: ALS <i>Positive/Indicative</i>		
8VWJK6	No Blood like stains could be observed Conclusion: Not Tested	Kind's reagent <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
8Y4DL2	Not tested due to observations: No visible signs of possible Blood. Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Not tested due to reported case circumstances and already identified Semen. Conclusion: Not Tested
8YL7UQ	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> p30-ABAcad <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Negative</i> Conclusion: Negative
8YVTV7	Rapid test for the detection of Human hemoglobin for use in legal medicine (Hexagon OBTI) <i>Negative</i> Conclusion: Negative	Sperm identification by Contrast Staining Christmas tree <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Seratec – Amylase Test <i>Negative</i> Conclusion: Negative
946QJ6	Conclusion: Not Tested	AP <i>Positive/Indicative</i> Spermatozoa Microscopic ID <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
963BVP	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas Negative Conclusion: Negative
98U8JJ	Micro Negative ALS Negative KM Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	ALS Positive/Indicative RSID Negative Conclusion: Negative
9GXGNN	Micro Negative KM Negative Conclusion: Negative	ALS Positive/Indicative FBB Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
9H2PM2	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
9H9FT6	Visual Negative O-tol Negative Conclusion: Negative	Visual Positive/Indicative ALS Positive/Indicative AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
9HBZPK	No RBS observed. Conclusion: Not Tested	ALS Positive/Indicative AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	[Laboratory] does not routinely test for Saliva. Conclusion: Not Tested
9KDWA9	Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
9KH9X6	Bluestar Obti Negative LGM Negative Conclusion: Negative	Seratec PSA Semiquant Positive/Indicative Conclusion: Positive/Indicative	Phadebas® Amylase Test Tablets Negative Conclusion: Negative
9T4NET	KM Negative Conclusion: Negative	Conclusion: Not Tested	Conclusion: Not Tested

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
9UXQ9H	NO RBS Conclusion: Not Tested	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Amylase PLATE Negative Conclusion: Negative
A4FUMV	KM Negative Conclusion: Negative	not in our scope Conclusion: Not Tested	not in our scope Conclusion: Not Tested
A4QTAY	Visual Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative RSID-Semen Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID-Saliva Negative Conclusion: Negative
AAEF3J	scenario Conclusion: Not Tested	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	not testing for Saliva Conclusion: Not Tested
AANRDT	No red/brown discoloration Conclusion: Not Tested	Acid Phosphatase Positive/Indicative Microscopic Sperm Search Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
<u>Other Tests Performed</u> Semen/Saliva - ALS Positive/Indicative			
AD6FYL	Not tested due to no red-brown staining present. Conclusion: Not Tested	Conclusion: Not Tested	Conclusion: Not Tested
ADC2PF	Visual Negative KM Negative Conclusion: Negative	ALS Positive/Indicative p30 Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
ADPV82	Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
AHKH3B	Conclusion: Not Tested	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
AJHZYK	not tested due to no visible RBS on tampon Conclusion: Not Tested	AP Positive/Indicative Microscopic exam Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
<u>Other Tests Performed</u> CrimeLite 82S Blue (with yellow goggles) Positive/Indicative			
APAZLT	No staining consistent with Blood was observed on the tampon, only yellow staining. Conclusion: Not Tested	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Scenario did not provide info to support testing for Saliva. Conclusion: Not Tested
AWH3H9	HemaTrace (hHb) Positive/Indicative Conclusion: Positive/Indicative	PSA Positive/Indicative Conclusion: Positive/Indicative	is not within the scope of this laboratory, we only perform identification of human blood and semen (P30 enzyme) Conclusion: Not Tested
BAYGU7	No red-brown staining Conclusion: Not Tested	AP Positive/Indicative p30 Positive/Indicative Microscopy Positive/Indicative ALS Positive/Indicative Conclusion: Positive/Indicative	Not tested per history Conclusion: Not Tested
BKPKDT	LMG Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
BPVRLH	Negative Conclusion: Negative	AP, PSA Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
BTGR6Q	Phenolphthalein Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
BXGUFK	No red coloring observed on tampon Conclusion: Not Tested	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Amylase Diffusion Negative Conclusion: Negative
BXPWFJ	Phenolphthalein Negative Conclusion: Negative	ALS Positive/Indicative BCIP Positive/Indicative p30 Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	ALS Positive/Indicative Phadebas Negative Conclusion: Negative

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
CFMNR6	Conclusion: Not Tested	ALS <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
CHEUWG	Conclusion: Not Tested	Conclusion: Not Tested	Conclusion: Not Tested
CLAAZR	FOB <i>Negative</i> Hemophan <i>Negative</i> Conclusion: Negative	PSA <i>Positive/Indicative</i> RSID Semen <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID Saliva <i>Negative</i> Conclusion: Negative
CLWA6M	LMG <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
CQ7C8H	Conclusion: Not Tested	ALS <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
CTEPNK	No dry red stain present Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> ALS <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Saliva not suspected to be on tampon based on scenario Conclusion: Not Tested
CU4ETT	Kastle Meyer <i>Negative</i> Hexagon OBTI <i>Negative</i> Conclusion: Negative	not in our scope Conclusion: Not Tested	not in our scope Conclusion: Not Tested
CXB6NT	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Negative</i> Conclusion: Negative
CZGMMP	Visual exam did not support chemical testing for Blood <i>Negative</i> Conclusion: Not Tested	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Scenario did not support testing for Saliva Conclusion: Not Tested
D7LP7V	no Blood testing conducted due to no observation of red/brown staining on the tampon Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	not testing for Saliva due to [Laboratory]'s policies and no allegations in this region Conclusion: Not Tested
	<u>Other Tests Performed</u> ALS (blue light) <i>Positive/Indicative</i>		

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
D8PPVU	Blood not tested based on case scenario and no stain observed Conclusion: Not Tested	AP Positive/Indicative RSID-Semen Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	Saliva not tested based on case scenario Conclusion: Not Tested
D9NLU6	Visual Examination Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative RSID Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	ALS Positive/Indicative RSID Negative Conclusion: Negative
D9VJDG	KM Negative Conclusion: Negative	Acid Phosphatase Positive/Indicative Sperm Microscopy Positive/Indicative Conclusion: Positive/Indicative	Phadebas Negative Conclusion: Negative
DGJNZM	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas Negative Conclusion: Negative
DHUHPV	Blood not tested based on case scenario and no stain observed. Conclusion: Not Tested	AP Positive/Indicative RSID-Semen Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	Saliva not tested based on case scenario. Conclusion: Not Tested
DKDWTV	Blood not tested based on case scenario and no stain observed. Conclusion: Not Tested	AP Positive/Indicative RSID-Semen Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	Saliva not tested based on case scenario. Conclusion: Not Tested
DRK8DB	Micro Negative KM Negative ALS Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative ALS Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
END392	Visual Negative O-tol Negative Conclusion: Negative	Visual Positive/Indicative ALS Positive/Indicative AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
EYWP69	KM Negative TMB Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	This analyst does not perform Saliva testing. Conclusion: Not Tested

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
F36YQ6	Conclusion: Not Tested	AP <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
	<u>Other Tests Performed</u> ALS <i>Positive/Indicative</i>		
FN7W9X	Conclusion: Not Tested	Christmas Tree Stain (Micro) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
FUZFGT	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Not in analyst's protocol to test for Saliva. Conclusion: Not Tested
FW4ZGL	not tested as no obvious rbs discoloration Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Blood - IR <i>Negative</i>		
FWLA6G	Blood testing not conducted due to no visible RBS Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> Crime-Lite 82S Blue ALS <i>Positive/Indicative</i>		
FYRKAW	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
FZGMQZ	HemDirect Hemoglobin (SERATEC) <i>Negative</i> Conclusion: Negative	Christmas Tree Stain Slide (Micro) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase Test (SERATEC) <i>Negative</i> Conclusion: Negative
G2DR8Y	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
G6XGVH	There did not appear to be any reddish staining on tampon Conclusion: Not Tested	AP <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Our Crime Scene unit does not test for Saliva Conclusion: Not Tested
GC99WW	Molecular Serology <i>Negative</i> Conclusion: Negative	Molecular Serology <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Molecular Serology <i>Negative</i> Conclusion: Negative

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
GCA8Q3	Phenol <i>Negative</i> Conclusion: Negative	PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
	<u>Other Tests Performed</u> ALS light - stains <i>Positive/Indicative</i>		
GD34A6	Conclusion: Not Tested	ALS <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
GXYAE	Synopsis doesn't indicate a reason to test for this Conclusion: Not Tested	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Synopsis doesn't indicate a reason to test for this Conclusion: Not Tested
H24RD4	Conclusion: Not Tested	ALS <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
H8YVMF	No DRS observed Conclusion: Not Tested	PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> ALS <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID-Saliva <i>Negative</i> Conclusion: Negative
H9DDT8	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
HJFLMB	Hemastix <i>Negative</i> Conclusion: Negative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Microscope slide <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Exhibit type and case circumstances Conclusion: Not Tested
HK3CJA	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
HQ68XW	SERATEC HemDirect <i>Negative</i> Conclusion: Negative	Christmas Tree <i>Positive/Indicative</i> Conclusion: Positive/Indicative	SERATEC Amylase Test <i>Negative</i> Conclusion: Negative

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
J74M6N	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas Amylase <i>Negative</i> Conclusion: Negative
J9AY8J	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
JBD4EC	Visual, no red color Conclusion: Not Tested	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase Diffusion <i>Negative</i> Conclusion: Negative
JF76CK	 Conclusion: Not Tested	ALS-Blue <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
JGETMY	 Conclusion: Not Tested	ALS <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
JKBG7R	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
JRX62W	No stains consistent with Blood were observed Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> ALS <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
JU8UHJ	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
JUFNLR	KM <i>Negative</i> ALS <i>Negative</i> Conclusion: Negative	PSA <i>Positive/Indicative</i> ALS <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
JWCGWJ	Blood not tested based on case scenario and no stain observed Conclusion: Not Tested	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Saliva not tested based on case scenario Conclusion: Not Tested

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
JXBCVT	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
JXJDAU	Molecular Serology <i>Negative</i> Conclusion: Negative	Molecular Serology <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Molecular Serology <i>Negative</i> Conclusion: Negative
JYAL9W	LMG <i>Positive/Indicative</i> HemaTrace <i>Positive/Indicative</i> Conclusion: Positive/Indicative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase Radial Diffusion <i>Negative</i> Conclusion: Negative
JZDH48	 Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> ALS <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
K62L7A	No brown stained observed. Conclusion: Not Tested	ALS and Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	The laboratory has no validated procedure for Saliva testing. Conclusion: Not Tested
K68HAF	 Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
K93MAD	No red-brown staining Conclusion: Not Tested	Acid Phosphatase <i>Positive/Indicative</i> Microscopic Sperm <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID Saliva <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> CrimeLite Blue light <i>Positive/Indicative</i>		
KDVHHD	no r/b staining Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Our lab doesn't test for Saliva, testing for Amylase is conducted Conclusion: Not Tested
	<u>Other Tests Performed</u> Amylase Diffusion <i>Negative</i> ALS <i>Positive/Indicative</i>		
KG2YYG	No visible staining Conclusion: Not Tested	PSA (p30) <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Possible presence of urine interferes with Saliva testing. Conclusion: Not Tested

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
KL4PPV	No Blood like staining observed Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative ALS Positive/Indicative Conclusion: Positive/Indicative	Not qualified Conclusion: Not Tested
KPFR9R	No red stains visible Conclusion: Not Tested	CS Proficiency for Blood screening only Conclusion: Not Tested	CS Proficiency for Blood screening only Conclusion: Not Tested
KPU8K8	Kastle Meyer Negative Conclusion: Negative	Acid Phosphatase Positive/Indicative Semen Microscopy Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
KUY28N	no red/brown staining observed, no Blood testing performed Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
KX88VB	Item did not have red/brown discoloration Conclusion: Not Tested	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Amylase Diffusion Negative Conclusion: Negative
KXNE3D	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	Phadebas Negative Conclusion: Negative
KYUEVT	Molecular Serology Negative Conclusion: Negative	Molecular Serology Positive/Indicative Conclusion: Positive/Indicative	Molecular Serology Negative Conclusion: Negative
KZMBDV	Hexagon Negative Conclusion: Negative	Micro Christmas tree Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	Seratec Negative Conclusion: Negative
LNCNP2	Conclusion: Not Tested	Acid Phosphatase Positive/Indicative PSA Positive/Indicative Microscopy Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
LUH87E	No visible staining observed, item white in color Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
LWYWA	No staining observed Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
	<u>Other Tests Performed</u> Blood - Infrared ALS Negative Semen/Saliva - ALS Positive/Indicative		
LXYP8V	Visual Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative RSID Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	ALS Positive/Indicative RSID Negative Conclusion: Negative
M3BX8F	LMG Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	 Conclusion: Not Tested
MEKR4Q	No Visual indications of Blood Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative ALS Positive/Indicative Conclusion: Positive/Indicative	RSID Saliva Negative Conclusion: Negative
MKKFX7	ABACard HemaTrace Negative Conclusion: Negative	ABACard p30 Positive/Indicative Conclusion: Positive/Indicative	RSID Saliva Negative Conclusion: Negative
MKTUPU	Visual Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative RSID Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
MRPREH	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phad (Micro tech) Negative Conclusion: Negative
MTZTZG	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas Amylase Negative Conclusion: Negative
MU7NDM	staining w/Visual appearance of Blood not observed Conclusion: Not Tested	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	scenario did not indicate need for testing Conclusion: Not Tested

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
N4E6QJ	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Negative</i> Conclusion: Negative
N7HDRH	ALS <i>Negative</i> Conclusion: Negative	AP, ALS <i>Positive/Indicative</i> Conclusion: Positive/Indicative	CSI does not test for Saliva Conclusion: Not Tested
NBP26J	Scenario does not indicate the need for this test. Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Scenario does not indicate the need for this test. Conclusion: Not Tested
	<u>Other Tests Performed</u> ALS <i>Positive/Indicative</i>		
NE3WGH	NT per case scenario. Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	NT per case scenario. Conclusion: Not Tested
NJPXAF	No red-brown staining observed on item. Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> ALS <i>Positive/Indicative</i>		
NKEQ4T	HemaTrace (hHb) <i>Positive/Indicative</i> Conclusion: Positive/Indicative	PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	is not within the scope of this laboratory, we only perform identification of human blood and semen (P30 enzyme) Conclusion: Not Tested
P3Y9WY	No apparent staining present. Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
P7VMXN	Visual exam- staining not consistent with possible Blood. ALS -staining not consistent with possible Blood. Conclusion: Not Tested	AP Micro ALS- staining fluoresced Visual Exam -staining consistent with possible Semen. Conclusion: Positive/Indicative	Based on case scenario not needed. Conclusion: Not Tested
PF6RCZ	ALS <i>Negative</i> Conclusion: Negative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Our lab does not test for Saliva Conclusion: Not Tested

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
PJ6CF6	Not tested due to the appearance Conclusion: Not Tested	ALS Positive/Indicative AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	The [Laboratory] does not perform tests for Saliva Conclusion: Not Tested
PJZK2N	No stains Visually consistent with Blood were observed. Conclusion: Not Tested	ALS Positive/Indicative AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
PKKRVR	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative ALS Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
PLPWZU	Conclusion: Not Tested	ALS Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
PP7CBL	HemDirect Hemoglobin (SERATEC) Negative Conclusion: Negative	Christmas Tree Coloring Positive/Indicative Immunochromatographic Test (Semen-RSID) Positive/Indicative Conclusion: Positive/Indicative	Amylase Test (SERATEC) Negative Conclusion: Negative
PQ6UUG	no red brown stains observed Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	not indicated by scenario Conclusion: Not Tested
PVZNRH	no reddish brown staining observed Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID-Saliva Negative Conclusion: Negative
	<u>Other Tests Performed</u> ALS Positive/Indicative		
PWAG6Q	Fibre Optic light Negative KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Not suggested by case scenario Conclusion: Not Tested
PWE34X	KM Negative Micro Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative ALS Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
PXZMHV	Conclusion: Not Tested	AP <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
	<u>Other Tests Performed</u> Micro <i>Positive/Indicative</i>		
QDBYT8	no staining consistent with Blood Conclusion: Not Tested	Acid Phosphatase <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID Saliva <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> ALS (Blue) <i>Positive/Indicative</i>		
QGEQ4T	no visible staining Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: No Response	Conclusion: Not Tested
	KM Test <i>Negative</i> Conclusion: Negative	Brentamine Test <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas Forensic Press Test <i>Negative</i> Conclusion: Negative
QKUQYU	Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
	No visible staining. Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	No Saliva testing on commercial PTs. Conclusion: Not Tested
QWM64W	KM <i>Negative</i> Conclusion: Negative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
	No brown/red staining present on the tampon. Conclusion: Not Tested	ABACard p30 <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Test not performed as method verification for RSID field kit is currently being conducted. Conclusion: Not Tested
R6CW3F	Not tested due to absence of red brown staining Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Microscopic Identification of Spermatozoa <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
RBLTBV	Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
RKF2NZ	Visual examination Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Phadebas Negative Conclusion: Negative
RLN9YE	no reddish-brown staining/possible Blood observed Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
	<u>Other Tests Performed</u> ALS (blue light) Positive/Indicative		
RPRUWD	no red staining indicating Blood Conclusion: Not Tested	AP Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	no indication of oral sex performed on the victim Conclusion: Not Tested
RTMQHX	KM Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
RWNZ6K	Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
T8MP4A	Visual exam Negative Conclusion: Not Tested	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	not supported by scenario Conclusion: Not Tested
TFHUWB	Visual Negative no staining consistent with Blood observed Conclusion: Not Tested	Visual Positive/Indicative ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	no documentation of oral contact Conclusion: Not Tested
THH6LJ	Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
TKQCQW	KM Negative Conclusion: Negative	AP Positive/Indicative PSA Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
TLYEQA	KM <i>Negative</i> Conclusion: Negative	Conclusion: Not Tested	Conclusion: Not Tested
TUKQGY	LMG <i>Negative</i> RSID-Blood <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID-Saliva <i>Inconclusive</i> Conclusion: Inconclusive
TY7HZB	Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
U77Y9V	lab policy, case scenario Conclusion: Not Tested	AP Conclusion: Positive/Indicative	lab policy, case scenario Conclusion: Not Tested
U8QUAN	No Bloodlike stains observed Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: No Response	Amylase paper <i>Negative</i> Conclusion: No Response
UBCY2C	Blood not tested based on case scenario and no stain observed. Conclusion: Not Tested	AP <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Saliva not tested based on case scenario. Conclusion: Not Tested
UETNJJ	LMG <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	This laboratory does not test for Saliva. Conclusion: Not Tested
UF93XK	Phenolphthalein <i>Negative</i> Visual Exam <i>Negative</i> Conclusion: No Response	PSA test <i>Positive/Indicative</i> Acid Phosphatase <i>Negative</i> Conclusion: No Response	Conclusion: Not Tested
	<u>Other Tests Performed</u> ALS CSS wavelength w/ orange filter- fluoresced <i>Positive/Indicative</i>		
UGXDND	LMG <i>Negative</i> HemaTrace <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> p30 <i>Positive/Indicative</i> Microscopic examination for Sperm cells <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase <i>Inconclusive</i> Conclusion: Inconclusive
	<u>Other Tests Performed</u> ALS <i>Positive/Indicative</i>		

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
UH277Z	LMG <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
UZEWY4	Not tested due to no dried red staining observed Conclusion: Not Tested	ALS <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Not applicable based on case scenario Conclusion: Not Tested
V3FG72	no RBS observed on tampon Conclusion: Not Tested <u>Other Tests Performed</u> ALS <i>Positive/Indicative</i>	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	RSID <i>Negative</i> Conclusion: Negative
V9QJFB	Not indicated Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Not indicated Conclusion: Not Tested
VHPP6X	LMG Conclusion: Negative <u>Other Tests Performed</u> H&E SLIDE: +4H/T <i>Positive/Indicative</i>	AP Conclusion: Positive/Indicative	Phadebas Conclusion: Negative
VKJRCT	no relevant staining Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested
VHLQB	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Negative</i> Conclusion: Negative
W238Q9	no apparent rbs observed Conclusion: Not Tested <u>Other Tests Performed</u> ALS <i>Positive/Indicative</i>	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Per case details and due to positive Semen testing, Saliva testing not performed. Conclusion: Not Tested
WUV8A7	KM <i>Negative</i> Conclusion: No Response	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Conclusion: Positive/Indicative	 Conclusion: Not Tested

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
WV79W6	ALS <i>Positive/Indicative</i> KM <i>Negative</i> RSID <i>Negative</i> Conclusion: Negative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> RSID <i>Positive/Indicative</i> MIC <i>Positive/Indicative</i> Conclusion: Positive/Indicative	ALS <i>Positive/Indicative</i> RSID <i>Negative</i> Conclusion: Negative
WVXB33	LMG <i>Negative</i> Conclusion: Negative	ALS <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
WWQCWW	LMG <i>Negative</i> Conclusion: Negative	ALS <i>Positive/Indicative</i> AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Amylase Diffusion <i>Negative</i> Conclusion: Negative
WZB9GG	Conclusion: Not Tested	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
X824XC	Phenolphthalein <i>Negative</i> Conclusion: Negative	Acid Phosphatase <i>Positive/Indicative</i> Microscopic Sperm search christmas tree staining <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Saliva RSID <i>Negative</i> Phadebas press test <i>Negative</i> Conclusion: Negative
	<u>Other Tests Performed</u> ALS <i>Positive/Indicative</i>		
XAEG39	KM <i>Negative</i> Conclusion: Negative	AP <i>Positive/Indicative</i> Micro <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Phadebas <i>Negative</i> Conclusion: Negative
XE3WZ9	No Blood like stains observed Conclusion: Not Tested	Conclusion: Not Tested	Conclusion: Not Tested
XLJL39	Blood not tested based on case scenario and no stain observed. Conclusion: Not Tested	AP <i>Positive/Indicative</i> PSA <i>Positive/Indicative</i> RSID-Semen <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Saliva not tested based on case scenario. Conclusion: Not Tested
XMEKXM	Conclusion: Not Tested	ALS <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Conclusion: Not Tested
XXDLD7	No obvious red/brown staining, therefore no KM conducted Conclusion: Not Tested	AP <i>Positive/Indicative</i> H&E slide <i>Positive/Indicative</i> Conclusion: Positive/Indicative	Not part of examination strategy, given case circumstances. Conclusion: Not Tested

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD	SEMEN	SALIVA
	Tests Performed Results	Tests Performed Results	Tests Performed Results
Y74ZXA	none done due to no red/brown staining present Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	RSID-Saliva Negative Conclusion: Negative
<u>Other Tests Performed</u> Body Fluid- Crime-lite 82S Blue ALS Positive/Indicative			
Y9D7F4	LMG Negative Conclusion: Negative	AP Positive/Indicative Micro Positive/Indicative PSA Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
YA8RBG	Visual Negative Conclusion: Negative	ALS Positive/Indicative AP Positive/Indicative RSID Positive/Indicative Microscopic Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
YBF4NC	Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
YD44NM	no visible stain Conclusion: Not Tested	AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Conclusion: Not Tested
<u>Other Tests Performed</u> FLS Positive/Indicative			
YFVAQD	KM Negative Conclusion: Negative	AP Positive/Indicative p30 Positive/Indicative Micro Positive/Indicative ALS Positive/Indicative Conclusion: Positive/Indicative	RSID Negative Conclusion: Negative
YQRBGG	LC HR/MS Negative Conclusion: Negative	LC HR/MS Positive/Indicative Conclusion: Positive/Indicative	LC HR/MS Negative Conclusion: Negative
Z6BHHU	No red/brown color observed Conclusion: Not Tested	ALS Positive/Indicative AP Positive/Indicative Micro Positive/Indicative Conclusion: Positive/Indicative	Amylase Diffusion Negative Conclusion: Negative

Table 1 - Item 6
Evidence - Tampon containing Semen

WebCode	BLOOD		SEMEN		SALIVA	
	Tests Performed	Results	Tests Performed	Results	Tests Performed	Results
ZC9797	Blood not tested based on case scenario and no stain observed		AP Positive/Indicative RSID-Semen Positive/Indicative PSA Positive/Indicative		Saliva not tested based on case scenario	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Not Tested	
ZD4X6C	No Visual indications of Blood, no red brown staining noted so item was not tested for Blood.		AP Positive/Indicative Micro Positive/Indicative ALS Positive/Indicative		RSID-Saliva Negative	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: Negative	
ZFMHRM	No visible R/B		ALS, 420-470 Positive/Indicative A/P Positive/Indicative Control Positive/Indicative		No Saliva test on screening PT worksheet	
	Conclusion: Not Tested		Conclusion: Positive/Indicative		Conclusion: No Response	
ZFR8GY	KM Negative		AP Positive/Indicative Micro Positive/Indicative		Phadebas Negative	
	Conclusion: Negative		Conclusion: Positive/Indicative		Conclusion: Negative	
ZUEEDZ	item not applicable for crime scene processing		item not applicable for crime scene processing			
	Conclusion: Not Tested		Conclusion: Not Tested		Conclusion: Not Tested	

Response Summary for Table 1 - Item 6: Tampon containing Semen

	Blood	Semen	Saliva
Positive/Indicative	4	210	0
Negative	104	0	94
Inconclusive	0	0	2
Not Tested (NT)	113	10	126
No Response	3	4	2
Total	224	224	224

Additional Comments

TABLE 2

WebCode	Additional Comments
2G2TAX	Item 4 swab was tested for urea. The test was inconclusive due to a false positive result. For urea tests using bromothymol blue-infused filter paper and a urease solution, our manual states: "A blue color change at this step (after adding one drop of deionized water onto the suspected urine stain and filter paper, but before adding one drop of the urease solution) indicates a false positive result."
2Q8KTZ	My laboratory does not perform screening tests for saliva.
2Y2QL2	My laboratory does not perform tests for saliva.
47AY3T	For item 2 - CTS paperwork description says Four replicate white handled swabs but item was actually four pink handled swabs.
4W3QGF	Where visual / odour tests have been carried and negative (for blood / urine / faeces), no further presumptive testing was carried out. Mini-taping and fibre taping was carried out on Item 5. Visual check for possible hairs & fibres for Items 3 & 4 were carried out - none were seen.
62AAYQ	Under Item 2 heading, it says: Four replicate white handled swabs.
6VNAP7	Not probative is based on scenario information provided. Saliva indicated is based on detecting levels of amylase equal or greater than the 1:200 saliva dilution standard.
7LHP24	Item 2 = Amylase activity detected halo (There were very low levels of amylase detected which are not indicative of the presence of saliva).
7N8RBX	FURTHER EXPLANATION OF RESULTS: i. Blood was identified on Q02, Q03, and Q05. ii. No blood was found on Q01, Q04, Q06, and Q07. iii. Semen was identified on Q01, Q02, Q03, Q04, Q05, and Q06. iv. No semen was found on Q07. v. Appropriate samples are being submitted for DNA testing. If a DNA profile is developed, known oral swabs are needed from the victim and any developed suspect(s). vi. RETAINED SAMPLES: 1. Item(s) Q01-Q07. 2. Cutting(s) from Q05 and Q06. 3. Tape lift(s) from Q05
82LQ3H	No test for saliva validated for Field Investigation Unit (FIU) use. Non-swab Items 5 and 6 not tested because this type of testing is performed in the laboratory and not by FIU field responders.
8ERLLK	The [Laboratory] does not conduct forensic serological tests for saliva. Therefore, items 1, 2, 3, 4, 5, and 6 did not have saliva tests performed on them.
8RNLK4	Swab packet (substrate control) noted, returned to evidence packet, not itemized or examined.
8YVTV7	There is no observation or suggestion of the case.
9H2PM2	Saliva was not tested because our Forensics Lab does not test for saliva.
9H9FT6	Laboratory does not confirm the presence of saliva.
9KDWA9	This test is being completed as part of a male screening protocol and only tests that may be performed by the analyst in casework were performed on this proficiency test.
9KH9X6	Item 4: Bluestar Obti positive result is a weak positive.
9T4NET	This laboratory does not perform semen testing. Any evidence items from a suspected sexual assault would be forwarded to the state lab for semen testing. This laboratory does not perform saliva testing. Any evidence that is suspected of having saliva on it would be forwarded for DNA testing.
A4QTAY	No hairs were observed on the Tube top or on the tampon.
AAEF3J	RESULTS AND CONCLUSIONS: No blood was found on Items 1 or 4. Items 2, 3, and 5 contain blood. Presumptive* testing for human origin was positive. No semen was found on Items 1 or 2. Items 4, 5, and 6 contain semen. No further analysis was performed. *presumptive tests do not confirm the presence of bodily fluids or human origin METHODS: Phenolphthalein Test – Items 1, 2, 3, 4, and 5. HemaTrace® Test – Items 2, 3, and 5. Acid Phosphatase Test – Items 1, 2, 4, 5, and 6 p30 Test – Items 1, 2, 4, 5, and 6. Microscopy – Items 1, 2, 4, 5, and 6. Results relate only to the items or sampling of items tested.
AD6FYL	Items not tested for semen or saliva. Individual is only authorized to perform blood testing.
ADC2PF	Our laboratory does not test for saliva.

TABLE 2

WebCode	Additional Comments
ADPV82	This test is being completed as part of a male screening protocol and only tests that may be performed by the analyst in casework were performed on this proficiency test.
AHKH3B	Items 2, 3, and 5 contain blood. Presumptive* testing for human origin was positive. No blood was found on Item 1. Items 4, 5, and 6 contain semen. No semen was found on Items 1 or 2. No further analysis was conducted. *presumptive tests do not confirm the presence of bodily fluids or human origin
AWH3H9	The search for biological saliva fluid was not performed, since it is outside the scope of accreditation of this laboratory. Within the scope of accreditation there is only the identification of human blood and semen (P30 enzyme)
BAYGU7	No blood, semen, seminal fluid, sperm, or saliva was detected on the tested control swab. Per our technical procedures a minimal or partial clearing on an amylase radial diffusion test is to be reported as "inconclusive" for amylase.
BKPKT	My laboratory does not perform screening tests for saliva.
BXGUFK	I recognize the financial practicality of using a tampon as an evidence item (Item 6- One tampon worn by the victim after the assault); however, I was disappointed in how the sample was prepared for this proficiency. A single fluorescent semen stain on one side of the tampon is not realistic. I understand that you want analysts to pass the proficiencies, but if that is the case bring back socks (white socks).
CHEUWG	No test for saliva validated for Field Investigation Unit (FIU) use. Non-swab Items 5 and 6 not tested because this type of testing is performed in the laboratory and not by FIU field responders
CLWA6M	My laboratory does not perform screening tests for saliva.
CQ7C8H	Our laboratory have no validated methods for saliva testing. Kastle-Meyer Test for Blood was not performed on Item 1A (original Item 1) and 1D (original Item 4), and Sample 1F1 (original Item 6) because they were not brown stained. Visual examination with an Alternate Light Source was performed on Items 1E (original Item 5) and 1F (original Item 6). Kastle-Meyer Test for Blood Presumptive testing detected the possible presence of blood on items 1B (original Item 2) and 1C (original Item 3), and Sample 1E1 (original Item 5). Microscopic Examination for Spermatozoa and p30 Testing Spermatozoa were observed on Item 1D (original Item 4), and Samples 1E1 (original Item 5) and 1F1 (original Item 6), which confirms the presence of semen. No spermatozoa was observed on Items 1A (original Item 1), 1B (original Item 2), and 1C (original Item 3).
D9NLU6	All microscopic examinations were performed after the Kernechtrot-Picroindigocarmine staining procedure. No hairs were observed on the item 5 "tube top" or the item 6 tampon.
END392	This laboratory does not confirm the presence of saliva.
FW4ZGL	Items 1, 4 and 6 not tested for blood as no red brown staining was noted on items.
FZGMQZ	- Item 1: Four (4) swabs are POSITIVE for saliva and NEGATIVE for semen and blood. - Item 2: Four (4) swabs are POSITIVE for blood and NEGATIVE for semen and saliva. - Item 3: Four (4) swabs are POSITIVE for blood and NEGATIVE for semen and saliva. - Item 4: Four (4) swabs are POSITIVE for semen and saliva and NEGATIVE for blood. - Item 5: One tube top worn by victim during the assault is POSITIVE for blood and semen and NEGATIVE for saliva. - Item 6: One tampon worn by victim after the assault is POSITIVE for semen and NEGATIVE for blood and saliva.
H9DDT8	Items 4, 5, and 6 contain semen. No semen was found on Items 1, 2, or 3. Items 2, 3, and 5 contain blood. Presumptive* testing for human origin was positive. No blood was found on Items 1, 4, or 6. No further analysis was conducted. *presumptive tests do not confirm the presence of bodily fluids or human origin
HK3CJA	Items 2, 3, and 5 contain blood. Presumptive* testing for human origin was positive. No blood was found on Items 1 or 6. Items 4, 5, and 6 contain semen. No semen was found on Items 1 or 2. No further analysis was conducted. Results relate only to the items or sampling of items tested. *presumptive tests do not confirm the presence of bodily fluids or human origin

TABLE 2

WebCode	Additional Comments
HQ68XW	Item 1: Four swab with white plastic handle is POSITIVE for saliva and Negative for semen and blood. Item 2: Four swab with pink plastic handle is POSITIVE for blood and Negative for Semen and Saliva. Item 3: Four swab with blue plastic handle is POSITIVE for blood and Negative for Semen and Saliva. Item 4: Four swab with yellow plastic handle is POSITIVE for Semen and Saliva and Negative for blood. Item 5: One tube top worn by victim is POSITIVE for Semen and Blood, and Negative for Saliva. Item 6: One tampon worn by victim is POSITIVE for Semen and Negative for Saliva and Blood.
JF76CK	This test is being completed as part of a male screening protocol and only tests that may be performed by the analyst in casework were performed on this proficiency test.
JKBG7R	Our lab does not test for the presence of saliva.
JUFNLR	None of the items were tested for saliva. The [Laboratory] does not presumptive test for saliva at this time.
K62L7A	Presumptive testing detected the possible presence of blood on Items 1B (Item 2), 1C (Item 3) and Sample 1E1 (Item 5). Spermatozoa were observed on Item 1D (Item 4), sample 1E1 (Item 5) and sample 1F1 (Item 6), which confirms the presence of semen. Semen was not detected on Items 1A (Item 1), 1B (Item 2) and 1C (Item 3). Visual examination with an alternate light source was performed on Items 1E (Item 5) and 1F (Item 6).
K93MAD	Item 5 had one additional area of very bright spot fluorescence that was tested with Acid Phosphatase, resulting negative.
KPU8K8	Laboratory records show that on 11 April 2023 a number of items were received for examination in relation to this case. I understand from accompanying documentation that the above items were submitted in relation to an investigation into an alleged sexual assault. The results and conclusions provided in this statement form my expert opinion, which is based on my scientific knowledge, experience and training. The results apply to the items as received and relate only to the items tested. I have received assistance from scientific support staff for the examinations completed in this case. The reporting of blood in this statement is based on its appearance and a positive result for a presumptive chemical test for blood. A presumptive test is used to screen for the presence of a substance, but it cannot confirm its presence. The items received were four sets of four swabs (items 1 – 4), a tube top (item 5) and a used tampon (item 6). Only one swab from the set of four swabs were tested for items 1, 3 and 4. All four swabs were tested in item 2. The swabs in items 2 and 3 were bloodstained. No bloodstaining was observed on swabs in item 1 and 4. A chemical test for the presence of semen was undertaken on the swabs in items 1, 2 and 3. No semen was detected on these swabs. A microscope slide was prepared from one of the swabs in item 4 as I understand from accompanying documentation that this swab was taken from the victim's stomach where the suspect allegedly ejaculated. This slide was examined for the presence of spermatozoa (sperm), the cellular constituent of semen. Sperm were observed on this slide indicating the presence of semen on the corresponding swab. The swabs in items 1, 2, 3 and 4 were examined for the presence of amylase. Amylase is a substance which occurs at high levels in saliva, however it can also be found at much lower levels in other body fluids such as vaginal secretions. Amylase was detected on the swabs in items 1 and 4, indicating that saliva is probably present on these swabs. No amylase, hence no saliva was detected on the swabs in items 2 and 3. The tube top was white and in an overall good condition. A bloodstain was detected on the front left side of the top. No other bloodstaining was observed on the top. The top was examined for the presence of semen. A semen stain was detected on the front left side of the top. This semen stain corresponded with the location of the bloodstain. The top was also examined for the presence of amylase. No amylase, hence no saliva was detected on the top. Pale yellow staining was observed on the outer surfaces of the tampon. No obvious bloodstaining was observed on this tampon. A small orange stain observed on the tampon was chemically tested for blood. No blood was detected on this testing. The outer surfaces of the tampon were examined for the presence of semen. The upper end of the tampon, opposite to the string end, was semen stained.
LNCNP2	Items 2, 3, and 5 contain blood. Presumptive* testing for human origin was positive. No blood was found on Item 1. Items 4, 5, and 6 contain semen. No semen was found on Items 1 or 2.
LVVYWA	ABAcad PSA test was used for item 4. Seratec PSA test was used for item 5.
LXYP8V	No hairs were observed on the tube top or tampon.

TABLE 2

WebCode	Additional Comments
M3BX8F	My laboratory does not perform screening tests for saliva.
MEKR4Q	Our laboratory tests for blood, seminal fluid, and saliva. I used all of the tests our laboratory employs during this test. Normally, much of our testing is scenario-based and therefore, there are some test that I performed on this test that I would not have performed during normal testing. I have recorded them in my notes kept at our laboratory. I will not document them all here.
MU7NDM	Amylase Radial Diffusion: A faint/partial clearing was observed in the well around item 02; therefore, the result is inconclusive. Amylase is present in body fluids including blood, which this item tested positive for.
NBP26J	Amylase Radial Diffusion: A faint/partial clearing was observed in the well around item 2; therefore, the result is inconclusive. Amylase is present in other body fluids, including blood, which this swab also tested positive for.
NKEQ4T	The search for biological saliva fluid was not performed, since this study is not within the scope of accreditation of this laboratory. Within the scope of accreditation there is only the identification of human blood and semen (P30 enzyme)
P3Y9WY	Items 2, 3, and 5 contain blood. Presumptive* testing for human origin was positive. Items 4, 5, and 6 contain semen. No blood was found on Items 1 or 4. No semen was found on Items 1, 2, or 3. Results relate only to the items or sampling of items tested. *presumptive tests do not confirm the presence of bodily fluids or human origin
P7VMXN	Item 2, based on the visual exam, the staining present was consistent with possible blood and not semen or saliva; without a case scenario to indicate the possible presence of semen or saliva, testing was not performed for these body fluids as that is outside of standard practice; Item 5 and Item 6, based on the case scenario, the bra band/tube top and the tampon were not applicable for possible saliva testing, as the chest and vagina were not locations that possible saliva was introduced.
PP7CBL	Item 1: Human saliva was detected. Human blood and semen were not detected. Item 2: Human blood was detected. Semen and Human saliva were not detected. Item 3: Human blood was detected. Semen and Human saliva were not detected. Item 4: Semen and Human saliva was detected. Human blood were not detected. Item 5: Human blood and semen was detected. Human saliva were not detected. Item 6: Semen was detected. Human blood and Human saliva were not detected.
QGEQ4T	RESULTS AND CONCLUSIONS: No blood was found on Item 1.1. Items 1.2, 1.3 and 1.5 contain blood. Presumptive* testing for human origin was positive. *Presumptive tests do not confirm the presence of bodily fluids or human origin. Items 1.4, 1.5 and 1.6 contain semen. No semen was found on Items 1.1, 1.2 or 1.3. No further analysis was performed. METHODS: Phenolphthalein Test – Items 1.1, 1.2, 1.3 and 1.5. HemaTrace® Test – Items 1.2, 1.3 and 1.5. Acid Phosphatase Test – Items 1.1, 1.2, 1.3, 1.4, 1.5 and 1.6 p30 Test – Items 1.1, 1.2, 1.3, 1.4, 1.5 and 1.6. Microscopy – Items 1.1, 1.2, 1.3, 1.4, 1.5 and 1.6. Results relate only to the items or sampling of items tested.
QKUQYU	Saliva testing is not performed in our lab; therefore, examination for saliva was not performed.
QLAK2U	No blood was found on Item 1. Items 2, 3 and 5 contain blood. Presumptive* testing for human origin was positive. No semen was found on Items 1 or 2. Items 4, 5 and 6 contain semen. No further analysis was performed. *presumptive tests do not confirm the presence of bodily fluids or human origin. METHODS: Phenolphthalein Test – Items 1, 2, 3 and 5 HemaTrace® Test – Items 2, 3 and 5. Acid Phosphatase Test – Items 1, 2, 4, 5 and 6 p30 Test – Items 1, 2, 4, 5 and 6. Microscopy - Items 1, 2, 4, 5 and 6. Results relate only to the items or sampling of items tested.
QWM64W	ALS inconclusive on item 5 - faint fluorescence observed.
R4TNA8	All examinations were conducted in accordance with my experience and qualifications as a crime scene examiner. Results obtained from most items were expected given the details that were provided in the scenario. Saliva tests were not performed as method verification for RSID saliva field kit are currently being carried out within my work place.
RBLTBV	This test is being completed as part of a male screening protocol and only tests that may be performed by the analyst in casework were performed on this proficiency test.
RTMQHX	"control swabs"- received, not analyzed.

TABLE 2

WebCode	Additional Comments
RWNZ6K	This test is being completed as part of a male screening protocol and only tests that may be performed by the analyst in casework were performed on this proficiency test.
TFHUWB	Too many items every year in this test. Overly complicated scenario with missing information. Too many scenario items- please limit to one set of swabs and one physical item.
THH6LJ	This test is being completed as a part of a male screening protocol and only tests that may be performed by the analyst in casework were performed on this proficiency test.
TKQCQW	Items 2, 3, and 5 contain blood. Presumptive* testing for human origin was positive. No blood was found on Items 1, 4, or 6. Items 4, 5, and 6 contain semen. No semen was found on Items 1, 2, or 3. No further analysis was conducted. Results relate only to the items or sampling of items tested. *presumptive tests do not confirm the presence of bodily fluids or human origin
TLYEQA	The laboratory does not perform semen testing, saliva testing, or DNA analysis. If testing is needed for any of the above, the item could be submitted to another laboratory.
TY7HZB	This test is being completed as part of a male screening protocol and only tests that may be performed by the analyst in casework were performed on this proficiency test.
UF93XK	Presumptive tests for the presence of blood, and semen were conducted with the following results: Item #1: negative for the presence of blood and semen. Item #2: positive for the presence of blood, negative for the presence of semen. Item #3: positive for the presence of blood, negative for the presence of semen. Item #4: positive for the presence of semen, negative for the presence of blood. Item #5: positive for the presence of blood, negative for the presence of semen. Item #6: positive for the presence of semen, negative for the presence of blood. Please submit items #1-6 to the [Laboratory] if you desire further DNA analysis.
UGXDND	The results of the amylase test for items 2, 3, 5, and 6 are inconclusive as a partial clearing was observed for each. Elevated levels of amylase are normally associated with the presence of saliva. However, amylase is present in other body fluids such as vaginal secretions, sweat, and fecal material.
UH277Z	My laboratory does not perform screening tests for saliva
V9QJFB	No testing was performed on control swabs
VKJRCT	Saliva testing: not necessary. This test is being completed as part of a male screening protocol and only tests that may be performed by the analyst in casework were performed on this proficiency test.
WVXB33	Visible stains were observed on Items #2, 3, and 5. LMG was positive on these items. Positive LMG is not a confirmation of the identity of blood. No fluorescence was observed with ALS on Items #1, 2, 3, and 5. However, the presence of red-brown stains on Items #2, 3, and 5 could mask the fluorescence in those stains. Fluorescence with ALS was observed in Items #4 and 6. A positive fluorescence is not a confirmation that semen is present.
WZB9GG	This test is being completed as part of a male screening protocol and only tests that may be performed by the analyst in casework were performed on this proficiency test.
XE3WZ9	We do not perform semen or saliva tests. LMG is our only testing method.
XMEKXM	No visible red-brown staining on Items 1, 4 and 6 - not tested with LMG. No fluorescent staining seen on Item 5, however presence of the red-brown questioned stain could mask fluorescence in that stain. Item 6 - fluorescence observed under ALS and reported as pos./ind., however no further tests for semen performed and this is not a confirmation that semen is present.
Y9D7F4	My laboratory does not perform screening tests for saliva.
YA8RBG	No hairs were observed on Item #5 tube top and Item #6 tampon
YBF4NC	This test is being completed as part of a male screening protocol and only tests that may be performed by the analyst in casework were performed on this proficiency test.
YD44NM	Lab does not test/screen for saliva amylase.

TABLE 2

WebCode	Additional Comments
YQRBGG	Body Fluids identification is based on LC HR/MS method using proteomic approach. The same MS/MS method is used to identify a specific peptides list for each bodyfluids. We have also looking for specific peptides of vaginal fluid but results are negatives on each item. We expected to found them on item 6 considering the case work. To conclude, results are : Item 1 : Saliva, Item 2 : Blood, Item 3: Blood, Item 4 : Semen and saliva (mix), Item 5: Semen and blood (mix), Item 6 : Semen.
ZUEEDZ	[Laboratory] does not presumptively test for saliva.

-End of Report-
(Appendix may follow)

(149)

Test No. 23-5781: Body Fluid Identification

DATA MUST BE SUBMITTED BY **May 22, 2023, 11:59 p.m. EDT** TO BE INCLUDED IN THE REPORT

Participant Code: U1234A

WebCode: L3E9F6

The Accreditation Release section can be accessed by using the "Continue to Final Submission" button above. This information can be entered at any time prior to submitting to CTS.

Scenario (this scenario applies only to Items 3-6):

Police are investigating a sexual assault involving a female victim and male suspect. The victim reports the male hit her several times, causing her to severely bleed. He then forced her to have vaginal sex. She reports the suspect licked her stomach and ejaculated on her stomach and chest. The victim reports that she inserted a tampon immediately after the assault. The suspect has been apprehended. The following items are being submitted to your laboratory for examination: One envelope containing four duplicate swabs (Item 3) taken from a stain on the suspect's hand. One envelope containing four duplicate swabs (Item 4) taken from the victim's stomach. One stained tube top (Item 5) worn by the victim during the assault. One tampon (Item 6) worn by the victim directly following the assault.

Items Submitted (Sample Pack BF11):

Item 1: Swab with white handle (provided as 4 replicates) - please examine this item for body fluids tested in normal laboratory casework; scenario not applicable

Item 2: Swab with pink handle (provided as 4 replicates) - please examine this item for body fluids tested in normal laboratory casework; scenario not applicable

Item 3: Swab with blue handle (provided as 4 replicates) - taken from stain on suspect's hand

Item 4: Swab with yellow handle (provided as 4 replicates) - taken from the victim's stomach

Item 5: One tube top worn by victim during the assault

Item 6: One tampon worn by victim after the assault

Control swabs (unused): One package to be used for substrate control testing, as necessary.

Instructions:

Please indicate the overall conclusion for each body fluid using the check boxes provided. List the tests performed and the results of each test. Item descriptions (as applicable) are provided at the top of each Item page.

For the purposes of this test, all swabs within an Item have been prepared as duplicates in order to provide sufficient sample quantity. A package of unused swabs is also provided for control purposes.

Please use the abbreviations listed in this response key to fill in the Screening Test tables. This is not an all inclusive list of tests, and should not be used to determine what tests should be performed. TESTS NOT ON THIS LIST MAY BE USED FOR SCREENING.

Test	Abbreviation	Test	Abbreviation
Acid Phosphatase	AP	Ortho-tolidine	O-tol
Alternate Light Source	ALS	Phenolphthalein-Tetramethyl benzidine	PTMB
Kastle Meyer	KM	Prostate Specific Antigen	PSA
Leucomalachite Green	LMG	Rapid Stain Identification	RSID
Microscopic	Micro	Tetramethyl benzidine	TMB

**Example
Blood Tests**

		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for blood?	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	HemaTrace	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	
	KM	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Pos./Ind. - Positive/Indicative; Neg. - Negative; Inc. - Inconclusive; NT - Not Tested

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Item 1: Cotton tipped swabs with white plastic handle.

Four replicate white handled swabs. Scenario is not applicable for this item. Utilizing your laboratory policies, please examine this item for body fluids commonly encountered in normal laboratory casework.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for blood?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for semen?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for saliva?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Pos./Ind.	Neg.	Inc.
Please report Other Tests as 'Body Fluid - Test Performed'				
<input type="text"/>		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

Item 2: Cotton tipped swabs with pink plastic handle.

Four replicate white handled swabs. Scenario is not applicable for this item. Utilizing your laboratory policies, please examine this item for body fluids commonly encountered in normal laboratory casework.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for blood?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for semen?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for saliva?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Pos./Ind.	Neg.	Inc.
	Please report Other Tests as 'Body Fluid - Test Performed'			
	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

Item 3: Cotton tipped swabs with blue plastic handle.

Four replicate blue handled swabs taken from the stain on the suspect's hand.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for blood?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for semen?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for saliva?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Pos./Ind.	Neg.	Inc.
Please report Other Tests as 'Body Fluid - Test Performed'				
<input type="text"/>		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

Item 4: Cotton tipped swabs with yellow plastic handle.

Four replicate yellow handled swabs taken from the victim's stomach.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for blood?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for semen?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
What is the overall conclusion for saliva?		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Pos./Ind.	Neg.	Inc.
Please report Other Tests as 'Body Fluid - Test Performed'				
<input type="text"/>		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

Item 5: One tube top worn by victim during assault.

This item of evidence is a tube top worn by the victim during the assault.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for blood?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for semen?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for saliva?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Pos./Ind.	Neg.	Inc.
	Please report Other Tests as 'Body Fluid - Test Performed'			
	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

Item 6: One tampon worn by victim after assault.

This item of evidence is a tampon worn by the victim after the assault.

Blood Tests		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for blood?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Semen Tests		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for semen?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

Saliva Tests		Pos./Ind.	Neg.	Inc.	NT*
	What is the overall conclusion for saliva?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Test(s) Performed	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	

*Should a body fluid conclusion be marked as NT (Not Tested), please document any Item specific reasons under the "Test(s) Performed" section. All other reasons should be documented in the Additional Comments section of this data sheet.

Other Tests		Pos./Ind.	Neg.	Inc.
	Please report Other Tests as 'Body Fluid - Test Performed'			
	<input type="text"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

Additional Comments

Please note: Any additional formatting applied in the free form space below will not transfer to the Summary Report and may cause your information to be illegible. This includes additional spacing and returns that present your responses in lists and tabular formats.

RELEASE OF DATA TO ACCREDITATION BODIES

The Accreditation Release is accessed by pressing the "Continue to Final Submission" button online and can be completed at any time prior to submission to CTS.

CTS submits external proficiency test data directly to ASCLD/LAB, ANAB, and/or A2LA. Please select one of the following statements to ensure your data is handled appropriately.

- This participant's data is intended for submission to ASCLD/LAB, ANAB, and/or A2LA. (Accreditation Release section below must be completed.)
- This participant's data is not intended for submission to ASCLD/LAB, ANAB, and/or A2LA.

Have the laboratory's designated individual complete the following steps **only if your laboratory is accredited in this testing/calibration discipline** by one or more of the following Accreditation Bodies.

Step 1: Provide the applicable Accreditation Certificate Number(s) for your laboratory.

ANAB Certificate No.
(Include ASCLD/LAB Certificate here)

A2LA Certificate No.

Step 2: Complete the Laboratory Identifying Information in its entirety.

Authorized Contact Person and Title

Laboratory Name

Location (City/State)